CLASS «B» STERILIZATION AUTOCLAVE
AUTOCLAVE ESTERILIZACIÓN CLASE «B»

AUTESTER ST DRY PV B 12L

4001737

Indice

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Index

1 Packing list	16
1.1 Standard package	16
2 Safety Information	17
2.1 Safety icons	17
2.2 Risks that the user can be subjected to	17
2.3 Qualified users	17
2.4 Modifications	17
3 Description	18
4 Installation	19
4.1 Suitable location	19
4.2 Connection to power supply	19
4.3 Filling water reservoir	19
4.4 Draining water	19
5 Technical specifications	20
6 Description and usage	20
6.1 Location of elements	20
6.2 Control panel	22
7 Operation	23
7.1 Start a program	23
7.2 Select a program	23
7.3 Step selection	23
7.4 Tests	23
8 Maintenance	24
8.1 List of failures	24
8.2 Door adjusting	25
8.3 Cleaning drainage filter	25
8.4 Replacing/cleaning the sealing ring	26
8.5 Replacing the power fuses	26
8.6 Cleaning cycle	26
8.7 Printer user guide	27
8.8 Function settings	27
8.9 Cleaning the water tank	27
<i>NOTICE TO CUSTOMERS</i>	28

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
1 Packing list

1.1 Standard package

Accessories	Qty	Fig.
• Tray holder + 3 trays	1	(1)
• Power cord	1	(2)
• Main board fuse	1	(3)
• Drain water pipe	1	(4)
• Try holder	1	(5)
• Fuse	2	(6)
• Door adjusting spanner	1	(7)
• Internal Filter	1	(8)

2 Safety Information

The equipment comes complete with safety features.

This manual indicates the areas of possible risk.

2.1 Safety icons

They identify risky situations and the safety measures to be taken.

The icons affect the paragraphs marked by a grey line.

Risk danger

There's a risk of danger.

Be careful, follow the instructions as indicated.

Electrical risk

There's a risk of electric shock when accessing zones or doing operations shown with this icon.

Be careful, follow the instructions as indicated.

Risk of burns in high temperature areas

The temperature in the zone indicated with this icon can exceed 60°C. Use thermal protective gloves to complete the described functions.

Be careful, follow the instructions as indicated.

Important information

- Important information in obtaining the best results or optimum performance of the equipment.
- Important information to extend the life of the equipment and to maintain its optimum performance.

2.2 Risks that the user can be subjected to

- Possibility of touching surfaces that are over 60°C.
- Possible escape of steam vapours.
- Electrical risk.

2.3 Qualified users

This equipment should only be used by personnel that are qualified users.

This equipment should only be used by the staff that has understood this manual or has previously read the instructions of the steam sterilizers use.

2.4 Modifications

The modification of the safety devices operation does not agree with the manufacturer, and it could put the user under unexpected risks.

DO NOT OPEN THE LID OF THE AUTOCLAVE IF THE CHAMBER IS UNDER PRESSURE.

3 Description

J.P. SELECTA, s.a.u. designs and makes pressurized steam autoclaves especially to meet the needs for instruments sterilization and devices in clinics, hospitals and laboratories.

Users will include doctors or his/her assistances under doctors supervisions.

The autoclave is computerized, user friendly and easy to operate.

Its real time display shows all data during the sterilization cycle.

The software program provides autodiagnose, protection against over heating or over pressurizing, and thus ensures reliability of the sterilization results.

Build-in condenser water collecting system prevents discharge of waste steam and keeps the machine clean and safe.

AUTESTER uses the advanced technology of 3 times pre-vacuum and vacuum drying.

It enables the hot steam to penetrate in thin tubes and porous material.

This product is suitable for dentists, ophthalmology, surgical theatres and sterilization laboratories.





Sterilization programs:

Step	134°C	134°C	121°C	121°C	134°C
	Wrapped	Unwrapped	Wrapped	Unwrapped	Universal B+
Sterilization temperature	134	134	121	121	134
Pressure	2.10	2.10	1.10	1.10	2.10
Preheating time	0-7	0-7	0-7	0-7	0-7
Vacuum time (estimated)	8 min	4 min	6 min	6 min	8 min
Sterilization time	6 min	4 min	20 min	20 min	18 min
Drying time	15 min	10 min	15 min	10 min	15 min
Total time	40 min	25 min	45 min	40 min	50 min

Qualification and cleaning programs:

Step	Bowie	Vacuum	Cleaning
	Dick	Test	
Sterilization temperature	134	-	-
Pressure	2.10	-	-
Preheating time	0-7	-	-
Vacuum time	10 min	4 min	20 min
Sterilization time	3.5 min	10 min	10 min
Drying time	1 min	-	-
Total time	20 min	17 min	30 min

4 Installation

Once the unit is unpacked, visually check the packing list (See 1.1) and that there are no parts damaged during transportation.

Save the packing for some days for a later sending.

4.1 Suitable location

Place the unit over a flat level and stable surface, suitable for the unit weight.

Use the leveling instrument to check the flat placement.

Note: The front of the machine can not be lower than the rear part. In this case the unit won't operate properly and could display E3 message (Dry not completed).

Place the equipment not far than 1.5m from a power supply plug.

No tap water and waste will be needed.

Autoclave dimensions are: heigh- 460mm, width- 460mm & depth- 600mm.

To work comfortably, let a free space of 50cm on each side and not least than 20cm on the rear part.

Prepare to refill the water reservoir from the top of the unit.

4.2 Connection to power supply

Locate a power supply socket (rated to 16A) near the unit.

Use the power cord supplied or another of similar specifications.

For safety reason, the socket must have ground terminal, since the unit must work grounded.

Before plugging on the socket, check that its voltage matches the voltage marked on the unit plate.

4.3 Filling water reservoir

Use distilled or deionized water with this autoclave.

Open the top cover to fill the reservoir tank.

If there is not enough water, the "water" indication led will blink and an acustic alarm will beep.

Pay attention with the maximum filling level, in the red mark. Do not exceed it.

4.4 Draining water

To empty the tank of used water, connect the output shown in the picture to the tube and open the valve by turning it counterclockwise and clockwise.

This tank should be frequently emptied. We recommend emptying it whenever you fill the tank with clean water.

To empty the clean water tank, connect the outlet indicated in the picture to the supplied pipe and open the valve by turning it counterclockwise.

This tank does not normally empty, only in cases of tank cleaning or equipment transportation.

5 Technical specifications

Chamber inner dimensions (diameter x depth):	200X358mm
Power supply:	according to nameplate
Sterilizing temperature/pressure:	121°C/(1.0-1.3) bar
Sterilizing temperature/pressure:	134°C/(2.0-2.4) bar
Power fuse:	F12A 250V
Water reservoir capacity:	3.2 litres
Operation temperature range:	0-40°C
Water minimum level:	1L
Environmental conditions:	+5°C to +40°C
Relative humidity:	35% to 75%

6 Description and usage

The autoclave ST Dry PV 12 litres is an equipment useful for sterilizing processes in clinics departments and hospitals, of instruments and devices wrapped and unwrapped, prions and drying. The software programme provides auto diagnose, protections against over heating or over pressurizing, and ensure reliability of the sterilization results. They have a condensed water collection box which keeps the equipment safe.

This autoclave uses the advanced technology of 3 times fractionated and pulsatory drying vacuum. It enables the hot steam to penetrate into the porous material.

6.1 Location of elements


- 1 Door knob
- 2 Control panel

- 3 Air inlet filter
- 4 Sealing ring
- 5 Waste water tank used
- 6 Water reservoir outlet
- 7 Main switch
- 8 Drainage filter


- 9 Reservoir tank water filling








- 10 Supply socket
- 11 Fuses
- 12 RS232 socket (printer)
- 13 Outlet of excess waste water
- 14 Safety valve

6.2 Control panel

The control panel allows function selection and shows the autoclave alarm status. Once the unit is plugged in, put it on, by switching on the main switch. The unit goes to stand by status (ON).

- 15 Pressure display
It displays chamber pressure in real-time during the cycle.
- 16 Temperature
It displays chamber temperature during a cycle. The units is per Celsius degrees (°C).
- 17 State displays/Error codes/Time
It displays all the state of autoclave and the error codes when alarm goes off. You may correct the errors accordingly the error list.
- 18 SELECT key
Programme selection key. It is used for selecting a temperature (121 or 134) or the textile mode (B+).
It can also be used to make a washing cycle.
- 19 PACKAGE key
It is used to select wrapped or unwrapped instrumental.
- 20 TEST key
You can select the Bowie Dick test or the vacuum test.
- 21 STEP key
This is a post cycle selection key. By pressing this key, the autoclave keeps the chamber warm, and thus it saves time in the following cycle.
If you want to cancel this operation, press this key again. The indicator will turn off.
- 22 Quick DRY key
This key short drying time (1 or 5 min) is selected. It is only recommended in processes with not wrapped material.
- 23 START/STOP key
This key is used to start programmes and abort them manually. If you press it for more that 5" during the operation, it aborts the cycle and starts a 2-minute drying cycle.

7 Operation

7.1 Start a program

- 1 Plug the autoclave into the mains.
- 2 Switch on the main switch (7).
The control panel does an automatic test. It indicates an error in case of fail. It indicates "Ld" if it passes the test.
- 3 Put the material to be sterilized into the chamber and close it by using the knob (1).
- 4 Start the autoclave by pressing "Start/Stop" key (23). "HE" message appears.
- 5 The programme stops automatically when it ends. You can abort it if necessary by pressing again the "Start/Stop" key (23). The message "Ed" is displayed. The "running" indicator goes off and the "ready" indicator lights up.

7.2 Select a program

- 1 Select a temperature. Use the "Select" key (18) to select 121°C, 134°C or textile (B+). The selection will be indicated by the corresponding led.
- 2 To select wrapped or unwrapped instruments, use the "Package" key (19). The selection will be indicated by the corresponding led.
- 3 Select "Quick" drying (22) if necessary. Only in unwrapped processes. The selection will be indicated by the corresponding led.

7.3 Step selection

This option allows a pre-warming of the chamber in order to shorten the time of the cycles. The selection will be indicated by the corresponding led.

7.4 Tests

This autoclave has two test programs:

- 1 Bowie Dick: Useful to check the sterilization. Select it by the "Test" key (20). The led "B&D" must light.
- 2 Vacuum test: Useful to found leakages. Select it by the "Test" key (20). The led "Leak" must light. Vacuum test is recommended to be done everyday before starting the batch sterilization.
Ensure that the vacuum gauge in the upper left corner of the screen reaches at least a value of -0.7 BAR.

Caution: When unloading material from autoclave. Hazard of burns if in contact with hot surfaces.

The temperature in this area can exceed 60°C.


Use the tray holder to load or unload the tray with instruments in it in order to avoid burning injuries.

8 Maintenance

8.1 List of failures

Error Code	Malfunction	Cause	Troubleshooting
E00	Program Exit	Manually exit during the operation	Restart the Program
E10	Pressure Sensor Malfunction	Pressure over detected, $\geq 300\text{KPA}$	Check the pressure sensor or main board
E11	Internal Temp. Sensor Malfunction	Internal Chamber temp. over detected, $\geq 200^\circ\text{C}$	Check the temperature sensor of the sterilization room or main board
E12	Steam Generator Sensor Malfunction	Steam generator temp. over detected, $\geq 250^\circ\text{C}$	Check the temperature sensor of the steam generator or main board
E13	External Temp. Sensor Malfunction	External chamber temp. over detected, $\geq 200^\circ\text{C}$	Check the temperature sensor of the chamber or main board
E41	St.Failure	Actual temp.exceeds sterilization temp. $+4^\circ\text{C}$	Check the temperature sensor or main board, Whether or not the pipe is blocked
E42	St.Failure	Actual temp.is lower than sterilization temp.	Check the temperature sensor or main board, door sealing ring, solenoid valve and related pipes to see if they are leaking
E02	Pressure Over	Actual pressure.exceeds sterilization temp.	Check the pressure sensor, main board, exhaust valve, drain valve and related pipes to see if they are blocked
E03	Chamber Temp. Over	Chamber temperature exceeding set value	Check the temperature sensor of the chamber or main board
E05	Steam Exhaust Failure	Pressure discharge rate in chamber is lower than preset value	Check whether exhaust valve, drain valve and related pipes are blocked
E06	Door Malufunction	Door switch loosened	Check whether the door switch, door switch wire and door body are adjusted in right place
E07	Steam generator Temp. Over	The temperature of the steam generator is over 210 degrees	Check the temperature sensor or main board, inlet valve, water pump
E08	Heating Time Over	Sterilization preheating time exceeds preset value	Check the temperature sensor, main board, temperature protector, connecting wires to see if they are falling off
E91	The First Time Vacuum Failure	During the first vacuuming, the vacuum degree in the chamber cannot reach the preset value within the prescribed time	Check the vacuum pump, main board, two positions three ways valve, drain valve, exhaust valve and the filter in the chamber to see if it's blocked
E92	The Second Time Vacuum Failure	During the second vacuuming, the vacuum degree in the chamber cannot reach the preset value within the prescribed time	Check the vacuum pump, main board, two positions three ways valve, drain valve, exhaust valve and the filter in the chamber to see if it's blocked
E93	The Third Time Vacuum Failure	During the third vacuuming, the vacuum degree in the chamber cannot reach the preset value within the prescribed time	Check the vacuum pump, main board, two positions three ways valve, drain valve, exhaust valve and the filter in the chamber to see if it's blocked
E94	LEAK Program Vacuum Failure	The vacuum degree in the chamber cannot reach the preset value within the prescribed time	Check the vacuum pump, main board, two positions three ways valve, drain valve, exhaust valve and the filter in the chamber to see if it's blocked
E95	The First Time Pressurization Over	For the first time pressure boost, the pressure is less than 100KPA in 30 minutes	Check water pump, inlet valve, steam generator, intake pipe or water tank to see if it's short of water
E96	The Second Time Pressurization Over	For the second time pressure boost, the pressure is less than 100KPA in 30 minutes	Check water pump, inlet valve, steam generator, intake pipe or water tank to see if it's short of water
E97	The Third Time Pressurization Over	For the third time pressure boost, the pressure is less than 100KPA in 30 minutes	Check water pump, inlet valve, steam generator, intake pipe or water tank to see if it's short of water
A01	Water Quality Alert	Water quality in water tank is unqualified	Please use distilled water and still check the main board
A02	Water Refill Alert	Water shortage in water tank	Please add water, still check water level switch, main board
A03	Wasted Water Tank Full Alert	Water full in used water tank	Discharge used water, still check water level switch, main board

INCIDENCE: When the two TEST key leds are blinking means that a washing cycle is necessary (every 100 cycles). With the door open, press SELECT key (18) for 5 seconds. Leds B and D (20) will flash. Close the door and press START.

Clean the water reservoir every week by using medical disinfectant. First drain the distilled water in the reservoir, pumping the medical disinfectant and then wait for an hour. Flush the reservoir with distilled water after medical disinfectant is drained.

Clean and disinfect the chamber every month by brushing the chamber inner surface with medical alcohol.

8.2 Door adjusting

If you note that there are steam leakages through the door, you can adjust it by doing the following:

- 1 Stop the autoclave and let it cold.
- 2 Insert the spanner into the gap beneath the plastic door cover.
- 3 Turn the door adjusting nut counter clockwise using the supplied spanner as showed in the pictures.

NEVER ADJUST THE DOOR WHILE THE AUTOCLAVE IS UNDER PRESSURE OR IF IT IS LOOKED BECAUSE YOU MAY SUFFER SCALDS AND BURNS.

8.3 Cleaning drainage filter

After a long time usage, some small impurity or particles could deposit on the draining filter (26) and may affect vacuum and drying processes. The impurity could come from grease, dust or dirty that the instruments to be sterilized or the water may contain.

Clean the filter with a soft brush every month. To do it, please do as follows:

1. Locate the water filter
2. Pull off the rubber filter from the connection
3. After you clean or change a new filter shim, plug in the water filter back on the back chamber connection.

1. Chamber water filter

2. Back-out the water filter

3. Take out the filter net, clean it by soft brush or replace a new one.

8.4 Replacing/cleaning the sealing ring

Extract the old ring. You will need a flat screwdriver to begin the extraction. Then, insert the new ring easily with the hands.


1. Hold the edge of the seal with one hand and insert the screwdriver beneath the seal with the other hand. Pull out the seal slowly.
2. Once a section of the seal is out, you may pull the whole seal out with both hands. Then, clean the seal groove and the sea. If it is broken or damaged, replace it with a new one.
3. After the groove and the seal are cleaned, mount the seal back in place.
Note: the seal should be mounted evenly inside the groove with first 4 positioning points in place. You may use your hands to press the seal in place evenly.
4. When inserting the seal, the inner edge might evaginate. You may use a screwdriver to press the edge down in the groove.

To avoid steam leakage, we recommend cleaning the support edge of the recommended gasket with a soft cloth dampened with distilled water.

8.5 Replacing the power fuses

To replace the power fuses, do the following:

1. Switch the power off. Unscrew the fuse holder with a screwdriver counter clockwise.
2. Pull the fuse holder out with your fingers.
3. Make sure to replace the fuse by a suitable one.
4. Put back the fuse holder and lock it by turning it clockwise with a screwdriver.

8.6 Cleaning cycle

Every 100 cycles the equipment will indicate us that it needs a washing cycle. So, the two TEST key leds will blink.

1. When the tank is empty, with no material inside, pulse SELECT for more than 5 seconds until "CLE" appears on the screen.
2. Press START and wait for 30 minutes, the equipment will make the cleaning cycle automatically. Once finished, the equipment will work normally again.

POWER THE MACHINE OFF AND LET IT COOL DOWN BEFORE CHANGING THE SEAL AND AVOID ANY BURNING INJURIES

Sterilization	Records
Date: Y11M05D11	
Cycle: 00001	
C.N.: 134 2.10	
-B-SP=6mDP=10m	
C.S.T.: 13:05:30	
Vp1: 13:11:27	0.00 048.9
Pp1: 13:12:57	-0.81 050.2
Vp2: 13:16:33	0.29 107.5
Pp2: 13:18:01	-0.81 072.1
Vp3: 13:20:36	0.30 112.4
Pp3: 13:22:14	-0.81 081.4
H.S.T.: 13:27:55	
Max: 2.14 134.6	
Min: 2.02 134.5	
Ave: 2.08 134.5	
H.E.T.: 13:34:06	
D.S.T.: 13:34:43	0.30 115.5
C.E.T.: 13:45:32	
Result: Passed	
Operator:	
Thanks	

8.7 Printer user guide

According to the requirement of customer, steam sterilizer can adopt prescriptive printer. Connect the printer and the steam sterilizer with prescriptive cord. Switch on the power, two lights on the printer panel will go on, and this means that the printer has been installed.

Caution: steam sterilizer can only select micro. Printer offered by us.

Print information explaining:

- Date: Sterilize date.
- CN: Name of the cycle.
- C.S.T.: The starting time of the cycle.
- Vp1.Vp2.Vp3: Vacuum pressure1, Vacuum pressure2, Vacuum pressure3.
- Pp1.Pp2.Pp3: The phase of pressure1, the phase of pressure2, the phase of pressure3.

If you have activated the printing option at the end of the cycle, when the cycle finishes, the display will show "Ed" and at the same time "Prt" will appear, which means that the information is sent to the printer.

8.8 Function settings

a) Press the "STEP" key for 20 seconds. Without any active program (message "Ld").

b) By means of the "SELECT" key you can select the alternative options:

Option 1	YEAR	"YER"
Option 2	MONTH	"MtM"
Option 3	DAY	"DAY"
Option 4	HOUR	"HOR"
Option 5	MINUTES	"MIN"
Option 6	printing	"Prt" (OFF, ON)
Option 7	EXIT	"Eit"

c) Press the "STEP" key (arrow up) or "DRY" key (arrow down) to select the desired value.

d) Confirm the value by pressing the "START/STOP" key. The equipment will emit a confirmation beep.

8.9 Cleaning the water tank

1) Empty the tanks water by draining it. Connect the water outlet using the one end of the tube without connector, and the other end of the tube to a drainage sink. Then turn the draining knob anti-clock wise to drain the water from the machine.

2) Open the tank cover board, you can use a screwdriver to loose the five screws.


3) Pull up the water tank cover by one hand.

4) Pull up the tank cover completely and then disassembly the cover.


5) Then you can start cleaning the tank. Use a cotton with some alcohol or disinfectant to clean the tank. Wash it with distilled water, then empty the tank and dry it.

6) After cleaning the tank, assembly the cover and tighten the screws.

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**Notice to customers:**

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.