AUTOCLAVE PARA ESTERILIZACIÓN A VAPOR
AUTOCLAVE FOR STEAM STERILIZATION

MICRO 8
MED 12
MED 20

4001745
4001756
4001757

Indice

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Index

1 Safety	15
1.1 Safety icons	15
1.2 Risks to which the operator is subjected to.....	15
1.3 Qualified users	15
2 General Information for Users	16
2.1 Legal requirements for installation, use and maintenance of this autoclave	16
2.2 Autoclave reception	16
2.3 Documentation.....	16
2.4 Other interesting information.....	16
2.5 Warranty period.....	16
3 Technical Specifications	17
4 Contents List	18
4.1 Standard contents	18
5 Introduction	18
5.1 Sterilization	18
5.2 Disinfection	18
5.3 Culture media.....	18
5.4 Purge (extracting the vessel air)	18
5.5 Saturated vapour.....	18
5.6 Atmospheric purge	18
5.7 Recommendations for a perfect sterilization.....	18
5.7.1 Loose instruments sterilization:	19
5.7.2 Containers sterilization:	19
5.7.3 Liquids sterilization:.....	19
6 Equipment description	19
6.1 Design features decisive for the equipment life	19
6.2 Control panel.....	20
6.3 Lateral part	20
6.4 General view	20
7 Installation	21
7.1 Location.....	21
8 Operation.....	21
8.1 Preliminary operations for switching on	21
8.2 Opening and closing the lid	21
8.3 Starting up	22
8.4 Sterilization end.....	22
9 Alarm indicators	23
10 Spare Parts	24
Notice to customers	24

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ENGLISH VERSION

1 Safety

The equipment includes the appropriate safety features.

Risky situations that must be respected are indicated in this manual.

1.1 Safety icons

They identify risky situations and safety measures that should be taken.

The icons make reference to the paragraph marked with a grey line.

Danger risk

It indicates danger risk.

Please, respect the indicated instructions to carry out the described operation.

Electrical risk

It indicates risk of electric shock when accessing zones shown by this icon or when taking actions indicated in this manual with this icon.

Please, respect the indicated instructions to carry out the described operation.

Burn risk due to high temperature areas contact

The temperature in the area indicated by this icon could exceed 60°C. Please use the thermal protective gloves to carry out the described operation.

Please, respect the indicated instructions to carry out the described operation.

Important information

- Important information to obtain good results or an optimum equipment performance.
- Important information to extend the equipment life or to avoid degradation of some of its components.

1.2 Risks to which the operator is subjected to

- Possibility of touching hot parts over 60°C.
- Possibility of steam exposure.
- Electrical risk.

1.3 Qualified users

This equipment can only be used by qualified staff.

This equipment can only be used by users who have read and understood these instructions or those who have been qualified for the use of this equipment.

2 General Information for Users

2.1 Legal requirements for installation, use and maintenance of this autoclave

Valid reference regulation over pressure equipments:

- Directive 2014/68/UE

Installation requirements are established in section 7 of this manual.

2.2 Autoclave reception

Handle the equipment with care in a suitable way according to its weight and considering it as a FRAGILE equipment.

Unpack and check that the contents coincide with the «packing list». If you see any damages or something missing, please notify your distributor as soon as possible.

Keep the original package for some days. Afterwards, throw away in the right way, by separating its components into cardboard, wood and plastic.

2.3 Documentation

The following documents are supplied with the autoclave:

- Instructions manual:

Do not install nor use this equipment before reading carefully this user's manual.

These instructions are an important part of this equipment and must be available to all users.

Any doubt should be asked to the autoclave distributor.

- CE Declaration of conformity:

This is a necessary document for authenticating the autoclave installation.

2.4 Other interesting information

Autoclaves sent to J.P. SELECTA, s.a.u. for repairing or maintenance should be disinfected.

Any modification, elimination or lack of maintenance of any autoclave device is forbidden.

It is forbidden also the use of this equipment in explosive environments or with steamed substances or to make explosive or inflammable mixtures.

If this autoclave is used in a way not specified by J.P. SELECTA, s.a.u. the guaranteed protection and its functionality can be committed.

2.5 Warranty period

Autoclaves are guaranteed for one year. The warranty does not cover damages caused by incorrect use or causes beyond the control of J.P. SELECTA, s.a.u.

Any manipulation of the machine by not authorized personnel of J.P. SELECTA, s.a.u., automatically cancels the guarantee.

3 Technical Specifications

Model	4001745	4001756	4001757
Autoclave volume	8 litres	12 litres	20 litres
Mains voltage:	I 230V	I 230V	I 230V
Electrical power: (W)	1000	1800	2500
Maximum intensity: (A)	4.5	7.8	10.8
Fuses: (fast. H type, dimensions: 10x38mm)	5	10	15
Network connection cable:	Phase, Neutral, Earth (section 0.75mm ²)		
Weight (Net Kg):	13	18	41
Useful measures (cm Ø x depth)	Ø24x15	Ø24x26	Ø24x41
Ext. measures: (cm Depth x Width x Height)	40x38x38	46x38x38	70x38x38
Volume: (Litre)	8	12	20
Water quality	Demineralized between 15 and 200 microS/cm; pH between 5 and 7		
Operation system:	Continuous with 20 minutes intervals between cycles.		
Environmental conditions:	Inner use. Altitude up to 2000m Room temperature between 5°C and 40°C. Max. relative humidity 80% for temperatures untill 31°C, reducing in a linear way up to a 50% of 40°C relative humidity.		
Overpressure level:	Category II		
Pollution degree:	2		
Vapour contact material:	Stainless steel, copper, teflon and brass		
Frame material:	Steel with cover		
Vessel material:	Stainless steel		
Pipes material:	Copper, brass and teflon		

4 Contents List

4.1 Standard contents

The standard contents consist of the following items:

Litres	8L	12L	20L
Model	4001745	4001756	4001757
Heater element cover	59106	59106	59106
Lid gasket	21230	21230	21230
Inox basket	1004771	1004772	1004773 (x2)
Instructions manual	80244	80244	80244

5 Introduction

Micro 8, Med 12 and Med 20 autoclaves are equipment suitable for a wide range of applications in several fields: health-care, industrial processes and quality control. These autoclaves allow solids without packaging sterilization processes and also of liquids and culture media.

5.1 Sterilization

Sterilization is the destruction or elimination of all forms of microbial life, including spores, which are present in inanimate objects.

5.2 Disinfection

Disinfection is the process of destruction of infectious agents. It is normally made at a lower temperature than sterilization. Disinfection only achieves the elimination of some vegetable forms.

5.3 Culture media

Culture media is a nutritional material in which microorganisms can be recovered, multiplied and isolated, as well as used in susceptibility tests. They are generally offered as fine or granular dry powder but also can be offered hydrated and prepared. They must not be used before being sterilized.

5.4 Purge (extracting the vessel air)

Purge is the cycle used to eliminate the air inside the vessel to achieve saturated steam vapour.

5.5 Saturated vapour

It is water vapour at a certain temperature corresponding to the boiling point of the original liquid.

5.6 Atmospheric purge

In the atmospheric purge, the air is expelled by vapour in a gravimetric way for a determined period of time.

5.7 Recommendations for a perfect sterilization

The material to be sterilized must be perfectly clean, free from residues. It is recommended to clean it with a good detergent and distilled water. Rinse later with plenty of water.

Do not excessively charge the trays, the racks or the baskets. Try always to leave a free space of 1 or 2 cm between them, in order to let the steam pass through and also to make drying more easily.

5.7.1 Loose instruments sterilization

Place the instruments over the water resisting paper in the open position and not in contact with each other.

Please, avoid setting in the same tray different kind of metal instruments.

5.7.2 Containers sterilization

Never set containers hermetically sealed.

Place the containers downwards to avoid water deposits.

5.7.3 Liquids sterilization

Please, set the liquids for sterilization in appropriate containers prepared to support sterilization temperatures. Place them over the trays, so it could be easier to pick possible spillage.

The containers should be filled about 2/3 of their capacity. DO NOT CLOSE THEM hermetically, just cover them with a cotton or some kind of stopper to let the air exit and not to form any pressure.

Avoid using narrow bottleneck containers.

At the end of the sterilization cycle, let the autoclave freely cool down (till pressure will be 0 kg/cm²).

DO NOT STERILIZE TUBES NOR WRAPPED MATERIAL

This kind of material is reserved to the autoclaves with a prevacuum purge system.

6 Equipment Description

Micro 8, Med 12 and Med 20 autoclaves are saturated vapour equipments which are distinguished thanks to the following features:

- They are not suitable for fixing sterilization units.
- Suitable for sterilizing:
 - Unpacked metal load: Clamps, Scalpels, ...
 - Liquid containers load.
- Water manual wastepipe.
- Automatic purge (MED12 and MED20) and manual purge (MICRO8).

6.1 Design features decisive for the equipment life


The maximum service temperature for these equipments is much lower than the one for which used materials show damage for slow fluids.

- Fatigue: These equipments have been designed to safety resist the maximum admissible pressure $PS = 2,4$ bar, as well as the fatigue in materials due to pressure cycles (from vacuum to maximum pressure).


The theoretical expected life for these equipments is about 20 years, with a work system of 4 operations / day during 365 days a year. However, the equipments real life is subjected to the periodical tests every 10 years and the annual checks.

- Corrosion: The use of genuine stainless steel in the autoclave vessel guarantees the resistance to the corrosion of water, steam and products to be sterilized.

6.2 Control panel

1. Pressure gauge
2. Overtemperature indicator
3. Open door indicator
4. Temperature selector 121°C/134°C
5. Timer knob
6. Lid release button
7. Operation at 121°C indicator
8. Operation at 134°C indicator
9. Sterilization phase indicator
10. End of cycle indicator

6.3 Lateral part

11. Safety valve
12. Purge, vacuum and closed manual valve
13. Purged steam and drainage outlet
17. Network connection cable
18. Fuses
19. Safety thermostat (reset button)

Important: The safety thermostat comes adjusted from the factory. It is advised not to handle it except to reset it by pressing its button.

6.4 General view

14. Control panel
15. Open /Close door handle
16. Main power switch

7 Installation

7.1 Location

Locate the autoclave near an appropriate power source for the machine.

According to the current legislation, the autoclave should be located in a way that the safety valve discharge outlet DO NOT point to any people or could not reach them in case of steam outlet.

Install the autoclave over a flat, stable horizontal surface, suitable for the machine weight (see chapter 3 «Technical Specifications»), leaving a free space of at least 20 cm around the equipment.

Fix a pressure hose to the outlet (13).

WARNING

If you use the machine for quick depressurization, please fix a pressure hose.

Fix the hose in such a way that cannot move by the effect of the pressure at the end of the cycle.

Do not hold the hose with the hands. It could reach temperatures higher than 60°C.

8 Operation

8.1 Preliminary operations for switching on

- Check that the lateral valve (12) is in the «CLOSED» position.
On model MICRO-8, in the «STEAM» position.
- Place the heater element cover tray at the bottom of the vessel.
- Fill the sterilization chamber with water until the heater element cover rack.
- Plug the autoclave to the mains.
- Press the main switch on (16).
- Place the material to be sterilized.
- Close the lid.

8.2 Opening and closing the lid


CLOSING


- To close, just press over the lid with one hand and turn the knob (15) with the other hand till the «CLOSED» position.

OPENING

- It is not possible to open the lid if the pressure is different to 0 bar.
- To open, press the button (6) and turn the knob (15) with the other hand till the «OPEN» position.

8.3 Starting up

Once the operations in chapter 8.1 are made, please do as follows:

- Make sure that the lateral valve (12) is in the «CLOSED» position. On MICRO-8 model, in the «STEAM» position.
- Select the sterilization temperature (121°C/134°C) by means of the selector (4).
- Select the sterilization time with the timer knob (5).

Recommended sterilization temperatures and times:

Material	Temp.	Time
Liquids	121°C	30'
Metallic tools and glass	134°C	10'

IMPORTANT

Is advised to use decalcified water. Distilled water and deionized water produce a chemical reaction which shorten life of the heating elements.


iWARNING!

The timer has a double scale:

- **Blue for equipments connected to a 60 Hz frequency.**
- **Red for equipments connected to a 50 Hz frequency.**

- The autoclave will begin the purge process by means of which the air is expelled through the exit (13). At the end of the purge process, the autoclave closes automatically the steam exit. On MICRO-8 model, monitor the purge outlet. Approximately 2 minutes after the steam has come out, change the valve position to «CLOSED».

- Once the valve is closed (12), the temperature and pressure increase till they arrive to 121°C/1 bar or 134°C/2 bar, according to what it has been selected.
- When the autoclave reaches the operating temperature, the green indicator lights on (9) and time begins to deduct from the sterilization time.
- When the sterilization time finishes, the indicator lights (10) and the autoclave finishes the sterilization cycle.

8.4 Sterilization end

Once the sterilization time ends, one of these two situations can be carried out:

- Solids sterilization: a quick devaporization can be made by setting the valve (12) to the «STEAM» position.

!WARNING!

The steam outlet (13) expels high temperature steam and could produce serious burns if you don't take precautions.

- Liquids sterilization: do not open the valve when it cools down till the pressure gauge (1) indicates 0 bar. Then, set the valve (12) to the «Steam» position.

In both cases, please open the autoclave lid when the pressure gauge indicates 0 bar, as explained in chapter 8.2.


When the autoclave is cooling with the valve closed, some vacuum can be produced inside the chamber. That can avoid the lid from opening.

In case of not being able to open the lid, break the vacuum by turning the valve (12) onto «steam» position.

Do not forget to turn the valve back on the «closed» position before the next sterilization.

9 Alarm indicators

Overtemperature indicator (2).

The safety thermostat has disconnected the heater:

Possible cause 1:

- The equipment has been switched on without or with not enough water.

Possible solution:

1. Add water until the heater element cover rack, as indicated in chapter 8.1.
2. Press the safety thermostat restart button (19) and begin the sterilization process once again.

Possible cause 2:

- The thermostatic valve did not get closed after the purge phase, and the water in the sterilization chamber has been evaporated.

Possible solution:

1. It must be solved by a technical service authorized by J.P. Selecta. s.a.u.


Open lid indicator (3)

Possible cause:

- The indicator is on while the lid is open and must be off when closing it. If it does not switch off or on when operating, depressurize carefully the autoclave and contact an authorized technical service of JP Selecta, s.a.u.

10 Spare Parts

07001	Network connection cable
13007	Restart and heater contactors
15318	Fuse carrier
15595	Fuse (MICRO 8)-5A
16078	Pressure gauge
16212	Pressure switch 2 bar (134°C)
16220	Pressure switch 1 bar (121°C)
20061	Temperature selector
20096	Main switch
21230	Lid gasket (MED 12/20)
22003	Amber indicator light
22004	Red indicator light
22063	Green indicator light
25089	Spring lid elevator
34057	Three ways valve
34200	Safety valve
37025	Timer
39295	Heater element 1000W (MICRO 8)
39296	Heater element 1800W (MED 12)
39297	Heater element 2500W (MED 20)
34050	Purger (MED 12/20)
15577	Fuse 15A (MED 20)
15551	Fuse 10A (MED 12)

Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.