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
## Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.

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# 1. Information for users

## 1.1 Autoclave incoming

Handle the equipment with care. Unpack and check that all items coincide with the delivery note. If you see any discrepancy or damage of any kind, notify your distributor as soon as possible.

Keep original package some days, just in case of good returns.

## 1.2 Documentation

This documentation comes with the package or could be get from the manufacturer:

### - Instruction manual:

Read this manual before using the equipment.

These instructions are an important part of the equipment and should be made available to all users.

If you are in any doubt or want to clarify the use of this instrument, contact your local or the service department of J.P. Selecta, s.a.u.

### - Factory test:

Record from manufacturer checkings.

## 1.3 Other important information

The autoclaves sent to J.P. SELECTA, s.a.u. for repairs or maintenance must be decontaminated.

Any unauthorized modification, lack of maintenance or deliverate damage to the machine its not allowed by law. And could inhibit the warranty.

Do not use the equipment for fluids that can produce vapours which can form inflammable or explosive mixtures.

Use it only for the intended use of the equipment as specified in the catalogue, or manuals, ...

# 2. Warranty

This product is guaranteed for one year. It does not cover damages caused by an incorrect use or causes beyond the control of J. P. SELECTA, s.a.u.

Any manipulation of the equipment by personnel not authorized by J.P. SELECTA, s.a.u. automatically cancels the guarantee.

### 3. Technical Specifications

|  |  |    |
|--|--|----|
| Autoclave model:                           | Autester ST DRY PVL III 25L.   |    |
| Mains voltage:                             | 230V I   |    |
| Electrical power: (W)                      | 3000   |    |
| Maximum intensity: (A)                     | 13   |    |
| Fuses: (fast. type H, dimensions: 10x38mm) | 16A  |    |
| Network connection cable:                  | Phase, Neutral, Earth (section 1.5 mm <sup>2</sup> )   |    |
| Weight (Net Kg):                           | 60   |    |
| Useful measures: (cm Ø x depth)            | Ø30x35   |    |
| Ext. measures: (cm Depth x Width x Height) | 58x46x59   |    |
| Volume: (Litre)                            | 25   |    |
| Maximum load: (Kg) (Metal)                 | 5  |    |
| Maximum load : (flasks)                    | 250ml flasks   | 12 |
|  | 500ml flasks   | 9  |
|  | 1000ml flasks  | 3  |
| Water volume per cycle: (Litres)           | 0.8  |    |
| Water quality:                             | Demineralized over 100 micreS/cm   |    |
| Clean water tank volume: (Litres)          | 7  |    |
| Used water tank volume: (Litres)           | 7  |    |
| Air inlet filter:                          | Porosity 0.20 micres   |    |
| Acoustic level:                            | Less than 70 dBA   |    |
| Environmental conditions:                  | Inner use.<br>Altitude up to 2000m.<br>Temperature between 5°C & 40°C.<br>Max. relative humidity 80% for temperatures till 31°C, reducing in a linear way up to 50% of 40°C relative humidity. |    |
| Overpressure level:                        | Category II  |    |
| Pollution degree:                          | 2  |    |
| Vapour contact material                    | Stainless steel, copper, teflon, brass   |    |
| Frame material:                            | Steel with cover   |    |
| Vessel material:                           | Stainless steel AISI 304L  |    |
| Vapour generator material:                 | Stainless steel AISI 304L  |    |
| Pipes material:                            | Copper, brass and teflon   |    |
| RS-232 outlet:                             | 9 pin male connector. Pin 2 & 3 crossed, common pin 5. Speed 9600, no parity, data 8 bits, stop 1 bit.   |    |

### 4. Packing List

Autoclave delivered together with:

- Instruction manual (M.80330)
- Internal grid (59027)
- Tray holder (9740)
- 3 x trays (30237)
- Hex socket key to adjust door (44474)
- Bag rack (9746)
- Tray extractor (9236)

## 5. Overview

The **Autester ST DRY PV** are versatile and suitable autoclaves for a wide range of applications, sanitation, industrial processes and quality control. These autoclaves can sterilize solids, with or without being wrapped, liquids, culture media. They also have disinfecting drying cycles.

### 5.1 Purge

The purging operation is performed by the autoclave to remove air inside the boiler to fill the vapor space. Where there is no air steam.

The purge system must be suitable to sterilize intended load.

### 5.2 Atmospheric purge

In atmospheric vented comes out gravimetrically, impulsed by steam for a certain time.

Suitable for solid unwrapped (Tweezers, needles, scalpels ...).

### 5.3 Pulsed vacuum purge

A vacuum pump is used to remove air from the boiler. During one or more extractions purged air is produced.

Suitable for wrapped loads, catheters, ...

### 5.4 Bowie Dick test and Vacuum test

These tests are pre-programmed cycles.

Bowie Dick test used to verify operation of the purging by vacuum fractional.

The Vacuum Test is designed to verify the operation of the vacuum pump and sealing the autoclave. This test must be carried out with a cold equipment (it is usual to do it before starting the work sequence).

### 5.5 Sterilization not wrapped material

Place the instruments on waterproof paper and they do not touch each other.

### 5.6 Sterilization of wrapped materials and pipes

Only one instrument per bag.

Placing the bags upright they do not touch one another nor the walls of the autoclave.

The tubes must be clean, rinsed and drained. Ensure that the tubes do not form folds that obstruct the inner diameter.

### 5.7 Sterilization of vessels

Never place sealed containers.

Place the containers upside down to prevent water tanks.

### 5.8 Sterilization of Liquids

Placing fluids in suitable containers at a temperature and above trays so that they can collect spills. Avoid using containers with narrow necks.

The containers were filled about 2/3 full. Not close tightly, they will be plugged with cotton or some sort of cap to facilitate exhaust of the pressure vessel and not be formed.

At the end of the sterilization cycle, the autoclave was let cool freely (until pressure is 0 kg / cm<sup>2</sup>).

### 5.9 Water to use

Type of water is of vital importance in the operation of the equipment and its subsequent maintenance.

It is recommended to use distilled water with a conductivity level higher than 100 uS / cm so that the steam generator water level sensor works correctly.


Mains water may contain calcareous and ferritic residues that can damage the passage through valves and heating components.

It is also not advisable to use pure or ultra pure water because it is very corrosive and can even damage the stainless steel of the chamber.

### 5.10 Preheating

When starting the equipment, to accelerate the sterilization cycle, it heats the internal walls of the chamber. The temperature parameter can be selected from the SETTINGS menu.

## 6. Equipment description

### 6.1 Cycle types and temperature range

The autoclave allows you to configure the sterilization cycle according to the characteristics of the load to be sterilized:

| Program type | Purge type         | With internal probe | With external probe | Material type          |
|--------------|--------------------|---------------------|---------------------|------------------------|
| SOL1         | Atmospheric        | 105°C-134°C         | 105°C-130°C         | Unwrapped instruments  |
| SOL2         | Vacuum 2 fractions | 115°C-134°C         |                     | Bagged instruments     |
| SOL3         | Vacuum 3 fractions | 115°C-134°C         |                     | Double bagging, porous |
| LIQ          | Atmospheric        | 105°C-134°C         | 105°C-130°C         | Liquids, culture media |

Sterilizing time range: From 3 to 60 minutes

Drying time range: From 0 to 60 minutes

### 6.2 Pre-set cycles

The autoclave has 10 preset cycles, which can be modified according to the user's needs, except the cleaning program.

| Nr | Name                  | Temperature<br>°C | Dry<br>(minutes) | Time | Purge            |
|----|-----------------------|-------------------|------------------|------|------------------|
| 0  | Unwrapped instruments | 134               | 4                | 10   | 2x vacuum cycles |
| 1  | Wrapped instruments   | 134               | 12               | 30   | 3x vacuum cycles |
| 2  | Prions                | 134               | 20               | 20   | 3x vacuum cycles |
| 3  | Delicate unwrapped    | 121               | 15               | 15   | 2x vacuum cycles |
| 4  | Delicate wrapped      | 121               | 30               | 30   | 3x vacuum cycles |
| 5  | Instruments           | 126               | 12               | 25   | 3x vacuum cycles |
| 6  | Disinfection          | 105               | 25               | 15   | Atmospheric      |
| 7  | Liquid                | 121               | 30               | -    | Atmospheric      |
| 8  | Solid                 | 121               | 15               | 30   | Atmospheric      |
| 9  | Cleaning              | 134               | 4                | 30   | Atmospheric      |

### 6.3 Test cycles

Autoclave Autester **ST DRY PV III** has two test programs

|                    |         |        |
|--------------------|---------|--------|
| Test de Bowie Dick | 134°C   | 3' 30" |
| Vacuum Test        | 0,7 bar | 5+ 10' |

## 7. Installation

Place the autoclave near a power outlet suitable for the machine.

The autoclave will sit on a flat, horizontal, stable surface, leaving a space of 30 cm around the machine.

Open the lid. Check that there are no packaging elements inside the boiler. Verify that the boiler is clean and it does not contain particles that could damage the passage through the valves.

In case of overpressure, the safety valve at the rear part of the equipment will act. Avoid exposing this area to people or objects, to avoid possible damage.

## 8. Operation

**The autoclave lid only opens with the autoclave connected to the mains (Main Switch lighted)**

**Do not force the opening lever**

Open / Close Lid selecting button

### 8.1 Load

Open the lid.


Load the previously cleaned and settled material trays.

Check that the seal is in place.

Close the lid.

**Lift the lid closure handle, located to the right of the door, and bring the door closer until it makes contact with the camera. Once this is done, release the handle (which will fall under its own weight) and click on the "close" icon on the screen. A slight vacuum will be produced in the chamber that helps to automatically close the door, dropping the handle completely and allowing the start of a program.**

**Similarly, when the cycle is finished and you want to open, press the same button on the screen, in this way it is possible to raise the closing handle and the door can be opened more easily.**


### 8.2 Select Program

Select the suitable program for sterilization and drying the load.

To set the program (cycle) see 8.6.

Start the cycle by pressing RUN.

Shows running program status


### 8.2.1 Application of the cycle

During the cycle (also called program), autoclave's screen will indicate the segment that is being made.


Cycle Segments are:

- WATER : Water inlet to the boiler.
- VACUUM : Removing the air.
- HEAT : Steam generation.
- STERIL. : Sterilization of the load.
- EMPTY : Draining the boiler steam and water.
- UNPRESS : Depressurization boiler.
- DRY : Moisture extraction and drying
- OK : Sterilization conditions have been achieved according to cycle parameters.
- FAIL : Not achieved all the conditions for the sterilization.


This is the "End of cycle" message.


### 8.3 Quick guide for start auto-clave cycle (program).


1 Touch over «PROGRAMMING».


2 Choose a cycle  


3 Select the program.

4 Touch over  to go to main menu.

5 Press over «TASK».



6 Press over «RUN» to start the cycle.

### 8.4 Main menu

The autoclave control panel its from the touch screen.  
Press, with your finger or with a pointer, in the center of each icon.


### 8.5 Settings

Setting operation parameters as: Language, units, control probe...

Point over «SETTINGS».

Press on the parameter to be established and once established, press on the different keys to configure them.


#### 8.5.1 Language:

Set language between:

- Spanish.
- English.
- French.


Press to return to main menu or press to modify other parameters.

#### 8.5.2 Date and time:

Press over icons to modify day, month, ...

Press to return to main menu or press to modify other parameters.


#### 8.5.3 Units:

Choose units between °C or °F.

Select thousands separator format and decimal point.

Select date format.

Press to return to main menu or press to modify other parameters.

#### 8.5.4 Log download:



Select this option to download the logs to a USB memory. Follow the instructions indicated on the screen (have a USB stick at hand).

#### 8.5.5 Parameters protected by a password:


The password protects the parameters that may affect the operation of the autoclave. It is only accessible by contacting to J.P. Selecta, s.a.u. technical service.

- **Calibration of measuring devices "password 0000":**

- HEART: Offset value for heart temperature probe.
- WATER: Offset value for autoclave water temperature probe.  
Offset value for pressure sensor.  
Temperature for strip heaters.

Access to calibration and configuration parameters.


**- Activates /disables the water entry**

With this parameter set to «**activate**» the autoclave will take water from an external network for the steam generation.


With this parameter set to «**disabled**», the autoclave will take water from the tank itself, for the steam generation.

This option is not available on this model, it should always be disabled.

**- Printer (accessory only factory installable):**

Enables the record on the printer.  
"Printer Feed" feeds printer paper.


**- Select internal or external (heart) probe:**

**Internal probe:** The autoclave controller takes the vessel water temperature to control the cycle.

**External probe (accessory only factory installable):** Also called «heart probe», the autoclave takes the sample temperature to control the cycle.

Press to return to main menu or press to modify other parameters.


**8.5.6 Test cycles (Vacumm test & Bowie Dick)**

The autoclave has the two EN13060 test cycles to facilitate the validation and commissioning of the autoclave:

- Vacuum test: Qualifies the autoclave tightness and its capability to get suitable vacuum according to EN13060.
- BOWIE DICK test: Qualifies the steam quality and steam capability to penetrate into wrapped loads.


This test cycles should be run periodically, and results must be stored as quality system records. Press TEST button and select the desired test to be done.

**8.6 Sterilization programs**

Up to 9 cycles could be stored on the autoclave.

The cycle sets the treatment to apply over samples, by means of the setting parameters.


**8.6.1 Modify & save a program**

- Select the program to modify by pointing over:

- The program sets the treatment to apply over samples, by means of the setting parameters.

- Purge type.
- Sterilization temperature.
- Sterilization time.
- Drying time.

**Note:** Pointing over this icon will use selected cycle as a default.


Purge type

- **Setting purge type.**

Press on the purge type to be selected:

- Solid 1: 1 vacuum cycle. For unwrapped loads.
- Solid 2: 2 vacuum cycles. For wrapped loads.
- Solid 3: 3 vacuum cycles. For double wrapped loads, catheters, hollow devices, ...
- Liquid Not vacuum. Culture media in demi-close bottles.


Sterilization temperature

- **Setting sterilization temperature.**

Press over temperature label to set open the range from 115°C ... 134°C.


(Note: if the heart probe is used, the temperature range is affected from 115°C ... 130°C to avoid overpressure)

Sterilization time

- **Setting sterilization time.**


Press over sterilization time to set time from 3 to 60 minutes.

Drying time


- **Setting drying time.**

Press over sterilization time to set time from 0 to 60 minutes.

Save program changes

- **Store changes.**

Quit from main menu

- **Return to main menu.**

```
*****
14/05/2015 09:52
=====
PROG: 0
MODE: SOLID2
TEMP: 134°C
TIME: 04'
DRY : 10'
=====
LOT: 000010
*****
```

```
TIME TEMP. PRESS. STATUS

023°C +0.0bar VACUUM
021°C -0.2bar VACUUM
021°C -0.4bar VACUUM
021°C -0.5bar VACUUM
022°C -0.5bar HEAT
023°C -0.5bar HEAT
023°C -0.5bar HEAT
025°C -0.5bar HEAT
060°C -0.3bar HEAT
071°C -0.1bar HEAT
080°C +0.0bar HEAT
089°C +0.2bar HEAT
095°C +0.3bar HEAT
100°C +0.5bar HEAT
095°C +0.0bar VACUUM
090°C -0.1bar VACUUM
086°C -0.3bar VACUUM
083°C -0.4bar VACUUM
081°C -0.4bar VACUUM
078°C -0.5bar VACUUM
076°C -0.5bar VACUUM
079°C -0.4bar HEAT
088°C -0.2bar HEAT
095°C +0.0bar HEAT
101°C +0.1bar HEAT
106°C +0.3bar HEAT
110°C +0.5bar HEAT
114°C +0.7bar HEAT
119°C +0.9bar HEAT
122°C +1.2bar HEAT
126°C +1.4bar HEAT
129°C +1.7bar HEAT
133°C +2.0bar HEAT
134°C +2.0bar HEAT
03' 134°C +2.0bar STERIL
02' 134°C +2.1bar STERIL
01' 134°C +2.1bar STERIL
00' 134°C +2.1bar STERIL
    132°C +1.9bar EMPTY
    124°C +1.1bar EMPTY
    112°C +0.4bar EMPTY
    104°C +0.1bar EMPTY
00' 099°C +0.0bar DRY
01' 093°C -0.2bar DRY
02' 089°C -0.3bar DRY
03' 085°C -0.4bar DRY
04' 082°C -0.4bar DRY
05' 079°C -0.5bar DRY
06' 076°C -0.5bar DRY
07' 073°C -0.6bar DRY
08' 071°C -0.6bar DRY
09' 069°C -0.6bar DRY
    068°C +0.0bar UNPRES
RESULT: PASS
```

## 9. USB data recorder

The equipment has an internal memory for data acquisition (Datalogger). Once the cycle is finished, it is possible to download this information to a USB memory through a front USB connector.

The device (USB pen drive) must be formatted to FAT 32.

### Operation

- By default, the logger is active with a 1 minute logging interval. Follow the instructions to modify the recording interval (see 8.5).

With the cycle completed:

- Connect the «USB pen drive» to the front USB connector.

- When you connect a “USB pen drive”, the USB symbol will appear in the lower bar of the screen. (If the USB memory is not detected, try turning the autoclave off and on with the USB memory already connected.)

- From the touch screen select:  
SETTINGS> DOWNLOAD LOGS

A download button will appear and by pressing it all logs will be downloaded and the internal memory emptied.

- The screen will indicate when the recording of the logs has finished. Disconnect the «USB pen drive» from the device.

- Connect the «USB pen drive» to a computer to download the data file.


**Important Note:** The «USB pen drive» must have a device name. It cannot be blank.

The created file is of text type, that can be opened with WordPad, WORD or EXCEL.

**Note:** If the data is not downloaded, when the device memory is at 80% capacity, a warning message will appear indicating that the memory needs to be downloaded.

If the message is ignored and the memory is completely filled, the oldest cycles will be erased.

Approximately the memory has capacity to store about 150 cycles (it will depend on the total time each cycle has).

## 10. Alarms

If a failure occurs during operation, a message will appear.

Please contact Technical Service or a distributor if you get one of these messages.

| ERROR | MESSAGE               | DESCRIPTION   |
|-------|-----------------------|---|
| 01    | PRESSURE SENSOR FAIL  | Pressure sensor error. Failure sensor or disconnected cable.  |
| 02    | INTERN UNDETEMP       | Low temperature. Internal temperature does not reach <b>8°C</b> .   |
| 03    | INTERNAL OVERTEMP     | Over temperature error. Internal temperature has exceeded maximum limit of <b>150°C</b> .   |
| 04    | EXTERNAL OVERTEMP     | External heating element overtemperature. External temperature has exceeded maximum limit of <b>150°C</b> .   |
| 05    | OVERPRESSURE          | The pressure has exceeded maximum admissible pressure of <b>+2.5bar</b> .   |
| 06    | WATER TANK EMPTY      | Clean water reservoir empty. Fill the equipment water tank.   |
| 07    | WATER TANK FULL       | The used water tank is full. Once the cycle is finished, connect the hose to the front and empty the tank.  |
| 08    | WATER ERROR           | Failed to get water. Blocked filter.  |
| 09    | LEVEL ERROR           | Steam generator sensor error or blocked tubes. More than 200s filling water.  |
| 10    | DOOR OPEN             | Door open. Close the door to start the computer.  |
| 11    | PRESSURE FAIL         | Pressure error during purge phase.<br>Insufficient water or water heater failure.   |
| 12    | SETPOINT UNREACHEABLE | Set value error not reached. No water in the system or heater failure.  |
| 13    | STERILIZATION FAIL    | Sterilization error. Set point margins overcome. There is air inside.   |
| 14    | VACUUM TEST FAIL      | Vacuum error. Vacuum leakage or pump failure.   |
| 16    | DEPRESS FAIL          | Maximum depressurization time exceeded error.   |
| 17    | PURGE ERROR           | Maximum purged time exceeded error.   |
| 18    | EMPTY TOO LONG        | Empty error, takes too long to empty the system. Valves or drain blocked or failure.  |
| 23    | PRINTER PAPER END     | No paper error, printer requires paper.<br>(ONLY THERMAL PRINTER RS-485)  |
| 24    | PRINTER ERROR         | Printer error. Possible paper jam. (ONLY THERMAL PRINTER RS-485)  |
| 25    | POWER FAIL            | Power error. Detected in the start-up.<br>The equipment has been switch off before the end of the cycle.  |
| 26    | COMMUNICATION FAIL    | Communication error between touch screen and autoclave. This error message will appear when, after several attempts, messages do not receive any reply. |
| 27    | PRINTER OFFLINE       | Printer not connected. (ONLY THERMAL PRINTER RS-485)  |

## **11. Safety Systems**

### **11.1 Safety pressure switch**

When an abnormal overpressure is produced in the sterilization chamber, the safety pressure switch acts in a way to stop the autoclave and it opens the aeration electrovalve.

### **11.2 Safety valve**

If an abnormal overpressure is produced, the safety valve is opened, expelling the steam.

This valve is calibrated to initiate the opening between 2.2 kg/cm<sup>2</sup> and 2.5 kg/cm<sup>2</sup>.

In case of safety valve reset, the machine must be prepared to advise a technical personnel.

### **11.3 Safety thermostat**

When an important increase or decrease of water level is produced the safety thermostat acts by stopping the heater.

It is necessary to reset it so that the autoclave returns to its normal operation. Press the white button on the back of the device.

### **11.4 Door blocked**

In case of open lid blocked by a safety system failure make as follows:

- 1 Stop the machine.
- 2 Disconnect the machine from the mains.
- 3 Ensure that the boiler has no pressure
- 4 Take the necessary precautions to avoid possible burns.
- 5 Remove the device cover unlock the safety piston. The access is in the right side of the equipment.

## 12. Maintenance

**Before removing the lid, disconnect the machine from the mains.**

**The control board should only be changed by a technical service or authorized personnel.**

### **CLEANING:**

- Wall panels: with a damp cloth.
- Limescale deposits: use home cleaning products.
- TFT Screen: use glasses cleaning products.
- Cover and plastic areas: use non-abrasive cloth dipped in alcohol.

### **Monthly operations:**

- Change water in the tanks (clean and used)
- Record results from vacuum test and Bowie dick test.

### **Three months operations:**

- Change air entry filter of 0.2 microns from the back.

### **Yearly operations:**

- Check and contrast temperature autoclave against traceable temperature standard.
- Check safety valve operation.

**ATTENTION:** this operation should be run for qualified personnel. Dangerous contact with hot steam.

- Check, visually the state of electrical wiring, ground and mains.