STERIL-PRES

Autoclave para esterilización
Autoclave for sterilization

STERIL-PRES USB 80L	code 4002442
STERIL-PRES USB 150L	code 4002443
STERIL-PRES GMP 80L	code 4002444
STERIL-PRES GMP 150L	code 4002445

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INDEX

1. General Information	46
2. Packing List	46
3. Technical Specification	46
3.1 Electrical installation	47
3.2 Water inlet installation	48
3.3 Installation of drain system	48
4. Introduction	48
4.1 Temperature control	49
5. Autoclave Description	50
6. Installation	52
6.1 Location	52
6.2 Water Supply & Connections	52
6.3 Connection to the drain	53
6.4 Connection to the electrical network	54
6.5 Protective earthing	54
7. Operation	54
7.1 Start-up	55
7.2 Cleaning the chamber	55
7.3 Load	55
7.4 HMI (Human-Machine Interface) Operation	55
7.4.1 Main Menu	55
7.4.2 Configuration Menu:	57
7.4.3 Recipe Menu:	63
7.4.4 Home/Trends Menu:	66
7.4.5 Hydraulics Menu:	67
7.4.6 User Administrator Menu:	68
8. Error Messages & Troubleshooting	68
9. Values for F ₀ :	72
10. Safety Thermostat:	73
10.1 Rearming the safety thermostat	73
10.2 Adjusting the safety thermostat	73
11. Emergency manual lid opening:	74
12. Recycling & dismantling of the autoclave:	75
13. Spare parts:	76
13.1 Spares for STERIL-PRES code 4002442 & 4002444 (80L)	76
13.2 Spares for STERIL-PRES code 4002443 & 4002445 (150L)	79
14. Maintenance:	82
14.1 Cleaning	82
14.2 Periodic maintenance tasks	82
15. Warranty:	83

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Action	Access level						
	Guest	Operator	Maintenance	Administrator	Selecta		
Visualization of screens and logging in and silence alarms	X	X	X	X	X		
Visualization of screens, logging off and change his password.		X	X	X	X		
Start and Stop.		X	X	X	X		
Alarms control		X	X	X	X		
Modify date and time			X	X	X		
Analog inputs adjustment.			X	X	X		
Modify general parameters				X	X		
Activate and deactivate PLC outputs			X	X	X		
Create, modify, unlock and delete users.				X	X		
Variables modify					X		

TABLE 8.1: ACCESS LEVELS SUMMARY

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ENGLISH

1. General Information

- Handle the package with care. Unpack it and check that the content matches what is indicated in the "Packing list" section. If any damaged component or the absence of any is observed, notify the distributor quickly.
- Do not install or use the equipment without first reading this instruction manual. These instructions are an inseparable part of the equipment and must be available to all its users.
- Consult any doubt to the technical service of J.P. SELECTA, s.a.u.
- Do not use the equipment in explosive, flammable or corrosive atmospheres.

2. Packing List

Item	150L	80L
	Cod. 4002443 / 4002445	Cod.4002442 / 4002444
Autoclave STERIL-PRES USB / GMP	(x1) cod. 4002443 / cod. 4002445	(x1) cod.4002442 / cod.4002444
Reinforced Hose 3/4" x 1500mm (Steam-Water-Drain)		(x3) cod. 46167
Instruction Manual		M.80461
EU Declaration of Conformity		

3. Technical Specification

Item	4002443 / 4002445	4002442 / 4002444
	(150L)	(80L)
Required Voltage (V):	3ph+N+T/T 400V	1ph+N+T/T 230V
Installed Power (W):	9000	5600
Max. Intensity (A):	23 A	26 A
Fuses (quick, type H, dimensions: 10x38mm):	32 A	32 A
Power Cord:	3f+N+T/T (4mm ²)	1f+N+T/T (4mm ²)
Net weight (Kg):	250	125
Useful Dimensions (cm Ø x prof.):	Ø50x70	Ø40x60
External Dimensions : (cm depth. x wide x height):	90x78x122	72x58x118
Nominal volume: (Liters):	150	80
Max. Load: (Kg) (Metal):	50	20
Total dimensions + lids:	Add 1.2 m in height	Add 40 cm in height
Min pressure water inlet (bar):	3	

Item	4002443 / 4002445 (150L)	4002442 / 4002444 (80L)
Volume of water per cycle (Liters)	25	20
Water quality	* Soft water (see clasification table)	
Dissipated energy to environment: (W / h)	600	525
Water Inlet connection	3/4"	
Accoustic level:	Inferior a 70 dBA	
Operating conditions of cycle:	Continuous with a breaktime of 20 minutes between each cycle.	
Environmental conditions:	Indoor use. max. altitude 2000m	
Room temperature	Temperature range between 5°C & 40°C.	
Max relative humidity	80% max.	
Overvoltage rate	Cat. II	
Pollution grade	2	
Steam contact materials	Stainless steel AISI 304, Copper, PTFE, brass	
Chamber material:	Stainless steel AISI 304	
Piping Material:	Copper & PTFE	
USB Connection	YES	YES
ETHERNET connection	4002445	4002444
Printer connection (RS-232)	YES	YES

3.1 Electrical installation

Check the electrical installation of the enclosure where the autoclave will be installed.

For your safety, we advise you to feed the autoclave with a breaker and a differential switch.

Size	Voltage	Power (W)	Breaker
80L	230V (1ph+N+PE)	5600	32A
150L	400V (3ph+N+PE)	9000	32A

3.2 Water inlet installation

In areas with “hard” water (conductivity above 100uS / cm) it is necessary to use water decalcification equipments.

It is recommended to work with soft water, to extend the equipment life

The water inlet must have a pressure between 3 and 4 bar (absolute pressure)

If it is necessary, install the following:

- A pressure group, if the pressure of the water network is less than 2 bar.
- A hand-operated device if the water network pressure is higher than 4 bar.

Water hardness classification table:

Water type	mg/l	*fH	*dH	*eH
Soft water	≤17	≤1.7	≤0.95	≤1.19
Slightly hard water	≤60	≤6.0	≤3.35	≤4.20
Moderately hard water	≤120	≤12.0	≤6.70	≤8.39
Hard water	≤180	≤18.0	≤10.05	≤12.59
Very hard water	>180	>18.0	>10.05	>12.59

mg/l: milligrams per liter of calcium carbonate water (CaCO₃)
***fH:** French grade (10.0 mg CaCO₃/l)
***dH:** German grade (17.8 mg CaCO₃/l)
***eH:** English grade (14.3 mg CaCO₃/l)

3.3 Installation of drain system

Located at ground level or below.

Material to resist temperature of 100°C.

(Gray PVC plumbing is suitable). Please see chapter 6.3 for details.


4. Introduction

The “STERIL-PRES” autoclave has been designed to perform thermal treatments of sterilization of containers (metal, glass, plastic) with any content.

The autoclave has the capacity to store 60 heat treatment programs with 8 steps each.

All programs can be programmed by the user.

The autoclave works by spraying overheated water to transmit the heat (temperature) on the containers subjected to heat treatment.

4.1 Temperature control

The autoclave uses two temperature probes:

- **Water temperature probe:**

It indicates the water temperature of the autoclave and it is located at the bottom of the autoclave, inside the autoclave chamber.

- **Product temperature or “heart” probe:**

It indicates the temperature of the product, this is required for the F_0 control.

The F_0 value is calculated by the autoclave software from the temperature recorded by the heart probe inserted in a container and the thermoresistance characteristics of four types of microorganism according to the thermal treatment to which the product is subjected.

In case of a sterilization process: *Geobacillus stearothermophilus* or user definable.

The F_0 value is therefore subjected to time and temperature to which the products are submitted during the respective thermal treatment (above 100°C in sterilization processes).

The total numerical value F_0 , obtained when completing the heat treatment, will indicate the degree of sterilization of the product. There are reference F_0 values depending on the type of product and the nature of the preserves to be obtained.

The product probe will give us the temperature readings inside the package in which it is inserted, these temperature values used to count the time that the product is subjected to as scheduled, as long as you are in product probe mode.

Remember:


For the F_0 calculation to be valid, the “heart” (or “product”) probe must be inside a sample container. Please see chapter 7.4.3.1 & chapter 9 for more info.

5. Autoclave Description




Steril-Pres autoclave belongs to a series of over-heated water circulation autoclaves that are distinguished for the following features:

1. Arm / joist locking system.
2. PLC controlled system, with HMI (human machine interface) to monitorize all systems and programming.
3. Safety locking system that only allows lid opening with pressure 0 bar and a temperature under 85°C.
4. Safety valve and maximum pressure limiter that disconnects the heating element in case of overpressure, max. 2.5 bar.
5. Limiter thermostat to protect the heating element, set at 150°C.

AUTOCLAVE CODE. 4002442 / 4002444 (80L)

AUTOCLAVE CODE. 4002443 & 4002444 (150L)

Excess steam Outlet
Water Inlet
Water Outlet


6. Installation

6.1 Location

Left a free space of 50 cm around the autoclave.

Provide a water tap and a drain box, following next picture, see additional info in chapter 6.2 & 6.3

Locking Wheels of 80L autoclave

Locking Wheels of 150L autoclave

Place the autoclave on a flat, horizontal and stable surface, suitable to the weight of the autoclave.

For an ergonomic use, leave a free space of about 50cm around the machine. To facilitate maintenance operations, it is advisable to have direct access to all sides of the autoclave.

Immobilize the autoclave, adjusting it to the ground, by using the stops located in the front part of the appliance, provided for that purpose.

Place the autoclave near an electrical outlet, water, and drain.

According to the legislation on safety, the autoclave must be located in such a way that the discharge outlet of the safety valve does NOT point to people or can reach them in the case of steam escaping.

In the autoclave, the safety valve outlet points to the rear.

6.2 Water Supply & Connections

Water type	mg/l	*fH	*dH	*eH
Soft water	≤17	≤1.7	≤0.95	≤1.19
Slightly hard water	≤60	≤6.0	≤3.35	≤4.20
Moderately hard water	≤120	≤12.0	≤6.70	≤8.39
Hard water	≤180	≤18.0	≤10.05	≤12.59
Very hard water	>180	>18.0	>10.05	>12.59

mg/l: milligrams per liter of calcium carbonate water (CaCO₃)
 *fH: French grade (10.0 mg CaCO₃/l)
 *dH: German grade (17.8 mg CaCO₃/l)
 *eH: English grade (14.3 mg CaCO₃/l)

It's required a water tap with G-3/4" thread (same as washing machines)

It's required a water pressure between 3 and 4 bar. If the water pressure is greater than 4 bar, it is advisable to install a reducer.

Must be connected one of the hoses between the tap and the "water inlet".

Example of installation for water inlet. It's required a G3/4" to connect the water inlet hose to the autoclave

6.3 Connection to the drain

ATTENTION

Fix the hose so that it cannot move due to the pressure at the end of the cycle.

Do not hold the hose with your hands. This hose can reach temperatures above 60°C.



The autoclave will eject excess steam, hot air and hot water during sterilization. They should be connected to a drain, catch basin..

There are required two hoses:

«**CONDENSED OUTPUT**»: This hose must **be firmly attached to the drain**, to prevent the force of the steam from removing the hose from the drain.

«**WATER OUTPUT**»: This hose must **be firmly attached to the drain**, to prevent the force of the steam from removing the hose from the drain.

BOTH HOSES SHOULD DISCHARGE TO A LOWER LEVEL THAN THE HOSES CONNECTIONS AT THE AUTOCLAVE

Only the hoses are supplied with the autoclave. Final customer must supply the water inlet tap & draining piping to attach the hoses (threads type G 3/4")

6.4 Connection to the electrical network

Plug into an electrical outlet, ready to supply 32A continuously.

In the model of 150L, the electrical cable doesn't include any connector or plug, must be directly connected by means of a terminal block or a suitable plug & socket. In the model of 80L, a CETAC 32A (1ph+N+Earth) is included.

It is advisable to install an ICP bipolar circuit breaker of about 32A located at an accessible point quickly and without obstacles.

Consult the type and level of the supply voltage and the maximum intensity value on the equipment's rating plate.

6.5 Protective earthing

The installation of electric power supply for the autoclave **MUST** have a protective earth conductor.

The protection against the electrical risks of the autoclave requires that the network cable, marked with a yellow-green color, to be connected to the ground conductor of the installation.

7. Operation

ATTENTION

**THE AUTOCLAVE IS SUPPLIED WITH THE CHAMBER LID CLOSED.
THE DOOR OPENING IS ONLY POSSIBLE BY ACTIVATING THE ELECTRIC POWER SUPPLY OF
THE AUTOCLAVE.**

7.1 Start-up

Main switch

1. Start up by using the main switch.
2. Make sure that the water tap is open.
3. Make sure that the drain is secure.

It is advisable to use water with little content of lime or deionized / demineralized water. DO NOT use distilled water.

7.2 Cleaning the chamber

To do easy the cleaning of the autoclave chamber, when the main switch is pressed on, the chamber drain is opened.

It's recommended to spray water with a small amount of cleaner (JP Selecta BIO-SEL cleaner recommended) the walls of the chamber. The water will go down the drain.

Run the cleaning program periodically.

7.3 Load

Place the **product probe** in a container with the same characteristics as the rest of the containers to be sterilized.


7.4 HMI (Human-Machine Interface) Operation

7.4.1 Main Menu

All control panel functions are selected via the touch screen.

Press with your finger or a pointer in the center of the touch screen.

The autoclave operation diagram is shown in the following diagram:

Below is what the start menu looks like when the HMI starts up, the entire circuit of the autoclave and its main components:

At the left of the screen, there is the button to start the recipe that has been programmed, (see more info in this chapter)

At the left of the screen, there is the manual opening button for the autoclave lid. It will only be activated if the safety conditions of the equipment allow it, in which case it will be visible the status, by means the padlock drawing.

7.4.2 Configuration Menu:

To access the configuration menu, press at the home menu, the button:

7.4.2.1 User profile:

By default, the system works with a "User-Operator" profile; to manage other users and permissions, they must be activated through the corresponding screen. To do this, click on the "users" button:

There are 5 predefined user levels, with the following permissions:


Action	Access level				
	Guest	Operator	Maintenance	Administrator	Selecta
Visualization of screens and logging in and silence alarms	X	X	X	X	X
Visualization of screens, logging off and change his password,		X	X	X	X
Start and Stop,		X	X	X	X
Alarms control		X	X	X	X
Modify date and time			X	X	X
Analog inputs adjustment,			X	X	X
Modify general parameters				X	X
Activate and deactivate PLC outputs			X	X	X
Create, modify, unlock and delete users,				X	X
Variables modify					X

TABLE 8.- ACCESS LEVELS SUMMARY

7.4.2.2 Date / Time settings:

Allows the selection of the date and time, to change, press directly on the values shown on the screen.


7.4.2.3 Language:

This menu, allows to choose the language of the menus and settings. ES / EN / FR.

7.4.2.4 Settings:

It allows the visualization of the temperature and pressure values of the different probes and measuring elements of the equipment. It allows the application of offsets to correct possible deviations that may occur in the measuring elements. Depending on the permissions available to the user, it may or may not be possible to make changes.

Prior to modify any internal parameter, please contact with JP Selecta Aftersales Dpt.

In order to modify any parameter will be done by means the "Internal parameters" menu.

- Water probe: Located inside the autoclave chamber at the bottom. It's the main probe and controls general temperature of the autoclave.
- Product Probe: Probe located in contact with the product to measure "real" temperature inside the jar.
- Pressure: Allows to check the pressure inside the chamber in real time.

By means these buttons allows to see more parameters in different screens.


This button allows to return to adjustments menu

This is the list of internal parameters and their values:

Parameter	Default Values		Units	Description
	80 L	150 L		
A	120	120	(s)	Time limit to reach maximum level
B	10800	10800	(s)	Maximum time to reach setpoint
C	85	85	(°C)	Maximum temperature for door opening
D	28	28	(s)	Delay in low-level switch detection
E	4	4	(s)	Hysteresis product probe to overtemp
F	30	40	(s)	Drain time at the end of the process
G	5	5	(s)	Filling time during cooling
H	30	30	(s)	Waiting time during cooling
I	2	2	(°C)	Setpoint offset of the second temperature probe


7.4.2.5 Maintenance:



Screen that allows you to monitor the behavior of any element of the autoclave and check that the autoclave alarms work correctly, manually forcing the signals of the different components or sensors.

In addition, the autoclave have a batch counter. When arriving to the defined “batch limit”, a message will appear on screen to notify to the user than a maintenance revision must be performed. The autoclave can continue on operation but it’s recommended the maintenance revision as soon as possible.

At the Steril-Press Check can be revised any existing alarm, and locate the problem easily.

All alarms are identified by means a number and a description. In case of existing any alarm, will be lightened to identify it clearly.

In some cases will be required to contact JP Selecta aftersales dpt. to solve the problem. In other cases can be solved by the user. (more info in chapter 8)

By means the “active alarms”, more info will be shown, and it’s possible to download a historical alarms report by means the USB port.

7.4.2.6 Communications (only for models including ETHERNET connection):

It’s required a ftp server to enable the autoclave communications service. The user must fill the parameters according the local net configuration.

Very Important: In order to store properly the data from the autoclave, the user must create into the ftp server, at the root, a folder named:

\sterilpres_serial number/

Where serial number is the effective serial number of the autoclave


By means this button it's possible to modify the parameters.

7.4.2.7 Report Backup:


It allows to send the stored internal memory information by means FTP or USB. *Please, see chapter 7.3.2.6, about storing folder to be implemented into the FTP server.

The files that can be transferred are:

- Process files *.pdf
- Alarm list *.pdf
- Process file *.csv (allowing to be post-processed into Ms-Excel or others programs.)


7.4.3 Recipe Menu:

Located at the home menu, can be activated by means the indicated button.

This menú, allows to generate a tailor-made process of sterilization, according the type of product to be treated.

7.4.3.1 F₀ Custom:

The equipment incorporates three lethality tables of three reference microorganisms and a self-calculated lethality table so that the client can incorporate the microorganism of their choice to calculate the total lethality value of their heating treatment.

Lethality (L): Time (minutes) which it is necessary to heat a product to a certain temperature (T) in order to obtain the same lethal effect as heating to the reference temperature (Tx) in 1 minute, singular value for each microorganism.

The lethality values sum (L) during the time (minutes) the heat treatment lasts will give us the total lethality value or the final F₀ sterilization value.

Default lethality values (F₀)

Default tables:

F₀ STERILE: Lethality table of microorganism *Geobacillus stearothermophilus*. This is the reference microorganism, for all thermal sterilization treatments.

Calculated lethality values (F₀ Calculated)

F₀ Calculated option allows the creation of a lethality table according to the reference microorganism that the client wishes, to obtain the total lethality value of the programmed heating treatment.

Note: When editing the program, the value of “F₀ Calculated” will be calculated automatically.

The partial lethality values are obtained by applying the formula:

$$L = 10^{((T-T_x) / z)}$$

By entering the **Z** and **T_x** values (unique to each microorganism), the custom lethality table will be created. The definition of these parameters is:

Z: Temperature to reduce a reference microbial population by a logarithmic power.


T_x: Reference temperature.

Most common reference values:

Reference bacteria	T _x (°C)	Z (°C)
<i>Geobacillus stearothermophilus</i>	121,1	10

7.4.3.2 Recipe Programming:

In this chapter, it's explained the process to configure a sterilization process, by means the recipe menu:

- Step 1: **PRE** (water preheating): This step heats the water in the autoclave to be hotter than the product temperature.
 - Step 2: **PUR** (air purge): In this step, it's generating a depressurization of the chamber along some time.
 - Step 3: Thermal treatment cycle: Can be programmed from 1 to 8 steps, to define the complete thermal treatment: heating stage; stabilization stage & cooling stage:
- a) **Heating stage:** Initial step in establishing the temperature, time and pressure to reach as a step prior to the maintenance step. This step can be ignored if there is no technological interest in carrying out previous steps before reaching the maintenance phase.
- b) **Stabilization stage:**

Step where the heat treatment (sterilization) itself is carried out. In this step it is required to establish time, temperature and pressure to which the product will be subjected (activated product probe) or the water in the autoclave (activated water probe) to reach the F₀ value of sterilization.

c) Cooling stage:

Stage in which once the desired degree of sterilization is reached, the product is cooled down. It is necessary to establish a cooling temperature and a pressure during determined times for each step. It is important to maintain an overpressure throughout this phase to avoid, especially in large containers, that product can escape from the containers due to excessive differences between the internal pressure of the product and the autoclave chamber.

The equipment during the cooling stage, follows the sequence of filling/emptying the chamber until the programmed temperature is reached.


The required time of this phase has the purpose of maintaining the inner temperature during the programmed time.

Note that cooling is performed by using the inlet water. It is not possible to cool below the inlet water temperature. If the programmed temperature is lower than the inlet water temperature, the autoclave will show an error because it will not be able to achieve this temperature.

During the operation, it's possible to modify the time values for each step in order to adjust the programming process to the product features.

It's possible to store a recipe for a determined product, and give a "customized name", by clicking inside recipe edit screen at the cell "name". A keyboard will be shown to define the new name.

7.4.3.3 Programming according samples type:

The composition characteristics of products are a key factor that determine the heat penetration capacity in the product internals during the thermal treatment. This fact is demanding the ability of the heat generated by the autoclave to heat the entire mass of products in a homogeneous way. According the thermal transfer capacity of products, they can be: not very conductive and also conductive products.

In case of non-conductive products, it is recommended to:

- The product probe must be placed centered into the container.
- Keep in mind that it will take longer to reach the maintenance temperature.

In case of conductive products, it is recommended to:

- Keep in mind the thermal inertia in the cooling phase, the product cools quickly and the final heat treatment may not reach the predetermined values.

Suggestion for counter-pressure for every temperature:

Temperature (°C)	Pressure
20 ...100°C	1 bar
100...115°C	1.6 bar
116...121°C	1.7 bar


7.4.4 Home/Trends Menu:

Located at the home menu, can be activated by means the indicated button.


This menú, allows to visualize graphically in real time the evolution of defined parameters:

- Temperature detected by water probe
- Temperature detected by product probe
- Pressure.
- Calculated F₀


By pressing this button it's possible to see the historical graphic in detail, including visualization tools, located at the top of the graphic.


7.4.5 Hydraulics Menu:

Located at the home menu, shows the PID diagram of the autoclave and the status of each component.

7.4.6 User Administrator Menu:

Located at the home menu, allows to get access to the user, with the defined rights (please, see chapter 7.4.2)

8. Error Messages & Troubleshooting

As defined in chapter 7.4.2.5, the Steril-Pres autoclave, can detect an abnormal situation, showing an alarm to identify the problem.

In case the cycle is interrupted, for example: due to power failure, lack of water, etc ... it must be considered that the load has not been sterilized and the sterilization must be repeated.

As a general rule: disconnect the autoclave, by the general switch, wait a few seconds, and start it up again. If the door is blocked, wait about 20 minutes for the temperature of the chamber to descend spontaneously.

ATTENTION

The process can be accelerated by opening the manual depressurization valve. Close the valve before restarting. Take precautions, do not stand in front of the outlet of the manual depressurization valve. Open the valve slowly.

The following table indicates the error messages, their meaning and the possible solution:

1	E2PROM Fail	2	Pressure sensor failure	10	Empty water tank
12	Door open	13	Overpressure	16	Overtemperature water probe
19	Unreachable setpoint	20	Drain time too long	23	Power failure
29	Safety thermostat	30	Level error	33	Product probe failure
34	Undertemperature water probe	36	Under-temperature sterilization	37	over-temperature sterilization

38 Vacuum drain pressure timer	39 Sterilization pressure failure	41 Preheating failure
42 Heating failure during purging	44 Short circuit in product probe	45 Comm. error product probe
48 Short circuit in water probe	49 Comm. error water probe	53 Pressure swith PSH010

Code	Name	Description
01	<p align="center">E2PROM FAIL</p> <p>Cause: Controller malfunction. Solution: Contact the technical service.</p>	(Internal error E2PROM)
02	<p align="center">PRES SENSOR FAIL</p> <p>Cause: Pressure sensor does not work properly Solution: Check the pressure sensor wiring.</p>	(Pressure sensor reading <0 bar)
10	<p align="center">MINIMUM WATER LEVEL</p> <p>Cause: There's not water enough inside chamber Solution: Revise water supply or check level control</p>	(Minimum water level)
12	<p align="center">DOOR OPEN</p> <p>Cause: "Open door" switch does not work properly. Solution: Revise the locking swith device and the activation</p>	(Door open.)
13	<p align="center">OVERPRESSURE</p> <p>Cause: The pressure has reached a dangerous value. Solution:</p> <ol style="list-style-type: none"> 5) Stop the autoclave by using the main switch. 6) Depressurize the autoclave with the manual valve.TAKE PRECAUTIONS for depressurization: <ol style="list-style-type: none"> a. Do not stand in front of the valve outlet. b. Open the valve slowly. c. Contact the technical service. 7) DO NOT USE THE AUTOCLAVE. 	(Press. sensor reading >3.5 bar)
16	<p align="center">INTERNAL OVERTEMP WATER PROBE</p> <p>Cause: The pressure has reached a dangerous value. Solution: The pressure has reached a dangerous value.</p> <ol style="list-style-type: none"> 1) Stop the autoclave by using the main switch. 2) Depressurize the autoclave with the manual valve.TAKE PRECAUTIONS for depressurization: <ol style="list-style-type: none"> a. Do not stand in front of the valve outlet. b. Open the valve slowly. 3) Check the connections of the water and heart probes. 4) If the connections are correct: <ol style="list-style-type: none"> a. Contact the technical service. 5) DO NOT USE THE AUTOCLAVE 	Temperature higher than 150°C

Code	Name	Description
19	<p>UNREACHABLE SETPOINT</p> <p>Cause: The pre-set temperature cannot be reached. Solution: Increase the temperature of the initial preheat phase. Jars too big. Product with little thermal conduction. Heating resistance(s) cut off. Controller failure SSR failure (solid state relay).</p>	<p>It has not been possible to reach the set point during the sterilization phase in less than 360 minutes.</p>
20	<p>DRAIN TIME TOO LONG</p> <p>Cause: The chamber has not been cleared after 30/40 minutes Solution: Check that the water temperature is low enough to cool the cans. Check drain valve. Check and clean the outlet filter. Check the connection of the drain hose.</p>	<p>Some of the phases of the emptying take more than 30/40 min</p>
23	<p>POWER FAILURE</p> <p>Cause: During an operation cycle, the power has failed. Solution: Restart the cycle.</p>	<p>(Power cut off with the equipment operating)</p>
29	<p>SAFETY THERMOSTAT</p> <p>Cause: Inner temperature has reached a too high value. Solution: Wait for a few minutes and reset the thermostat Increase parameter "B" in 5s.</p>	<p>(Safety performance of the boiler thermostat)</p>
30	<p>LEVEL ERROR</p> <p>Cause: No water enters the chamber. Solution: Check the water inlet tap Increase parameter "B" , more water enters.</p>	<p>(Water level is not detected)</p>
33	<p>PRODUCT PROBE FAILURE</p> <p>Cause: Abnormal temperature reading. Solution: Check probe status & connection and replace if required.</p>	<p>(Product Probe failure)</p>
34	<p>UNDERTEMP. WATER PROBE</p> <p>Cause: If this happens at the beginning of the cycle, it is likely that the water in the network is very cold (winter). Solution: In the autoclave loaded with cans, pour half a bucket of hot water into the boiler and then start the autoclave.</p>	<p>(Reading water probe temperature lower than 8°C for more than 5 minutes during the HEAT phase or in any moment during any other phase when the program is operating).</p>
36	<p>STERILIZ. UNDERTEMP</p> <p>Cause: In the sterilization phase, the temperature has dropped 2°C below the set point for 30s. Solution: Check the temperature probe and the location inside the jar. Make sure there has been no power interruptions.</p>	<p>(Temperature drop)</p>
37	<p>STERILIZ. OVERTEMP</p> <p>Cause: In the sterilization phase, the temperature has exceeded the set point at 4°C for 5s. Solution: Revise cycle parameters.</p>	<p>Excessive temperature</p>

Code	Name	Description
38	<p>VACUUM DRAIN PRESSURE TIMER</p> <p>Cause: The temperature has been above the cooling temperature and the pressure has been below the cooling pressure for 60 minutes.</p> <p>Solution: Revise cycle parameters.</p>	(Timeout)
39	<p>STERILIZATION PRESSURE FAILURE</p> <p>Cause: During the sterilization phase, have not reach the purge pressure in 30 seconds.</p> <p>Solution: Check that the manual release pressure valve is not open, that the equipment is not leaking, that no any solenoid valve is opened</p>	
41	<p>PREHEATING FAILURE</p> <p>Cause: The preheat temperature is not reached before 60'.</p> <p>Solution: Check (Reset) the safety thermostat. Check the heating elements</p>	Excessive water pre-heating time
42	<p>HEATING FAILURE DURING PURGING</p> <p>Solution: Check (Reset) the safety thermostat. The product is very cold, increase the temperature of the PRE phase.</p>	The purge temperature is not reached before 60'
44	<p>SHORT CIRCUIT IN PRODUCT PROBE</p> <p>Solution: Revise Product probe & connections to locate the damage</p>	
45	<p>COMM. ERROR PRODUCT PROBE</p> <p>Solution: Revise Product probe & connections to locate the damage</p>	
48	<p>SHORT CIRCUIT IN WATER PROBE</p> <p>Solution: Revise water probe & connections to locate the damage</p>	
49	<p>COMM. ERROR WATER PROBE</p> <p>Solution: Revise water probe & connections to locate the damage</p>	
53	<p>PRESSURE SWITCH PSH010</p> <p>Cause: The reading from the pressure switch is abnormal</p> <p>Solution: Check connections & replace if required</p>	

9. Values for F₀:

It is a numerical value (minutes) that parameterizes the time required to destroy a known number of reference pathogenic bacteria or spores (sterilization), at a determining temperature. The final F₀ value is the result of a combination of times and temperatures that occur inside the container (at the coldest point) during the heat treatment. A certain F₀ value (sterilization value) can be reached with different combinations of times and temperatures (see table of Sterilization). A higher temperature requires less time and vice versa, to obtain the same degree of sterilization (F₀) as the case may be.



Temp.	F ₀
100	0,007
101	0,009
102	0,012
103	0,015
104	0,019
105	0,024
106	0,030
107	0,038
108	0,048
109	0,061
110	0,077
111	0,097
112	0,122
113	0,154
114	0,194
115	0,244
116	0,303
117	0,388
118	0,488
119	0,615
120	0,774
121	0,974
122	1,227
123	1,544
124	1,944
125	2,448
126	3,082
127	3,830
128	4,885
129	6,150
130	7,745

Table 3: Sterilization F₀ values

* tables according to Reichert, 1977

10. Safety Thermostat:

As indicated in chapter 5, this safety device is located at the rear right hand of the autoclave, and limits the temperature inside the autoclave chamber.

8. Safety thermostat regulation screw (grey).
9. Safety thermostat reset button (white).

10.1 Rearming the safety thermostat

Safety thermostat is an element to prevent the temperature from exceeding a value. If by a fortuitous event the thermostat is triggered, that is, it opens the power circuit, the heating element will stop working. Usually there is an indicator light, or an error message, that indicates that the thermostat is in this position.

If the thermostat is triggered, check that the equipment looks good and it is not overheating at any point.

To change its status, simply press the white button (9), and you will hear a “click” that indicates that the thermostat has been reset. If the equipment has not cooled down sufficiently, the thermostat cannot be reset, and you must wait for the temperature to drop before it can be reset.

10.2 Adjusting the safety thermostat

In some equipment such as ovens and baths and with certain processes it may be useful to adjust the thermostat trigger value, but this operation should NEVER be performed to the rest of equipments.

From factory, the thermostat is adjusted to its trigger temperature, set at 150°C.

Stop the equipment and leave it to cool with the lid open for 15 minutes.

11. Emergency manual lid opening:

In case the door opening is blocked due to a failure of the security system, follow the following protocol:

- 1 Stop the machine.
- 2 Unplug it from the network.
- 3 Put a nipple and a tube on the aeration valve.
- 4 **Open the valve carefully so that the steam escapes. Take the necessary precautions to avoid burns.**
- 5 Wait for the steam to completely empty.
- 6 Remove the upper casing where the display is located to move it aside (it is not necessary to disconnect any of the cables). Does not use screws, pull firmly.
- 7 Remove the plastic cap (PULL) and pull lightly until you check that the lid closing handle is released.

ATTENTION

IMPORTANT:

MANUALLY OPENING THE AUTOCLAVE LID WHILE PRESSURE AND TEMPERATURE EXIST, MAY CAUSE AN ACCIDENT DUE TO BURNS PRODUCED BY THE HOT WATER AND EXISTING STEAM, PLEASE CHECK AND FOLLOW ALL THE INDICATED STEPS.

12. Recycling & dismantling of the autoclave:

Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.

13. Spare parts:

13.1 Spares for STERIL-PRES code 4002442 & 4002444 (80L)

		
16070 Solenoid Valve G1/4"	16067 Water Valve	43121 Pressure Sensor
		
39370 Water Heater	16138 Pressure Switch	43172 Temp. Water Probe
		
43173 Product Probe L=100 mm Cable ø 4mm	43173 Product Probe L=145 mm Cable ø 6mm	43174 Product Probe L=200 mm Cable ø 6mm
		
15712 Product Probe gland for cable ø6mm	21309 Gasket for probe 43173 L=100mm ø 4mm	25094 door gas strut

		
<p>20134 Water Level Sensor</p>	<p>20096 Main Switch</p>	<p>16170 3-Way Valve G1/4"</p>
		
<p>13019 Contactor</p>	<p>36028 Solid State Relay</p>	<p>10054 Door Locking Switch</p>
		
<p>36030 Relay</p>	<p>24361 Compressor</p>	<p>34336 Check Valve</p>
		
<p>34358 Full Cone Nozzle</p>	<p>34400 Safety Valve</p>	<p>34556 2-Way Valve</p>
		
<p>43168 Thermostat</p>	<p>46167 Water hose G3/4"</p>	<p>7488 32 Amp. Connector</p>

		
<p>24353 Water Pump</p>	<p>15409 Power Fuse 32A</p>	<p>21271 Door Gasket</p>
		
<p>40039 Autoclave Wheel</p>	<p>16144 Solenoid Valve 1/2" NC</p>	<p>16145 Solenoid Valve 1/2" NO</p>
		
<p>85283 Temperat. Transmitter</p>	<p>55303 Filter Mesh</p>	

13.2 Spares for STERIL-PRES code 4002443 & 4002445 (150L)

		
<p>40026 Autoclave Wheel</p>	<p>16067 Water Valve</p>	<p>43121 Pressure Sensor</p>
		
<p>39372 Water Heater</p>	<p>16138 Pressure Switch</p>	<p>43172 Temp. Water Probe</p>
		
<p>43173 Product Probe L=100 mm Cable ø 4mm</p>	<p>43173 Product Probe L=145 mm Cable ø 6mm</p>	<p>43174 Product Probe L=200 mm Cable ø 6mm</p>
		
<p>15712 Product Probe gland for cable ø6mm</p>	<p>21309 Gasket for probe 43173 L=100mm ø 4mm</p>	<p>25080 door gas strut</p>

		
20134 Water Level Sensor	20096 Main Switch	85283 Temperat. Transmitter
		
13019 Contactor	36028 Solid State Relay	10054 Door Locking Switch
		
36045 Relay	24362 Compressor	34412 Check Valve
		
34358 Full Cone Nozzle	34400 Safety Valve	21206 Door Gasket
		
43168 Thermostat	46167 Water hose G3/4"	55303 Filter Mesh

		
<p>24363 Water Pump</p>	<p>15410 Power Fuse 20A</p>	<p>16079 3-Way Valve G1/4"</p>
		
<p>16145 Solenoid Valve 1/2" NO</p>	<p>16146 Solenoid Valve 3/4" NO</p>	<p>16147 Solenoid Valve 3/4" NC</p>

14. Maintenance:

14.1 Cleaning

To avoid damaging both the wear and the appearance of the working chamber in the long term, it is important to periodically clean and prevent the formation of remains.

You have to avoid working with products that may get rusty. The oxide residues cause the stainless steel to become contaminated. If the usual work is contaminated, they should be cleaned and polished immediately.

14.2 Periodic maintenance tasks

Daily

Check the condition of the cover gasket and its placement to prevent closing the lid being pinched or cut.

Weekly

Check the condition of the product probe, that the cable does not have cuts in the cover, or excessive folds that could damage the internal cables.

Check that the low level control buoy (at the bottom of the chamber) is free to move.

Carry out a cycle with only clean water to eliminate possible product residues. You can add a few drops of anti-degreasing, dishwasher-like, when doing this cycle, as it will help to clean the pipes (recommended JP Selecta BIO-SEL).

Monthly

Clean the drain filter (on the front of the equipment).

Check that the consumption of the resistance is adequate. In the three-phase model, check that the consumption is the same in each phase.

Check the water diffuser nozzles on the top of the boiler. Disassemble and ensure the correct operation.

Annually

Remove the side covers and check that there is no leak. In the event of a leak, repair the union of the fitting.

Check the operation of the water diffusers. If necessary, disassemble and clean.

Clean the compressor inlet filter.

If the use of the equipment is high (more than 1 cycle per day) it is advisable to replace the door seal.

NOTE: The Steril-Pres is a “pressure equipment”, and it has been designed according the European Pressure Equipment Directive 2014/68/UE.

In all european countries, exists a regulation about the periodical inspections to be performed on this kind of devices according the category of the autoclave.

It's a responsibility of the autoclave owner, to follow the national regulations about the installation and inspections.

15. Warranty:

The autoclave has a one-year warranty.

The warranty does not cover damage caused by improper use or due to causes beyond the control of J.P. SELECTA, s.a.u.

Any manipulation of the device by personnel not authorized by J.P. SELECTA, s.a.u. automatically voids the benefits of the warranty.