Autoclaves con ciclos de vacío y secado Dry & Vacuum cycles autoclaves

Autester ST Dry-S USB 50L code 4002521
Autester ST Dry-S USB 80L code 4002522
Autester ST Dry-S USB 150L code 4002523

Autester ST Dry GMP-S 50L code 4002524
Autester ST Dry GMP-S 80L code 4002525
Autester ST Dry GMP-S 150L code 4002526

INDICE

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INDEX

1. Qualified installation of the autoclave	45
2. General information for the user	45
2.1 Receiving the autoclave	45
2.2 Documentation	45
2.3 Other information of interest	46
2.4 Warranty period	46
3. Technical Specifications	47
4. Packing List	48
5. Introduction and recommendations	48
5.1 Purging	48
5.2 Atmospheric purging	48
5.3 Fractional vacuum purging	49
5.4 Bowie Dick Test and Vacuum Test	49
5.5 Sterilization of unwrapped material	49
5.6 Sterilization of wrapped material and tubes	49
5.7 Containers Sterilization	49
5.8 Liquid sterilization	49
5.9 Water to use	50
5.10 Preheating	50
6. Equipment Description	51
6.1 Front View (mod. 50 & 80L)	51
6.2 Rear View (mod. 50 & 80L)	52
6.3 Front View (mod. 150L)	53
6.4 Rear View (mod. 150L)	54
7. Installation	55
7.1 Outside water ingress	55
7.2 Drainage	55
7.3 Measurements of the external connections of the autoclave	56
8. Operation	57
8.1 Main Menu	57
8.2 Settings Menu:	58
8.2.1 User Management	58
8.2.2 Date and time	59
8.2.3 Language	60
8.2.4 Settings	60
8.2.5 Maintenance	61
8.3 Recipe Menu:	64

8.3.1 Recipe Programming	65
8.4 Start-Trend Menu:.....	66
8.5 User Change Menu:.....	67
8.6 FTP Start-Transfer Menu:.....	67
9. Protection Systems – Security	68
9.1 Safety pressure switch.....	68
9.2 Safety Valve	68
9.3. Unlocking the lid	68
10. Maintenance Plan	69
11. Troubleshooting.....	70
12. Spare parts	75
12.1 50/80L Model	75
12.2 150L Model	78

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AUTOCLAVE MANUFACTURED IN CONFORMITY WITH THE PRESSURE EQUIPMENT DIRECTIVE 2014/68/UE

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Action	Access level					
	Guest	Operator	Maintenance	Administrator	Selecta	
Visualization of screens and logging in and silence alarms	X	X	X	X	X	
Visualization of screens, logging off and change his password.		X	X	X	X	
Start and Stop the AUTESTER		X	X	X	X	
Alarms control		X	X	X	X	
Modify date and time			X	X	X	
Analog inputs adjustment.			X	X	X	
Modify general parameters				X	X	
Activate and deactivate PLC outputs			X	X	X	
Create, modify, unlock and delete users.				X	X	
Variables modify					X	

TABLE 8.: ACCESS LEVELS SUMMARY

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ENGLISH

1. Qualified installation of the autoclave

In the European Community, Directive 2014/68/EU establishes the mandatory qualified installation of this autoclave.

The owner of the equipment is responsible for complying with the provisions of Royal Decree 809/2021 on pressure equipment, regarding installation, legalization, and periodic inspections. (applicable to Spain)

This installation must be performed by a "Qualified Pressure Vessel Installer."

For installation in Spain, you can request information from J.P. SELECTA, S.A.U.'s technical service for a list of qualified installers in your autonomous community.

The documents required for processing the qualified installation are as follows:

- The EC Declaration (document that comes with the equipment).

The following documents can also be requested from **calidad@jpselecta.es**:

- Functional test.
- EC Declaration of the safety valve.

The qualified installer will also inform you about periodic inspections.

2. General information for the user

2.1 Receiving the autoclave

Handle the package appropriately for its weight and consider it FRAGILE.

Unpack and check the "Packing List." If you notice any damaged or missing components, notify the distributor.

Keep the original packaging for a period of time. After this time, dispose of it properly, separating its components: cardboard, wood, and plastic.

2.2 Documentation

The following documents are supplied with the autoclave:

- Instruction manual:

Do not install or use the equipment without first reading this supplied instruction manual.

These instructions are an integral part of the equipment and must be available to all users.

Any questions should be consulted with the autoclave supplier.

- EU Declaration of Conformity.

2.3 Other information of interest

Autoclaves sent to J.P. SELECTA, s.a.u. for repair or maintenance must be decontaminated. Failure to comply with this standard will result in the equipment being returned to the sender.

Any modification, removal, or lack of maintenance of any autoclave device is prohibited.

The autoclave may not be used in explosive environments or with substances that may release vapors or form explosive or flammable mixtures.

If this autoclave is used in a manner not specified by J.P. SELECTA, s.a.u., the protection provided and its functionality may be compromised.

2.4 Warranty period

The autoclave has a one-year warranty.

The warranty does not cover damage caused by improper use or causes beyond the control of J.P. SELECTA, s.a.u.

Any tampering with the appliance by personnel not authorized by J.P. SELECTA, s.a.u. automatically voids the warranty.

3. Technical Specifications

Autoclave code	4002521 / 4002524	4002522 / 4002525	4002523 / 4002526
Autoclave volume	50L	80L	150L
Net capacity	48L	78L	140L
Mains Voltage: (other under request)	Monophase 230V-50Hz		Triphase 400V-50 Hz
Supplied connector	CETAC 32A – 1PH+N+T/T		None
Electrical Power: (W)	3600	4850	14000
Max. Intensity: (A)	16	21	22.5
Fuses: (quick type. type H, dimensions)	20A	32A	32A
Cable section	1Ph+N+Earth (sec 2,5 mm ²)	1Ph+N+Earth (sec 4 mm ²)	3 Phase+N+Earth (sec 6 mm ²)
Net Weight (Kg):	108	136	250
Useful Dimensions (cm Ø x depth)	Ø30x63	Ø40x60	Ø50x70
Ext dimensions: (cm depth x width x height)	117X48X62	119x58x72	118x80x90
Required space + lid	Add 0,7m to the height		Add 1m to the height
Water supply max. pressure (bar)	3		
Water supply min. pressure (bar)	0.5		
Required water per cycle (Liters)	0,75	1,5	5,0
Water quality	Demineralized, conductivity > 100 microS/cm		
Thermal energy transferred to room (W/h)	750		850
Water volume of tank (Liters)	10		15
External water connection	3/4"		
Air filter type	Porosity 0.20 microns cod. 15645		
Accoustic Level	Lower than 70 dBA		
Operation running level	Continuous, 20 minutes of intermittent stops between cycles		
Room conditions	Indoor use. Max. height 2000m Room temperatura between 5°C to 40°C. Relative humidity max. 80%		
Overvoltage Level	Category II		
Pollution Level	2		
Materials in contact with steam	Stainless steel, Copper, PTFE, Brass		
Cabinet material	Galvanized Steel		

Autoclave code	4002521 / 4002524	4002522 / 4002525	4002523 / 4002526
Chamber material	Stainless Steel AISI 316		
Steam generator material	Stainless Steel AISI 304		
USB Connection	YES		
ETHERNET Connection	YES		
Printer Port (RS-232)	YES		
GMP – Audit Trail	Only GMP models: 4002521 / 4002522 / 4002523		

4. Packing List

The autoclave is delivered with:

- Instruction manual
- EU Declaration of Conformity.
- Inner grate at the bottom of the boiler.

5. Introduction and recommendations

These autoclaves allow sterilization processes of solids with and without packaging, liquids, culture media and also allow disinfection as well as drying cycles.

The material must be clean, free of any type of incrustation or residue.

Do not overload the trays, racks or baskets and leave a space between them of 1 or 2 cm, in order to allow the passage of steam and facilitate drying.

5.1 Purging

Purging is the operation performed by the autoclave to remove the air from inside the boiler and thus be able to fill the entire space with steam.

The purging system must be suitable for the load to be sterilized.

5.2 Atmospheric purging

In atmospheric purging, the air exits the outside in a gravimetric method, driven by steam for a certain time.

Suitable for solid unwrapped material (tweezers, needles, scalpels...).

5.3 Fractional vacuum purging

A vacuum pump is used to extract air from the boiler. During purging, one or more air extractions occur. This system is much more effective than atmospheric purging.

Suitable for wrapped loads, catheters, ...

5.4 Bowie Dick Test and Vacuum Test

The Bowie Dick and Vacuum tests are pre-programmed cycles.

The Bowie Dick test is used to verify the performance of fractional vacuum purging.

The Vacuum Test is intended to verify the operation of the vacuum pump and the tightness of the autoclave. This test should be carried out with cold equipment (it is usually carried out before starting the work sequence).

5.5 Sterilization of unwrapped material

The instruments must be placed on water-repellent paper and the pieces must not touch each other.

5.6 Sterilization of wrapped material and tubes

Only one instrument per bag must be placed.

Place the bags in an upright position that does not touch each other or the walls of the autoclave.

The tubes should be cleaned, rinsed and drained. Make sure that the tubes do not form bends that obstruct the inner diameter.

5.7 Containers Sterilization

Never place hermetically sealed containers.

Place the containers upside down to avoid water deposits.

5.8 Liquid sterilization

Place liquids in containers suitable for working at sterilization temperature and place them on top of trays, so that possible spills can be collected.

The containers will be filled about 2/3 of their capacity, not exceeding this value. DO NOT CLOSE them hermetically, they will be covered with cotton or some type of stopper that facilitates the exit of air from the container and does not form pressure inside them.

Avoid using narrow-necked containers.

At the end of the sterilization cycle, allow the autoclave to cool freely (until the pressure is 0 kg/cm²).

5.9 Water to use

The type of water is of vital importance in the operation of the equipment and in its subsequent maintenance.

It is recommended to use distilled water with a conductivity degree greater than 100 uS/cm for the steam generator water level sensor to work properly.

Mains water may contain limescale and ferritic residues that can damage the passage through valves and heated components. Use soft water (see table) to avoid fouling.

It is also not advisable to use pure or ultra-pure water because it is very corrosive and can even damage the stainless steel of the camera.

Water Type	mg/L	*fH	*dH	*Yes
Soft water	≤17	≤1.7	≤0.95	≤1.19
Slightly hard water	≤60	≤6.0	≤3.35	≤4.20
Moderately hard water	≤120	≤12.0	≤6.70	≤8.39
Hard water	≤180	≤18.0	≤10.05	≤12.59
Very hard water	>180	>18.0	>10.05	>12.59


mg/l: milligrams per liter of calcium carbonate (CaCO₃) water
***fH:** French grade (10.0 mg CaCO₃/l)
***dH:** German grade (17.8 mg CaCO₃/l)
***eH:** English grade (14.3 mg CaCO₃/l)

5.10 Preheating

When starting the equipment, to speed up the sterilization cycle, heaters are located at the walls of the chamber.

6. Equipment Description



6.1 Front View (mod. 50 & 80L)


Item	Description
1	Touch screen
2	Lid opening and closing lever
3	General switch
4	USB connection
5	Lid gasket
6	Water tank filler cap
7	External probe access port
8	Outdoor water intake 3/4"G
9	Safety valve
10	Manual steam drain valve
11	Mains plug
12	Air inlet filter
13	Condensate emptying socket ø10mm
14	ETHERNET socket
15	RS-232 Printer Connection
16	Wheel
17	Brake leg
18	Water tank emptying valve ø10mm
19	Temperature probe


6.2 Rear View (mod. 50 & 80L)



6.3 Front View (mod. 150L)

Item	Description
1	Touch screen
2	Lid opening and closing lever
3	General switch
4	USB connection
5	Water tank filler cap
6	Lid gasket
7	External probe access port
8	Outdoor water intake 3/4"G
9	Safety valve
10	Manual steam drain valve
11	Mains cable
12	Air inlet filter
13	Condensate drain socket 3/4"G
14	ETHERNET socket
15	RS-232 Printer Connection
16	Wheel
17	Brake leg
18	Water tank emptying valve \varnothing 10mm
19	Temperature probe
20	Drain Mesh Filter

6.4 Rear View (mod. 150L)

7. Installation

Place the autoclave near a power outlet suitable for the machine's consumption.

The autoclave will sit on a flat, horizontal and stable surface, leaving a free space of 30cm around the machine.

Place the autoclave, adjusting it to the ground by means of the stops located on the front of the equipment, provided for this purpose.

Open the lid. Place the rack at the bottom of the boiler. Check that there is no element of the packaging inside. Check that the boiler is clean and that it does not contain particles that could damage the passage through the valves.

Attach a tube to the "CONDENSATE OUTLET" spigot to collect the condensate that is formed during the phases in which the vacuum pump operates. Do not place the collection container at a level higher than that of the teat.

7.1 Outside water ingress

The autoclave can be filled via the outside network using a 3/4" socket. **The supply pressure should never exceed 3bar.**

This form of filling must be configured from the configuration panel.

7.2 Drainage


ATENCIÓN

Steam output in case of overpressure.

Risk of burns.

Position this equipment in the appropriate way to prevent this risk from occurring.


ATTENTION

Do not use the external emptying option without first connecting a pressure hose. Fix the hose so that it cannot move due to the effect of pressure at the end of the cycle. Do not hold the hose with your hands, you may suffer burns.

7.3 Measurements of the external connections of the autoclave

Model	50L	80L	150L
Water Inlet	Thread 3/4" G	Thread 3/4" G	Thread 3/4" G
Emptying clean water tank	Spigot ø10mm	Spigot ø10mm	Spigot ø10mm

8. Operation

The autoclave lid is only opened with the autoclave plugged into the mains (illuminated switch)

Do not FORCE the opening lever


Plug the cable into the base of the mains socket.

Connect the device by operating the main switch (4) and the main screen will appear.


8.1 Main Menu

All functions on the control panel are selected by the touchscreen.

Tap with your finger or a pointer in the center of the touchscreen. The operating diagram of the autoclave is shown in the following diagram:

Below is what the start menu looks like when you start the autoclave, in it, you can see the entire circuit of the computer and its main components:

On the left of the screen, there is the button to start the recipe that has been programmed, as described in point 8.3 of this manual.

On the left of the screen, there is the manual opening button of the autoclave lid, it will only be activated if the safety conditions of the equipment allow it, in which case it will be visible in green

8.2 Settings Menu:


To access the settings menu, click on the icon:

8.2.1 User Management

By default, the system works with a "User-Operator" profile, to manage other users and permissions, it must be activated through the corresponding screen, to do this, click on the user management icon:


There are 5 predefined user levels, with the following permissions:

Action	Access level						
	Guest	Operator	Maintenance	Administrator	Selecta		
Visualization of screens and logging in and silence alarms	X	X	X	X	X		
Visualization of screens, logging off and change his password.		X	X	X	X		
Start and Stop the AUTESTER		X	X	X	X		
Alarms control		X	X	X	X		
Modify date and time			X	X	X		
Analog inputs adjustment.			X	X	X		
Modify general parameters				X	X		
Activate and deactivate PLC outputs			X	X	X		
Create, modify, unlock and delete users.				X	X		
Variables modify					X		

TABLE B.: ACCESS LEVELS SUMMARY


8.2.2 Date and time

It allows the selection of the date and time, to change, to click directly on the values shown on the screen.


8.2.3 Language

Allows choice of language from menus and settings. English or Spanish.


8.2.4 Settings

It allows the visualization of the temperature and pressure values of the different probes and measuring elements of the equipment. It allows the application of offsets to correct possible deviations that may occur in the measurement elements. Depending on the permissions available to the user, it will be possible or not to make changes.

- Internal Probe: Water probe installed at the bottom of the equipment basin.
- External Probe: Probe installed in the heating jacket of the equipment bucket.
- Product Probe: Temperature probe to monitor the temperature of the product to be treated.
- PSG: Steam Generator temperature probe.
- Pressure: Allows monitoring of the pressure value recorded in the autoclave tray.
- Internal parameters: Allows you to modify certain values of the internal operation of the autoclave. **(consult with our technical department before making any modification of parameters)**

Using this key it is possible to advance through all the screens of internal parameters to locate the one you want to view and/or modify.


Using this key it is possible to return to the previous parameter screen.


This key allows you to return to the settings screen.


8.2.5 Maintenance


Screen that allows you to monitor the behavior of any autoclave element and check that the autoclave alarms are working correctly, manually forcing the signals of the different components or sensors.

With the "Check autester" function, a synoptic is displayed, which allows you to change its status by clicking on each component (valves, etc.) to verify that it is working properly. When you click on the component, it changes color.

With the "Check alarms" function, you can access a list with all the defined alarms, and if there is an active alarm, the color of the alarm changes, making it possible to immediately identify the problem


By pressing the active alarms button, you can see in detail the problem and the history.

By pressing this button you can view the alarm history of the device.

Detail of the screen with the history of alarms generated by the equipment. The search for alarms can be filtered by different methods and it is also possible to generate a PDF report with the entire alarm history that will be sent via FTP.

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

It will be essential to have an ftp server to enable communications. The user must enter the parameters based on their network configuration.

This key provides access to the configuration screen, to enable the autoclave HMI to access the customer's network.

8.3 Recipe Menu:

To access the recipes menu, from the start menu, you must press the indicated key.

The system allows you to edit recipes defined by the customer or those predefined on the computer. In the same way, there are defined vacuum and Bowie-Dick tests.


8.3.1 Recipe Programming

If you click on the type of product to be processed, a drop-down menu will appear, where you can select the type of program, with the factory presets that can be modified:

Program	Purge	Temperature	Time	Drying
Non-wrapped tools	3	121	20	3
Wrapped tools	2	121	20	3
Prions	3	134	10	0
Disinfection	0	115	20	3
Liquid	0	121	20	0
Agar	0	121	20	0
Liquid (product probe)	0	121	20	0
Cleaning	0	134	20	0
Bowie & Dick 121	3	121	15	0
Bowie & Dick 134	3	134	4	0



NOTE: Both the Bowie Dick tests and the vacuum test cannot be modified

To program a recipe, one of the predefined types will initially be loaded and the time and temperature values will be modified at will, using the "Edit recipe" key. The "load recipe" key will then be pressed and the programmed values will appear in the recipe on the right side of the screen, ready to be executed in the autoclave

8.4 Start-Trend Menu:

It allows the visualization in real time of the parameters of the water probe temperature and product and pressure of the autoclave. This menu is accessed, from the start menu, by means of the indicated key

If desired, it is possible to view the graph in its entirety, accessing it by means of the indicated key

8.5 User Change Menu:

Allows User change. It is accessed from the autoclave's home menu, using the indicated key:

User Name:

Password:

Show password

The user must have registered in the system via the configuration menu, as described in point 8.2

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9. Protection Systems – Security

9.1 Safety pressure switch

When there is an abnormal overpressure in the sterilization chamber, the safety pressure switch acts to stop the autoclave and open the aeration solenoid valve.

9.2 Safety Valve

If an abnormal overpressure occurs, the safety valve will be opened and the steam will be expelled.

ATTENTION: the steam outlet must be oriented in a direction that cannot reach the personnel.

This valve is calibrated to 2.3 kg/cm²


In the event of a trip by the safety valve, the machine must be stopped and the service service must be notified.

9.3. Unlocking the lid

In the event of a blockage of the door opening due to a breakdown of the security system, follow the following protocol:

- Stop the machine.
- Unplug it from the mains.
- Put a teat and tube in the aeration valve.
- Open the valve carefully to allow the steam to escape.
- Take the necessary precautions to avoid burns.wait for the total emptying of the steam.
- Remove the upper case where the display is to move it to one side (it is not necessary to disconnect any of the cables). Do not use screws, pull tightly.
- Remove the plastic cap and pull lightly until the handle is released.

10. Maintenance Plan

Before removing the cover from the machine, disconnect it from the mains.

The control plate should only be handled by duly authorized maintenance personnel.

Cleaning:

- Outer Enclosure: Clean with stainless steel polish.
- Lime incrustations. Use household anti-limescale cleaner.
- Touch screen: use a cloth and a glasses cleaning liquid.
- Plastic caps: use a cloth and a glasses-cleaning liquid.
- Perform a **Monthly Cleaning Program:**
- Clean the drain filter located at the bottom front.
- Changing the water in the tank. Empty the water tank using valve 18 (in the 50 & 80L autoclave)
- Record the result of the Bowie Dick Test and the Vacuum Test.

Quarterly:

- Change the air intake filter.

Annual:

- Contrast the temperature and pressure of the autoclave against equipment traceable to recognized standards.
- Check the operation of the safety valve.
- **ATTENTION:** this operation must be carried out by qualified personnel. There is a danger of burns from contact with the steam.
- Visually check the status of the electrical connections, earthing and connection cable to the electrical network.

- ***Note - Cleaning schedule:**

- Make sure the boiler is empty.
- Select the cleaning schedule as indicated in point 8.3
- Fill with water manually (about 2 liters of water) and pour through the boiler door.
- Close the lid.
- Run the program normally.

It is recommended to execute this program after an intense day of work or once a week to avoid deposits in the pipes.

11. Troubleshooting

The AUTESTER GMP/USB autoclave can detect an abnormal operating situation, showing the user an alarm and thus identifying what the problem is.

In case the cycle has been interrupted, for example: power failure, lack of water, etc.. The load should be considered not to have been sterilized and therefore the heat treatment cycle should be repeated.

As a general rule: disconnect the autoclave, using the start switch, wait about 20 seconds and start again. If the door is locked, wait about 20 minutes for the temperature of the chamber to drop spontaneously.

The process can be streamlined by opening the manual depressurization valve. Close this valve again before restarting any process. Take precautions, do not stand in front of the manual depressurization valve. Open this valve slowly, high temperature steam may come out.

The following table indicates the error messages set, their meaning, and their possible solution:

1 E2PROM Fail	2 Pressure sensor fail	10 Max filling time
12 Door open	13 Overpressure	14 Vacuum fail
15 Pressure fail	16 Internal probe overtemp	17 External probe overtemp
18 Sterilization fail	19 Setpoint unreachable	20 Too long emptying
23 Power failure	29 Safety thermostat	30 Level error
31 Purge error	33 Product probe overtemperature	34 Internal probe undertemperature
36 Sterilization undertemp	37 Sterilization overtemp	44 Product probe shortcircuit
45 Product probe comm. error	46 External probe shortcircuit	47 External probe cut wire
48 Internal probe shortcircuit	49 Internal probe cut wire	50 PSG probe shortcircuit
51 PSG probe cut wire	52 PSG probe overtemp	53 Pressure switch PSH010
54 Pressure switch PSH024		

Cod.	Name	Description
01	E2PROM FAIL Cause: Controller failure. Solution: Contact Technical Service	(Internal error E2PROM)
02	PRESSURE SENSOR FAILURE Cause: Pressure sensor not working properly Solution: Check the sensor wiring.	(Pressure sensor reading <0 bar)
10	MAXIMUM FILLING TIME Cause: The predefined time has been exceeded Solution: Check water supply or level sensor status	(Tank filling time)
12	DOOR OPEN Cause: Door Hold Actuator Fails Solution: Check the door latch or its relay	(Lid is open)
13	OVERPRESSURE Cause: The pressure has reached a dangerous level. Solution: 1. Stop the autoclave with the main switch. 2. Depressurize the autoclave with the manual depressurization valve. 3. Take precautions: - Do not stand in front of the depressurization valve (high steam temp.) - Open the valve slowly. - Contact technical service. 4. DO NOT USE THE AUTOCLAVE.	(Pressure sensor >3.5 bar)
14	VACUUM FAILURE Cause: Chamber vacuum is not maintained Solution: Check the fit of the door and vacuum pump hoses. The machine must be cold, if the autoclave is started, as it preheats up, the volume increases and can generate the error	(Vacuum chamber failure)
15	PRESSURE FAILURE Cause: The pressure detected is abnormal and does not reach the set set Solution: Check for Leaks in Hose	(Pressure failure)
16	INTERNAL PROBE OVERTEMPERATURE Cause: The temperature has reached a dangerous level. Solution: 1. Stop the autoclave with the main switch. 2. Depressurize the autoclave with the manual depressurization valve. 3. Take precautions: - Do not stand in front of the depressurization valve (high steam temp.) - Open the valve slowly. - Contact technical service. 4. DO NOT USE THE AUTOCLAVE.	Temperature > 150°C

Cod.	Name	Description
17	<p>EXTERNAL PROBE OVERTEMPERATURE</p> <p>Cause: The temperature has reached a dangerous level.</p> <p>Solution:</p> <ol style="list-style-type: none"> 1. Stop the autoclave with the main switch. 2. Depressurize the autoclave with the manual depressurization valve. 3. Take precautions: <ul style="list-style-type: none"> - Do not stand in front of the depressurization valve (high steam temp.) - Open the valve slowly. - Contact technical service. 4. DO NOT USE THE AUTOCLAVE. 	Temperature > 150°C
18	<p>STERILIZATION FAILURE</p> <p>Cause: The temperature during the sterilization phase has reached temperatures below the setpoint temperature or has exceeded the permissible temperature according to EN 13060</p> <p>Solution: Check Temperature Probe</p>	(Sterilization failure)
19	<p>SETPOINT UNREACHABLE</p> <p>Cause: The selected temperature is not reached.</p> <p>Solution:</p> <ul style="list-style-type: none"> Extra-large jars Product with low driving value Heating Resistor Fails Driver failure Failure of the solid-state relay. 	Impossible to reach temperature during the sterilization phase in less than 360 min.
20	<p>EMPTYING TOO LONG</p> <p>Cause: The chamber has not been emptied in 30-40 minutes</p> <p>Solution: Check that the cooling temperature is low enough</p> <ul style="list-style-type: none"> Check Drain Valve Check and clean filter. Check hose connections 	The emptying phase exceeds 30-40 minutes.
23	<p>POWER FAILURE</p> <p>Cause: During a cycle the power has failed</p> <p>Solution: Restart cycle</p>	(Power failure with equipment in operation)
29	<p>SAFETY THERMOSTAT (steam generator)</p> <p>Cause: Excessively high temperature.</p> <p>Solution: Wait a few minutes and reset (modify parameter value B)</p>	(High temp caldera thermostat)
30	<p>LEVEL ERROR</p> <p>Cause: No water enters the chamber.</p> <p>Solution: Check the water supply faucet (wall)</p> <ul style="list-style-type: none"> Increase parameter "B", more water enters 	(Insufficient water level)
31	<p>PURGE ERROR</p> <p>Cause: Air is not properly removed from the chamber</p> <p>Solution: Check pump and purge solenoid valve</p>	(Purge Error)

Cod.	Name	Description
33	PRODUCT PROBE OVERTEMPERATURE Cause: Abnormal reading of the product probe. Solution: Check the status of the probe and its connections.	(Product probe failure)
34	UNDERTEMPERATURE WATER PROBE Cause: The temperature of the inlet water is excessively low (winter) Solution: In the autoclave, pour a bucket of WARM water into the tank	(Water probe reading below 8°C for more than 5 minutes during the warm-up phase or any other time or phase)
36	UNDERTEMPERATURE STERILIZATION Cause: During the sterilization phase, the temperature has dropped more than 2°C from the setpoint value for more than 30 seconds. Solution: Check the placement of the product probe inside the bottle and there has not been a power failure.	(Temperature drop)
37	OVERTEMPERATURE STERILIZATION Cause: During sterilization, the temperature has exceeded 4°C the setpoint temperature for more than 5 seconds. Solution: Review the recipe schedule	(Excess temperature)
44	SHORT CIRCUIT PROBE PRODUCT Solution: Check product probe and its connections	
45	PRODUCT PROBE COMMUNICATION ERROR Solution: Check product probe and its connections	
46	SHORT CIRCUIT EXTERNAL PROBE Solution: Check Probe of the external resistance and its connections	
47	EXTERNAL PROBE CUT CABLE Solution: Check the external resistance probe and its connections	
48	SHORT CIRCUIT WATER PROBE Solution: Check water probe and its connections	
49	WATER PROBE CUT CABLE Solution: Check water probe and its connections	
50	SHORT CIRCUIT STEAM GENERATOR PROBE Solution: Check Probe and Its Connections	
51	STEAM GENERATOR PROBE CUT CABLE Solution: Check Probe and Its Connections	
52	STEAM GENERATOR PROBE - OVERTEMPERATURE Solution: Check the probe's operation and replace if appropriate	

Cod.	Name	Description
53	PRESSURE SWITCH PSH010 Cause: Pressure switch has reached the setting pressure, and therefore the pressure is close to the opening of the safety valve (operating limit) Solution: Review Operating Conditions and Component Status	
54	PRESSURE SWITCH PSH024 Cause: Pressure switch of the steam generator has reached the setting pressure (operating limit) Solution: Review Operating Conditions and Component Status	

12. Spare parts

12.1 50/80L Model

		
40039 Autoclave Wheel	15645 Air microfilter	43121 Pressure Transmitter
		
32051 Steam generator	16138 Pressure Switch 2,3bar	43172 Water temperature probe
		
16173 Pressure Switch Generat	43173 Probe Product L=145 mm Cable ø 6mm	43174 Probe Product L=200 mm Cable ø 6mm
		
15712 Product Probe gland cable ø6mm	21271 Lid Gasket 80L 21193 Lid Gasket 50L	25106 Gas Strut 80L 25094 Gas Strut 50L

		
<p>20134 Water level sensor generator</p>	<p>20096 General Switch</p>	<p>15409 Fuse 32A</p>
		
<p>13019 Contactor</p>	<p>36028 Solid State Relay</p>	<p>10054 Lid closure</p>
		
<p>36030 Relay</p>	<p>24360 Vacuum pump</p>	<p>34336 Check Valve</p>
		
<p>24359 Diaphragm pump</p>	<p>34400 Safety valve</p>	<p>34271 2-way manual valve</p>
		
<p>28036 Tank cap</p>	<p>46167 Hose G3/4"</p>	<p>7488 32 Amp. Plug</p>

		
<p>16057 Solenoid Valve NO ø3-24DC</p>	<p>16067 Solenoid Valve Water Inlet</p>	<p>16070 Solenoid Valve 1/4" 24VDC</p>
		
<p>16079 Solenoid Valve 3 ways</p>	<p>16170 Solenoid Valve. 3 ways G1/4"</p>	<p>16213 Solenoid Valve NO 122K8363</p>
		
<p>16227 Solenoid Valve 139V08</p>	<p>85289 Temperature Transmitter</p>	<p>16143 Power supply</p>
		
<p>36034 Tank Level sensor vert</p>	<p>36035 Tank Level sensor horiz</p>	<p>39286 Autoclave ext. heater</p>
		
<p>43071 Probe PT-100 ext. heater</p>	<p>85279 Autoclave HMI Display</p>	

12.2 150L Model

		
<p>40026 Autoclave wheel</p>	<p>15645 Air microfilter</p>	<p>43121 Pressure Transmitter</p>
		
<p>43168 Thermostat</p>	<p>16138 Pressure Switch 2.3bar</p>	<p>43172 Water te</p>
		
<p>16173 Steam gen. pres. switch</p>	<p>43173 Product Probe L=145 mm Cable ø 6mm</p>	<p>43174 Product Probe L=200 mm Cable ø 6mm</p>
		
<p>15712 Product Probe gland cable ø6mm</p>	<p>21206 Lid gasket 150L</p>	<p>25080 Gas Strut 150L</p>


		
<p>20134 Steam Gen. level sensor</p>	<p>20096 General Switch</p>	<p>15409 Fuse 32A</p>
		
<p>13019 Contactor</p>	<p>36028 Solid State Relay</p>	<p>10054 Lid closure</p>
		
<p>36030 Relay</p>	<p>24362 Vacuum pump</p>	<p>28036 Tank cap</p>
		
<p>24364 Diaphragm pump</p>	<p>34400 Safety valve</p>	<p>34271 2-way manual valve</p>
		
<p>46167 Hose G3/4"</p>	<p>46167 Hose G3/4"</p>	<p>46167 Hose G3/4"</p>

		
<p>16057 Sol. Valve NO ø3-24VDC</p>	<p>16067 Solenoid Valve water Inlet</p>	<p>16148 Solenoid Valve 2w-160-NC</p>
		
<p>16079 3-way solenoid valve</p>	<p>16149 Solenoid Valve 2w-15bk NO</p>	<p>16150 Solenoid Valve 2w-15B NC</p>
		
<p>36060 Three-phase solid-state relay</p>	<p>85289 Temperature Transmitter</p>	<p>16143 Power supply</p>
		
<p>36034 Tank water level vertical</p>	<p>36035 Tank water level horizontal</p>	<p>85279 Autoclave HMI Display</p>

Notice to customers about recycling:

The product consists of various components and various materials that must be recycled or, failing that, deposited at the appropriate debris disposal sites when the life of the product has been completed or when it is otherwise necessary to dispose of it. To do this, the fine user who acquires the product must know the current regulations of each municipality and/or locality depending on electrical and electronic waste. The user purchasing this product must be aware of and responsible for the possible effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen reach if it is a previous dismantling and knowledge of the components incorporated. If you do not know the procedure to follow, please consult with the city council for more information.