SOFTWARE V.WAVE BASIC

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/

**USER MANUAL
(FOR SPECTROPHOTOMETER UV-2005)**

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INDEX

1. Functions	12
1.1 Main functions.....	12
2. Setup	13
2.1 System requirements for M.Wave Basic	13
2.2 Setup M.Wave Basic to PC	13
2.3 RemoveM.Wave Basic from PC	14
2.4 Run M.Wave Basic	14
2.5 Set communication port	14
2.6 Connect/Realease spectrophotometer	14
3. Introduction	15
3.1 Main screen.....	15
3.2 Menur bar and toolbar.....	16
4. Operation	17
4.1 Quantitative Test	17
4.2 Kinetics	20
4.3 Additional Function	22
4.4 File Operation	22

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1. Functions

This section introduces the functions of the M.wave Basic.

1.1 Main Functions

Single wavelength photometric measurement

- 2 methods to set up the regression curve.
 - Up to 20 standards to set up the regression curve. The M.wave Basic will calculate the working curve using a linear equation that fits the data. Enter factor values to generate regression curve.
- 2 methods for curve fit.
 - Linear fit and Quadratic fit.

Kinetics

- Allow user to set scan Interval (0.5, 1.0, 2.0, 5.0, 10.0, 30.0 and 60.0s).
- Spectrum display mode can be changed (Time-100%Trans. or Time-Abs.).

2 Setup

This section introduces how to setup the M.Wave Basic to PC.

2.1 System Requirements for the M.Wave Basic

- Pentium processor-based personal computer
- CD-ROM driver
- USB port
- 8 MB of RAM (16 MB or more recommended)
- 20 MB of available hard disk space
- Microsoft Windows 95, Windows 98/Me, Windows 2000 or Windows XP

2.2 Setup the M.Wave Basic to PC

- Insert the CD-ROM of the M.Wave Basic into CD-ROM driver.
- Open the directory of CD-ROM.
- Double-click the icon **Setup.exe** to start to setup (Fig. 2-1). Click **Next**.
- Input user's information (Fig. 2-2). Click **Next**.
- Select installation folder (Fig. 2-3). Click **Next**.
- Select shortcut folder (Fig. 2-4). Click **Next** and the installation information will show (Fig. 2-5). Click **Next** to starting copy files to computer.
- After finished copy all files, click **Finish** to finish setup(Fig. 2-6).


Fig. 2-1 Welcome


Fig. 2-2 User Information


Fig. 2-3 Installation Folder

Fig. 2-4 Shortcut Folder


Fig. 2-5 Show Installation Information

Fig. 2-6 Finished Installation

2.3 Remove the M.Wave Basic from PC

- Start→All Program→Mapada→Uninstall M.Wave Basic 1.0, Uninstallation information will show (Fig. 2-7). Click **Next** to uninstall.
- After remove all files from PC (Fig. 2-8), click **Finish** to finish uninstall.


Fig. 2-7 Uninstallation Information


Fig. 2-8 Finished Uninstallation

2.4 Run M.Wave Basic

There are two ways to start the M.Wave Basic.

1. Double-click shortcut icon  on the desktop.
2. Start→All Program→M.Wave Basic1.0.

1) Set Communication Port

Start the M.Wave Basic, on the **View** menu, click **Options** appears the following box (Fig. 2-9), select the Port type (Serial Port), Select the Comm. Port No. and Baud Rate (19200bps), click **OK**. Tick the **Connect to Spectrophotometer when start** will automatically connect to the Spectrophotometer when starting M.Wave Basic next time.


Fig. 2-9 Communication Setup

2) Connect/Release the Spectrophotometer

Click  on the toolbar to connect/release the Spectrophotometer.

3. Introduction

This chapter introduces the interface of the M.Wave Basic.

3.1 Main Interface

After running the M.Wave Basic, the Main Form appears (Fig. 3-1).


Fig. 3-1

3.2 Menu Bar and Toolbar

Menu bar and **Toolbar** are both provided in the software offering you two ways to select a desired function.

- On the menu bar, use your keypad or mouse to select the desired function.
- Almost all the functions listed in the menu bar can be reached by clicking a corresponding button in the toolbar.

Main Menu	Sub Menu-1	Sub Menu-2	Tool	Function
File	New	Quantative Analysis		New a Quantative Analysis Measurement
		Kinetics		New a Kinetics Measurement
	Open...			Open a Measurement File
	Close			Close Current Measurement
	Save			Save Current Measurement
	Print...			Print Current Measurement
	Exit			Exit M.Wave
View	ToolsBar	Standard		Show/Hide Standard ToolsBar
		Spectrum		Show/Hide Spectrum ToolsBar
	StatusBar			Show/Hide Statusbar
	Option			Setup System Parameter
UV-Spectro-photometer	Connect/Release			Connect/Release Instrument
	Turn On/Off W Lamp			Turn On/Off W Lamp
	Turn On/Off D2 Lamp			Turn On/Off D2 Lamp
	Go to wavelength...			Set and Goto a result Wavelength
Operation	Start			Start operate
	Cancel			Cancel Current Operate
	Modify			Modify Current Result
	Delete			Delete Current Result
	Zero/Full			Zero and Full
	Setting...			Setup Measurement Parameter
Spectrum	Display Mode	Abs.		Display as Abs. Mode
		Trans.		Display as Trans. Mode
Window	Close			Close Current Measurement
	Cascade			Multi Windows Display in a Cascade
	Tile Vertically			Multi Windows Display in a Tile Vertically
	Tile Horizotally			Multi Windows Display in a Tile Horizotally
Help	About...			Information about M.Wave

Tab. 3-1

4. Operation

This chapter introduces the operation of the M.Wave Basic.

4.1 Quantitative Test

M.Wave Professional use Standard Curve to test the sample's concentration at a fixed wavelength point.

4.1.1 Measure

1. Click the icon  in the Tools Bar to create a new Quantitative Test. (Fig. 4-1)


Fig. 4-1 the Main Interface of Quantitative Test

2. Click the icon  in the Tools Bar to set the parameters of Quantitative Test. (Fig. 4-2)


Fig. 4-2 Parameters Setting Interface of Quantitative Test

3. Select Test Method in **Area A** (Fig.4-2)Pls. Ref. **Appendix 1**4. Input the wavelength you want in **Area B** (Fig.4-2).5. Choose the Concentration Unit in **Area C** (Fig.4-2).6. Choose curve fitting method in **Area D** (Fig.4-2)

7. Create Standard Curve.

Two methods are under your choice to create a new curve.


1) Method 1: Coefficient1) Click **Coefficient** option in **Area E**. (Fig.4-2)2) Input the curve equation's coefficients in corresponding frames in **Area E**. (Fig.4-2).3) Click **OK** to finish setting.**2) Method 2: Use Standard Samples to create a new curve**1) Click **Standard Sampling** option in **Area E** (Fig.4-2).2) Choose the quantity (≤ 20) of the standard samples in **Area E**. (Fig.4-2)3) Input the standard samples' concentration in corresponding frame in **Area F**. (Fig.4-2)4) Click **OK** to finish setting.5) Pull the Reference into the light path and click the icon **Z1**, then the system will move the wavelength to the one you just set, it blanks then.6) Pull 1# standard sample into the light path and move the cursor on the first frame of **Abs.** in **Area C**. (Fig.4-3) The sample's Absorbance value will appear in this frame after double click.7) Measure other standard samples in the same way as **step (6)** shows. The curve displays automatically in **Area B** (Fig.4-3) when finish.

Fig. 4-3 Create Standard Curve and Measure

Note: You can modify the curve's coordinates by clicking the icon  (Fig. 4-4) .

You can also save the curve by clicking  , and you can reload it by clicking  next time.


Fig. 4-4 Rescale the Curve's Coordinates

- 8) Pull the reference into the light path, click the icon **Zt** again to set Zero and Full scale. Pull the unknown concentration sample into the light path, click the icon  in the tools bar to test. And the test result will display in the data sheet in **Area A**. (Fig.4-3)
- 9) Finish all the test following the way as **step 8** does.

4.1.2 Modify a Test Result

- 1) Move the cursor on the data you want to amend and click.
- 2) Pull the sample into the light path.
- 3) Click the icon  in the Tools Bar to re-measure the sample, and the new value will replace the old one.

4.1.3 Delete a Test Result

- 1) Move the cursor on the data you want to delete and click.
- 2) Click the button  in the Tools Menu to delete the data.

4.1.4 Name a Sample

- 1) Move the cursor on the sample line you want to name.
- 2) Input the sample name after double click the frame, then press **Enter** to confirm.

4.2 Kinetics

We'll introduce you how to measure a sample's change of Absorbance or Transmittance over a selected period of time.

4.2.1 Measure a sample

- 1) Click the icon  to begin a new test of Kinetics. (Fig.4-5) .


Fig.4-5 Main Interface of Kinetics

- 2) Click the button of  to set the parameters (Fig.4-6)


Fig.4-6 Set Parameters of Kinetics

- 3) Choose Photometric **Display Mode**, input the **wavelength**, **Delay Time**, **Total Time**, **Intervals** and coordinates values
- 4) Click **OK** to finish and cancel setting.
- 5) Pull the reference into the light path and click icon  in the Tools Bar to set zero and full scale.
- 6) Pull the sample into the light path and click the button of  to begin the test. (Fig. 4-7), you can click icon  in the Tools Bar to cancel at any time.


Fig. 4-7 Kinetics

4.2.2 Change Display Mode

Click the icon **A** or **T** in the Tools Bar to choose the photometric mode you want. (Abs.-Time, Tran.-Time)

4.2.3 Magnify the Selected Zoon.

Click the icon  in the Tools Bar, choose the area you want to magnify, then you'll see the magnified one, and the Mouse's cursor turns into the form of "+".
Click the icon  again to cancel the current status.

4.2.4 Rescale the coordinates

Click the icon  in the Tools Bar, input the new value you want in the new frame (Fig. 4-8) and click **OK** to confirm. Then the curve with new coordinates will display on the screen.


Fig.4-8 Rescale the Coordinate

4.2.5 Switch the Display Status

Click the icon  or  to display the curve's upper or latter status, click  to recover the original status.

4.2.6 Search Peaks

Click the icon  in the Tools Bar to search the peaks, they will be labeled automatically on the curve, and the peaks' height will be listed in the Right Frame. (Fig.4-7)

Click the icon  in the Tools Bar. You should input the peaks' height that you want to find in the new frame and click **OK**. (Fig.4-9) Then the peaks whose height is lower than the inputted value will be automatically labeled.


Fig.4-9 Set peak's height

4.3 Additional Function

This chapter will introduce some other function of M.Wave Basic.

4.3.1 Set and go to a wavelength

Click  on the toolbar to set a wavelength (Fig. 4-10). Key in the wavelength value in the **Wavelength** box, then click button **Goto...** to goto the result wavelength, after do Zero/Full automatically the Abs. value and Trans. will display in real time. Click button **zero/Full** to get bright current and dark current. If you tick the **Resampling Dark Current** check box, it will resample the dark current when do Zero/full.


Fig. 4-10

4.3.2 Turn on/off W lamp

Click  on the toolbar to turn on/off W lamp.

4.3.3 Turn on/off D2 lamp

Click  on the toolbar to turn on/off D2 lamp.

4.4 File Operation

4.4.1 Save a File

Click , a **Save** dialog box will be displayed. Type in a file name, click **Save**.

4.4.2 Open a File

Click , a **Open** dialog box will be displayed.. Select a folder and filename. Click **Open** to open the selected file.

4.4.3 Print a Measurement

Click , a **Print** dialog box will be displayed. Setup the printer parameter. Click **Print** to print the measurement.