V-1100

4120025

**Espectrofotómetro de rango visible
Spectrophotometer visible range**

Indice

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Index

1 INTRODUCTION	16
1.1 Package contents	16
1.2 Appearance	16
1.3 Operating Panel	17
1.4 Descripción del Teclado	18
2 SAFETY	18
3 OPERATING PRINCIPLE	19
4 INSTALLATION	19
5 OPERATION	20
5.1 Local Software Structure	20
5.2 Basic Operation	20
6 MEASUREMENT	21
6.1 Measure the Absorbance	21
6.2 Measure the Transmittance	22
6.3 Two-point Method	23
6.4 Coefficient Method	23
7 SPECIFICATIONS	25
8 WARRANTY	25
9 TROUBLESHOOTING	26
10 CE CONFORMITY DECLARATION / DECLARACIÓN DE CONFORMIDAD	27

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ENGLISH VERSION

Contact with dangerous substances

Electrical accident

1 INTRODUCTION

Spectrophotometers V-1100 described in this manual are designed to be used by properly trained personnel in a suitable equipped laboratory. For the correct and safety use of this equipment, it is essential that laboratory personnel follow generally accepted safe procedures in addition to the safety precautions called for in this manual.

The covers on this instrument may be removed for servicing. However, the inside of the power supply unit is a hazardous area and its cover should not be removed under any circumstances. There are no serviceable components inside this power supply unit.

1.1 Package contents

Carefully unpack the contents and check the materials against the following packing list to ensure that you have received everything in good condition.

Packing List

Description	Quantity
Spectrophotometer	1
Mains Lead	1
Glass cuvettes	4
Manual	1


1.2 Appearance

Please see Fig. 3-1:

- 1 Lid of the compartment
- 2 Cell Holder
- 3 Pole
- 4 Operating Panel
- 5 Wavelength Knob
- 6 USB Port
- 7 Print Port
- 8 Power Socket 85-265VAC
- 9 Power Switch

Fig. 3-1

1.3 Operating Panel

Please see the panel of V-1100 (Fig. 3-2):

Fig. 3-2

1 LCD (128 X 64)

2 Key Button

Switch the Photometric Mode

Confirm/Print

Decrease Number/Set Zero

Increase Number/Set 100%T

2 SAFETY

Spectrophotometer VR-1100 is an electrical test instrument widely used in laboratories.

This equipment has been designed to be used mainly in laboratories.

It is important that the laboratory staff know the safety general procedures, together with the safety precautions explained in this manual.

Some of the chemical products used in this spectrophotometer are corrosive and/or inflammable and the samples can be radioactive, toxic or potentially infectious.

Be careful to follow the safety standardized procedures in laboratories for samples and chemical products manipulation.

DANGER:
contact with
dangerous
substances

DANGER:
Electrical accident

Electrical connection:

Power requirements of V-1100 are from 85V to 265VAC. Make sure that the local power supply is within this range.

The power cord shall be inserted in a socket provided with a protective earth contact.

Warning:

Any interruption of the protective conductor inside or outside the equipment is likely to make it dangerous.


NEVER touch or handle the power supply on V-1100 due to the high voltage in the inside, even when it is disconnected from the mains.

3 OPERATING PRINCIPLE

V-1100 equipments are composed of five parts:

- An halogen lamp to supply the light.
- A monochromator to isolate the wavelength of interest and to eliminate unwanted second order radiations.
- A compartment for sample cuvettes.
- A photodetector to receive light transmittance and to convert it to electrical sign.
- A digital display where absorbance or transmittance.

Block diagram for the spectrophotometer.

4 INSTALLATION

After carefully unpacking the contents, check the materials with the packing list to ensure that you have received everything in good condition.

Keep the packaging for some days. If there is something missing, please contact your distributor as soon as possible.

Place the instrument in a stable level surface, suitable for the equipment weight and away from direct sunlight.

Set the unit up in an area that is free of dust, corrosive gases and strong vibrations.


Leave a free space of 30cm around the equipment for a good ventilation.

Turn on the spectrophotometer V-1100, it will begin to self test.

After a 20 minutes' pre-warm, you could take any readings.

5 OPERATION

5.1 Local Software Structure

5.2 Basic Operation

1) Select Test Mode

Press  to select the measure mode.

2) Set Wavelength

Turn the Wavelength Setting Knob to select the wavelength you want, and the wavelength value can be displayed on the screen in real-time.

3) Input the Coefficient

When the system prompts you to input C, K or B, press  or  button to change the value till it displays the one you want, then press  to confirm.

Note: The system will always memorize your last input, until you input another value.

4) Calibrate Zero

Close the lid of the compartment and press  to calibrate Zero T%.

5) Calibrate 100%T

Pull the Reference in the light path, press  to calibrate 100%T0.

6) Print the result

Press the  button to print the test result.

6 MEASUREMENT

6.1 Measure the Absorbance


Fig. 4-1

- 1) Press  to choose «A» mode (Fig.4-1).
- 2) Turn the wavelength knob till it displays the wavelength value you want.


Fig. 4-2

- 3) Pull the reference in the light path, press  to set 100%T (Fig.4-2).
- 4) Pull the unknown concentration sample in the light path and then record the displayed value (Fig. 4-3).

Fig. 4-3

- 5) Press  to print the test result.
- 6) Repeat steps 4) and 5) to test other unknown concentration samples.

6.2 Measure the Transmittance


Fig. 4-4

- 1) Press  to choose the «T» mode (Fig. 4-4).
- 2) Turn the wavelength knob to set the wavelength at the point you want.


Fig. 4-5

- 3) Pull the Reference in the light path to calibrate 100%T (Fig. 4-5).
- 4) Pull the unknown concentration sample in the light path, then the value displayed on the screen is the one you need (Fig. 4-60).


Fig. 4-6

- 5) Press  to print the test result.
- 6) Repeat steps 4) and 5) to test other samples.


Fig. 4-7






Fig. 4-8

6.3 Two-points Method

This method allows measuring a concentration from two pattern concentrations.

- 1) Choose «A» or «T» mode and the wavelength.
- 2) Insert the «blank» cuvette, press  to get 100%T.
- 3) Insert the «pattern» cuvette and press  to choose «C» mode.
- 4) Press  or  to input the standard concentration, and press  to confirm. (Fig. 4-7).
- 5) Repeat to insert the second pattern.
- 6) Set the unknown concentration sample in the cell holder, then its concentration will display on the screen. (Fig.4-8).
- 7) Press  to print the test result.
- 8) Repeat step 6) and 7) to test other samples.

6.4 Coefficient Method

- 1) Set the wavelength.
- 2) Press  to choose «F» Mode.
- 3) Press  or  to set the coefficient value of «K» and «B», followed with  pressed to confirm. (Fig. 4-9). From now on, you can test the samples.
- 4) Press  to set 100%T transmittance or zero absorbance (Fig. 4-10).
- 5) Set the samples in the cell holder. Its concentration will be displayed on the screen automatically (Fig.4-11).
- 6) Press  to print the test result.

7 SPECIFICATIONS

Model	V-1100D
Wavelength Range	325-1000nm
Spectral Bandwidth	4nm
Optical System	Single Beam, Grating 1200lines/mm
Wavelength Accuracy	±2nm
Wavelength Repeatability	1nm
Wavelength Resolution	0.5nm
Photometric Accuracy	±0.5%T
Photometric Repeatability	±0.3%T
Photometric Range	-0.3-3A, 0-200%T
Stray Light	0.3%T@360nm
Stability	±0.004A/h @500nm
Display	Graphic LCD (128X64 dpi)
Keyboard	4 Button keypad
Photometric Mode	T, A, C, F
Detector	Photodiode
Sample Compartment	Standard 10mm path length cuvette holder (with optional holder)
Light Source	Tungsten lamp
Output	USB Port Parallel Port (printer)
Power Requirement	85V ~ 265VAC
Dimensions (W x D x H)	480 x 360 x 160mm
Weight	10 kg

8 WARRANTY

This product is guaranteed for one year from the date of purchase and it covers defective components and assembly.

Please keep the purchase ticket as it is the only valid document for the purchase date.

This guarantee does not cover damage caused by an incorrect use of the equipment. This equipment must be properly validated before its use in the value sample analysis for the user.

Any return without the authorization of J.P. SELECTA, s.a.u. is not accepted.

9 TROUBLESHOOTING

PROBLEM	Possible Causes	Solution
No functioning after power is on	-Power cord not connected to outlet	-Plug instrument in
	-Dead Power outlet	-Change to a different outlet
	-Internal fuse melted or defective electronic component	-Change fuse or call an authorized service engineer.
Instrument cannot set 100%T (0A)	-Light beam blocked: Holder misaligned	-Check sample holder
	-Lamp is old or defective	-Replace lamp
	-Lamp is off alignment	-Readjust the lamp
	-Defective electronic component	-Call an authorized service engineer
Unstable display	-Insufficient warm up time	-Warm up at least 20 minutes
	-Poor grounding	-Check ground connection
	-Excessive vibration, strong air current near light source, or strong external light.	-Improve working conditions
	-Lamp old or defective	-Replace with a new lamp
	-Unstable power supply	-Equip a manostat to make the voltage stable
	-Defective or dirty detector or defective electronic component	-Call an authorized service engineer
Incorrect readings obtained	-Insufficient sample volume	-Fill cuvette with more samples
	-Wrong wavelength setting Failed to blank (0A/100%T)	-Check analytical procedure and wavelength setting.
	-Stray sample preparation vapours	-Prepare sample away from instrument. Use proper ventilation
	-Bubbles or particles in solution	-Check sample preparation and analytical procedure
	-Instrument out of electronic calibration	-Call an authorized service engineer