VR 2000

4120026

**ESPECTROFOTÓMETRO
SPECTROPHOTOMETER**

Indice

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Index

1 INTRODUCTION	18
1.1 Package contents	18
1.2 Appearance	18
2 SAFETY	19
2.1 Risks to which the operator is exposed	20
2.2 Staff qualification	20
2.3 Equipment use	20
2.4 Modifications	20
3 OPERATING PRINCIPLE	21
3.1 Functions	21
4 INSTALLATION	22
4.1 Location	22
4.2 Power supply connection	22
5 OPERATION	23
5.1 Keyboard description	23
5.2 Starting up	23
5.3 Basic Operation	23
6 ANALYZE SAMPLE	26
6.1 Photometry	26
6.2 Quantitative	26
7 MAINTENANCE	31
7.1 Lamps replacing	31
7.2 Battery replacement	32
8 WARRANTY	32
9 TECHNICAL SPECIFICATIONS	32

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
1 INTRODUCTION

VR-2000 spectrophotometer has a wide range wavelength, high sensitivity, powerful function, ease of use, and a simple structure. Besides this, its large LCD, high precise A/D and easy to store RAM makes the instrument much more superior than other originals. It is widely used in chemistry, pharmaceuticals, biochemical, metallurgy, light industry, textile and medical and educational environments. It is one of the most important instruments in quality control and it is essential in all laboratories.

1.1 Package contents

- Spectrophotometer
- Power connection cable
- Instructions manual
- Glass cuvettes (4)

1.2 Appearance

2 SAFETY

Spectrophotometer VR-2000 is an electrical test instrument widely used in laboratories.

This equipment has been designed to be used by technicians who have been qualified for its use. For the correct and safety use of this equipment it is important that the laboratory staff know the safety general procedures, together with the safety precautions explained in this manual.

Some of the chemical products used in this spectrophotometer are corrosive and/or inflammable and the samples can be radioactive, toxic or potentially infectious. Be careful to follow the safety standardized procedures in laboratories for samples and chemical products manipulation.

The risky situations are indicated along this manual by means of the safety icons and together with some instructions that must be respected.

2.1 Safety icons

These symbols identify risky situations and the safety measures to be taken.

The icons make reference to the paragraph marked in grey.

Danger risk

Risk of danger.

Please, respect the instructions indicated to do the described operation.

Electric risk

Risk of electric accident when accessing the areas labelled with this signal or when doing the operations indicated with this icon in the present manual.

Please, respect the instructions indicated to do the described operation.

Burn risk due to contact with high temperature areas

The temperature in the area indicated by this icon is very hot. Please, use anti thermally gloves to do the described operation.

Please, respect the instructions indicated to do the described operation.

Gloves use compulsory

It is compulsory the use of gloves to do the described operation.

A recommended model is indicated in the «Accessories» section.

Eyes protection use compulsory

It is compulsory the use of safety glasses or complete facial protection to do the described operation.

A recommended model is indicated in the «Accessories» section.

Important information

- Important information to obtain good results or for an excellent operation of the equipment.
- Important information to lengthen the equipment life or to avoid the degradation of any of its components.

2.1 Risks to which the operator is exposed

The working process with the described equipment exposes the operator to the following risky situations:

- Corrosive and / or inflammable substances manipulation.
- Sample manipulation which can be radioactive, toxic or potentially infectious.
- Glass pieces manipulation.
- Electric risk.

2.2 Staff qualification

This equipment can only be used by staff qualified in knowing dangers exposed in chemical analysis laboratories.

This equipment can only be used by conveniently trained staff who have read the instruction manual or have been qualified in the use of this equipment.

2.3 Equipment use

This equipment is expected to be used in clinical, pharmaceutical and biochemical laboratories.

The use of this equipment in a way not indicated in this manual can commit the protection guaranteed by the equipment.

This equipment is NOT prepared to work in explosive atmospheres.


2.4 Modifications

The modification of the equipment safety systems operation or manipulation by personnel not authorized by the manufacturer can expose the operator to risks not foreseen in this manual.

3 OPERATING PRINCIPLE

VR-2000 equipment is composed of five parts:

- An halogen lamp to supply the light.
- A monochromator to isolate the wavelength of interest and to eliminate unwanted second order radiations.
- A compartment for sample cuvettes.
- A photodetector to receive light transmittance and to convert it to electrical sign.
- A digital display where absorbance, transmittance and concentration values are shown.

The halogen lamp light, focus on the monochromator entrance where a collimated mirror points the light beam towards the diffraction grating.

The diffraction grating diffuses the light beam to produce an spectrum. A part of it centers in the monochromator output slot by means of another collimated mirror.

The wavelength light beam selected passes through a filter which helps eliminating second order radiations from the diffraction grating.

This light beam goes through the sample, and then the photodetector produces an electrical sign used to make the transmittance or absorbance calculations, shown in the digital display.

3.1 Functions

VR-2000 has 3 functions:

- Photometry:



Test mode can be switched between ABS, %T and Energy. Test results can be stored (up to 200 groups) and print out.

- Quantitative:

Two methods are included in the menu: Coefficient Method and Standard Curve Method. The regression equation (up to 20 groups) and the test results (up to 200 groups) can be stored in the RAM and printed out.

- Utility:

User can set their favourite items to keep the equipment work in the best conditions. You can calibrate Dark Current, Calibrate Wavelength, Set Light Source Switch Point and Load Default Setting.

4 INSTALLATION

Once the equipment is unpacked, check if it is all right. Preserve the packing for some days. Any damage or anything missing, please contact your dealer as soon as possible.

4.1 Location

Place the equipment over a flat, levelled and stable surface, suitable for the equipment weight and not exposed to direct sunlight.

With the aim of having the best equipment performance, keep it isolated, as much as possible, from any magnetic or electrical field source or any other electrical device that could generate high frequency fields.

Install the unit in an area free of dust, corrosive gases and hard vibrations.


Place the equipment at a distance lower than a 1.5m from a power socket.

The power required is 220±22V@50±1Hz or 110±11V@60±1Hz.

To work in a comfortable way, left a free space of 50cm around the equipment and a minimum of 20cm behind.

The best working temperature range is 16-35°C and the humidity is 45-80%.

Remove obstructions or materials that could obstruct the air flow around and underneath the equipment.

4.2 Power supply connection

Please, select a socket near the equipment and suitable for its power.


Make sure the printer's power is off. Link the printer's data cable to the equipment's parallel port.

Make sure the machine's power switch is in the OFF condition, link the power cord to the equipment and insert the other end in the socket provided with a protective earth contact. Use only the supplied cord or similar ones to connect the equipment.

For your safety, the equipment must always be earthed.

Before switching on, make sure that the voltage and the power frequency corresponds to the indicated in the equipment, if not a severe damage may occur.

Switch ON the power. The equipment will begin to self-test. After a self-test and a 20 minutes' pre-warm, the machine will be ready to work. To skip the 20 minutes pre-warm, press any key.

Note: DO NOT OPEN SAMPLE COMPARTMENT LID DURING SELF-TESTING.

Attention:

After transportation or a period of use, the working condition of the instrument may change because of environment changing or any other cause, such as the attenuation of power light, vibrations, etc. Users need to re-calibrate the system to keep it in its best working condition.

5 OPERATION

The control panel allows the equipment functions selection and informs about the current operating status.

Keyboard

5.1 Keyboard description

[SET]	To program parameters
[GOTO]	To establish the wavelength
[ZERO]	To make blank, 100%T/0Abs
[PRINT]	To print out results
[-]	«-» symbol that functions as indicated in the display
[>] [<]	Cursors for selecting and scroll

5.2 Starting up

Switch the spectrophotometer on by pressing the main switch. The equipment will start.

First of all, do the self-testing.

Note:

Do not open the sample compartment during the self-testing.

Then, the equipment will enter the waiting mode for 20 minutes to achieve the working pattern.

To leave the waiting mode, press any key.

Note:

If you do not respect the waiting time, it is possible that a little drift will be shown in the equipment readings.

5.3 Basic Operation

The display screen is user-friendly, always you can complete all the operation according to the clues on the screen.

5.3.1 Set Wavelength

Press **[GOTO]** into wavelength setting screen. Then press **[>] [<]** to select the wavelength you want, and press **OK** to confirm. After the wavelength has been changed, it returns to the main menu automatically (Fig. 2-1).


Fig. 2-1


Fig.2-2 Blank

5.3.2 Blank

Pull the reference cuvette into the light path, press **[ZERO]** to blank (Fig.2-2).

5.3.3 Utility

Press **[SET]** to go into system functions screen, use **[>]** **[<]** to select, then press **OK** to go into the corresponding item (Fig. 2-3).


Fig.2-3 go into utility setting interface

5.3.3.1 Test Mode

Move the cursor on to «Test Mode», and press **OK**, then it will go into the mode selection screen (Fig. 2-4).


Fig. 2-4 select test mode

5.3.3.2 D2 Lamp Management

Move the cursor on «D2 Lamp On», press **OK**, it will go into another screen.

Choose the D2 lamp condition you want, press **OK** to confirm, then press **Cancel** to return to upper screen.


Fig. 2-5 Deuterium lamp management


Fig. 2-6 | W lamp manage

5.3.3.3 W Lamp Management

Move the cursor on «**W Lamp On**», press **ok**, it will go into another screen.

Choose the W lamp condition you want, press **ok** to confirm, then press **Cancel** to return to upper screen.


Fig. 2-7 Get dark current

5.3.3.4 Dark Current

Move the cursor on «**Dark Current**», followed with **ok** pressed to get dark current. Once the environment changed, such as temperature, instrument's station, etc, you should go to this operation.



Fig. 2-8 Calibrate Wavelength

5.3.3.5 Reset WL

Move the cursor on «**Reset WL.**», then press **OK**, system begins to reset wavelength.

5.3.3.6 See the version

Move the cursor on «**Versión**», then press **OK**, you'll find the version of this instrument's hardware and software. Press any key to return.



Fig 3-1 go into Basic test Interface



Fig 3-2 Basic test



Fig 3-3 print or clear



Fig 3-4 go into quantitative Interface



Fig 3-5 go into method choice Interface

6 ANALYZE SAMPLE

Before measurement, make a blank reference solution by filling a clean cuvette half full with distilled or de-ionized water or other specified solvent. Wipe the cuvette with tissue to remove fingerprints and droplets of liquid. Fit the blank cuvette into the 4-cell linear changer and place the cuvette into the nearest slot. Push the changer so that the cuvette is in the light path (push the rod in). Close the lid.

6.1 Photometry

Firstly, set the wavelength you need.

6.1.1 Basic mode

Pull the blank cuvette into the light path, then press **Basic**. After automatic blanking, it will go into Basic test interface.

6.1.2 Begin to test

Pull the sample into the light path, then press **Test** to measure. The test result is displayed on the screen. And it will be numbered automatically.

Repeat this step to measure other samples, the results will be numbered orderly. It can be saved to a total of 200 groups of test data in the RAM, 3 groups per screen. Use [>] [<] keys to scroll.

Press **Cancel** to return to the main menu.

6.1.3 Clear or Print out the data

Press the [PRINT] key as shown.

Move the cursor on the item you want to select and then press **Ok** to confirm.

Item 1 means all the data will be deleted after printing.

Item 2 means all the data will be deleted from screen and RAM without printing.

Item 3 means returning to upper interface without any choice, also you can press **Cancel** to return.

Note: You can also **Set wavelength** and **Blank** by pressing [GOTO] and [ZERO] individually in the Basic Test interface.

6.2 Quantitative

Press **Quantitative** to go into quantitative test main menu. Three items are under your choice. You can create, load and delete curves.

6.2.1 Create curve

Use [>] [<] to move the cursor on the first item, then press **OK** to confirm. You have two methods to choose.

If you know a curve's equation, you can choose «Coefficient» method. Input all the coefficients to complete the equation, then to measure unknown samples.

If you don't know the equation, you'll have to create a standard curve by standard samples. At most 9 samples can be used.


Fig.3-6

6.2.1.1 Coefficient

1) Move the cursor on to «Coefficient», then press **ok** and it will go into the wavelength setting interface automatically.

2) Choose the wavelength you need by pressing [**>>**] [**<<**] and then press **ok**. It will go into the coefficient setting interface.


Fig.3-7

3) Input the K coefficient

There are five zeros on the screen, the cursor is on the first one at the beginning. Press [**>>**] [**<<**], the first zero will vary between «0~9». Choose the number or symbol you need, then press **ok** to confirm, the cursor then moves to the next zero automatically. The second to the fifth bits vary between «0~9» and decimal point (it's between 9 and 8). Complete other bits of K Coefficient in the same way.


Fig.3-8

4) Input the B coefficient

When the last bit of K value is confirmed and **Ok** pressed, the interface goes into the B Coefficient setting interface.

The first bit is a symbol bit, and varies from «+» to «-». The other bits are number bits. The way to Input B is the same as you have done with K.

After the last bit of B is confirmed, the curve will be displayed on the screen.


Fig.3-9

5) Begin to test

When the curve displays on the screen, you can use it to test the unknown samples' concentration. Just Press **Test** to do it.


Fig.3-10

6) Print

Press **[PRINT]** in the concentration measure interface, and the test results will be printed.


Fig.3-11

6.2.1.2 Standard Curve

1) Insert Blank Reference

Move the cursor on to «Standard Curve», then press **ok**. The system will ask you to insert the blank reference.

- a. Pull the blank cuvette in the light path, then close the compartment lid.
- b. Press **[GOTO]** to set wavelength.
- c. Press **Ok** to blank.

2) Input the number of standards

After blanking, the system will ask you to input the number of the standards. Use **[>]** **[<]** to choose the number you need, then press **Ok** to confirm.

Note: At most 9 standards are permitted. Other number is invalid.


Fig.3-12

3) Input the concentration of standards

When you have finished the last step, the system will ask you to input the standards' concentration one by one.

You should input the No.1 standard sample's concentration first as indicated in the screen.


Fig. 3-13

4) See the curve.

When the last one is finished, for example if a total of 6 standards are used, the curve will be drawn on the screen, and the curve equation will also be displayed on the screen. The curve and the curve equation will be saved in the RAM automatically. A total of 200 curves can be saved, and they'll be numbered orderly.


Fig.3-14

5) Test

Insert the unknown concentration sample in the holder, then close the compartment lid and pull the sample in the light path. Press **Test** to measure. The results will be displayed on the screen.


Fig. 3-15

6) Print

You can print the test results and the curve. For printing, just press **[PRINT]**.


Fig.3-16 go into curve choice interface

6.2.2 Load Curve

Use **[>>]** **[<<]** to move the cursor on to the second item, then press **OK** to confirm and you'll go into the curve choice interface.

All the curves you create will be saved in the RAM automatically. You can select the one you want by using **[>>]** **[<<]** and then pressing **Ok**.

Note: A total of 200 curves can be saved, you can find them by pressing **[>>]** **[<<]**.


Fig.3-17



Fig 3-18 Go into curve choice interface

6.2.3 Delete Curve

Use [>>] [<<] to move the cursor on to the third item, then press **OK** to confirm and you'll go into the curve deleting interface.

Move the cursor on the equation you want to remove, then press **Ok**. The curve will be deleted from the screen and from the RAM.

7 MAINTENANCE

7.1 Lamps replacing

Tungsten-halogen lamp

Switch off and unplug the equipment.

Remove the four screws from the spectrophotometer sides.

Remove the **sample selection bar** unscrewing the counter clock.

Remove the equipment cover very carefully and place it to the right.

BE AWARE of not PULLING OUT the cable.


Unscrew the two lamp chamber screws and remove the cover.

Attention

The cover could be very hot, remember to use gloves!

Unplug and remove the lamp from the ceramic base (black connector). Insert the new lamp, pushing it as far as possible.

CAUTION:

- 1) Do not touch the lamp with bare fingers. Use a cloth to handle it.
- 2) There's no polarity difference between both tungsten lamp pins.

Connect the equipment and adjust the lamp in order that the light beam focus on the monochromator inlet slot.

Unplug the equipment, cover the lamp chamber again and install the cover. Pay attention not to pinch the cables during the process.

7.2 Battery replacement

VR-2000 can save many data in the RAM memory, and so a battery is needed in the circuit board.

You will have to replace the battery when it is exhausted.

Switch off and unplug the equipment
(VERY IMPORTANT: high voltage).

Battery

Unscrew the board's 13 screws and remove it.
Remove the old battery and place the new one.
Place again the board and the 13 screws.

8 WARRANTY

This product is guaranteed for one year from the date of purchase and it covers defective components and assembly.

Please keep the purchase ticket as it is the only valid document for the purchase date.

This guarantee does not cover damage caused by an incorrect use of the equipment. This equipment must be properly validated before its use in the value sample analysis for the user.

Any return without the authorization of J.P. SELECTA, s.a. is not accepted.

9 TECHNICAL SPECIFICATIONS

Model	VR 2000
Optical system	one single beam
Wavelength range	325-1100nm
Stray light	0,2%T @ 220nm & 360nm
Band width:	4nm
Photometric range:	0-200%T, -0,3-3.0A
Wavelength accuracy:	±1.0nm
Wavelength repeatability:	≤0.5nm
Photometric accuracy:	±0,5%T
Stability:	0,002A/h @ 500nm
Display:	128x64 dots matrix LCD
Data output:	USB, parallel port
Weight (Net):	12 Kg
Dimensions: Depth x Width x Height	490 x 360 x 210 mm