PH-METRO DIGITAL PH-3004

4120403

Contenido

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Formato de informe:

```
***** Report Title
-----
Measure Time:2021/01/19 12:27:28
Operator: Operator 1 Model: PH-3004 pH Meter Serial Number:
SW Version: Ver 1.00
-----
.....MATCHED INFO
Stored Num:                28
Matched Num:               1
Stored No.:                15
-----
.....CALIB INFO
Calib Operator:           REX Team Calib Time: 2020/05/13 08:30:00
Calib Num:                3
.....CALIB RESULT STD 1: 4.00pH    177.3mV 25.0c
STD 2: 6.86pH            8.0mV 25.0c STD 3: 9.18pH -129.1mV 25.0c
pH Slope 1:              100.00%
pH E0 1:                  0.0mV
pH Slope 2:              100.00%
pH E0 2:                  0.0mV
-----
.....BRIEF INFO
Reading Mode:             Timed Reading
Stable Type:              Medium
Temp Comp Type:           ATC
-----
.....SAMPLE INFO
Sample ID:                Sample 1
-----
.....RESULT
Result:                   7.00pH
Signal Value:             -0.0mV
Temp Value:               25.0c
-----
-----
Signature:
```

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1 Introduction

1.1 Introduction

PH-3004 portable pH meter is a newly designed functional instrument that can be widely used in universities, environmental protection, medicine, food, sanitation, geological prospecting, metallurgy, ocean exploration and other fields measurements for acid rain detection, industrial wastewater, surface water, drinking water, beverages, daily chemical products, textiles and other related industries.

General Features

- LCD display screen, 3.5 inches.
- Multi-reading feature allows auto-read, timed-read and continuous-read.
- Automatic/Manual temperature compensation ensures accurate results.
- Auto-hold feature senses and locks the measurement endpoint.
- Data capacity of up to 500 sets for each parameter (GLP-compliant).
- Support for USB communication.
- Auto-power off feature effectively extends the battery service life.
- Reset feature automatically resumes all settings back to factory default options.
- IP65 waterproof.
- 1-5 points calibration with standard recognition.
- Selectable pH buffer groups, including NIST, DIN, GB.
- Automatic electrode diagnosis with pH slope and offset display.

1.2 Technical Specification

Table 1-1 Instrument Specifications

Model		PH-3004
pH level		0.01 pH
mV	Range	(-2000.0~2000.0)mV
	Minimum resolution	0.1mV
	Electronic unit indication error	±0.1% or ±0.3mV
	Electronic unit repeatability	1mV
	Electronic unit input current	≤1×10 ⁻¹² A
	Electronic unit input impedance	≥1×10 ¹² Ω
pH	Range	(-2.00~20.00)pH
	Minimum resolution	0.01pH
	Electronic unit indication error	±0.01pH
	Electronic unit repeatability	0.005pH
	Instrument indication error	±0.02pH
	Instrument repeatability	0.01pH
Temperature	Range	(-5.0~110.0)°C
	Minimum resolution	0.1°C
	Electronic unit indication error	±0.2°C
	Instrument indication error	±0.4°C(0.0°C-60.0°C) ±1.0°C(Other range)
Use environment		Environment temperature: (0~40)°C Relative humidity: not more than 85%
Dimensions(l×b×h), Weight(kg)		80mm ×225mm×35mm, About 0.4kg
Power supply		Rechargeable lithium battery, power adapter (Input AC 100~240V; Output DC 5V)

● Equipped

1.3 Function Introduction

Table 1-2 Functions Specification

Features		Explanation
Basic Function	Backlight adjustment	●
	Automatic diagnostics	●
	Factory reset	●
	Default parameter	●
	Prompt Sound	●
	Time settings	●
	Power failure protection	●
	Firmware upgrade	●
	Anti-interference automatic recovery	●
	Automatic shutdown	●
	Protection	IP65
Reading Function	Reading balance settings	●
	Auto-lock reading	●
	Reading Mode	●
	Sample ID	●
Data Management	Storage	500 sets of measurement parameters each
	View	●
	Delete	●
	GLP	●
Communications and external devices	Printer connection	●
	Content and format customization	●
	Connect to a PC for instrument control	●
pH/mV Measurement	pH electrode status/performance	Slope, Electrode status
	Multi-point calibration	5 points
	Automatic standard solutions recognition	4 groups
	Standards customization	●
	Standard groups customization	1 group
	Automatic temperature compensation	●
	Manual temperature compensation	●
	pH electrode diagnostics	●
Temperature Measurement	Temperature units	°C, °F

● Equipped

2 Safety Notices

Please read the entire contents of this manual carefully before use, and please keep this manual properly. The user **MUST** use the instrument following this manual to avoid damage to the user and equipment.

Before using the meter, **READ** the following notes:

- **DO NOT DISASSEMBLE** the device for inspection or repair.
- To prevent electric shock or damage to the device, **DO NOT** place cables and connectors in any liquid, wet or corrosive environment.
- Please use the defaulted power adapter, **DO NOT** use it if the power cord is damaged (the wire is exposed or broken).
- **DO NOT** use in flammable and explosive environments.
- **DO NOT** use if the user finds any abnormalities such as damage or deformation of the device.

The following identifiers will be used in this manual.

3 Terms Explanation

- **pH/pX Slope:** The amount of potential change generated by each 1 pH/pX change, expressed in mV/pH or by 100% Theoretical Slope (PTS). $pX = -\log[X]$, where $[X]$ means molar concentration (mol/L) of X ions.
- **E0 of pH:** Also known as "zero potential", it usually refers to the potential value at a pH of 7.
- **One-point calibration:** Calibration with a standard solution.
- **Two-point calibration:** Calibration with two standard solutions.
- **Multi-point calibration:** Calibration with more than two standard solutions.

4 Overview and Installation

4.1 Overview




Figure 4-1 Overview-Front View

- 6 Power sockets
- 7 Wristbands
- 8 Flip Tilting Stand
- 9 Electrode Holder
- 10 Electrodes

Figure 4-2 Overview- Back View

- 11 4-pin aviation connectors
- 12 pH electrode connectors
- 13 Electrode protection cap
- 14 pH electrodes

Figure 4-3 Electrodes and connector Table 4-1 Connector Specifications

Electrode type	Connector specifications	Electrodes Connection
pH electrode	BNC(Q9)	pH electrode, ORP electrode, ISE probe
Temp electrode	4-pin aviation	ATC probe

The PH310F adopted a pH combination Electrode with a four-pin aviation connector. For replacement purchasing, please choose electrodes with the proper connector.

4.2 Instrument Installation

4.2.1 Electrodes Stand Installation

Figure 4-4 Wristband Installation

Installation:

- 1) Pull the wristband switch at the head of the wrist so that its closed round hole opens.
- 2) Snap the opened round hole card slot into the shaft of the corresponding instrument housing and close the round hole.
- 3) The lower wristband card slot is also operated as above.

4.2.2 Electrodes Connection

Push the pH Electrode into the electrode holder. Remove the socket protector cap of the pH electrode. Connect the pH electrode into the right socket. If the ATC probe is applied, or ATC has been integrated into the pH probe, please connect the ATC probe onto the DO/T electrode socket.

5 Instrument Operation

5.1 Switch On/Off

Press and release  to switch on the meter. The startup screen shows software version and other related information. After the self-test program, the screen turns to the homepage and the meter are ready to measure. If the meter does not turn on, charge the meter for 15 minutes. Otherwise, please contact the manufacturer for further assistance.

The meter equipped with 7 function keys. Users press and hold  the key for more than 3 seconds and release to shut down.

5.2 Screen Icons


Figure 5-1 Screen icons explanation

1 System time. 2 Reading mode. 3 Measurement parameters. 4 Measurement information. 5 Calibration information. 6 Sample ID. 7 User ID. 8 Operation tips. 9 Power information. 10 Reading status. 11 Soft function keys.

The instrument displays symbol identification that has the following functional implications:

Table 5-1 Symbol Explanation

No.	Symbol	Explanation
1		Reading status, display the measurement status of reading, stable, locked each indicates that the processing, stable, and reading completed.
2		The percentage slope of the pH electrode calibration data
3		Standard buffer solution for calibration
4		Time to Calibrate Electrodes
5		Standard buffer solution used to calibrate the electrode
6		Automatic temperature compensation
7		Manual temperature compensation
8		User ID
9		Sample ID
10		Operation Tips

5.3 Function Key



Figure 5-2 Function keys explanation

Table 5-2 Key Function Explanation

No.	Key	Explanation	Note
1		Power	Press to switch on/off
2		Setting	Set the parameters and settings
3		Cancel	Cancel the operation
4		Enter	Confirm the option
5		F1	Function key, Corresponds to the function options on the screen
6		F2	Function key, Corresponds to the function options on the screen
7		F3	Function key, Corresponds to the function options on the screen

5.4 Parameter Settings

When the meter is in the standby status or measuring status, press "Setting" to reset the instrument setting.

5.4.1 Tutorial settings

For the first use, please follow the guide to settings the measurement parameters. After all the settings, press the "Enter" to return to the previous page.

5.4.2 Reading Mode Settings

The meter provides three reading modes, including continuous reading, auto reading, and timed reading.

- Continuous reading: The instrument displays real-time measurement results. User can end the measurement at any time and save the last result.
- Auto-reading: The measurement reached the balance, and the meter locked the reading result. The meter offers "Fast", "Medium", "Strict" and "Custom" four options for endpoint detection conditions.
- Time reading method: Timed Reading contains two kinds of timed reading methods: "Interval measurement" and "Timed measurement".

"Interval Measurement" provide measurement results at interval time and timed reading provide measurement result after a set time.

Table 5-3 Reading Parameters Settings

Stability Type		pH
Fast	Stable time	4s
	Fluctuation	1mV
Medium	Stable time	6s
	Fluctuation	0.5mV
Strict	Stable time	8s
	Fluctuation	0.1mV
Custom (Recommended value)	Stable time	1 to 30s
	Fluctuation	0.1~1mV

5.4.3 pH Parameter Settings

5.4.3.1 pH standard groups

The meter provides various Standards Group including DIN, NIST, GB and USA. And allows the user to prepare the customized Standard groups.

Table 5-4 Standard Solution Groups

Groups	Contents
NIST	1.68pH, 4.01pH, 6.86pH, 7.00pH, 7.42pH, 10.01pH, 12.47pH
USA	1.680pH, 4.010pH, 7.000pH, 10.010pH
DIN	1.68pH, 2.00pH, 3.56pH, 3.78pH, 4.01pH, 6.87pH, 7.00pH, 7.42pH, 9.18pH, 10.01pH, 12.45pH
GB	1.68pH, 3.56pH, 4.00pH, 6.86pH, 7.41pH, 9.18pH, 12.46pH

The meter supports up to five-points calibrations. Neighboring standards (pH gap<2) choice in the group may be frozen for accurate calibration. For neighboring standards, please choose the customization to perform calibration.

If the selected standard solution group is different from the pH standard buffer solution used, it will lead to wrong calibration results.


Figure 5-3 Selection of standard groups and standard solution

5.4.3.2 Recognition

Auto Mode and Manual Mode.

In some special cases, it is necessary to use some non-standard pH buffer solutions, or use two very close pH standard buffer solutions for electrode calibration. In this case, the manual standard solution identification function can be used. When set to "Manual Mode", the pH value of the current standard solution can be input during and used for electrode calibration.

5.4.3.3 Resolution settings

The pH measurement resolution of the instrument is adjustable. pH resolution: 0.01pH and 0.1pH.
mV resolution: 0.1 mV and 1 mV.

5.4.4 Temperature Parameter Settings

The temperature unit of the meter is selectable in °C and °F. Temperature compensation mode: ATC and MTC.

ATC means automatic compensation.

MTC means manual compensation. It allows user to input the temperature.

ATC input selection.

For single-parameter measurement, please follow the above instruction.

When conductivity compensation electrode dissolved oxygen electrode is applied. The default ATC input is from DO electrode.

If the ATC input selection is incorrect. Try resetting or manually entering the temperature.

5.4.5 Data Management Settings

5.4.5.1 Sample ID type

The instrument supports three setting methods of Sample ID: number order, time order, and manual.

- Number order: The sample ID No. is increasing with series number.
- Time order: The sample ID No. is increasing with sample measuring time. Format: Year/Y, Month/M, Day/D, Hour/H, Minutes/M, Second/S
- Manual: Manually set the sample ID No. It allows samples to manually enter the sample ID when saving or printing data.

5.4.5.2 Result Autosave

When this function is enabled, the meter saves the results when the reading is stable in the auto-reading and interval timed reading mode.

5.4.5.3 Data Overwrite

The meter provides 1000 sets of measurement results storage space. When this function is enabled, the results data that exceeds capacity will overwrite the old results data.

5.4.6 Output option

The data formats are GLP, STD Format, and Custom. It could select one data format to output the result.

5.4.7 User ID Settings

Set the user ID.

5.4.8 System Parameter Settings

5.4.8.1 System Date & Time Settings of system date and time.

5.4.8.2 Buzzer setting

Users can set the key sound by this setting.

5.4.8.3 Brightness setting

Users can adjust the screen brightness by this setting.

5.4.8.4 Auto Power off

The meter provides auto shutdown function. When the meter is not using, the meter switches off automatically.

5.4.8.5 Restore Default

The meter supports "Restore Default" and "Restore Parameters". "Restoring Default" will restore all meter parameters to the factory state. "Restoring parameters" will restore the measurement parameters to the factory state.

5.4.8.6 Software version

Users can find the software version information on the general setting page.

5.5 pH Measurement

5.5.1 Calibration Preparation

The electrode slope and zero potential of pH electrodes drift slightly over time. To accurately measure pH, it is recommended to calibrate the Ph electrode before use, the instrument supports 1-5 point calibration.

One point calibration is a calibration process with a single standard solution, commonly applied in a quick test. The calibration slope is 100% in here.

Two-point calibration is to use two pH standard buffer solutions to calibrate the electrode, and construct a linear calibration curve through two points. Two-point calibration is the most commonly used calibration method, and it is usually recommended that the pH value of the solution to be measured lies between the two standard buffer solutions. Two-point calibration can improve pH measurement accuracy.

Multi-point calibration is a calibration process with more than one standard solution. It is recommended to calibrate between two standard buffer solutions at the pH of the solution to be tested. Multi-point calibration covers a wider measurement range for accurate pH measurement. Before starting calibration, please prepare one or more pH standard buffer solutions.

5.5.2 Standards group selection

Before starting calibration, please prepare one or more pH standard buffer solutions. The meter has standards recognition function. Please set the Standard Group before the measurement.

You can also set the identification type to "Manual Mode" and manually enter the nominal value during the calibration process.

5.5.3 pH Calibration

Go to pH electrode calibration by "pH Calibration". Connect the ATC probe or enter the reference temperature manually.


Figure 5-4 pH electrode calibration information

The calibration process is as follows:

1. Setting.
 - 1) Set the parameters (e.g. pH).
 - 2) Select NIST standard solution group, and check pH 4.01, pH 7.00 and pH 10.01 three standard solutions.
 - 3) Set the Auto Mode recognition.
2. Connect the ATC probe or enter the temperature manually.
3. Press the F2 "Calibrate" for one-parameter measurements or press the F2 "Calibrate"- "pH Calibration" for multi-parameter measurements.
4. Put the cleaned electrode into pH 4.01 standard solution.
5. Wait for the instrument to display "Auto Mode Matched".
6. If only 1-point calibration is required, after 1-point calibration is completed, press the "Enter" key to complete the calibration.
7. If multi-point calibration is required, please replace the pH7.01 and pH10.01 standard buffer solutions. After cleaning the electrode, put the electrode into the standard solution. After the instrument recognizes it successfully, the instrument reads stably, press the F3 "Next Point" to complete the calibration.
8. After completing the calibration, press the "Enter" key to complete the calibration, save the calibration results and end the calibration, directly enter the start interface. If the checked standard solution group is 5, automatically end the calibration after five points of calibration.

The meter will not save the date if the calibration results are not confirmed.

5.5.4 pH Measurement

The measurement process is as follows:

1. Setting.
 - 1) Set the parameters (e.g. pH).
 - 2) Set the reading mode (e.g. continuous reading, auto-reading, or timed format).
2. Connect the ATC probe or enter the temperature manually.
3. Rinse the pH electrode with DI water, dry out.
4. Put the electrode into test solution under test.
5. Put the measurement end of the electrode into the sample solution.
6. In the measurement status, press "Enter" to enter into measurement status.
7. When the reading is stable, read the results.

8. Press the "Save" to save the measurement results and print the result.
9. Between measurements, stored pH electrode in distilled or deionized water.
10. After measurement, rinse the pH electrode with deionized water thoroughly.


Figure 5-5 pH measurement information

The measurement end of the electrode should well be immersed into the sample solution.

For high accuracy measurement, make sure the measurement is carried out at the lab with constant temperature and pressure.

5.6 Data Management

Press "Data" to view the detail of results.

The meter stores the measurement results independently according to the measured parameters. The meter provides data Storage 500 sets for each parameter (pH/mV).

The user can press "Delete" into the delete menu. It allows users to select the parameter data or all data to delete.

The user can view the data filter by parameter, locate No. or stored date.

By the filter setting, press "Enter" to look up the data. The filter data shows in a graph. Press "←" and "→" to choose data. User can choose one and press "Enter" key to see the detail result. Users can press "Operate" to the operation menu. In the operation menu, it allows to delete current result, output the current result, output matched result and output all result.


Figure 5-6 Results setting view

Figure 5-7 Results data graph


Figure 5-8 Measurement result output settings

Please connect the meter with printer when the meter is switch off.

The Brad rate is 9600bps for printer connecting. The default setting is 8 data bits, one start bit, one stop bit, and no parity.

The output format is followed:

```
***** Report Title
-----
Measure Time:2021/01/19 12:27:28
Operator: Operator 1 Model: PH-3004 pH Meter Serial Number:
SW Version: Ver 1.00
-----
.....MATCHED INFO
Stored Num:                28
Matched Num:               1
Stored No.:                15
-----
.....CALIB INFO
Calib Operator:           REX Team Calib Time: 2020/05/13 08:30:00
Calib Num:                3
.....CALIB RESULT STD 1: 4.00pH    177.3mV  25.0c
STD 2: 6.86pH            8.0mV  25.0c STD 3: 9.18pH -129.1mV  25.0c
pH Slope 1:              100.00%
pH E0 1:                  0.0mV
pH Slope 2:              100.00%
pH E0 2:                  0.0mV
-----
.....BRIEF INFO
Reading Mode:             Timed Reading
Stable Type:              Medium
Temp Comp Type:          ATC
-----
.....SAMPLE INFO
Sample ID:                Sample 1
-----
.....RESULT
Result:                   7.00pH
Signal Value:             -0.0mV
Temp Value:               25.0c
-----
-----
Signature:
```

6 Maintenance/Troubleshooting

6.1 Meter Maintenance

The correct use and maintenance of the instrument can ensure the accurate and reliable performance of the instrument. Maintenance plays an important role in pH/pX measurements, which have high input impedance and are susceptible to damage from static electricity and electromagnetic field disturbances. Additionally, exposure to chemicals or harsh use environments can affect performance.

The pH/pX electrode socket has a protective plug, when the meter is not in use, please insert the protective plug into the pH/pX socket.

The protective cover can avoid surface damage.

Keep the meter and accessories clean and away from acids, alkalis, and any corrosive solutions/gases. Users can clean the meter surface with poor waters and detergent.

The electrode socket of the instrument must be kept clean and dry, and should not be in contact with acid, alkali, and salt solutions.

6.2 Electrodes Maintenance

For more detailed information, please refer to the electrode instruction manual.

6.3 Battery Maintenance

The instrument uses a built-in lithium battery. The battery takes 6 hours to fully charge. Follow the below instruction to maintain battery:

- When the meter is turned on and not charged, the power light flashes four times when the power flashes once when the power is 25% and is 100%.
- When the meter switches off, the power light keeps the lighting for a full charge.
- Precautions for the use of lithium batteries:
 - Keep the meter away from the high-temperature environment.
 - Keep away from water.
 - Keep away from flammable and explosive substances.
 - Please charge for 15 minutes if not in use for a long time.
 - Please stop charging in time after fully charged.
 - Continuous charging may result in decreased battery life.
 - Please charge and discharge it every 3 months to extend the life of the lithium battery.

6.4 Troubleshooting

Table 6-1 Troubleshooting

Phenomenon	Probable reasons	Solutions
1. No Display	Damage to the meter.	Replace or repair as required.
2. Incorrect mV measurement is	1. The electrode is out of service life. 2. The electrode plug is in poor contact.	1. Replace the electrodes. 2. Connect the protection plug, if the potential is not 0mV, please contact the after-sales service.
3. Incorrect pH measurement	1. Refer to as 2.2. 2. Refer to as 2.2. 3. The electrodes are not calibrated or are calibrated incorrectly.	1. Refer to as 2.2. 2. Refer to as 2.2. 3. Recalibrate the electrode or replace the standard solution.

If the meter still does not work, please contact your local dealer for further assistance.

7 Technical Supports

7.1 Accessories

Please refer to the accessories table for purchasing recommendations.

Table 7-1 Meter accessories


Name	Description
E-301-QC 3 in 1pH composite electrode	pH Measurement Probe
Standard pH buffer solution 4.01/7.00/10.01	Standard solution

8 Appendixes

8.1 pH-Temperature Relationship Table of pH Standard Solutions

Temperature (°C)	1.68	4.01	7.00	10.01
5	1.67	4.00	7.09	10.25
10	1.67	4.00	7.06	10.18
15	1.67	4.00	7.04	10.12
20	1.68	4.00	7.02	10.06
25	1.68	4.01	7.00	10.01
30	1.68	4.01	6.99	9.97
35	1.69	4.02	6.98	9.93
40	1.69	4.03	6.97	9.89
45	1.7	4.04	6.97	9.86
50	1.71	4.06	6.97	9.83

Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.