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```
*****  
Report Title  
-----  
Measure Time:2021/02/06 13:31:51  
Operator: Admin  
Model: PH-3006 pH Meter  
Serial Number:  
SW Version: Ver 1.00  
-----  
.....SAMPLE INFO  
Sample ID: Sample 1  
.....ELECTRODE INFO  
pH EC ID: pH Electrode  
.....CALIB INFO  
Calib Operator: REX Team  
Calib Time: 2020/06/18 12:13:10  
Calib Num: 3  
.....CALIB RESULT  
STD 1: 4.01pH 176.8mV 25.0c  
STD 2: 7.00pH 0.0mV 25.0c  
STD 3: 10.01pH -178.1mV 25.0c  
pH Slope 1: 100.00%  
pH EO 1: 0.0mV  
pH Slope 2: 100.00%  
pH EO 2: 0.0mV  
.....BRIEF INFO  
Reading Mode: Continuous Reading  
Stable Type: Medium  
Temp Comp Type: ATC  
.....RESULT  
Result: 7.00pH  
Signal Value: -0.2mV  
Temp Value: 25.0c  
-----  
Signature:
```

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1 Introduction

1.1 Introduction

PH-3006 pH-meter is a newly designed laboratory high-precision intelligent analytical instrument that can be widely used in universities, environmental protection, medicine, food, sanitation, geological prospecting, metallurgy, ocean exploration and other fields measurements for acid rain detection, industrial wastewater, surface water, drinking water, beverages, daily chemical products, textiles and other related industries.

General Features

- Colour High contrast LCD touchscreen, 7.0 inches.
- Intelligent operation system provides features including user management, calibration management, method management, data management and etc.
- Multi-reading feature allows auto-read, timed-read and continuous-read.
- Automatic/Manual temperature compensation ensures accurate results.
- Data Storage 1000 sets for each parameter (GLP-compliant).
- Data analysis feature helps user review, compare and recalculate results.
- Support for USB or RS-232 communication.
- Reset feature automatically resumes all settings back to factory default options.
- IP54 waterproof.
- 1-6 points calibration with Standard Recognition.
- Selectable pH buffer groups, including DIN, NIST, USA, MERK, JIS, GB.
- Support for pH standard group customization.
- Automatic electrode diagnosis with pH slope and offset display.

1.2 Technical Specification

Table 1-1 Instrument Specifications

Model		PH-3006
pH level		0.001pH
mV	Range	(-2000.00~2000.00)mV
	Resolution	0.01mV
	Accuracy	±0.03% or ±0.1mV
	Repeatability	0.1mV
	Input Current	≤1×10 ⁻¹² A
	Input Impedance	≥3×10 ¹² Ω
pH	Range	(-2.000~20.000)pH
	Resolution	0.001pH
	Accuracy	±0.002pH
	Repeatability	0.001pH
	Measurement Accuracy	±0.01pH
	Measurement Repeatability	0.005pH
Temperature	Range	(-10.0~135.0) °C/(14.0~275.0) °F
	Resolution	0.1 °C/0.1 °F
	Accuracy	±0.1 °C
	Instrument indication error	±0.3°C (0°C~60°C), ±1.0 °C (Other range)
Work environment		Environment temperature: (0~40)°C Relative humidity: not more than 85%
Dimensions (L×B×H), weight (kg)		220mm ×195mm×68mm, about 0.95kg
Power supply		AC Adapter,100-240VAC input, DC24V output

1.3 Function Introduction

Table 1-2 Functions Specification

Function	Explanation	
Basic Function	Backlight adjustment	•
	Automatic diagnostics	•
	Reset settings	•
	Default parameter	•
	Prompt Sound	•
	Time settings	•
	Firmware upgrade	•
	Password	•
	Anti-interference automatic recovery	•
	Power failure protection	•
	Automatic shutdown	•
	Protection	IP54
Reading Function	Reading stability settings	•
	Auto-lock reading	•
	Reading Mode	Continuous, Auto, Timed
	Sample ID	•
	Alarm	•
Data Management	Storage	1000 sets of measurement parameters each. 100 sets of verification parameters each. 100 sets of calibration parameters each.
Data Management	Filter	Storage number, time, operator, method name, sample ID and electrode ID.
	View	Sheet and Graph
	Delete	•
	Alarm	•
	GLP	•
Communications and external devices	Connect to Printer	RS232
	U Disk	•
	Content and format customization	GLP, Standard, custom
	PC Output	•
	Connect to a PC for instrument control	•
	Scanner	•
	Auto sampler	•
	pH electrode status/performance	Slope, Electrode status (Excellent, good, bad)

Measurement	Multipoint calibration	6 points
	Automatic standard solutions recognition	6 groups
	Standards customization	•
	Standard groups customization	1 group
	Temperature compensation	Automatic/ Manual
	pH electrode diagnostics	•
	pH electrode calibration reminder	•
Measurement	Mandatory calibration of pH electrodes	•
	pH standard solution verification	•
	pH Mandatory Verification	•
	ORP measurement	•
	ORP electrode calibration reminder	•
	Mandatory calibration of ORP electrodes	•
Temperature Function	Temperature Unit	°C , °F
	Temperature calibration	•
Management	GMP Mode	•
	System & Running log	•
	User management	3 Level permissions, 8 users
	Method management	•
	Sample list	•

• Equipped

2 Safety Notices

Please read the entire contents of this manual carefully before use, and please keep this manual properly. The user **MUST** use the instrument following this manual to avoid damage to the user and equipment.

Before using the meter, **READ** the following notes:

- **DO NOT DISASSEMBLE** the device for inspection or repair.
- To prevent electric shock or damage to the device, **DO NOT** place cables and connectors in any liquid, wet or corrosive environment.
- Please use the defaulted power adapter, **DO NOT** use it if the power cord is damaged (the wire is exposed or broken).
- **DO NOT** use in flammable and explosive environments.
- **DO NOT** use if the user finds any abnormalities such as damage or deformation of the device.

The following identifier will be used in this manual.

3 Terms Explanation

pH

- **pH Slope:** The amount of potential change generated by each 1 pH change, expressed in mV/pH or by 100% Theoretical Slope (PTS). pH
 - $= -\log[H]$, where [H] means molar concentration (mol/L) of H ions.
- **E0 of pH:** Also known as "zero potential", it usually refers to the potential value at a pH of 7.
- **One-point calibration:** Calibration with a standard solution.
- **Two-point calibration:** Calibration with two standard solutions.
- **Multi-point calibration:** Calibration with more than two standard solutions.

Redox potential (ORP)

- **Redox potential:** Also known as ORP value, it refers to the potential difference between the potential difference of a measuring battery relative to a standard hydrogen electrode consisting of an indicator electrode, a reference electrode, and a solution under test. The symbol for ORP measurement is "RmV", and the unit is mV.
- **One-point calibration:** Calibration with an ORP standard solution.
- **Offset:** $= ERmV - E$ (ERmV potential relative to the hydrogen electrode, E measured potential). At a certain temperature, when the ORP calibration solution is used for calibration, the potential value "RmV" and the potential difference between the measured system.
- **Calibration reminder (Recommended):** "Electrode calibration" prompts on the homepage, user can start measurement without last calibration.
- **Calibration reminder (Mandatory):** "Electrode calibration" prompts on the homepage, calibration is required for accurate measurement.
- **Verification reminder (Recommended):** "Electrode verification" prompts on the homepage, user can start measurement without last verification.
- **Upper limit:** The upper limit value is monitored for the measured or calibrated data, and the upper limit value must not be lower than the lower limit value.
- **Lower limit:** Lower limit monitoring of measured or calibrated data.

4 Overview and Installation

4.1 Overview


Figure 4-2 Overview- Back View


Figure 4-3 Electrodes and connectors

Table 4-1 Connector Specifications

Connector specifications	Electrode type
BNC(Q9)	pH electrode, ORP electrode
4-pin aviation	Temperature electrode
Banana	Reference electrode

4.2 Instrument Installation

4.2.1 Electrodes Stand Installation

- 1 Fixed screw
- 2 Electrode stand
- 3 Electrode base
- 4 Fixing hole


Figure 4-4 Electrode Stand Installation

- 1) Put the electrode base on the table.
- 2) Screw the fixed screw into the fixing hole.
- 3) Place the electrode holder from top to bottom to fix at the vertical shaft of the fixing screw.

4.2.2 Electrodes Connection

Push the pH electrode into the electrode holder. Remove the protector cap of the pH electrode. Connect the pH electrode into the right socket. If the ATC probe is applied, or ATC has been integrated into the pH probe, please connect the ATC probe onto the Temperature electrode socket.

When using a 3 in 1 composite pH electrode, the Q9 connector can be inserted into the pH electrode socket, and at the same time, the 4-pin connector can be inserted into the temperature electrode socket.

5 Instrument Operation

5.1 Switch On/Off

Press and  to switch on the meter. The start-up screen shows the software version and other related information.

After the self-test program, the screen turns to the homepage and the meter is ready to measure. Users press and hold  key for more than 3 seconds and release to shut down.

5.2 Screen Icons


Figure 5-1 Screen icons explanation

1 User ID. 2 Measurement parameters. 3 Reading states. 4 Main measurement box. 5 Method information. 6 Tips information. 7 System time. 8 Calibration information. 9 Sample ID. 10 Function buttons.

The instrument displays symbol identification that has the following functional implications:

Table 5-1 Symbol Explanation

No.	Symbol	Explanation
1		Reading status, display the measurement status of reading, stable, locked, each indicates that the processing, stable, and reading completed.
2	PTS	The percentage slope of the pH electrode calibration data.
3	ATC	Automatic temperature compensation
4	MTC	Manual temperature compensation
5	ORP	Redox potential value, in mV
6	Offset	Offset potential, in mV
7	STDCHK	Standard Check
8	CALTime	Time to Calibrate Electrodes
9	BUFF	Standard buffer solution for calibration
10		Standard buffer solution used to calibrate the electrode
11		Measurement method management, display the current method information
12		User ID
13		Sample ID
14	Auto Mode	Auto-recognition of Standards
15	Manual Mode	Manual-recognition of Standards

5.3 Methods Management

The meter provides a library of built-in methods. Select the proper method to measure, system loads method's preset information, including the method name, brief overview, creator, creation date, and measurement parameters.

In the method settings, it supports the navigation setting including the method information, parameters, reading mode, pH parameters setting, ORP parameters settings, temperature setting and data management setting. For the first use, please follow the guide to settings the measurement parameters.

Figure 5-2 Methods Information

If a modified built-in method is not saved to library, the “Method management” logo will turn to red. When the user returns to the method management page, the system will prompt.

5.3.1 Update Method

When the method is revised by the setting, there is a red mark in the method management logo to remind user that the current method is different from the method in the library.

When press the method in the method library, there is a system tip to give user a choose to update or create a new method.

5.3.2 Create a new Method

There are three modes to create a new method: create method step by step in the system, save the setting as a new method in the system, and save the method with a new name, description and location. The location can be the system and U disk.


Figure 5-3 Method setting view

5.3.3 Delete a Method

There are 3 methods saved in the meter as the default setting in the library method. The No.001-No.003 method cannot be deleted, and the others can be deleted.

5.3.4 Save and Load a Method

It is support to store up to 50 methods in the meter, and save and load the methods in the U Disk.

5.4 Parameter Settings

Press "Setting"- to the setting menu. It includes tutorial setting, method information, select parameters, reading mode, pH parameters, ORP parameters, Temperature parameters and data management, output option, User management, system parameters and GMP Mode.

Figure 5-4 Parameters setting view

5.4.1 Tutorial settings

For the first use, please follow the guide to settings the measurement parameters. After all the settings, press

 Parameter Setting " to return to the previous page.

5.4.2 Select parameters

It could select one measurement parameter from pH and ORP every test.


Figure 5-5 Select parameters

5.4.3 Reading Mode Settings

The meter provides three reading modes, including continuous reading, auto reading, and timed reading.

- Continuous reading: The instrument displays real-time measurement results. User can end the measurement at any time and save the last result.
- Auto-reading: The measurement reached the balance, and the meter locked the reading result. The meter offers "Fast", "Medium", "Strict" and "Custom" four options for endpoint detection conditions.
- Time reading method: Timed Reading contains two kinds of timed reading methods: "Interval Measurement" and "Timed Measurement". "Interval Measurement" provide measurement results at interval time and "Timed Measurement" provide measurement result after a set time.

Table 5-2 Reading Parameters Settings

Stability Type		pH
Fast	Stable time	4s
	Fluctuation	0.6mV
Medium	Stable time	6s
	Fluctuation	0.1mV
Strict	Stable time	8s
	Fluctuation	0.03mV
Custom (Recommended)	Stable time	1 to 30s
	Fluctuation	0.03~1mV

5.4.4 pH Parameter Settings

5.4.4.1 pH Electrode management

The electrode information includes the name, serial No., user ID, registration time, service life, calibration message, calibration number, last calibration time. It allows user to choose one electrode to calibrate and measure samples.

The meter has one created default pH electrode, which is convenient for users. When using the built-in default electrode, the meter automatically loads the electrode serial number and the last calibration result, including calibration points, slope value, zero-point potential value, etc., and the calibration time and calibrator will also be loaded.

The screenshot shows the 'Parameter Setting/pH Parameter' screen. On the left is a navigation menu with options: pH Electrode (selected), pH Buffer Group, Resolution, pH Alarm Setting, Calibrate Reminder, and pH STD Check. The main area displays the following information:

- Name:** pH Electrode
- Serial No.:** 12345678
- Recognition:** Auto Mode
- CAL Result:** A table showing calibration results for four buffer groups (BUF1-BUF2, BUF2-BUF3, BUF3-BUF4, BUF4-BUF5). Each row shows 100.00% slope, 0.00mV EI, and Perfect comment.
- CAL Time:** 2022/05/08 16:23:43
- User ID:** REX Team

At the bottom right, there is an 'Enter' button.

Figure 5-6 EC electrode information

5.4.4.2 pH standard groups

The meter provides various Standards Group including GB, DIN, NIST, USA, MERK, and JIS. And allows the user to prepare the customized Standard groups.

Table 5-3 Standard Solution Groups

Groups	Contents
NIST	1.677pH, 4.008pH, 6.864pH, 7.000pH , 7.416pH, 10.014pH, 12.469pH
USA	1.680pH, 4.010pH, 7.000pH, 10.010pH
DIN	1.680pH, 2.000pH, 3.557pH, 3.775pH, 4.008pH, 6.865pH, 7.000pH, 7.416pH , 9.184pH, 10.014pH, 12.454pH
GB	1.680pH, 3.559pH, 4.003pH, 6.864pH, 7.409pH, 9.182pH, 12.460pH
MERK	2.000pH, 4.000pH, 7.000pH, 9.000pH, 12.000pH
JIS	1.680pH, 4.008pH, 6.865pH, 7.413pH, 9.180pH, 10.010pH

The meter supports up to eight-points calibrations. Neighbouring standards (pH gap<2) choice in the group may be frozen for accurate calibration. For neighbouring standards, please choose the customization to perform calibration.

Figure 5-7 Selection of standard groups and standard solution

5.4.4.3 Recognition

Auto Mode and Manual Mode.

In some special cases, it is necessary to use some non-standard pH buffer solutions, or use two very close pH standard buffer solutions for electrode calibration. In this case, the manual standard solution identification function can be used. When set to "Manual Mode", the pH value of the current standard solution can be input during and used for electrode calibration.

5.4.4.4 Resolution settings

The pH measurement resolution of the instrument is adjustable. pH resolution: 0.01pH and 0.1pH.
mV resolution: 0.1 mV and 1 mV.

5.4.4.5 Alarm setting

The meter support pH measurement and calibration result monitoring alarm. The setting includes slope limit, potential value limit and monitoring options. When the measurement results out of the selected range, the meter alarms.

The screenshot shows the 'Parameter Setting\pH Parameter' screen. The 'pH Alarm Setting' option is selected in the left-hand menu. The settings are as follows:

pH Electrode	pH CAL Alarm:	On
pH Buffer Group	Slope Max(%):	120.00%
Resolution	Slope Min(%):	80.00%
pH Alarm Setting	E0 Max:	30.00mV
Calibrate Reminder	E0 Min:	-30.00mV
pH STD Check	pH Alarm:	On
	pH Max:	10.000pH
	pH Min:	4.000pH

An 'Enter' button is located at the bottom right of the screen.

Figure 5-8 Alarm setting information

5.4.4.6 Calibration reminder settings

The meter provides calibration prompts function. From the setting, the user can select calibration reminder (recommended) and calibration reminder (mandatory) for future calibration.

Calibration reminder (recommended): When enable the calibration reminder, the instrument verifies whether the calibration data of pH electrode is within the validity period. If calibration information expires, a striking calibration reminder will appear at "Calibrate" on the main interface, but it will not affect the measurement and data saved.

Calibration reminder (mandatory): When enable the mandatory calibration, the instrument verifies whether the calibration data of pH electrode is within the validity period. If calibration information expires, a striking calibration reminder will appear at "Calibrate" on the main interface, the measurement data cannot be saved and output.

The screenshot shows the 'Parameter Setting\pH Parameter' screen. The 'Calibrate Reminder' option is selected in the left-hand menu. The settings are as follows:

pH Electrode	pH CAL Message:	On
pH Buffer Group	CAL Interval:	10
Resolution	Time Unit:	Hour
pH Alarm Setting	pH Compulsive CAL:	On
Calibrate Reminder	Validity:	10
pH STD Check	Time Unit:	Hour

An 'Enter' button is located at the bottom right of the screen.

Figure 5-9 Calibration reminder setting information

5.4.4.7 Verification Reminder settings

The meter provides verification prompts function. The electrode slope and zero potential of pH electrodes drift slightly over time, and this change will have a direct impact on the measurement results. So the verification of standards is important for accurate results.

Verification reminder (mandatory): When enable the verification reminder, the instrument verifies whether the verification data of pH electrode is within the validity period. If verification information expires, a striking verification reminder

will appear at "Calibrate" on the main interface, the measurement data cannot be saved and output. It allows user to input the pH standard value and the allowable deviation. If the deviation amount is set to 0, it means that no verification is required.


Figure 5-10 pH verification information

5.4.5 ORP Parameter Settings

5.4.5.1 ORP Electrode management

In the meter setting, a set of ORP Electrode info is defaulted. At the measurement, the meter loads the defaulted pH electrode's info including serial number, last calibration result.

When replacing a new ORP Electrode replacement, please created information set for the electrode.

5.4.5.2 Alarm setting

The meter support ORP measurement and calibration result monitoring alarm. The setting includes potential value limit and monitoring options. When the measurement results out of the selected range, the meter alarms.

5.4.5.3 Calibration reminder settings

The meter provides calibration prompts function. From the setting, the user can select calibration reminder (recommended) and calibration reminder (mandatory)for future calibration.

5.4.6 Temperature Parameter Settings

The temperature unit of the meter is selectable in °C and °F. Temperature compensation mode: ATC and MTC.

ATC means automatic compensation.

MTC means manual compensation. It allows user to input the temperature.

5.4.7 Data Management Settings

5.4.7.1 Sample ID type

The instrument supports three setting methods of Sample ID: number order, time order, and manual.

Number order: The sample ID No. is increasing with series number.

Time order: The sample ID No. is increasing with sample measuring time. Format: Year/Y, Month/M, Day/D, Hour/H, Minutes/M, Second/S

Manual: Manually set the sample ID No. I allows samples to manually enter the sample ID when saving or printing data.

5.4.7.2 Result Autosave

When this function is enabled, the meter saves the results when the reading is stable in the auto-reading and interval timed reading mode.

5.4.7.3 Data Overwrite

The meter provides 1000 sets of measurement results storage space. When this function is enabled, the results data that exceeds capacity will overwrite the old results data.

5.4.8 Output option

The data format is GLP, STD Format, and Custom. It could select one data format to output the result.

5.4.9 User ID management

The meter supports hierarchical user authority management and password security management. It allows user to revise the user ID, user type and password, and create a new account. A maximum of eight users can be created, and user type can be system admin, method admin, and operator. The default user is system admin and the initial password is blank. Only system admin has the permission to create all type of users.

5.4.10 System Parameter Settings

5.4.10.1 System Date & Time

Settings of system date and time.

5.4.10.2 Buzzer setting

Users can set the key sound by this setting.

5.4.10.3 Brightness setting

Users can adjust the screen brightness by this setting.

5.4.10.4 Auto Power off

The meter provides auto shutdown function. When the meter is not using, the meter switches off automatically.

5.4.10.5 Restore Default

The meter supports "Restore Default" and "Restore Parameters". "Restoring Default" will restore all meter parameters to the factory state. "Restoring parameters" will restore the measurement parameters to the factory state.

5.4.10.6 Software version

Users can find the software version information on the general setting page.

5.4.11 GMP Mode

The meter supports the GMP mode, in order to achieve strict control of the measurement process and data traceability. In the mode, the meter enables the login password protection, three-level user rights, data security protection, system log and run log and other security control.

Press the "Setting"- "GMP Mode" to set the mode.


Figure 5-11 GMP mode setting

- In the GMP Mode, please check the information as follows:
- It needs to set a login password for system admin because the default admin password is blank.
- It needs to identify the system time.
- It needs to create a method admin and an operator.

5.4.11.1 System Log

Access the System Log to view the log, such as GMP mode open log, GMP mode close log, system time modification log and creation log.


Figure 5-12 System log in GMP mode

5.4.11.2 Running Log

Access the Operation Log to view the log, such as the data log and calibration log, print the data log, and print the calibration log, log in or log out of an existing account.

In the GMP mode, the data cannot be deleted because of the data protection.


Figure 5-13 Operation log in GMP mode

5.5 pH Measurement

5.5.1 Calibration Preparation

The electrode slope and zero potential of pH electrodes drift slightly over time. To accurately measure pH, it is recommended to calibrate the pH electrode before use, the instrument supports 1-6 point calibration.

One-point calibration is a calibration process with a single standard solution, commonly applied in a quick test. The calibration slope is 100% in here.

Two-point calibration is to use two pH standard buffer solutions to calibrate the electrode, and construct a linear calibration curve through two points. Two-point calibration is the most commonly used calibration method, and it is usually recommended that the pH value of the solution to be measured lies between the two standard buffer solutions. Two-point calibration can improve pH measurement accuracy.

Multi-point calibration is a calibration process with more than one standard solution. It is recommended to calibrate between two standard buffer solutions at the pH of the solution to be tested. Multi-point calibration covers a wider measurement range for accurate pH measurement. Before starting calibration, please prepare one or more pH standard buffer solutions.

5.5.2 Standards group selection

Before starting calibration, please prepare one or more pH standard buffer solutions. The meter has standards recognition function. Please set the Standard Group before the measurement.

You can also set the identification type to "Manual Mode" and manually enter the nominal value during the calibration process.

5.5.3 pH Calibration


Figure 5-14 pH electrode calibration information

The calibration process is as follows:

1. Select a method. The method includes the parameter (e.g. pH), NIST standard solution group, pH 4.01, pH 7.00 and pH 10.01, auto Mode recognition.
2. Connect the ATC probe or enter the temperature manually.
3. Press the "Calibrate"- "pH Calibration".
4. Put the cleaned electrode into pH 4.01 standard solution.
5. Wait for the instrument to display "Auto Mode Matched".
6. When the pH and temperature reading are stable, press the "Start".
7. If only 1-point calibration is required, after 1-point calibration is completed, press the "End" key to complete the calibration.

- If multi-point calibration is required, please replace the pH7.00 and pH10.01 standard buffer solutions. After cleaning the electrode, put the electrode into the standard solution. After the instrument recognizes it successfully, the instrument reads stably, press the "Next Point" to complete the calibration.
- After completing the calibration, press the "Enter" key to complete the calibration, save the calibration results and end the calibration, directly enter the start interface. If the checked standard solution group is 8, automatically end the calibration after eight points of calibration.
The meter will not save the date if the calibration results are not confirmed.

5.5.4 pH Verification

The verification process is as follows:

- Select a method. The method includes the parameter (e.g. pH),
- Set the verification pH standard value (e.g. pH 4.01) and auto Mode recognition.
- Connect the ATC probe or enter the temperature manually.
- Press the "Calibrate"->"pH Verification".
- Put the cleaned electrode into pH 4.01 standard solution.
- Wait for the instrument to display "Auto Mode Matched".
- When the pH and temperature reading are stable, press the "Start".
- If the pH standard value is in the allowable range, the verification result show pass. If the pH standard value is not in the allowable range, the verification result shows no pass.
- Press the "Save" key to save the verification results. Press the "End" key complete the verification.


Figure 5-15 pH verification information

5.5.5 pH Measurement

The measurement process is as follows:

- Select a method. The method includes the parameter (e.g. pH), reading mode (e.g. continuous reading, auto-reading, or timed format).
- Connect the ATC probe or enter the temperature manually.
- Rinse the pH electrode with DI water, dry out.
- Put the electrode into test solution under test.
- Put the measurement end of the electrode into the sample solution.
- Press "Measure" to enter into measurement status.
- When the reading is stable, read the results.

8. Press the "Save" to save the measurement results and press the "Output" to print the result.
9. Between measurements, stored pH electrode in distilled or deionized water.
10. After measurement, rinse the pH electrode with deionized water thoroughly.


Figure 5-16 pH measurement information

- The measurement end of the electrode should well be immersed into the sample solution.
- For high accuracy measurement, make sure the measurement is carried out at the lab with constant temperature and pressure.

5.6 ORP Measurements

5.6.1 Calibration Preparation

In general, ORP electrodes need few calibrations. When the electrode is used for the first time or has not been used for a long time, a calibration is needed.

The instrument supports ORP 1-point calibration, the electrode can be calibrated using the ORP calibration solution. The instrument automatically calculates the offset and compensates for the measurement.

Before the calibration, prepare the ORP calibration standard solution.

5.6.2 ORP Calibration

The calibration process is as follows:

1. Select a method. The method includes the parameter (e.g. ORP).
2. Connect the ATC probe or enter the temperature manually.
3. Press the "Calibrate"- "ORP Calibration".
4. Put the cleaned electrode into ORP standard solution (e.g., 462mV ORP standard).
5. Press the "STD value" to input the ORP standard value.
6. When the reading is stable, press the "Start" to complete the first point calibration, and the instrument displays and stores the calibration results.
7. Press the "End" key to complete the calibration.

If the calibration process is not completed, the calibration will be ended, and the setting parameters will not be saved.


Figure 5-17 ORP electrode calibration information

5.6.3 ORP Measurement

The measurement process is as follows:

1. Select a method. The method includes the parameter (e.g. ORP), reading mode (e.g. continuous reading, auto-reading, or timed format).
2. Connect the ATC probe or enter the temperature manually.
3. Rinse the pH electrode with DI water, dry out.
4. Put the electrode into test solution under test.
5. Put the measurement end of the electrode into the sample solution.
6. Press "Measure" to enter into measurement status.
7. When the reading is stable, read the results.
8. Press the "Save" to save the measurement results and press the "Output" to print the result.
9. Between measurements, stored ORP electrode in distilled or deionized water.
10. After measurement, rinse the ORP electrode with deionized water thoroughly.


Figure 5-18 ORP measurement information

5.7 Batch Sample Measurement

Batch sample measurement is suitable for multiple sample measurements to improve the work efficiency. This function needs to match the automatic sampler to achieve multi-sample measurement.

Enter the batch sample measurement function by "Sample List".

5.8 Data Management

Press "Data" to view the detail of results.

The meter stores the measurement results independently according to the measured parameters. The meter provides data Storage 1000 sets for each parameter (pH/mV/ORP/Temp),100 sets of calibration data and 1000 sets of verification data.

5.8.1 Filter

The data type includes the measurement data and verification data, and calibration data. The parameter can be selected by pH and ORP. By the lookup filter, the user can view the result data by storage number, time, operator, method name, sample ID and electrode Name. The filter data shows in a graph or curve. By the filter setting, press "Start Search" to look up the data.

Figure 5-19 Data setting view

5.8.2 Result

In the result page, press "Page Up", "Page Down", "<<" and ">>" to choose data. Users can press "Delete" to delete the current result.

5.8.2.1 Measurement Data

In the measurement data, the result is displayed as a graph or a sheet, and the coordinates correspond to the relevant information.

No.	Time	Sample ID	Signal	Temp	Result
0001	2022/05/08 16:58:58	Sample ID	0.01mV	25.0°C	7.000pH
0002	2022/05/08 16:59:02	Sample ID	0.01mV	25.0°C	7.000pH
0003	2022/05/08 16:59:08	Sample ID	0.02mV	25.0°C	7.000pH
0004	2022/05/08 16:59:12	Sample ID	0.02mV	25.0°C	7.000pH

Figure 5-20 Matched measurement data in a sheet


Figure 5-21 Matched result data results in a graph

5.8.2.2 Verification Data

In the verification data, the pH verification result is displayed as a sheet, the coordinates correspond to the relevant information.

The screenshot shows the 'Database' interface with the 'Results' tab selected. The parameter is 'pH STD Check'. The table below shows the following data:

No.	Time	Sample ID	STD Value	Meas Value	Conclusion
0001	2022/05/08 16:57:33	Sample ID	7.000pH	7.000pH	Pass
0002	2022/05/08 16:57:36	Sample ID	7.000pH	7.000pH	Pass
0003	2022/05/08 16:57:40	Sample ID	7.000pH	7.000pH	Pass

Figure 5-22 Matched verification data results in a sheet

5.8.2.3 Calibration Data

In the calibration data, the pH calibration result or ORP calibration result is displayed as a sheet, and the coordinates correspond to the relevant information.

The screenshot shows the 'Database' interface with the 'Results' tab selected. The parameter is 'pH CAL Log'. The table below shows the following data:

No.	Time	Operator	Comment
0001	2022/05/08 16:24:27	Admin	PTS:100.00%

Figure 5-23 Matched calibration data results in a sheet

5.8.3 Report

The report shows relevant information of the data in detail. It includes the result, calibration result, and basic information.


Figure 5-24 Measurement result report

5.8.4 Statistic

The meter supports the statistical analysis function. By pressing "Statistic", the meter calculates basic statistical results, including maximum, minimum, average, standard deviation, relative standard deviation, and related statistical information of the results.


Figure 5-25 Statistical analysis

5.8.5 Output Settings

It is the same to the output setting.

The output device can be printer or U disk.

Output the measurement results by selecting the data format. Format: Standard, GLP and custom.

Figure 5-26 Output settings

- Please connect the meter with printer when the meter is switch off.
- The Brad rate is 9600bps for printer connecting. The default setting is 8 data bits, one start bit, one stop bit, and no parity.

The output format is as follows:

```

*****
Report Title
-----
Measure Time:2021/02/06 13:31:51
Operator: Admin
Model: PH-3006 pH Meter
Serial Number:
SW Version: Ver 1.00
-----
.....SAMPLE INFO
Sample ID: Sample 1
.....ELECTRODE INFO
pH EC ID: pH Electrode
.....CALIB INFO
Calib Operator: REX Team
Calib Time: 2020/06/18 12:13:10
Calib Num: 3
.....CALIB RESULT
STD 1: 4.01pH 176.8mV 25.0c
STD 2: 7.00pH 0.0mV 25.0c
STD 3: 10.01pH -178.1mV 25.0c
pH Slope 1: 100.00%
pH EO 1: 0.0mV
pH Slope 2: 100.00%
pH EO 2: 0.0mV
.....BRIEF INFO
Reading Mode: Continuous Reading
Stable Type: Medium
Temp Comp Type: ATC
.....RESULT
Result: 7.00pH
Signal Value: -0.2mV
Temp Value: 25.0c
-----
Signature:

```

6 Maintenance/Troubleshooting

6.1 Meter Maintenance

The correct use and maintenance of the instrument can ensure the accurate and reliable performance of the instrument. Additionally, exposure to chemicals or harsh use environments can affect performance.

- If the meter is not used for a long time, please disconnect the power supply.
- The electrode socket of the instrument must be kept clean and dry, and should not be in contact with acid, alkali, and salt solutions.
- Keep the meter and accessories clean and away from acids, alkalis, and any corrosive solutions/gases.
- Users can clean the meter surface with clean waters and detergent.
- When the meter is transported, please follow the instructions:
 - please remove all connected cables.
 - Please remove the electrode holder.
 - Please use original packaging in the long distance transport to avoid damage.

6.2 Electrodes Maintenance

For more detailed information, please refer to the electrode instruction manual.

6.3 Troubleshooting

Table 6-1 Troubleshooting

Phenomenon	Probable reasons	Solutions
1. No Display	Damage to the meter.	Connect the adapter and press the power key to turn it on, replace or repair as required.
2. Incorrect mV measurement	1. The electrode is out of service life 2. The electrode plug is in poor contact	1. Replace the electrodes 2. Connect the protection plug, if the potential is not 0mV, please contact the after-sales service.
3. Incorrect pH measurement	1. Refer to as 2.1 2. Refer to as 2.2 3. The electrodes are not calibrated or are calibrated incorrectly	1. Refer to as 2.1 2. Refer to as 2.2 3. Recalibrate the electrode or replace the standard solution

If the meter still does not work, please contact your local dealer for further assistance.

7 Technical Supports

Accessories

Please refer to the accessories table for purchasing recommendations.

Table 7-1 Meter accessories

Name	Description
E-301-QC 3 in 1 pH composite electrode	pH Measurement Probe
pH 4.01/7.00/10.01 standard pH buffer solutions	To prepare the pH standard
501 ORP composite electrode	ORP measurement probe
462mV/296mV /419mV ORP calibration reagent	Calibrate ORP electrode
REX-6 Electrode stand	Hold the electrodes during measurement

8 Appendixes

Appendix 1

pH-Temperature Relationship Table of pH Standard Solutions

Temperature (°C)	1.68	4.01	7.00	10.01
5	1.67	4.00	7.09	10.25
10	1.67	4.00	7.06	10.18
15	1.67	4.00	7.04	10.12
20	1.68	4.00	7.02	10.06
25	1.68	4.01	7.00	10.01
30	1.68	4.01	6.99	9.97
35	1.69	4.02	6.98	9.93
40	1.69	4.03	6.97	9.89
45	1.7	4.04	6.97	9.86
50	1.71	4.06	6.97	9.83

Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.