MACROFRIGER-BLT GR-50T
380V 50Hz III (ESPECIAL)

7001491

CENTRÍFUGA REFRIGERADA PARA BANCOS DE SANGRE
BLOOD BANK REFRIGERATED CENTRIFUGE

Indice

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Index

1 Safety	19
1.1 Risks to which the operator is subject	19
1.2 Personnel qualification	19
1.3 Use with hazardous substances	19
2 General information	20
3 Packing list	20
4 Technical Specifications	21
5 Installation	22
5.1 Transport	22
5.2 Location.....	22
5.3 Initial setup.....	23
5.4 Drain connection	24
5.5 Head location	24
5.6 Power supply connection.....	25
6 Operation.....	25
6.1 General conditions during operation	25
6.1.1 General information	25
6.1.2 Maximum density of liquids to be centrifuged.....	25
6.1.3 Filling and location of the buckets	26
6.1.4 Emergency stop.....	26
6.1.5 Safety measures for the staff	26
7 Operating menus	26
7.1 Settings	27
7.2 Programming.....	27
7.3 Run a program (WORKING).....	28
7.4 Pre-programmed programs.....	29
<i>Notice to customers.....</i>	<i>29</i>
8 Maintenance	30
8.1 Visual inspection every 3 months	30
8.2 Software maintenance	30
8.3 Cleaning	31
8.4 Periodic official review	31
Review Form.....	32

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Fecha: _____ Fecha: _____

1 Safety

The equipment includes appropriate security measures. This protection is compromised if the centrifuge is used in a way not specified in this manual.

Risky situations that must be respected are indicated along this manual and must be respected.

Danger Risk

Risk of danger.

Respect the indicated instructions to perform the operation.

Electrical risks

Risk of electrical accident when accessing the areas indicated by this sign or when performing the operations indicated by this icon in this manual.

Respect the indicated instructions to perform the operation.

1.1 Risks to which the operator is subject

- Risk of infection due to sample handling.
- Risk of parts expelled.
- Electric risk.

1.2 Personnel qualification

This equipment can only be used by personnel who have been properly qualified. This equipment can only be used by personnel who have read and understood these instructions or who has been properly qualified in the operation of this equipment.

Handle of biological specimens requires special knowledge.

1.3 Use with hazardous substances

When you know that the materials used in the centrifuge are toxic, radioactive or contaminated by pathogenic microorganisms, the following safety measures must be applied:

- Containers specified by the manufactures will be used, to ensure sealing with those materials.
- When materials of Risk Group II are handled according to the identification «Biological Safety for Laboratories Manual» edited by the World Health Organization, act accordingly to the instructions.

- In case of materials belonging to a higher group, it requires a higher level of protection.
- It is forbidden to centrifuge the following materials in these centrifuges:

- Flammable or explosive materials.
- Materials that chemically react with enough violence to cause a hazard.

2 General information

Specially designed for blood bank applications, where repeatability, precision and reliability of: speed, time, temperature, acceleration and deceleration are significant factors in obtaining perfect separations of the distinct blood fractions.

Serology investigations, molecular biology. Serum separation, pharmaceutical applications.

Lid with double safety lock.

Digital electronic circuit controlled by microprocessor for controlling all operating parameters:

- Speed
- Acceleration
- Brake
- RCF
- Temperature
- Time
- Possibility to set 10 programs in memory
- Alarm of operation

Drive unit: high frequency induction motor maintenance-free (no brushes nor collector).

Rotors, sealed buckets and adapters properly identify.

Windbreak hood with detector which reduces friction of the inner buckets and improves cooling.

Inner insulation which absorbs noise and vibrations.

Elastic suspension of the motor which guarantees a smooth operation with radial dampers.

Hermetically sealed compressor fixed on anti-vibration mounts with a recirculating evaporator unit around the interior of the centrifuge bowl.

Emergency stop brake: totally independent circuit that stops the equipment in case of emergency.

3 Packing list

The standard equipment consists of the following components:

- Centrifuge
- Rotor fixing wrench (1x)
- Allen wrench (1x)
- User manual (1x)
- Drain nipple
- Drain silicon tube $\varnothing 8 \times 14$ (2 mts, cod. 46041)

4 Technical Specifications

Code	7001491					
Power supply ($\pm 10\%$)	400V III + N 50Hz					
Motor + cooling consumption	2900 + 1600 VA Total 4500VA					
Maximum volume	200 vacutainer tubes (11ml) $\varnothing 16 \times 110 = 2200\text{ml}$					
Maximum admissible density	1.2 Kg/dm ³					
Speed range	300...3000 rpm					
Stability	± 5 rpm					
Resolution	1 rpm					
RCF range	29-2900					
RCF stability	± 5 g					
Resolution	4 g					
Temperature range	0-40°C					
Stability	$\pm 2^\circ\text{C}$					
Resolution	1°C					
Time range	99h 59 min					
Acceleration ramps	nº	1	2	3	4	5
rpm/minute:		4	6	12	17	24
Braking ramps	0 Free (unbraked)					
	1-7 arithmetic line					
	8-15 quadratics					
Measurements (Height x Width x Depth)	1050 x 960 x 960					
Fuses	Power: 20A					
	Control: 4 A					
	Type: gl-gG					
Weight	429 Kg.					
Sound level:	Less than 70 dBA					
Operating range:	Continuous intervals of 10 minutes between cycles.					
Environmental conditions:	(EN61010-1)					
	Indoor use					
	Altitude up to 2000m					
	Room temperature between 5°C and 40°C.					
	Maximum relative humidity of 80% for temperatures up to 31°C, linearly decreasing to 50% of relative humidity at 40°C					
Over voltage level:	Category II					
Pollution degree:	2					

5 Installation

WARNING !

The installation instructions indicated in this manual must be applied to guarantee the machine safety.

5.1 Transport

The centrifuge comes with four lifting supports, for positioning the unit. Note that the unit weights 429 Kg.

5.2 Location

Place the centrifuge on an horizontal flat and stable surface.

The surface must allow adherence of the centrifuge legs.


This equipment can cause vibrations affecting other nearby equipments.

Extreme conditions of temperature, humidity and contamination could negatively affect the machine security. Avoid direct exhibition to sunlight.

For normal operation, ambient working temperature should be between 10°C and 40°C, with a maximum humidity level of 80%.

Do not install the centrifuge in areas where flammable liquids are stored or special protection areas.

Leave a free space of 30 cm around the centrifuge to permit refrigeration of the equipment. This area is called **SAFETY ZONE**

SAFETY ZONE

Leave a free space of 30cm around the centrifuge.

With the centrifuge operating:

- Do not place hazardous substances in the safety zone.
- Personnel must be outside the safety zone.

DO NOT INSTALL THE EQUIPMENT IN EXPLOSIVE ATMOSPHERES OR WHERE SOLVENTS, OR OTHER PRODUCTS THAT MAY FORM EXPLOSIVE OR FLAMMABLE MIXTURES, ARE USED OR STORED.


Do not use the centrifuge:


If the speed reading (RPM) never reaches stabilization.
(in this case, actual speed could be very different from that indicated)

Contact with J.P. SELECTA's technical service.

5.3 Initial setup


Once the unit is installed, remove the lifting supports and put the hole taps supplied into place.

Adjust the 4 levelling rubber feet supports so that the unit is raised and do not rest on the wheels, and adjust until the machine is levelled.

Remove the caps to open manually.

Open the door by gently pulling the cables.

Remove the packaging protections from inside the chamber. Check that there is no packing material remaining in the unit.

5.4 Drain connection


To prevent water from accumulating condensation inside the tank, there is a back drain.

Connect the supplied nipple by screwing it onto the record. Connect the silicone tube to the drain of your installation.

5.5 Head location

Before placing the swing out rotor, check that the locking ring wedge is in the right position.


Place the rotor on the spindle, ensuring that the keys and their lodgings are perfectly aligned and that it fits perpendicular to the axis. The rotor should slide easily, **it is not necessary to use a mallet or hammer.**

Place and screw the locking nut in place and tighten using the spanner provided.

Tighten the Allen screws on the locking nut block.

Periodically examine the status of the rotor and the buckets to check and detect possible corrosion, micro fractures and mechanical deformations.

Verify that the adapters and the extractor buckets freely rotate and can perform the reciprocal movement without impediments.

5.6 Power supply connection

Installation of the mains must be performed by qualified personnel.

Connect the equipment to the main power, checking that the power voltage is as indicated on the machine's specification plate.

Install an emergency disconnect switch (3P+N three-phase 20Amp) to cut power supply in case of malfunction.

This switch must be installed at a minimum distance of 10m from the centrifuge (preferably outside the room where the centrifuge is). In some installations it may be advisable to install a 3P + N and 40A differential circuit breaker. This switch must be at a minimum distance of 10m from the centrifuge (preferably outside the room where the centrifuge is located).

Sections and colours of the connection cable:

Grey cable 6mm ² :	Phase
Brown cable 6mm ² :	Phase
Black cable 6mm ² :	Phase
Blue cable 6mm ² :	Neutral
Yellow/green cable 6mm ² :	Earth

It is supplied with a CETAC 16A 3P + N + T plug. Make sure that the installation has a CETAC base of 16A 3P + N + T.

Do not use the equipment if it is not connected to earth. Make sure that the rotation direction of the motor is counterclockwise.

6 Operation

6.1 General conditions during operation

Consider the following precautions for a safe use of the centrifuge.

6.1.1 General information

- Do not use the machine to centrifuge liquids capable of forming explosive, flammable mixtures, with hazardous products or materials that chemically react with virulence.
- Do not move the centrifuge in operation.
- Never open manually the lid while the centrifuge is running.

6.1.2 Maximum density of liquids to be centrifuged

The maximum speed recorded on the rotor is calculated to centrifuge fluids at a maximum density of 1.2gr/cm³.

To centrifuge liquids with a higher density than 1.2gr/cm³, the rotor's maximum speed must be reduced according to the following formula:

$$\text{Rpm} = \text{Max. Rpm} \times \sqrt{1/\text{density}}$$




where:

Rpm = Max speed of fluids higher than 1.2gr/cm³

Max Rpm = Maximum speed indicated on the rotor

Density = gr/cm³.

OPERATING TIME:	00004h 10m
ROTOR TIME:	00004h 10m
NORMAL BUCKET TIME:	00004h 10m
DOUBLE BUCKET TIME:	00004h 10m

Accumulated working times of components with limited service life

Main menu

6.1.3 Filling and location of the buckets

Once the buckets are filled, and before placing them in the rotor, please check that their weight do not differ in more than 15 grams.

Do not start the centrifuge if the buckets are not placed.

And under no circumstances start the centrifuge with the rotor loaded in an unbalanced way.

6.1.4 Emergency stop

If you have to stop the machine, quickly activate the front panel button (18).

Note: Do not use this button as an usual form of stoppage. Please only activate it in case of emergency.

6.1.5 Safety measures for the staff

Do not lean on the centrifuge.

Do not stay in the safety free space more than the strictly necessary.

Do not leave potentially hazardous materials inside the security free space.

7 Operating menus

Press the ON/Off switch.

The hour counter of the centrifugal elements appears with limited service life.

Rotor

Normal bucket.

Double bucket.

After a few seconds, the main menu screen appears.

You can select from the main menu, by pressing on the touch screen:


SETTINGS: To update date, time, language, data download, software updates, life counters initialization.

PROGRAMMING: To establish programs for samples treatment:

- Times
- Speed
- Temperature.

WORKING: To select one of the stored programs.

Return to main menu

Access authorization for Firmware update.

7.1 Settings

Set on this screen the following settings:

LANGUAGE: Spanish, English & French.

DATE & TIME.

TEMPERATURE UNITS in °F or °C.

ROTOR: Enters the rotor used in the current operation, essential for braking and acceleration ramps application.

ENTER / CHANGE PASSWORD: Password access for software updates or restart life counters.

DOWNLOAD LOGS: Download Time/ temperature/ RPM/RCF logs of any executed cycle.

7.2 Programming

In this menu the program parameters are set.

- Rotor: Enter the rotor identification. This limits the programming to avoid dangerous situations.
- Temperature: It sets the working temperature.

A program is composed of «phases» of:

- Acceleration (slope)

Approximate acceleration values:

SLOPE 0:	04 rpm x seconds
SLOPE 1:	06 rpm x seconds
SLOPE 2:	12 rpm x seconds
SLOPE 3:	17 rpm x seconds
SLOPE 4:	24 rpm x seconds

- Speed (rpm) and time (minutes, seconds) and:


- A final braking phase.
Select the program and press «EDIT».

Program

Mark as selected

Select the program

Add phase

Delete phase

Approximate linear braking ramps.

0 - Free braking (can take 30' to stop)

1 - 28 rpm x seconds

2 - 19 rpm x seconds

3 - 14 rpm x seconds

4 - 11 rpm x seconds

5 - 09 rpm x seconds

6 - 08 rpm x seconds

7 - 07 rpm x seconds

Approximate quadratic braking ramps

8 - 160 rpm x seconds


9 - 28 rpm x seconds

10- 19 rpm x seconds

11- 14 rpm x seconds

12- 11 rpm x seconds

13- 09 rpm x seconds

**7.3 Run a program (WORKING)**

From the main menu, press on «WORKING».

Program selected in paragraph 7 will be running.

Press on «RUN» to start the program.

The program starts and displays the following information:

- Current speed
- Time elapsed (minutes, seconds)
- Temperature in the bucket.

Press on «RPM» so that the display shows «RCF».

RPM: Revolutions per minute

RCF: Number of «g» (acceleration of gravity) to which the sample is subjected to.

The program alerts when it arrives to the end.

7.4 Pre-programmed programs

Pre-programmed program:


Attention: ALL programs are modified by the user.

Program 1:

Speed:	2650 r.p.m.
RCF:	1996 g
Acceleration:	5
Braking:	5
Time:	12 minutes
Temperature:	22 °C

Program 2:


Speed:	2550 r.p.m.
RCF:	1778 g
Acceleration:	9
Braking:	5
Time:	6 minutes
Temperature:	22 °C

Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.

8 Maintenance

Before removing the centrifuge cover to handle inside, disconnect it from the mains.


Manipulation of internal electronic circuits of the centrifuge by unauthorized personnel could cause irreparable damages. Take the equipment to one of J.P. SELECTA's technical service.

Please check on a monthly basis the fixing of the swing-out rotor. For inspection, proceed as follows:

1. Remove the protective cover of the rotor.
2. Using the Allen key supplied, loosen the screws (see picture).
3. Using the spanner provided, check the tightening of the nut by pressing in the clockwise direction (see picture).
4. Periodically check the rotor and the buckets status in order to detect possible signs of corrosion or mechanical deformations.

8.1 Visual inspection every 3 months


- Centrifuge lids locking status.
- That the screws of the machine continue very tighten.
- Deformations and visible fractures.
- Side walls.
- Cover gasket status.
- Continuity of earth connection.
- Levelling: legs are not loosening.

8.2 Software maintenance

Software maintenance functions.

- Insert new rotors or accessories.
- Update Firmware.
- Initialize rotors and buckets service life.
- Calibration and logging programs.
- Add programs for specific treatments.

These functions are only available under the supervision of the technical service.

8.3 Cleaning

For cleaning of different parts of the equipment, we recommend to use the following products:

- Cleaning of stainless steel: alcohol.
- Cleaning of covers and plastics: alcohol with cotton or a non-abrasive cloth.
- Clean with a brush or a condenser fins vacuum located on the back of the equipment.

Cleaning and decontamination of spills inside the centrifuge:

The user is responsible to make the appropriate decontamination if dangerous materials have spilled on or inside the equipment.

In case of samples spillage in the centrifuge chamber, clean it with a suitable biological disinfectant.

Maintenance operations:

The following components are the operator's protection against the impact of projected parts in case of accident, thus requiring a periodical inspection (quarterly):

Door hinges, closing and chassis: visually examine the appearance of cracks, deformations, rust spots, etc ... that may indicate decreased resistance of the material.

Rubber feet: verify that the material (rubber) maintains its flexibility and adherence to the table features.

8.4 Periodic official review

It is recommended to carry out an inspection by an official technical service before reaching 1000 operating hours, to avoid possible more serious breakdowns.

J.P. Selecta, s.a.u. is not responsible if these instructions are not followed periodically.

Doc: Q.00031.00-C

MACROFIBER BILT BLOOD CENTRIFUGE FUNCTIONAL REVIEW

Brand: **Machine Code:**
Model: **Work Center:**
Model Number: **Location:**

Preliminary equipment inspection

	Correct	No correct	Action
- Integrity of the network cable			
- Integrity of the control panel			
- Verification of acoustic and light alarms			
- General condition of the equipment			
- Machine operating hours			
- Head operating hours			
- Adapter cups operating hours			

Checkpoints

	Good	Bad
1. Visual inspection of:		
1.1- HEAD: There are no visible deformations in the vessel fixing part.		
1.2- HEAD: There are no visible cracks.		
1.3- BUCKETS: It has no bumps, deformations or wear.		
1.4- HEAD FAIRING: It does not have bumps, deformations or wear.		
2. Visual inspection of the IMHUBER:		
2.1- There are no cracks, they do not wear, there are no deformations.		
2.2- Check that the door scraper cleans and cannot be opened during operation.		
2.3- Check Mags (rotating cones, chimes and handles).		
3. Internal inspection:		
3.1- Disconnect the side cover.		
3.2- Check the resistance between metal parts of the enclosure and the grounding (from WGround).		
3.3- Visually check the state of the "brake-block" and the counterweight. (Brake-blocks: braked, loose, barely attached, rubber deformation, counterweight barely attached).		
3.4- Check that the centrifuge is level. Level with the legs, if necessary.		
4. Operating tests:		
4.1- Fill the centrifuge by using the most exact program, or the one with the highest speed. (With cold component activation)		
5. Measure and register:		
5.1- Supply voltage (V) VR= VS= VT=		
5.2- Consumption (A) IR= IS= IT=		
5.3- Check that rotation direction is well-characterized.		
5.4- Put the IM on.		
6. Check that there are no vibrations, rashes, out of the ordinary. MEM value measurement:		

Patterns used, models: _____

TEST RESULT

VALID TEST

NOT VALID TEST

REMARKS: _____

Technician: _____

Date: _____