**2000368**

**R-3L**

**2000369**

**R-8L**

**HORNO de MUFLA**

**MUFFLE FURNACE**

# Indice

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# Index

<b>1. General Information.....</b>	<b>12</b>
<b>2. Packing list .....</b>	<b>12</b>
<b>3. Technical specification .....</b>	<b>12</b>
<b>4. Information about safety in the use of the equipment.....</b>	<b>13</b>
4.1 Risks the operator is subjected to .....	13
4.2 Staff qualification.....	13
4.3 Use with hazardous substances .....	13
4.4. Maintenance.....	13
4.5. Warranty.....	13
<b>5. Equipment Description .....</b>	<b>14</b>
<b>6. Installation .....</b>	<b>14</b>
6.1 Location.....	14
6.2 Connection to the mains .....	14
<b>7. Operation.....</b>	<b>14</b>
7.1 Inicio rápido para trabajar a una temperatura .....	15
7.2 Setting screen .....	16
7.3 TASK screen .....	18
7.4 PROGRAMMING screen .....	19

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


## 1. General Information


- Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributors immediately.
- Do not install or use the equipment without reading this handbook first. This instruction manual must be attached to the equipment and must be available for all the users.
- If you have any doubts or enquiries, contact with J.P. SELECTA's technical service.
- Do not use the equipment in explosive, flammable or corrosive atmospheres.

## 2. Packing list

Product code	2000368	2000369
Size	3L	8L
Thermal shelf:	x1	x1

## 3. Technical specification

	2000368	2000369
Usable space:	3L	8L
Chamber dimensions: High/wide/depth (mm)	115/125/200	130/200/300
External dimensions: High/wide/depth (mm)	430/340/470	510/440/560
Power supply:	230 VAC 50Hz	
Electrical power (W):	1700	1800
Weight: (Kg)	18	28
Fuse:	Type: 10x5mm	
Temperature range:	From ... T <sup>a</sup> amb + 5°C ... to ... 1100°C	
Temperature resolution:	1°C	
Timed operation:	week days, hh.mm	
Temperature profiles:	Up to 10 profiles 6 segments. Segment from 1minute to 99hour.	
Data recording:	USB memory device	

### Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.

## 4. Information about safety in the use of the equipment

**Hazard risk**

**Electric risk**

**Risk of high temperature contact**

### 4.1 Risks the user is subjected to

- Risk of burns
- Electrical risk

### 4.2 Staff qualification

This equipment may only be used by properly trained personnel who have read and understood these instructions or have been properly trained in this equipment operation.

### 4.3 Use with hazardous substances

Do not treat flammable or explosive materials or those which chemically react with violence.

## ATTENTION!

### Burning risk

Take proper precautions  
(use of thermal gloves, clamps,...)  
to extract samples from inside the furnace.

**It can be at a high temperature.**

### 4.4. Maintenance

Unplug the equipment from the mains before removing the furnace's cover to handle inside.

**The manipulation of the equipment internal electronic circuits by unauthorized personnel can cause irreparable damage. Take it to one of J.P.SELECTA's authorized technical services.**

It is recommended to clean the accessible parts of the furnace with a cloth and a cleaner, or a similar solution. Do not use aggressive solvents.

### 4.5. Warranty

This product is guaranteed for one year. It does not cover damages caused by an incorrect use or causes beyond the control of J. P. SELECTA, s.a.u.

Any manipulation of the equipment by personnel not authorized by J.P. SELECTA, s.a.u. automatically cancels the guarantee.


Fig 1. Temperature profile

## 5. Equipment Description

Furnace muffles are used in sample thermal treatments, for its desiccation, cured and sterilization.

Fan assisted circulation gets a great temperature uniformity and a quick temperature rise.

## 6. Installation

Unpack the equipment and visually check that the contents are in good conditions. Keep the package for some days.

### 6.1 Location

Place the equipment on a flat level stable surface, suitable for the equipment weight.

To work in a comfortable way, please leave a free space of 50cm around the equipment and a minimum of 20cm at the rear part.

### 6.2 Connection to the mains

Plug the equipment to a near socket suitable for the equipment power. Use the connection cable supplied or a similar one.

For your safety, the equipment must be earthed.

Before plugging the machine, check that the voltage supply and the frequency corresponds to the one indicated in the equipment characteristics plate.

## 7. Operation

All the equipment functions are set with the touch screen. Gently press with your finger in the center of each icon.

When starting the computer, the starting screen will be shown (Fig 7.1)

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Fig 7.1. Initial screen

- **SETTINGS:**

- Time and date: menu to set the date and time.
- PID and Autotuning: menu to set the parameters of temperature control.
- Downloading recorded data.
- Temperature correction.
- Temperature alarm limits.

- **PROGRAMMING**


- Allows you to program the temperature profile, by defining temperatures, times and ramps.

- **TASK**


- Allows you to select and launch the operating temperature or a temperature profile.


### 7.1 Quick start to work at one temperature


**1** Press «TASK».

**2** Press over-temperature panel.

**3** Press   to set temperature.


**4** Press on  to start.

## 7.2 Setting screen

Icons meaning:

### 7.2.1 Set time and date

Press on the time/date icon. Press successively on the box to select the day, month, year, minute & second adjustment.

Adjust by pressing on the icons Increase and Decrease.

### 7.2.2 Set parameters of temperature control and «Autotuning»

For an optimal use of the furnace, P, I, D parameters of temperature control should be the appropriate value.

With a value of the appropriate P, I, D parameters, the furnace will reach the temperature in the shortest possible time and without exceeding the set temperature. However, those values can be quite different at high and low temperatures.

The furnace comes with the following values: P=10; I=100; D=25 which are appropriate for most jobs.

These values can be modified, but it's advisable to perform an «Autotuning» when you need a better accuracy.

The «Autotuning» automatically adjusts the optimal parameters for a given temperature and a load.

Press on the PID/Autotuning icon. Press successively on the box to select the P, I, D values adjustment and modify the value by pressing on the Increase and Decrease icons.

**Attention:** It is not recommended to change these parameters. Use the «Autotuning» option to automatically adjust these values.

Press on the PID/Autotuning icon.

Select the «Autotuning» temperature with:    
(for example, the temperature most used)

Start «Autotuning» by pressing on the icon: 

If for some reason you want to stop, press on: 


Fig 7.2. Communication

### 7.2.3 Communication type

The equipment can store readings from time to time and send them to a computer or store them in a USB memory pendrive.

- «**PORT**» Set the readings storing mode:
  - RS232: Send the temperature readings to a computer by a serial cable (DB9 type).
  - USB: Store the readings in a USB memory pen drive.
- «**INTERVAL**» Select with the register time interval and the data sending.
- «**MODE**» Select the data to be stored:
  - Start: Save all temperature readings.
  - Set point: Save temperature readings, since the set point temperature has been reached.
  - PC: Send data, it is not stored.
- «**SEND**» Number of stored data, pending from sending.


Fig 7.2. Temperature alarms

### 7.2.4 Set the temperature and the temperature alarms

It is recommended to verify calibration of temperature with some regularity (for example annually) by using for instance a thermometer or a temperature pattern.

If it is necessary to correct the reading, as a result of calibration, one must enter a correction.

To do this, press on this icon:

A box indicating the current correction value appears.

Adjust with keys until you get the indicated value on the temperature pattern.

Applied correction will appear in brackets.

#### 7.2.4.1 Temperature alarms

By pressing successively on the temperature correction box, one may enter an over-temperature and low-temperature limit. Both alarms activate a message on screen.

### 7.2.5 Temperature units

Press on the temperature units menu. Press successively to change between Celsius or Farenheit degrees. It is also possible to change the ramp units de °C/min or °C/hour.

### 7.2.6 Set language

Press on the language icon. Press successively on the box to select among Spanish, French or English.

Adjust by pressing on the Increase and Decrease icons.

### 7.3 TASK screen

There are two operation modes:

- Working at a unique temperature.
- Working according to a pre-programmed temperature profile.

Press on  to change the operation mode.


Fig 7.3.1. Run screen operating at one temperature.

#### 7.3.1 Operation at a temperature

- Select the set point temperature by pressing on the temperature indicator.
- Select a temperature increase ramp if you want to increase it in a slowly way.  
The value 0°C/min makes that the ramps increases as fast as possible.
- Select an operating time. This time starts to count down when the temperature reaches the set point.  
The 0°C value sets a continuous operation.

- Start by pressing on this icon: 

#### 7.3.2 Operation with a temperature profile

- Temperature profiles must be previously programmed from the «PROGRAMMATION» display.

- Select a profile (from 0 to 9) by pressing on keys   .

- Start by pressing on the icon:  .

- In the temperature box (top-left corner) the current temperature and the set point are indicated.


Fig 7.3. Run screen operating with temperature profile.



Fig 7.3.2 Working screen.


### 7.4 PROGRAMMING screen

Go to PROGRAMMING screen to:

- Set the temperature profiles.
- Set the timed operation for a profile.  
(For instance, you can set a time an a week day to start a profile. As: start, profile n<sup>er</sup> 1 on Sunday at 19h00')

#### 7.4.1 Programming a temperature profile

A temperature profile is a sequence of steps, each step could be a ramp or a dell.

To go from the end of a «dell» to the following «dell» could be done in two ways:

- With a «temperature ramp»  
(go from 100°C to 200°C following a ramp at 2°C per minute).
- As fast as possible.

#### Profile parameters input:

- Select the profile number between 0 and 10 with .
- Press on the parameter to be modified. (Temp, ramp, length ...)
- Press on these icons: to enter the value.
- Modify another parameter.

#### Special values:

- If you enter value **OFF** in **temperature** (1st column) (you have to press Increase key) the next steps are disabled and not executed.
- If you enter **INF** in **Step length**, the furnace remains permanently at this temperature.
- If you enter **OFF** in the **ramp**, the furnace makes the temperature rise at the fastest speed possible.

#### Starting of temperature profile

- To set the profile starting time, press on «←:→».
- Enter time by pressing in **hh:mm** format.
- Press on the days of the week when the profile should be run.
- Accept the programming by pressing on .