DRY-BIG	2002961	DRY-BIG	2002972
DRY-BIG	2002962	DRY-BIG	2003721
DRY-BIG	2002971	DRY-BIG	2003741

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1. General Information

- Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributors immediately.
- Do not install or use the equipment without reading first this handbook. This instruction manual must be attached to the equipment and must be available for all users.
- If you have any doubts or enquiries, contact with J.P. SELECTA's technical service.
- Do not use the equipment in explosive, flammable or corrosive atmospheres.

2. Packing List

Models:	2002961	2002971			
	2002962	2002972	2003721	2003741	2003741 (1 door)
Inner volume:	216L.	288L.	400L.	720L.	360L.
Shelf: (2x Oven)	2000062	2000062	2000063	2000064	475x600mm

Note: More shelves can be purchased separately.

3. Technical Specification

Models:	2002961	2002971			
	2002962	2002972	2003721	2003741	2003741(1 door)
Inner volume:	216L.	288L.	400L.	720L.	360L.
Inner dimensions: height/width/depth(mm)	600/600/600	800/600/600	1000/800/500	1200/1000/600	1200/500/600
External dimensions: height/width/depth (mm)	870/1120/840	1070/1120/840	1280/1320/740	1500/1520/800	1500/1000/800
Voltage:	230/400 V-3ph / 230 V-1ph AC 50Hz				
Electric power (W):	4000	5000	5250	6000	6000
External fuse:	Removable (Type: 10x5mm)				
Temperature range:	From 40°C to 250°C				
Display resolution:	1°C				
Temperature stability:	±0.4°C to 100°C				
Temperature uniformity:	± 2°C to 100°C				
Set point error:	± 2.5%				
Programmable operating time:	From 0h 1m to 99h				
Programmable timeout:	From 0h to 99h				
Data record output:	Serial RS-232				
Oven's remote control:	Serial RS-232				

4. Accesories

Digital microprocessor TFT:	2000010
-----------------------------	---------

5. Equipment Description

DRY-BIG ovens are used in sample thermal treatments, for its desiccation, cured and sterilization. Especially when samples can be treated by fan assisted circulation.

Fan assisted circulation gets a great temperature uniformity and a quick temperature rise.

The temperature controller includes the following functions:

- Temperature selection up to 250°C in increments of 1°C.
- Timeout selection, after which the oven will start.
- Operating time selection, after which the oven will stop.

The oven's remote control from a PC allows temperature profiles creation, data recording and control panel lock, by means of a RS-232 serial interface communication and a temperature recording in a USB pen-drive.

6. Installation

Unpack the equipment and visually check that the contents are in good conditions. Keep the package for some days.

¡ATENCIÓN!
THREE PHASE POWER SUPPLY.

1. ENSURE THAT THE DEVICE IS RATED AT MAIN POWER SUPPLY.
2. CHECK MOTOR FAN TURNING WAY. IF NOT MATCH THE INDICATED WAY SWAP TWO PHASES.

6.1 Location

Place the equipment on a flat level stable surface, suitable for the equipment weight.

To work in a comfortable way, please leave a free space of 50cm around the equipment and a minimum of 20cm at the rear part.

6.2 Connection to the mains

Plug the equipment to a near socket suitable for the equipment power. Use the connection cable supplied or a similar one.

For your safety, the equipment must be earthed.

Before plugging the machine, check that the voltage supply and the frequency corresponds to the one indicated in the equipment characteristics plate. Also, check that motor turns on the indicated way.

7. Operation

Unpack the equipment and visually check that the contents are in good conditions. Keep the package for some days.

7.1 Introduction

DRY-BIG basic operation is very simple: select the operating temperature (See 7.2) and press the «START» key.

The oven will start and a few minutes later, the inner temperature will be stabilized at the operating temperature.

To get the most out of the oven, it will be interesting you know all its additional features:

- Adding a timeout to the starting and an operating time.
- Adjusting the temperature measure in a precise way.
- Recording temperature readings in a USB pen-drive.
- Controlling the oven from a PC via a RS-232 serial.

ATTENTION
Burning risk

Take proper precautions (use of thermal gloves, clamps,...) to extract samples from the inside the oven.

It can be at a high temperature.


Fig 1. Display:

- a. : Heater On / Off
- b. : Operation time
- c. : Timeout
- d. : Oven running

7.2 Setting working temperature

- a. The «RUN» led must be off (Press «START/STOP»).
- b. Press and without releasing, press or until the display shows the desired operating temperature.

7.3 Setting operation time

This time begins to run from the moment the oven reaches the operating temperature. After this time, heater will stop. But ventilation won't stop.

It is not necessary to set an operating time. If this time is «0», the oven will run continuously (This value is set by default).

- a. The «RUN» led must be off (Press «START/STOP»).
- b. Press and without releasing, press or until the display shows the desired operating time (in hours and tenths of minute).
- c. After this time, and after starting, the oven will reach the operating temperature. At this time, the led will be lit indicating there's a programmed operation time.

7.4 Setting timeout

This time will begin counting from the moment the oven reaches the operating temperature. (See Fig. 2)

After this time, the heater will stop. However, the ventilation won't stop.

It is not necessary to set a timeout. If this time is «0», the oven will run when pressing «START» (This value is set on default).

- a. The «RUN» led must be off (press «START/STOP»).
- b. Press and, again, and without releasing it, press or until the display shows the desired operating time (in hours).
- c. When starting the oven , «TIMER» led is lit, indicating a timeout before it starts.


Fig 2. Timeout and operation time.

7.5 Calibrating (adjusting) the oven

DRY-BIG ovens come factory adjusted. However, if there is a deviation seen of the indicated temperature in the display regarding the measured temperature inside a reference thermometer, you can adjust the temperature controller.

Initially, try to correct the error by adding an «Offset» to set the temperature. See 7.5.1

For example:

The oven's display indicates 100°C.

Reference thermometer indicates 99°C.




An offset of 1.0 must be **added**.

7.5.1 Setting «Offset»

The «Offset» allows adding a quantity (positive or negative) to the measured temperature value, moving the reading.

An «offset» between -9.9 and 9.9 can be set.

- a. The «RUN» led must be off (press «START/STOP»).

- b. Press  and without releasing, press  for 5 seconds.
- c. Release  and .
- d. Press  once.
- e. Press  or  until the display shows the desired «Offset».
- f. Press  four times.

7.5.2 Setting the «Temperature coefficient»


The temperature coefficient adds a factor that multiplies the controller original reading.

A coefficient between 800 and 1000 can be set (The value 1000 is shown in the display as 000)

The coefficient default value is 1000.

This coefficient applies according to the formula:

$$T^a_{display} = T^a_{probe} \times (\text{coefficient} / 1000).$$

- a. The «RUN» led must be off (Press «START/STOP»).
- b. Press  and without releasing, press  for 5 seconds.
- c. Release  and .
- d. Press  twice.
- e. Press  or  until the display shows the desired «coefficient».
- f. Press  three times.

7.6 Setting the safety thermostat

DRY-BIG ovens are equipped with a safety thermostat (See Fig. 4). This device interrupts the oven's heater if it exceeds the set value.

It is intended as an additional safety device, in case the temperature electronic controller stops operating properly. In this case, its action prevents that the temperature could rise to excessive values.

Another use of the thermostat is to protect valuable samples, by setting it between 5 and 10 °C above the operating temperature.

The thermostat comes factory set to maximum: 100% (at full clockwise).

To set the safety thermostat shutdown temperature, follow the next instructions:

- a. Set the thermostat at full clockwise (100%).
- b. Be sure that the safety thermostat is reset by pressing «RESET». The red lamp must be off (See Fig. 4).



Fig 1. Display:

**Safety thermostat:**

1. «RESET» button.
2. Thermostat lock.
3. Release temperature setting button.
4. Safety thermostat lamp:

Off: Reset.

On: Released (Oven's heater cut)

Fig 4. Safety thermostat

c. Set the safety thermostat release temperature in the oven's control panel.

This temperature must be at least 5°C higher than the operating temperature.

d. Press «START» to start the oven.

e. When the oven's temperature is stabilized, slowly turn the safety thermostat knob clockwise until it releases (you can hear a «click») and the red lamp lights. The oven's heater switches off permanently.

Lock the safety thermostat knob.

f. Stop the oven by pressing the «START/STOP» key. Open the door and let it cool.

g. Press the «RESET» button (1) and the red lamp (4) must be off (See Fig. 4).

h. At this moment, the oven is ready to work.

NOTE: The safety thermostat setting process is long, that's the reason why it is only practical for users who usually work at the same temperature.

7.7 Serial output type selection

DRY-BIG are equipped with a serial outlet DB-9 of 9 pins. This outlet could be for logging data and controlling the oven from a computer.

To select the serial output type:


a. The «RUN» led must be off (Press «START/STOP»).

b. Press  and without releasing, press  for 5 seconds.

c. Release  and .

d. Press  three times.

e. Press  or  until the display shows «USB» or «PC».

f. Press  twice.

7.8 Recording interval selection

The recording interval is the time between readings registered in the «USB» mode:

a. The «RUN» led must be off (Press «START/STOP»).

b. Press  and without releasing, press  for 5 seconds.

c. Release  and .

d. Press  four times.

e. Press  or  to reduce / increase (interval minutes).

f. Press  once.

7.9 Serial identification «ID number» selection

The «ID number» identifies the equipment within a serial network.

You can control several ovens from a computer. Everyone has a different «ID number».

All the ovens have by default «ID number» = 1.

a. The «RUN» led must be off (Press «START/STOP»).

b. Press  and without releasing, press  for 5 seconds.

c. Release  and .

d. Press  once.

e. Press  or  to reduce or increase the «ID number».

f. Press  four times.

7.10 Error messages

Error messages shown in the display:

rtd The temperature probe reading is not correct. Check the probe connection.

e2p Memory reading or writing error in the electronic board. Press  to restart.

8. Temperature Record

DRY-BIG ovens include a temperature data recording system from the oven's inside.

9 pin female cable; 9 pin female:


Fig. 8.1: Null- modem cable

9 pin female cable; 25 pin male:


Fig. 8.2: Null- modem cable

Record via RS-232 with software you can request from JP Selecta s.a.u.

The material needed is:

- a. Communications software RS-232 from the oven.
- b. Serial communications cable, not supplied with the oven.
- c. PC computer with WINDOWS XP and RS-232 connection.

Installation:

a. You can acquire the components to build the serial communication cable it in any specialized shop. Please follow diagrams in figures 8.1 and 8.2 to build the cable.

Do not hesitate to contact J.P. SELECTA, s.a.u if you cannot get the cable.

- b. Install the requested software file «**DRY-BIG**».
- c. Connect the serial cable in a free port of the computer (normally COM1 or COM2).

9. Oven's Remote Control

To control the oven from a remote computer, you can build a software, from the serial commands sent to the temperature controller (**Note:** advanced programming knowledge is needed). Temperature profiles of different steps, ramps, etc. can be programmed in this way.

Please contact J.P. Selecta,s.a.u. for more specific details.

10. Maintenance

Unplug the equipment from the mains before removing the oven's cover to handle inside.

The manipulation of the equipment internal electronic circuits by unauthorized personnel can cause irreparable damage. Take it to one of J.P.SELECTA's authorized technical services.

It is recommended to clean the accessible parts of the oven with a cloth and a cleaner, or a similar solution. Do not use aggressive solvents.

List of more usual spare parts:

	2002961	2002971		
	2002962	2002972	2003721	2003741
Regulation plate	29091	29091	29091	29091
Contactors 7.1	13006	13006	13006	13006
Contactors LS	13008	13005	13005	13005
Condenser	15408	15408	15408	15408
Fuse 20Amp.	15410	15410	15410	15410
Connector	15503	15503	15503	15503
Fuse 4Amp.	15555	15555	15555	15555
Opening Selector	20040	20040	20040	20040
Silicone gasket	21027	21027	21027	21027
Indicator green lamp	22017	22017	22017	22017
Indicator red lamp	22021	22021	22021	22021
Motor	24010	24010	24010	24010
Fan	30334	30335	30335	30335
Thermal relay	36009	36009	36009	36009
Heater element	39046	39046	39046	39046
Pt100 probe	43076	43076	43076	43076
Triac	45055	45055	45055	45055

11. Warranty

This product is guaranteed for one year. Warranty does not cover damages caused by an incorrect use or causes beyond the control of J. P. SELECTA, s.a.u.

Any manipulation of the equipment by personnel not authorized by J.P. SELECTA, s.a.u. automatically cancels the guarantee.