HIGHTEMP

2001406

**ESTUFAS PARA ALTA TEMPERATURA
HIGH TEMPERATURE OVEN**

Indice

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Index

1. Safety	10
1.1 Hazards to which the operator is exposed	10
1.2 Staff qualification.....	10
1.3 Proper use of the equipment	10
1.4 Modifications	10
2. Warranty	10
3. Packing List	10
4. Equipment Description	11
5. Technical Features	11
6. Installation	11
6.1 Location.....	11
6.2 Power Supply Connection	11
7. Operation	12
7.1 Control Panel.....	12
7.2 Starting Up.....	13
7.3 Temperature Controllers	13
7.4 Motor Thermal Protection	13
8. Maintenance	14
8.1 Electrical Panel Components List with its Spare Part Numbers	14
8.2 Electrical Diagram for Oven 2001406 (Diagram Number: E-20093.10)	15
<i>Aviso a clientes / Notice to customers</i>	16

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


1. Safety

HIGHTEMP ovens have the suitable measures for its use in laboratories. However, the user must be informed of some of the risks.

1.1 Hazards to which the operator is exposed

It identifies risky situations and the safety measures to be taken.

Electrical shock risk

There's a risk of electrical accident. Please respect the indicated instructions to carry out the described operation.

Burn risk due to contact with high temperature surfaces

Temperature in the heater surface could exceed 70°C. Please use the specified tools to perform the described operation.

1.2 Staff qualification

This equipment must only be operated by staff suitable qualified to know the hazards to which the operator is exposed in a laboratory and who have read and understood this manual or have been qualified on this equipment operation.

1.3 Proper use of the equipment

This equipment is for laboratory use.

The improper use of this unit (by not following the manual's instructions) can endanger the protection guaranteed for the equipment.

This equipment is NOT designed to work in explosive environments areas except the ones foreseen in this manual.

1.4 Modifications

The modification of the operation or the manipulation of the safety systems unauthorized by the manufacturer, could expose the operator to hazards not indicated in this manual.

2. Warranty

This product is guaranteed for one year from the date of purchase. The warranty does not cover damages caused by an incorrect use of the equipment or causes beyond J.P. SELECTA, s.a.u.

Any manipulation by staff not authorised by J.P. SELECTA, s.a.u. automatically cancels the warranty.

3. Packing List

The standard equipment consists of the following components:

Part	Part No	Quantity
Oven	2001406	1
Shelf	2000071	2
Instructions manual	80028	1
Timer Instructions manual (optional)	80431	1
CE Certificate		1
Electrical board key		2

4. Equipment Description

High temperature ovens **HIGHTEMP** have been designed with fan assisted air circulation, that allow adjustable temperature from 60°C up to 400°C.

Temperature is controlled by a digital regulator type PID.

As a safety element, a second temperature regulator acts as a temperature limiter.

The motor can be activated during the cooling cycle, regardless the heater.

5. Technical Features

(See also the oven's features plate.)

• Maximum temperature:	400°C
• Minimum temperature:	60°C
• Volume:	80 litres
• Fuse size 10x38 qG: (Supply 400V III)	Power: 3x 16A Control: 1x 2A
• Fuse size 10x38 qG: (Supply 230V III)	Power: 3x 20A Control: 1x 2A
• Max. Power:	4600W
• Power Supply:	230 / 400V III 50 Hz
• Shelf guides:	4
• Noise level:	70 dB - 80dB
• Inner chamber materials:	Stainless steel
• Weight:	75 Kg

6. Installation

6.1 Location

You must choose a suitable location for the equipment in order to avoid accidents and troubles to other laboratory activities, bearing in mind that this machine produces an unpleasant noise level, and that the oven generates fumes when operating.


It is important that the location meets the following conditions:

- It should be wide enough to operate with safety.
- There should be a ventilated area, either by windows or extractor fan, just in case there's a fumes leakage.
- There should be a three-phase + neutral wire to the mains.
- The equipment should be placed over a flat solid surface, suitable for the equipment weight.
- For a better operation of the equipment, do not expose it near sunny places nor heating sources.

6.2 Power Supply Connection

Take care when connecting this oven to the mains.

As it has a three-phase supply, the oven motor would turn one-way or the other depending on how the phases are connected.

You must try different combinations between wires Black, Brown and Grey until you find that the motor turns in the right way (See figures 6.1 and 6.2).

Make sure that the motor turns in the suitable way when switching on the oven.

If the motor turns the wrong way, this will excessively warm and its life will shorten.

Ground wire

The plug must have a ground wire.

The oven must be connected to an electrical supply provided with an earth wire.

Do not use nor start up the oven if its ground wire is not grounded to the installation.


Fig. 6.1 Motor turn way

7. Operation

7.1 Control Panel

1. Main and motor switch: It is used as main switch and also to start the fan motor.

2. Heating switch: It switches the heater on.

3. Motor starting up light (green): It indicates that the motor is ON.

4. Motor overtemperature alarm light (red): It indicates that the motor has stopped due to overtemperature.

5. Heater starting up light (green): It indicates that the oven's heater is ON.

6. Heater overtemperature alarm light: It indicates that the heater has stopped due to overtemperature.

7. Temperature controller: It establishes the operating temperature.

8. Safety temperature controller: It establishes the safety temperature.

9. Timer (optional): To program an operation program (see manual 80431)

10. End of Cycle light (optional): It's activated when cycle is finished (see manual 80431).

11. One cycle /Multi-cycle switch (optional): The programmed cycle at the timer (9) only is activated 1 time or indefinitely. (see manual 80431)



Fig. 6.2 Wire identification for mains connection


Figure 7.3
Temperature control panel

7.2 Starting Up

- Press the main switch/selector (1) to the ON position to work in no-timer mode.
 - The light (3) will come on and the controllers (7) and (8) and the motor start.
 - Press the heater switch (2).
 - Light (5) will light up and the heater will start up.
 - Select the working temperature in the controller (7).
 - Select the safety temperature in the controller (8).
- (See annexed instructions)

7.3 Temperature Controllers

HIGHTEMP uses two temperature controllers:

• **The working temperature controller:**

Once the controller has finished its starting routine, it allows the selection of the working temperature (SV) by means of the «Up» and «Down» keys (See figure 7.3).

Once the working temperature is selected (SV), the oven’s current temperature (PV) will increase up to the selected value (SV).

• **The safety temperature controller:**

This controller works in a similar way as the former one.

By means of the «Up» and «Down» keys, a safety temperature is selected.

The safety temperature must be minimum 5°C over the working temperature (PV).

If, in any case, the oven’s current temperature (PV) exceeds the safety temperature, the oven’s heater automatically disconnects.


The red light (6) indicates this failure.

7.4 Motor Thermal Protection

The oven has a thermal relay protective device for the motor which, in case of a failure in the power supply, motor tightening or excessive temperature, etc ... will switch off the oven.

The red light (4) indicates this failure.

If the oven is switched off by means of the thermal relay, you should disconnect it from the mains, check the motor and once the failure has been repaired, press the RESET key. Then the oven is ready to work again.

8. Maintenance

HIGHTEMP oven has no periodical maintenance operation.

However, it is important to check the oven every 12 months:

Check the electrical board connections status:

You have to:

- Disconnect the oven from the mains.
- Open the electrical board closet.
- Visually check the connections status. If you notice that there are black areas, please replace the cable or tighten the connection.
- Occasionally, tighten the power circuit board connections.

Check the mains connection status:

- Visually check the connections (plug) to the mains status, specially the protection ground cable connection. If you notice that there are black areas, please replace the cable or tighten the connection.

The manipulation of the oven's electrical board by unauthorized personnel could cause irreparable damage. Please, make sure to take the equipment to one of the authorized J.P. SELECTA's technical service.




CLEANING:


For the oven's cleaning, either the internal and external parts, do not use organic solvents which could damage the equipment. It is recommended to use alcohol or window cleaner.

8.1 Electrical Panel Components List with its Spare Part Numbers


Spare part:	Code:	Diagram No:
Commutator with light	20061	1
Green neon	22017	2 , 8
Red neon	22021	3 , 11
Motor contactor	13006	4
Motor thermal relay	36029	5
Switch	20049	6
Heater contactor	13008	7
Alarm reset timer	37008	9
Heater alarm relay	13007	10
Temperature regulator	16123	12 , 13
Solid relay	36028	14 (3 units)
Motor	24024	15
Electrical heater element	39047	16 (6 units)

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**Notice to customers:**

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.