**Incudigit-TFT**

2001261	19L
2001262	36L
2001263	52L
2001264	80L
2001265	150L

Digitheat-TFT

2001251	19L
2001252	36L
2001253	52L
2001254	80L
2001255	150L

Digitronic-TFT

2005163	33L
2005165	47L
2005167	76L
2005169	145L

(con puerta vidrio)

2005164	33L
2005166	47L
2005168	76L
2005170	145L

Incubig-TFT

2000238	288L	2000239	400L	2000240	720L
---------	------	---------	------	---------	------

ESTUFAS DE CULTIVO / ESTUFAS DESECACIÓN
CULTURE OVENS / DESICCATION OVENS

Indice

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Index

1. General Information.....	16
2. Packing list	16
3. Technical specification	16
4. Information about safety in the use of the equipment.....	18
4.1 Risks the operator is subjected to	18
4.2 Staff qualification.....	18
4.3 Use with hazardous substances	18
4.4 Maintenance.....	18
4.5 Warranty.....	18
5. Equipment Description	19
6. Installation	19
6.1 Location.....	19
6.2 Connection to the mains	19
7. Operation.....	19
7.1 Quick start to work at a fixed temperature.....	20
7.2 Setting screen	21
7.3 OPERATION screen	23
7.4 PROGRAMMING screen	24
8. USB log	25
9. Maintenance	25
10. Safety thermostat	26
10.1 Rearming the safety thermostat.....	26
10.2 Adjusting the safety thermostat	26
Aviso a los clientes	27
Notice to customers	27

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ENGLISH VERSION

1. General Information

- Handle the parcel with care. Unpack and check that the contents coincide with the packing list. If any part is damaged or missing, please advise the distributors immediately.
- Do not install or use the equipment without reading this handbook first. This instruction manual must be attached to the equipment and must be available for all users.
- If you have any doubts or enquiries, contact with J.P. SELECTA's technical service.
- Do not use the equipment in explosive, flammable or corrosive atmospheres.

2. Packing list

Model:	2001261	2001262	2001263	2001264	2001265
Shelves guides:	2x 2000011	2x 2000012	2x 2000012	2x 2000013	2x 2000015
Shelves:	2x 2000021	2x 2000022	2x 2000024	2x 2000023	2x 2000025

Model:	2000238	2000239	2000240
Shelves guides:	2x 2002371	--	--
Shelves:	2x 2002372	2x 2000063	2x 2000064

Model:	2001251	200125	2001253	2001254	2001255
Shelves guides:	2x 2000011	2x 2000012	2x 2000012	2x 2000013	2x 2000015
Shelves:	2x 2000021	2x 2000022	2x 2000024	2x 2000023	2x 2000025

Model:	2005163	2005165	2005167	2005169
(glass door)	2005164	2005166	2005168	2005170
Shelves guides:	2x 2000012	2x 2000033	2x 2000013	2x 2000015
Shelves:	2x 2000072	2x 2000073	2x 2000074	2x 2000075

3. Technical specification

Type:	CULTURE OVEN				
Model:	2001261	2001262	2001263	2001264	2001265
Chamber space: (L)	19	36	52	80	150
Interior: High/wide/depth (cm)	30/25/25	40/30/30	33/47/33	50/40/40	50/60/50
External : High/wide/depth (cm)	51/57/49	60/62/54	53/79/57	70/72/64	70/92/74
Power supply:	230 VAC 50Hz/60Hz (according to characteristics label)				
Electrical power (W):	170	225	250	300	525
Weigth: (Kg)	26	36	46	54	75
Fuse:	Replaceable (Type: 10x5mm)				
Temperature range:	From ambient + 5°C ... up to ... 80°C				
Display resolution:	0.1°C				
Stability	+ - 0.1°C at 37°C				
Homogeneity	+ - 0.5°C at 37°C				
Timed operation:	week days, hh.mm			max 23h:59min	
Temperature Profiles:	Up to 10 profiles 6 segments. Segment from 1minute to 99hour.				
Data recording:	USB memory device				

Type:	CULTURE OVEN				
Model:	2000238	2000239		2000240	
Chamber space: (L)	288	400		720	
Interior: High/wide/depth (mm)	800/600/600	1000/800/500		1200/1000/600	
External: High/wide/depth (mm)	970/910/760	1300/1140/750		1520/1840/850	
Power supply:	230 VAC 50Hz/60Hz (according to characteristics label)				
Electrical power (W):	570	1100		1600	
Weigth: (Kg)	87	160		225	
Fuse:	Reemplazable (Type: 10x5mm)				
Temperature range:	From ambient + 5°C ... up to ... 80°C				
Display resolution:	0,1°C				
Stability	+ - 0.2°C at 37°C				
Homogeneity	+ - 0.7°C at 37°C				
Timed operation:	Week, hour, minute		max 23h:59min		
Temperature Profiles:	Up to 10 profiles 6 segments. Segment from 1minute to 99hour.				
Data recording:	USB memory device				
Type:	DRYING OVEN (NATURAL CONVECTION)				
Model:	2001251	2001252	2001253	2001254	2001255
Chamber space: (L)	19	36	52	80	150
Interior: High/wide/depth (cm)	30/25/25	40/30/30	33/47/33	50/40/40	50/60/50
External: High/wide/depth (cm)	50/60/44	60/65/49	53/82/52	70/74/59	70/95/68
Power supply:	230 VAC 50Hz/60Hz (according to characteristics label)				
Electrical power (W):	600	900	1000	1200	2100
Weigth: (Kg)	24	35	44	59	73
Fuse:	Replaceable (Type: 10x5mm)				
Temperature range:	From ambient + 5°C ... up to ... 250°C (*special model up to 300°C upon request)				
Display resolution:	1°C				
Stability	+ - 0.3°C at 150°C				
Homogeneity	+ - 1°C at 150°C				
Timed operation:	Week days, hh.mm		max 23h:59min		
Temperature Profiles:	Up to 10 profiles 6 segments. Segment from 1minute to 99hour.				
Data recording:	USB memory device				
Type:	DRYING OVEN (FORCED CONVECTION)				
Modelo:	2005163	2005165	2005167	2005169	
(glass door)	2005164	2005166	2005168	2005170	
Chamber space: (L)	33	47	76	145	
Interior: High/wide/depth (cm)	40/28/30	33/45/32	50/38/40	50/58/50	
External: High/wide/depth (cm)	60/65/55	53/81/58	70/75/65	70/95/72	
Power supply:	230 VAC 50Hz/60Hz (according to characteristics label)				
Electrical power (W):	1200	1200	1600	2000	
Weight (Kg):	38/40	46/50	58/64	74/79	
Fuse:	Replaceable (Tipo: 10x5mm)				
Temperature range:	From ambient + 5°C ... up to ... 250°C (*special model up to 300°C upon request)				
Display resolution:	1°C				
Stability	+ - 0.3°C at 100°C				
Homogeneity	+ - 1°C at 100°C				
Timed operation:	Week, hour, minute		max 23h:59min		
Temperature Profiles:	Up to 10 profiles 6 segments. Segment from 1minute to 99hour.				
Data recording:	USB memory device				

4. Information about safety in the use of the equipment

Hazard risk

Electric risk

Risk of high temperature contact

4.1 Risks the operator is subjected to

- Risk of burns
- Electrical risk

4.2 Staff qualification

This equipment may only be used by properly trained personnel who have read and understood these instructions or have been properly trained in this equipment operation.

4.3 Use with hazardous substances

Do not treat flammable or explosive materials or those which chemically react with violence.

ATTENTION!

Burning risk

Take proper precautions
(use of thermal gloves, clamps,...)
to extract samples from inside the oven.

It can be at a high temperature.

4.4 Maintenance

Unplug the equipment from the mains before removing the oven's cover to handle inside.

The manipulation of the equipment internal electronic circuits by unauthorized personnel can cause irreparable damage. Take it to one of J.P.SELECTA's authorized technical services.

It is recommended to clean the accessible parts of the oven with a cloth and a cleaner, or a similar solution. Do not use aggressive solvents.

4.5 Warranty

This product is guaranteed for one year. It does not cover damages caused by an incorrect use or causes beyond the control of J. P. SELECTA, s.a.u.

Any manipulation of the equipment by personnel not authorized by J.P. SELECTA, s.a.u. automatically cancels the guarantee.


Fig 1. Temperature profile

5. Equipment Description

Ovens are used in sample thermal treatments, for its desiccation, cured and culture incubation.

Giving the heat treatment possibility by following a "temperature profile" allows avoiding cracks, breaks, color changes, etc.

A temperature profile is a succession of "ramps" and "dells."

A ramp is specified by a slope of °C/minute or °C/hour.

A dell is specified by a temperature and a time.

The device has a temperature log in an USB memory drive.

6. Installation

Unpack the equipment and visually check that the contents are in good conditions. Keep the package for some days.

6.1 Location

Place the equipment on a flat level stable surface, suitable for the equipment weight.

To work in a comfortable way, please leave a free space of 50cm around the equipment and a minimum of 20cm at the rear part.

6.2 Connection to the mains

Plug the equipment to a near socket suitable for the equipment power. Use the connection cable supplied or a similar one.

For your safety, the equipment must be earthed.

Before plugging the machine, check that the voltage supply and the frequency corresponds to the one indicated in the equipment characteristics plate.

7. Operation

All the equipment functions are set with the touch screen. Gently press with your finger in the center of each icon.

When starting the computer, the starting screen will be shown (Fig 7.1).

Fig 7.1. Initial screen

• SETTINGS:

- Time and date: menu to set the date and time.
- PID and Autotuning: menu to set the parameters of temperature control.
- Downloading recorded data.
- Temperature correction.
- Temperature alarm limits.

• PROGRAMMING

- Allows you to define the temperature profile, by defining heating temperatures, times and ramps.


• OPERATION

- Allows you to select and launch the operating temperature or a temperature profile.


7.1 Quick start to work at a fixed temperature


1 Press on «OPERATION».


2 Press over temperature panel. The setpoint indicator will change to blue.

3 Press to select the setpoint temperature. Wait a few seconds to confirm the new value.

4 Press on to start.


5 Press on to end the cycle.

NOTE: at the end of the cycle, a message appears saying that the cycle has been completed by the user. To delete it, click on the message.

NOTE: If the oven is turned off with the ON/OFF switch and without pressing STOP at the end of the cycle, a power failure error message will appear each time the equipment is started up. By clicking on the warning, the message will disappear and you could work normally.

7.2 Setting screen

Icons meaning:

Note: it is important to keep the equipment on time so that alarms and logs could indicate the exact moment when they occurred.

7.2.1 Set time and date

Press on time/date. Press successively on the box to select the day, month, year, minute and second adjustment.

Adjust the values by pressing on the icons Increase and Decrease.

7.2.2 Control

Through this menu you can access the general parameters of the equipment, the model and the heating parameters.

With a value of the appropriate P, I, D parameters, the oven will reach the temperature in the shortest possible time and without exceeding the set temperature. However, those values can be different at high and low temperatures.

The oven comes with the values which are appropriate for most jobs.

These values can be modified, but mathematical knowledge is needed to calculate them. It is much more advisable to perform an "Autotuning" when you need more accuracy.

"Autotuning" automatically adjusts the optimal parameters for a given temperature and load. This process is long because it performs several heating / cooling phases to know the functioning of the environment.

Press on the Control menu. Press successively on the box to select the P, I, D values adjustment and modify the value by pressing on the Increase and Decrease icons.

Attention: It is not recommended to change these parameters. Use the «Autotuning» option to automatically adjust these values.

Select «Autotuning» 40 or 100 with:   .
(it is recommended to set the most used temperature). There are two calculation methods that are identified with 100 or with 40. The 100 method is more accurate but it requires more calculations. For most applications, 40 method is more than enough.

Start «Autotuning» by pressing on: 

If for some reason you want to stop, press on: 

Note: «Autotuning» process has a long duration, it is important that the conditions do not vary, therefore do not open the door or handle the oven.

NOTE: If you use the equipment in various temperatures, perform the «Autotuning» function at the lowest temperature.


Fig 7.2. Communication


Fig 7.2. Temperature alarms

7.2.3 Log

The equipment can store readings from time to time and send them to a computer or store them in an USB memory pen drive.

«PORT» Set the readings storing mode:

- RS232: Send the temperature readings to a computer by a serial cable (DB9 connector).
- USB: Store the readings in the internal memory and download later to a USB

«INTERVAL» Select with the register time interval or the data sending.

«MODE» Select the data to be stored:

- Start: Save all temperature readings.
- Set point: Save temperature readings, since the set point temperature has been reached.

«USB» Update firmware.

Note: The system only saves the last process. To save the process to an USB memory, simply connect the USB to the connector, and data will send automatically to the USB memory. A "sent OK" message will appear. If a process is being performed, this action will be performed after pressing STOP.

7.2.4 Temperature

It is recommended to verify calibration of temperature with some regularity (for example, annually) by using for instance a thermometer or a temperature pattern.

If it is necessary to correct the reading as a result of a calibration, one must apply a correction.

To do this, press on the temperature menu.

A box indicating the current correction value appears.

Adjust with keys until you get the indicated value on the temperature pattern.

Applied correction will appear in brackets.

7.2.4.1 Two-point calibration

This function allows you to calibrate the equipment at two different points. For example, at 40°C and at 100°C. In this way, a better linearity is achieved in calibration than if only offset is adjusted.

7.2.4.2 Temperature alarms

By pressing successively on the temperature correction box, one may enter an over-temperature and a low-temperature limit. Both alarms activate only a message on screen, if the temperature is outside the limits.

7.2.5 Units

Press on the temperature units menu. Press successively to change between Celsius or Farenheit degrees. It is also possible to change the ramp units de °C/min or °C/hour.

7.2.6 Language

Press on the language menu. Press successively on the box to select among Spanish, French or English.

Adjust by pressing on the Increase and Decrease icons.

7.3 OPERATION screen

There are two operation modes:

- Working at a fixed temperature.
- Working according to a pre-programmed temperature profile.

Press on  to commute the operation mode (fixed/program)


Fig 7.3.1. Run screen operating at one temperature.

7.3.1 Operation at a fixed temperature

- Select the set point temperature by pressing on the temperature indicator.
- Select a temperature increase ramp if you want to heat slowly. The value "—" heats as fast as possible.
- Select an operating time. This time starts to count down when the temperature reaches the set point. The INF value sets a continuous operation.

- Start by pressing on: 

7.3.2 Operation with a temperature profile

- Temperature profiles must be previously programmed from the «PROGRAMMATION» display (see 7.4).
- Select a profile (from 0 to 9) by pressing on keys:  
- Start by pressing on: 
- In the temperature box (top-left corner) the current temperature and the set point are indicated.

Fig 7.3. Run screen operating with temperature profile.

Fig 7.3.2 Working screen data.

Step target temperature	*	Program:	Current step duration:
Current temperature		Nº step on running.	Last ramp run value.
			Total elapsed time of the current step:


- * Point 1 (top) indicates the status of the oven:
 - . if it flashes red, the heating works.
 - . if it is yellow (only in programs) it indicates that the oven is in the resting phase (waiting to reach the heating phase).
- * Point 2 (lower) indicates the data log:
 - . if it flashes green, it indicates that data logging is activated.
 - . if it flashes white, it indicates the frequency of data recording, or what is the same, the recording interval.

7.4 PROGRAMMING screen

Go to PROGRAMMING screen to:

- Set the temperature profiles.
- Set the timed operation for a profile.
(For instance, you can set a time and a week day to start a profile. As start, profile nº 1 on Sunday at 19h00')

7.4.1 Programming a temperature profile:

A temperature profile is a sequence of steps, also called "dells".

To go from the end of a «dell» to the following «dell», it could be done in two ways:

- With a «temperature ramp».
(go from 100°C to 200°C following a ramp at 2°C per minute).
- As fast as possible, "--" indicator (maximum power).

Profile parameters input:

- Select the profile number between 0 and 10 with .
- Press on the parameter to be modified. (Temp, ramp, length ...)
- Press on these icons: to enter the value.
- Modify another parameter.


Special values:

- If you enter value **OFF** in **temperature** (you have to press Increase key) and the next steps are disabled and not executed.
- If you enter **INF** in **Step length**, the oven remains permanently at this temperature.
- If you enter -- in the **ramp**, the oven makes the temperature rise at the fastest speed possible.

Starting of temperature profile

- To set the profile starting time, press on «—:—» (09:00 in the picture).
- Enter time by pressing in hh:mm format.
- Press on the days of the week when the profile should be run.
- Accept the programming by pressing on .

Note: If you run a timed program, do not forget to press the RUN key to start the program. During the timeout, it will indicate that it is in Standby mode.

8. USB log

The device has an internal memory for data acquisition (Datalogger).

Once the cycle is finished, it is possible to download this information to an USB memory via an USB connector. The device (USB memory) must be formatted to FAT 32.

Functioning:

- By default, the registrar is not active.

Follow the instructions to change the log interval (see settings>log). Adjust this interval with the conditions of your program.

With the cycle finished:

- Connect the «USB memory» to the USB connector on the front of the device.
- When an «USB memory» is connected, the registration data will be automatically downloaded. A message appears indicating that a connected memory is found and the message "send ok" appears at the end of the download.
- This operation can be performed several times, as long as another cycle has not started (press RUN). When pressing RUN, internal memory of the device empties. The device only stores the last process performed.

Recover data from PC:

- Connect the «USB memory» to a computer to download the data file.
- Important note: The "USB memory" must have a device name. It cannot be blank.

The created file is text type, which can be opened with WordPad, WORD or EXCEL.

Nº lot.	00025	Batch number
Date	2018/03/21	Batch data
Start time	12:00	Registration start date and time
Interval	30 sr	Interval between logs
Nº Dat.	32	Number of logs, amount of data recorded
Temp.	11	Target temperature
	19	
	19	
	20	
	21	
	23	
	28	
	33	
	40	
	47	
	56	
	64	
	73	

Example of text file downloaded via USB

9. Maintenance

Before removing the equipment cover, disconnect it from the mains. The control board must only be handled by duly authorized maintenance personnel.


Cleaning:

- Inner enclosure: clean with polish for stainless steel.
- Outer enclosure: clean with a damp cloth or a mild detergent.
- Touch screen: use a cloth and a glass-cleaning liquid.

Annually:

- Contrast the temperature against recognized standart traceable equipments.
- Visually check the status of electrical connections, grounding and connection cable to the mains.

10. Safety thermostat

8. Safety thermostat regulation screw (grey).

9. Safety thermostat reset button (white) - only DIGITHEAT & DIGITRONIC models.

10.1 Rearming the safety thermostat

(Only for DIGITHEAT & DIGITRONIC models)

Safety thermostat is an element to prevent the temperature from exceeding a value. If by a fortuitous event the thermostat is triggered, that is, it opens the power circuit, the heating element will stop working. Usually there is an indicator light, or an error message, that indicates that the thermostat is in this position.

If the thermostat is triggered, check that the equipment looks good and it is not overheating at any point.

To change its status, simply press the white button (9), and you will hear a "click" that indicates that the thermostat has been reset. If the equipment has not cooled down sufficiently, the thermostat cannot be reset, and you must wait for the temperature to drop before it can be reset.

10.2 Adjusting the safety thermostat

From factory, the thermostat is adjusted to its maximum temperature. To adjust the safety thermostat temperature, do as follows:

1. With the equipment connected to the mains, turn the safety thermostat screw (8) clockwise to the end of the scale, with the help of a flat screwdriver.
2. Select the desired safety temperature using the control keys. This temperature has to be minimum at 3°C above the working temperature.
3. Turn on the equipment with the Start key.
4. When the equipment is stabilized at the above selected temperature, turn the safety thermostat screw (8) slowly to the left until it clicks, and the red indicator lamp lights up.
5. Stop the equipment and leave it to cool with the lid open for 15 minutes.
6. Select the working temperature by using the temperature selection keys (minimum at 3°C below the desired safety temperature).
7. Turn the equipment on again with the Start key. Once the temperature is stabilized, the machine is ready to work.

11. Spare Parts List

These devices are in constant development and depending on the age of the equipment, the spares could be different. Please consult with our After Sales Dept. to confirm the proper replacement part.

Model	Code	Description
DIGITRONIC-INCUDIGIT TFT	21062	Porous Silicone gasket for stove (1m)
DIGITRONIC TFT	43040	Temperature Limiter L-240 (stoves up to 250°C)
DIGITRONIC TFT	43176	Temperature Limiter L-300 (stoves up to 300°C)
INCUDIGIT TFT	43045	Temperature Limiter L-110C
DIGITRONIC TFT	16139	TOKY Temperature Controller
DIGITRONIC-INCUDIGIT TFT	16099	OMRON Temperature Controller
DIGITRONIC-INCUDIGIT TFT	29600	Stove Cabling Board 230V
DIGITRONIC-INCUDIGIT TFT	29602	Stove Cabling Board 115V
DIGITRONIC TFT	43152	Thermostat 300°C
INCUDIGIT TFT	43153	Thermostat 120°C
DIGITRONIC-INCUDIGIT TFT	43074	PT100 Sensor
DIGITRONIC TFT	24008	Fan motor
DIGITRONIC TFT	29589	Double Relay Board (only with motorized air valve)
DIGITRONIC TFT	16172	Motorized air valve.
DIGITRONIC-INCUDIGIT TFT	29560	CPU-TFT Stoves Board
DIGITRONIC-INCUDIGIT TFT	85167	Display-Digitizer TFT 4.3"

Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.