HOTCOLD-UCS

2101619

ESTUFA REFRIGERADA DE PRECISI3N
REFRIGERATED PRECISION CABINET

INDICE

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INDEX


GENERAL INFORMATION	16
PACKING LIST	16
ACCESSORIES	16
TECHNICAL FEATURES	17
EQUIPMENT DESCRIPTION	17
INSTALLATION	18
OPERATION	21
USB LOG	26
MAINTENANCE	27
SPARE PARTS	27
GUARANTEE	27
SAFETY THERMOSTAT	28

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Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.

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GENERAL INFORMATION

- 1) Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributor immediately.
- 2) Do not install or use the equipment without reading this handbook before.
- 3) This handbook must always be attached to the equipment and it must be available for all users.
- 4) If you have any doubts, please inquiries contacting with your supplier or J.P. Selecta's technical service.
- 5) **IMPORTANT! J.P. SELECTA WILL NOT ACCEPT ANY EQUIPMENT TO BE REPAIRED IF IT IS NOT DULY CLEANED.**
- 6) If any modification, elimination or lacking in maintenance of any device of the equipment by the user transgress the directive 2009/104/CE, the manufacturer is not responsible for the damage that can occur.
- 7) Do not use the equipment with liquids which can give off vapours capable of making explosive mixtures.

PACKING LIST

The standard equipment consist of the following components:

Description	Code	Quantity
Hotcold UC	2101515	1
Shelves	1001858	4
Instruction Manual	80460	1
EU Declaration of Conformity		1

ACCESSORIES

Hotcold UCS equipments allow connect stirrers and shakers without heating elements, which can be connected on inner sockets or through lateral hole.

TECHNICAL FEATURES

Voltage supply 115-230V 50/60 Hz according to the characteristics plate indications.

Features		
Code		2101619
Temperature range		-5°C to 65 °C
Capacity (litres)		235
Inner useful dimensions (cm) *with / without shelf	Height	109
	Width	51
	Depth	55
Overall dimensions (cm)	Height	156
	Width	61
	Depth	62
Consumption (W)		1000
Weight (kg)		90

EQUIPMENT DESCRIPTION

Hotcold UC refrigerated cabinets are built in stainless steel AISI 304. They have two electric sockets.

Fan air circulation which allows:

Stability*: $\pm 0.5^{\circ}\text{C}$
 Homogeneity: $\pm 2.5^{\circ}\text{C}$
 Resolution: $\pm 0.1^{\circ}\text{C}$

*Note: tests carried out in environmental conditions 20°C and 50% RH

They can work at temperatures between -5°C to 65°C. These equipments are designed for working under laboratory conditions of:

Temperature 18-27°C
 Relative humidity 25-75%

Working outside these margins can limit their benefits.

Temperature profile

These equipments are used in sample thermal treatments, for its desiccation, curing and others.

The possibility of giving thermal treatment by following a “temperature profile” allows to avoid cracks, breaks, colour changes, etc.

A temperature profile is a succession of “plains” and times.

A plain is specified by a temperature and a time.

The device has a temperature log on a USB memory pen drive.

INSTALLATION

Place the cabinet on a flat horizontal, level surface, leaving a space of about 10 cm at the back and sides of the equipment.

For a maximum efficiency, do not place the machine in direct sunlight or other sources of heat.

Connect the equipment to the mains. Be sure that the voltage supply is the same as the one indicated on the characteristics plate. Do not use the equipment if it is NOT earthed.

Note: Before connecting it to the mains, it is advisable to leave the machine at rest for 24h to stabilize oil and gas cooling unit.

If you have to change the plug cap, bear in mind the following:


- Blue cable: Neutral.
- Brown cable: Phase.
- Yellow/green cable: Earth.

Check that the circuit breaker on the electrical panel (4) is powered. See “Electrical panel” section on the “Operation” chapter.

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Do not use the cabinet for thermal treatment of products which can give off vapours capable of making an explosive mixture.

Inside view of the cabinet


INNER CHAMBER:

- 1. Internal light.
- 2. Air circulating fan
- 3. Side hole for connections inlet.

ELECTRICAL PANEL

- 4. Magnetothermic control, cold - heating group and plugs.
- 5. Safety thermostat

Initial screen

OPERATION

All the equipment functions are set with the touch screen. Gently press with your finger in the center of each icon.

When starting the computer, the starting screen will be shown.

• SETTINGS:

- Time and date: menu to set the date and time.
- PID and Autotuning: menu to set the parameters of temperature control.
- Downloading recorded data.
- Temperature correction.
- Temperature alarm limits.

• PROGRAMMING

- Allows you to program the temperature profile, by defining temperatures and times.

• TASK

- Allows you to select and launch the operating temperature or a temperature profile.

Fast temperature setting and start

1 Press on «PROGRAMMING».

2 Press over temperature panel.


3 Press   to set temperature

4 Press on  to start

Setting screen

Icons meaning:

Set time and date

Press on the time/date. Press successively on the box to select the days, months, years, minutes, seconds adjustment.

Adjust by pressing on the icons to increase and decrease.

Control

With a value of the appropriate P, I, D parameters, the cabinet will reach the temperature in the shortest possible time and without exceeding the set temperature. However, those values can be quite different at high and low temperatures.

These values can be modified, but it's advisable to perform an «Autotuning» when you need a better accuracy.

«Autotuning» automatically adjusts the optimal parameters for a given temperature and load. Do not manipulate the equipment during the process.

Press on the Control menu. Press successively on the box to select the P, I, D values adjustment and modify the value by pressing on the increase and decrease icons.


Attention: It is recommended to use «Autotuning» option to automatically adjust these values, instead of modifying them manually.

Press on the Control menu.


Select the «Autotuning» temperature with:   (it is recommended to set the most used temperature). There are two calculation methods that are identified with 100 or with 40. The 100 method has a slower duration, but it is more accurate. For most applications, 40 method is more than enough.

Start «Autotuning» by pressing on the icon 

If for some reason you want to stop, press on: 

LOG MENU

TEMPERATURE SETTINGS MENU

Log

The equipment can store readings from time to time and store them lately in a USB memory pen.

- «PORT» Set the readings storing mode:
 - USB: Store the readings in a USB memory pen drive.
- «INTERVAL» Select the log time interval of a data, with the keys
- «MODE» Select the data to be stored:
 - Start: Save all temperature readings.
 - Set point: Save temperature readings, since the set point temperature has been reached.
- «USB» Update the firmware from a file that is in a USB.

Temperature

It is recommended to verify calibration of temperature with some regularity (for example, annually), by using for instance a thermometer or a temperature pattern.

If it is necessary to correct the reading, as a result of calibration, one must enter a correction. To do this, select CORRECT T.

A box indicating the current correction value appears.

Adjust with keys until the you get the indicated value on the temperature pattern. Applied correction will appear in brackets.

Temperature alarms

By pressing successively on the temperature correction box, one may enter an over-temperature (**Alm sobre T**) and low-temperature (**Alm bajo T**) limit. Both alarms activate a message on screen.

T.S. settings

This parameter sets the maximum value of the electronic safety thermostat, it is recommended to leave it at its maximum value. It activates a message on the screen.

Comp. ON

Adjust the operating time of the compressor.

Comp. OFF

Adjust the compressor shutdown time.

RET.VET

Adjust the delay of entry into operation of the indoor fan.

Tª INI COMP

Adjust the temperature from which the compressor stops working.

Units

Press on the Units menu. Press successively to change between Celsius or Farenheit degrees.

Language

Press on the Language menu. Press successively on the box to select among Spanish, French or English.


Adjust by pressing on the increase and decrease icons.

OPERATION screen

There are two operation modes:

- Working at a fixed temperature and a time.
- Working according to a pre-programmed temperature profile.

Press on  to commute the operation mode.

Run screen operating at one temperature.

Run screen operating with temperature profile.

Operation at a temperature

- Select the set point temperature by pressing on the temperature indicator.
- Select an operating time. This time starts to count down when the temperature reaches the set point.
- The INF value sets a continuous operation.

Start by pressing on this icon: .

Operation with a temperature profile

• Temperature profiles must be previously programmed from the «PROGRAMMING» display (see next chapter).

• Select a profile (from 0 to 9) by pressing on keys  .

• Start by pressing on the icon: .

• In the temperature box (top-left corner), the current temperature and the set point are indicated.

• By pressing the graph icon you can switch the way of displaying the temperature evolution data either in graphical form or in a table format.

*

Step target temperature	*	Program:	Current step duration:
Current temperature		Nº step on running.	Total elapsed time of the current step:

Fig 7.3.2 Working screen data.

* Point 1 (top) indicates the status of the oven:
 . if it flashes red, the heating works.
 . if it is yellow (only in programs) it indicates that the stove is in the resting phase (waiting to reach the heating phase).


* Point 2 (lower) indicates the data log:
 . if it flashes green, it indicates that data logging is activated.
 . if it flashes white, it indicates the frequency of data recording, or what is the same, the recording interval.

PROGRAMMING

Go to PROGRAMMING screen to:

- Set the temperature profiles.
- Set the timed operation for a profile.
(For instance, you can set a time and a week day to start a profile. F.i.: start, profile nº 1 on Sunday at 19h00')

Programming a temperature profile

A temperature profile is a sequence of a maximum of 6 steps, also called dells.

Profile parameters input:

- Select the profile number with between 0 and 10.
- Press on the parameter to be modified (Temp, length ...)
- Press on these icons: to enter the value.
- Modify another parameter.

Special values:

- If you enter **INF** in **Step length**, the cabinet remains permanently at this temperature.

Starting of temperature profile:

- To set the profile starting time, press on «—:—»
- Enter time by pressing in **hh:mm** format.
- Press on the week days when the profile should be run.
- Accept the programming by pressing on .

STARTING UP

1. Check that the installation protocol has been correctly followed.
2. Switch on the main switch (see picture). The graphic screen of the temperature controller will light up.

NOTE: Maximum power you should connect on inner sockets is 500 W.

3. Configure the clock (real time clock) only the first time or when the summer / winter time change occurs. See "Clock configuration".

Selection offset value

"Offset" option allows increasing or reducing temperature reading in order to adjust it at other readings made with a standard temperature.

- 1 Check that the equipment is in RUN mode (in operation)
- 2 Follow the adjustment instructions to calibrate.

Temperature alarm / electronic safety thermostat

Control alarms can be set to show an alarm message as a warning, telling that temperature values reached are out of programmed range.

The equipment will continue normally operating despite the error condition. During the error condition it will be displayed on the screen and when the temperature is restored within the limits, the "alarm" status will disappear from the screen.

Note: Be aware when programming the alarm values, that opening the cabinet door may cause the temperature to reach values outside the programmed range and then the display will show the alarm message.

USB LOG

The device has an internal memory for data acquisition (Datalogger).

Once the cycle is finished, it is possible to download this information to an USB memory via an USB connector. The device (USB memory) must be formatted to FAT 32.

Functioning:

- By default, the registrar is not active.

Follow the instructions to change the log interval (see settings>log). Adjust this interval with the conditions of your program.

With the cycle finished:

- Connect the «USB memory» to the USB connector on the front of the device.
- When an «USB memory» is connected, the registration data will

be automatically downloaded. A message appears indicating that a connected memory is found and the message "send ok" appears at the end of the download.

- This operation can be performed several times, as long as another cycle has not started (press RUN). When pressing RUN, internal memory of the device empties. The device only stores the last process performed.

Recover data from PC:

- Connect the «USB memory» to a computer to download the data file.

Important note: The "USB memory" must have a device name. It cannot be blank.

The created file is text type, which can be opened with WordPad, WORD or EXCEL.

MAINTENANCE

Before removing the case, disconnect the equipment from the mains.

The manipulation of the internal electronic circuits of the cabinet by unauthorized personnel can cause irreparable damages. Take it to one of J.P.SELECTA's authorized technical services.

CLEANING:

For cleaning different parts, we recommend to use the following products:

Stainless steel cleaning: Alcohol

Plastic cleaning: Alcohol with cotton duster.

SPARE PARTS


16099 Temperature regulator
20075 Magnetothermal
29553 Control Circuit
36027 Solid relay 25A
39057 Electrical Resistance
43031 Safety thermostat
43071 Temperature probe
85167 4.3" TFT screen

GUARANTEE

This product is guaranteed for one year. The guarantee does not cover damages caused by an incorrect use or causes beyond the control of J.P. SELECTA, s.a.u.

Any manipulation of the equipment by unauthorized personnel by J.P. SELECTA, s.a.u. cancels the guarantee automatically.

SAFETY THERMOSTAT

8. Safety thermostat regulation screw (grey).
9. Safety thermostat reset button (white).

Rearming the safety thermostat

Safety thermostat is an element to prevent the temperature from exceeding a value. If by a fortuitous event the thermostat is triggered, that is, it opens the power circuit, the heating element will stop working. Usually there is an indicator light, or an error message, that indicates that the thermostat is in this position.

If the thermostat is triggered, check that the equipment looks good and it is not overheating at any point.

To change its status, simply press the white button (9), and you will hear a "click" that indicates that the thermostat has been reset. If the equipment has not cooled down sufficiently, the thermostat cannot be reset, and you must wait for the temperature to drop before it can be reset.

Adjusting the safety thermostat

From factory, the thermostat is adjusted to its maximum temperature. To adjust the safety thermostat temperature, do as follows:

1. With the equipment connected to the mains, turn the safety thermostat screw (8) clockwise to the end of the scale, with the help of a flat screwdriver.
2. Select the desired safety temperature using the control keys. This temperature has to be minimum at 3°C above the working temperature.
3. Turn on the equipment with the Start key.
4. When the equipment is stabilized at the above selected temperature, turn the safety thermostat screw (8) slowly to the left until it clicks, and the red indicator lamp lights up.
5. Stop the equipment and leave it to cool with the lid open for 15 minutes.
6. Select the working temperature by using the temperature selection keys (minimum at 3°C below the desired safety temperature).
7. Turn the equipment on again with the Start key. Once the temperature is stabilized, the machine is ready to work.