VACIOTEM-GMP

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ENGLISH

1. General info

Handle the package carefully. Unpack it and check that the contents match those indicated in the "Packaging Contents" section. If you notice any damaged or missing components, notify the distributor immediately.

Do not install or use the equipment without first reading this instruction manual.

These instructions are an integral part of the equipment and must be available to all users.

Any questions can be clarified by contacting the technical service of J.P. SELECTA, s.a.u.

ATTENTION! NO MACHINE WILL BE ACCEPTED FOR REPAIR THAT IS NOT PROPERLY CLEANED AND DISINFECTED.

Any modification, removal, or lack of maintenance of any device on the machine violates the User Directive 2009/104/EC, and the manufacturer is not responsible for any resulting damage.

DO NOT USE THE EQUIPMENT WITH FLUIDS THAT MAY RELEASE VAPORS THAT ARE CAPABLE OF FORMING EXPLOSIVE MIXTURES.

2. Packaging Contents

Item	
Vacuum Drying Oven VACIOTEM-GMP	(x1) cod. 4001491
Hose PRESCORD 20ATM Ø12/22 (1.5 mts)	cod. 46147
Hose Clamp INOX ø12/22	(x2) cod.44288
Silicone tube 10x6mm (1.5 mts)	cod. 46042
Instructions Manual	M.80462
EU Declaration of Conformity	

3. Technical Specification

Item	
Required Voltage (V):	(230V / 50Hz) otros valores bajo demanda
Installed Power (W):	800 (vacuum pump not included)
Máxima intensidad (A):	4 A
Fuses (quick type, type H, dimensions: 5x20mm):	8 A (x2)
Supplied Power Cord:	1ph+N+T/T (1.5mm ²) x 2 mts. EU type Schuko 16A (other connections on demand)
Net Weight (Kg):	75

Dimensions inside (cm Ø x depth.):	Ø34x52
External Dimensions.: (cm depth. x width x height):	70x76x54
Volume: (Liters):	47
Trays Nr.	2
Max height over tray (mm)	95
Tray dimensions (cm. Depth. x wide)	45 x 29
Temperature range (°C)	20..200 (setting temp. must be +5°C over room temperature)
Max vacuum (mbar / Pa)	5 / 500
Vacuum pump electrical socket?	YES
Hose connection for vacuum	1/2" (Ø12mm)
Thermal Homogeneity	+/- 1.5°C a 100°C
Thermal stability	+/-0.2°C a 100°C
Setpoint error	+/-1°C
Room conditions	Max. 40°C HR 70%
USB Connection	SÍ
ETHERNET Connection	SÍ

3.1 Electrical connection

Check the electrical installation in the location where the vacuum oven will be installed.

For your safety, we recommend to protect the device with a circuit breaker and a differential switch.

Please, consider additionally, the power of the required vacuum pump.

Model	Voltage	Power (W)	Circuit breaker
4001491	230V (1f+N+TT)	800	10A

3.2 Accessories

For the proper functioning of the vacuum oven, a vacuum system is required.

If a vacuum network is not available, JP Selecta recommends the following equipment:

Vacuum Pump cod. JP Selecta: 5900622

4. Equipment description

Vacuum ovens are equipment designed to perform thermal processes under controlled conditions, where air and moisture are removed to create a vacuum environment. This allows for precise results and improved quality of the treated products, especially in applications sensitive to oxidation or air degradation.

Examples of use:

a) Electronic components industry:

Drying electronic components after cleaning.

Removing moisture from printed circuit boards and PCBs before encapsulation.

b) Chemical products industry:

Drying of chemical and biological samples avoiding modification of their properties.

Controlled evaporation of solvents to obtain purified products.

c) Research & Development:

Material testing under specific vacuum and temperature conditions.

Heat treatment of new product prototypes.

The use of vacuum allows thermal processes to be more efficient, faster, and safer, especially in applications where it is crucial to avoid reactions with oxygen or preserve the integrity of the treated materials.


VACIOTEM-GMP vacuum ovens are designed to allow drying at adjustable temperatures from 20°C to 200°C with a stability of $\pm 0.2^\circ\text{C}$ and a homogeneity of $\pm 1.5^\circ\text{C}$. The equipment incorporates digital control of the temperature and operating modes. The heating elements are located beneath the trays for optimal temperature control. The VACIOTEM-GMP model's digital control also allows vacuum level control and the ability to communicate via Ethernet with FTP servers.

4.1 Placement

Place the oven on a flat, horizontal, and levelled surface, ensuring a 10 cm clearance at the back and sides of the oven.

Connect the vacuum pump, or if applicable, the vacuum net, and open the rear valve.

Connect the vacuum pump plug to the power socket on the back of the oven.


Connect the oven to the mains. Ensure that the equipment is connected to a mains voltage that matches the voltage indicated on the nameplate.

If the product to be used with the vacuum oven can emit vapors, only water jet pumps or pneumatic hydro-rotary pumps should be used directly connected to the oven. If an air ballast pump is used, a water-cooled condenser must be installed between the pump and the oven to collect vapors and prevent damage to the pump. In addition, condensation collectors must be added on both sides of the pump (suction and discharge).

For best stove performance, do not place it in direct sunlight or near significant heat or cold sources.

On the front of the unit are the USB port (1), manual vent valve (2), touchscreen (3), door opening lever (4), and stove power switch (5).

On the left side is the cover that protects all the unit's electronics, and on the back are the power cable connections (6), the vacuum inlet (7), the Ethernet connection (8), and the electrical connection for the vacuum pump (9). A couple of hoses have been included with the vacuum oven, one is for the vacuum pump connection (hose 12/22) for connection(7). The silicone tube (10/6), for the manual vent (2), in order to add any additional compound required in the process.

IN ORDER TO START THE OVEN, PRESS THE MAIN SWITCH AND WAIT APPROX. 30 SECONDS TILL THE SCREEN SHOWS THE JP SELECTA LOGO.

5. Operation

5.1. Home Menu.

The unit centralizes all its functions on the touch screen located on the front of the oven. The main functions are described below.

Back to the Home Page

Disconnect or change User

Allows to monitorize historical trends

Setting menu access

Recipe Menu

Allow to verify any alarm or historical.

Start / Stop thermal cycle

Manual activation of the venting valve to eliminate the vacuum inside chamber

In the Home Menu, the two trays and their temperatures are displayed. Physically can be disconnected the upper one (tray 2) to use the oven with larger containers.

To do this, remove the tray support, disconnect the corresponding connector, and unscrew the tray from the support.

In the Start Menu, click the corresponding tray 2 (if active, it will light up green).

5.2. Settings.

Once the device is in place and powered both electrically and pneumatically, can be activated by pressing the power button. By default, the device displays the Home Menu. By means the setting menu access, can be adjusted some parameters:


- **Language.**

From Home menu:

Settings → Language.

Choose the language at your convenience.


• **Date & time.**

From Home menu,

Settings → Date & time.

Please, modify date &/or time if required.

• **Users.**

From Home menu:

Settings → Users.


Within Users, can be created, modified, deleted, and managed the users who can log in to the device. This provides additional information in process reports, events, and alarms that occur.

There are different privileges according the user type.

	Unauthorized	Guest	Operator	Maintenance	Admin
Login, logout, modify date &time, visualize alarms & trends	x	x	x	x	x
Check alarms, modify offset, to do checks				x	x
Manage Users, modify internal parameters, backups					x

5.3. Communications (Ethernet).

The device can be connected to the local network to send service reports to a designated FTP server. To use this option, follow these steps:


Starting from the main menu, go to:

Configuration → Communications.


There, you must fill in the address, username, password, and port sections.

It is important to check the Ethernet configuration to ensure the device is on the network with access to the server described. To do so, within the same communications screen, press the Ethernet button to configure its parameters.

5.4. Calibration.

If required, can be applied an offset to the pressure and temperature sensors readings, as well as the ranges adjustment, by using the following key:


Configuration → Settings.

Within Settings, both the pressure and temperature sensors can be selected & modified

6. Operation Modes

Once the equipment is configured, with the sensors calibrated and connected to the network via Ethernet, the operation process can start.

6.1. Basic Mode


In the basic mode, enter the setting of temperature (between 20 and 200°C, always at least 5°C above ambient temperature), the required vacuum pressure (above 5 mbar), and the time to maintain the conditions indicated at the bottom of the main screen. The unit will start up immediately, aiming to reach the set values as quickly as possible, first ensuring the pressure and then the temperature. Once conditions are stable, the unit begins counting the stipulated time.

6.2. Custom Mode

In custom mode, the recipe structure will be used to guide the steps the oven should follow during the warm-up phase. A maximum of 6 setpoints are available, and a total of 50 recipes can be stored.


7. End of cycle

At the end of the heating cycle, the device automatically generates a report with the data obtained of the process.

Upon completion, it will wait for the user's notification and begin a new cycle.

Process reports are stored by default in the internal memory, with a maximum of 100 reports. Once the maximum is reached, it is necessary to extract the reports, as described in detail on chapter 8.

8. Memory backup

Reports generated and stored in internal memory can always be retrieved from the section:

Settings → Backup.

There are three different ways to retrieve the report:

a) USB Connection.

- Connect a USB flash drive to the unit's USB port and wait until it's recognized (FAT 32 format).

- Start downloading reports by pressing the indicated button and wait until it indicates that it is complete.

- Remove the USB flash drive and proceed with more copies or erase the stove's internal memory.

b) FTP.

- If the FTP server is configured, the address will appear on the screen. If it is correct, it's possible to send the information by pressing the indicated button.

- The process will start and, after finished, a notification will be generated.

- Can be done more copies or erasing the stove's internal memory.

c) Single process report.

It's possible to extract a report from a single process by following these steps:

- Before starting the process to be reported, connect the USB flash drive.
- Run the cycle normally and wait for the oven, to indicate that it has generated the report.
- Upon completion, will appear a notification that the process has ended and wait approximately 3 seconds.
- The USB flash drive can be removed, containing the process report for further analysis. You can extract a report from a single process by following these steps:

It can be done as many backup copies as required with no limits in any of the available options

The single process report doesn't eliminate the copy into the oven memory. It must be deleted manually.

8.1. Internal memory limits.

After reaching 100 reports or running out of free space, the device will prompt on screen to free up space in the internal memory. To do so, follow these steps:

After completing the necessary backups, press the Recycle Bin button.

Once finished, it's possible to return to the main menu and continue using the device.


Once deleted the internal memory content, the existing reports cannot be recovered

9. Maintenance.

Tampering the stove's internal electronic circuits by unauthorized personnel can cause damage that is difficult to repair. Be sure to take the equipment to one of J.P. SELECTA, S.A.U.'s authorized service centers.

9.1 Cleaning

To avoid long-term damage to both the wear and appearance of the working chamber, regular cleaning is important to avoid dust inside.

Avoid working with products that can cause rust. Rust particles can contaminate stainless steel. If contamination occurs through routine work, it should be cleaned and polished immediately.

To clean stainless steel and aluminum, use a cleaner that does not contain aggressive products (a universal household cleaner, for example).

To clean plastics, glass, and the front, use a household cleaner and a microfiber cloth.

9.2 Programmed maintenance

Daily

Check the condition of the lid gasket and its placement to prevent it from being pinched or cut when closing the lid.

Weekly

Check the condition of the internal wiring of the trays, ensuring that the cable does not have cuts in the sheath or excessive bends that could damage the internal cables.

Monthly

Check that the heater's power consumption is adequate.


Check that the vacuum hose and its connections and fasteners are in good condition.

Yearly

It is advisable to replace the door gasket if the equipment usage is high (more than 1 cycle per day).

10. Safety thermostat:

This safety device is located at the rear right hand of the vacuum oven, and limits the temperature inside the autoclave chamber.

8. Safety thermostat regulation screw (grey).
9. Safety thermostat reset button (white).

10.1 Rearming the safety thermostat

Safety thermostat is an element to prevent the temperature from exceeding a value. If by a fortuitous event the thermostat is triggered, that is, it opens the power circuit, the heating element will stop working. Usually there is an indicator light, or an error message, that indicates that the thermostat is in this position.

If the thermostat is triggered, check that the equipment looks good and it is not overheating at any point.

To change its status, simply press the white button (9), and you will hear a “click” that indicates that the thermostat has been reset. If the equipment has not cooled down sufficiently, the thermostat cannot be reset, and you must wait for the temperature to drop before it can be reset.

10.2 Adjusting the safety thermostat

In some equipment such as ovens and baths and with certain processes it may be useful to adjust the thermostat trigger value, but this operation should NEVER be performed to the rest of equipments.

From factory, the thermostat is adjusted to its trigger temperature, set over 200°C.

Stop the equipment and leave it to cool with the lid open for 15 minutes.

11. Recycling & equipment dismantling:

Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.

12 . Spares:

		
13007 Contactor CWC09	15551 Fuse IFO 520 10A	16215 Electrov. Vacuum NC 3/4" L153
		
16227 Electrov. 139V08	21065 Vaciotem Lid gasket	21066 Vaciotem chamber gasket
		
34135 Inlet Valve 1/2"	43017 Thermostat	36028 Solid State Relay 50A
		
16143 Power Supply	34086 Front manual valve (vacuum extraction)	15457 Front valve knob

		
<p>85284 PLC Screen</p>	<p>20017 Main Switch</p>	<p>59976 Oven tray</p>
		
<p>85286 PLC</p>	<p>85281 Converter A/D</p>	<p>85288 Temperature Transmitter</p>
		
<p>39374 Tray heater 400W</p>	<p>43121 Pressure Transmitter</p>	<p>43178 Thermostat Limit. 220°C</p>
		
<p>15743 Connector Female Heater. 7 pin</p>	<p>15744 Connector Male Heater 7 pin</p>	

13 . Error list & alarms:

Code	Name	Description
01	<p align="center">E2PROM FAIL</p> <p><u>Reason:</u> Controller Failure <u>Solution:</u> Contact aftersales dpt. JP Selecta</p>	(E2PROM_FAIL)
02	<p align="center">PRESSURE SENSOR FAIL</p> <p><u>Reason:</u> The pressure sensor signal is not as expected <u>Solution:</u> Revise sensor and wiring</p>	(PRES SENSOR FAIL)
14	<p align="center">VACUUM FAIL</p> <p><u>Reason:</u> It's not possible to achieve the vacuum <u>Solution:</u> Revise the vacuum system (centralized or pump), hoses and valves (open / close)</p>	(VACUUM FAIL)
19	<p align="center">UNREACHABLE SETTING</p> <p><u>Reason:</u> The oven cannot achieve the temperature setting after defined time <u>Solution:</u> Revise the power consumption from trays (400W each), and the solid state relays status.</p>	(SET_UNREACHABLE)
23	<p align="center">POWER FAIL</p> <p><u>Reason:</u> The equipment is recovering from a power fail. <u>Solution:</u> Revise the electrical connection (plug and cable), and the electrical protections (thermal breaker and differential switch)</p>	(POWER FAIL)
29	<p align="center">SAFETY THERMOSTAT</p> <p><u>Reason:</u> Safety thermostat is activated <u>Solution:</u> Revise the temperature setting and compare with a tested probe to verify the reading is correct.</p>	(SAFETY THERM)
57 & 60	<p align="center">HEATER TEMPERATURE PROBE - SHORTED</p> <p><u>Reason:</u> The reading from the trays heaters temperature probes are not as expected <u>Solution:</u> Trays heaters to be replaced</p>	(PROBE SHORT CIRCUIT)
58 & 61	<p align="center">HEATER TEMPERATURE PROBE - CUT</p> <p><u>Reason:</u> The reading from the trays heaters temperature probes are not as expected <u>Solution:</u> Trays heaters to be replaced</p>	(PROBE CUT WIRE)
59 & 62	<p align="center">OVERTEMPERATURE</p> <p><u>Reason:</u> The reading from the trays heaters temperature probes are not as expected (exceeding 250°C) <u>Solution:</u> Trays heaters to be replaced</p>	(OVERTEMP)

14. Warranty:

This VACIOTEM-GMP vacuum oven has a one-year warranty.

The warranty does not cover damage caused by improper use or causes beyond the control of JP Selecta.

Any tampering with components or the equipment circuit by personnel not authorized by JP Selecta automatically voids the warranty.