| | |
|-----------------------|----------------|
| UNIVEBA 110°C | 3000400 |
| UNIVEBA 200°C | 3000401 |
| TERMOBAT 110°C | 3000587 |
| TERMOBAT 200°C | 3000588 |
| AQUATERM | 3000916 |
| BAHER | 7000540 |
| TERMOFIN | 3000459 |

BAÑOS TERMOSTÁTICOS
THERMOSTATIC BATHS

INFORMACIÓN GENERAL**GENERAL INFORMATION**

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The standard equipment consists of the following components:

| DESCRIPCIÓN / DESCRIPTION | CÓDIGO / CODE | Cant/Qty. |
|---|---------------|-----------|
| UNIVEBA 110°C | 3000400 | 1 |
| UNIVEBA 200°C | 3000401 | 1 |
| TERMOBAT 110°C | 3000587 | 1 |
| TERMOBAT 200°C | 3000588 | 1 |
| AQUATERM | 3000916 | 1 |
| BAHER | 7000540 | 1 |
| TERMOFIN | 3000459 | 1 |
| Termómetro (solo TERMOFIN) / Thermometer (only TERMOFIN) | 1000120 | 1 |
| Tapa en aluminio teflonada (solo TERMOFIN) / Cover with P.T.F.E. coated aluminium (only TERMOFIN) | -- | 1 |
| Manual de instrucciones / Instruction manual | 80005 | 1 |

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Voltage supply 115-230V 50/60 Hz according to the characteristics plate indications.

| Código/ Code | Temperatura máxima °C/ Maximum temperature °C | Capacidad litros/ Capacity litres | Medidas útiles / Tank dimensions cm | | Medidas exteriores / Overall dimensions cm | | Consumo Consump. (W) | Estabilidad* Stability (K) | Peso Weight (Kg.) |
|-----------------|--|--|---|------------|--|-------------|----------------------------|----------------------------------|-------------------------|
| | | | Ø | Alto/Depth | Ø | Alto/Height | | | |
| 3000400 | 110 | 2,5 | 18 | 13 | 25 | 21 | 400 | ±1 | 2,5 |
| 3000401 | 220 | | | | | | 700 | | |
| 3000459 | 80 | 4 | 25 | 10 | 34 | 14 | 400 | --- | 2,3 |
| 3000587 | 110 | 11 | 29 | 17 | 29 | 25 | 1500 | ±2 | 5,5 |
| 3000588 | 200 | | | | | | 1700 | ±3 | |
| 3000916 | 110 | 1,5 | 20 | 11,5 | 22 | 20,5 | 700 | ±2 | 1,3 |
| 3000540 | 110 | 1 | 14 | 11 | 14 | 19 | 1500 | ±3 | 1 |

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EQUIPMENT DESCRIPTION

UNIVEBA baths are designed for test with rotatory evaporator and water and oil bath, for adjustable temperatures from ambient +5°C up to 110°C (Stability $\pm 1^\circ\text{C}$) or up to 200°C (Stability $\pm 1^\circ\text{C}$), depending on model.


TERMOBAT baths are designed for testing with rotatory evaporator (for flasks up to 10 litres) and water and oil baths, for adjustable temperatures from ambient +5°C up to 110°C (Stability $\pm 2^\circ\text{C}$) or up to 200°C (Stability $\pm 3^\circ\text{C}$), depending on model.

AQUATERM water bath, for adjustable temperatures from ambient +30°C up to 110°C (Stability $\pm 2^\circ\text{C}$).

BAHER boiling bath, with rapid boiling type (6 to 7 minutes), for adjustable temperatures from ambient +30°C up to 110°C (Stability $\pm 3^\circ\text{C}$).

TERMOFIN histology bath for paraffin, for melting and keeping paraffin in a liquid state (to allow the viewing of tissue sections), for adjustable temperatures from ambient +40°C up to 80°C.

To use the baths, fill the inner tank with liquid until the level covers at least the heater element. The liquid or mixture to be employed must be adequated to the bath model and the working temperature. As a guidance, here are some liquids recommended:

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INSTALLATION

Place the bath on a flat horizontal, level surface, leaving a space of 10 cm at the back and on each sides of the equipment.

Do not use the equipment if it is not earthed.

Be sure that the voltage supply is the same that the one indicated on the characteristics plate.

If you change the plug, bear in mind the following:

Blue cable: **Neutral.**
Brown cable: **Phase.**
Yellow/green cable: **Earth.**

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It is recommended the use of decalcified water. DO NOT USE DEIONIZED OR DEMINERALIZED WATER as it has a great tendency to rust. To prolong the life of the bath and avoid rusting of the parts, dissolve 1 gr of CO₃Na₂ (Sodium Carbonate) per litre of water used in the bath.

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

OPERATION

CONTROL PANEL UNIVEBA:

1. Main switch.
2. Temperature control thermostat.
3. Heater operation indicator lamp (amber).
4. Safety thermostat operation indicator lamp (red).

REAR PART:

5. Safety thermostat control screw.
6. Safety thermostat manual reset button (grey).

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
CONTROL PANEL TERMOBAT:


1. *Main switch.*
2. *Temperature control thermostat.*
3. *Heater operation indicator lamp (amber).*
5. *Mains indicator lamp (green).*

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CONTROL PANEL AQUATERM AND BAHER:


2. *Temperature control thermostat.*
3. *Heater operation indicator lamp (amber).*

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COMMON INDICATIONS FOR ALL BATHS:**STARTING UP:**

1. Fill the bath with liquid.
2. Connect the bath into the mains .
3. Turn on the main switch (1). The green indicator lamp (5) will light up.
4. Turn on the temperature control thermostat (2) from low to high according to the working temperature. The heater operation indicator lamp (3) (amber) will light up.
5. After some time, the heater operation indicator lamp (3) will begin to operate intermittently. When the cycle of these intermitences becomes uniform, that will indicate that temperature of the liquid is stabilized.
6. The thermostat scale graduation is only useful as a reference, because the real temperature is the one shown by a thermometer introduced into the liquid.
7. Be careful with liquid evaporation. The heater element must ALWAYS be covered with liquid. Replenish when necessary.

ADJUSTMENT FOR THE SAFETY THERMOSTAT:
(Only for the baths equipped with it)

The low limit of working temperature for the safety thermostat is nearly 50°C. The desired temperature of the safety thermostat must be at least 5°C above the working temperature. The adjusting process of the safety thermostat is long. For that reason, the safety thermostat is useful when the users normally work at the same temperature.

1. Turn the safety thermostat screw (5) clockwise to the maximum.
2. Select with the temperature control knob (2) the desired temperature of the safety thermostat.
3. Turn on the bath and leave it to stabilize.
4. When the bath is stabilized, slowly turn the safety thermostat screw (5) anticlockwise until it clicks and the red indicator lamp (4) lights up.
5. Wait until the bath cools to reset the safety thermostat with (6) button.
6. Select the working temperature of the bath with the temperature control knob (2).

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To guarantee the safety of the equipment, all spare parts must be purchased from J.P. SELECTA, s.a.u.

| Descripción / Description | Código / Code | | | | | | |
|--|---------------|---------|---------|---------|---------|---------|---------|
| | 3000400 | 3000401 | 3000587 | 3000588 | 3000916 | 7000540 | 3000459 |
| Modelo / Model | 3000400 | 3000401 | 3000587 | 3000588 | 3000916 | 7000540 | 3000459 |
| Resistencia Heating element | 39063 | 39049 | 39318 | 39318 | 39049 | 39049 | 39007 |
| Termostato de regulación Control thermostat | 43006 | 43005 | 43006 | 43005 | 43006 | 43006 | 43002 |
| Termostato de seguridad Safety thermostat | 43031 | 43017 | NO | | | | |
| Racord termostato 110°C Thermostat racord 110°C | 34140 | NO | 34140 | NO | | | |
| Racord termostato 200°C Thermostat racord 200°C | 34141 | | NO | 34141 | | NO | |
| Tuerca racord termostato Thermostat racord nut | 34142 | | | | | | NO |
| Alargo Schuko Schuko cable | 7002 | | 7001 | | 7002 | | |

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Before removing the cover, disconnect the equipment from the mains. The control and heating elements must only be manipulated by authorized personnel.

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The distilled or demineralized water has a great tendency to rust. To prolong the life of the bath and avoid rusting of the parts, dissolve 1 gr of CO_3Na_2 (sodium carbonate) per litre of water used in the bath.

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CLEANING:

If normal tapwater is used, calcareous deposits will appear in the bath and in the heating element. A heater covered with lime should not be used as it can lose heating power and deteriorate more quickly. Periodical cleaning of the inside of the tank and the heating element with a solution of 10% of hydrochloric acid (28% chlorhydric acid) or a solution of water and vinegar (acetic acid) is advised. Rinse with clean water afterwards.

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For the cleaning of the different parts use the following products:

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Cleaning of stainless steel: alcohol

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Cleaning of plastic: alcohol with cotton duster.

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GUARANTEE


This product is guaranteed for one year. The guarantee does not cover damages caused by an incorrect use or causes beyond the control of J.P. SELECTA, s.a.u.

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
Any manipulation of the equipment by unauthorized personnel by J.P. SELECTA, s.a.u. cancels the guarantee automatically.

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