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Index

| | |
|--|-----------|
| 1. General Information | 12 |
| 2. Packing list | 12 |
| 3. Technical specification | 12 |
| 4. Information about safety in the use of the equipment | 13 |
| 4.1 Risks the operator is subjected to | 13 |
| 4.2 Staff qualification..... | 13 |
| 4.3 Use with hazardous substances | 13 |
| 4.4 Maintenance..... | 13 |
| 4.5 Warranty..... | 13 |
| 5. Equipment Description | 14 |
| 6. Installation | 14 |
| 6.1 Bucket/Tank | 14 |
| 6.2 Connection to the mains | 14 |
| 7. Operation | 14 |
| 7.1 Quick operation | 16 |
| 7.2 Setting screen | 17 |
| 7.3 Programming screen..... | 17 |
| 7.4 Selecting a program for operation | 18 |
| 7.5 Operation screen | 18 |
| 7.6 Circulation pump | 18 |
| 7.7 External probe..... | 19 |
| 7.8 Communications | 19 |
| 7.9 Alarms | 19 |

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ENGLISH VERSION

1. General Information

- Handle the parcel with care.
- Check that contents coincide with the packing-list. If any part is damaged or missing, please advise the distributors immediately.
- Do not install or use the equipment without reading this handbook first. This instruction manual must be attached to the equipment and must be available for all users.
- If you have any doubts or enquiries, contact with J.P. SELECTA's technical service.
- Do not use the equipment in explosive, flammable or corrosive atmospheres.

2. Packing list

Units are supplied with their own manual.

3. Technical Specification

| | 3000545 | 3000554 |
|-------------------------------------|--|---|
| Dimensions: Height/width/depth (mm) | 360/280/360 | 360/280/360 |
| Power supply voltage: | 230 VAC 50Hz | 230 VAC 50Hz |
| Electrical power (W): | 2060 | 2060 |
| Weight: (Kg) | 9,5 | 9,5 |
| Fuse: | Replaceable (Type: 10x5mm) | |
| Temperature range: | From... ambient T ^a + 5°C ... up to ... 200°C | |
| | WITH EXTERNAL REFRIGERATION from -20°C up to 200°C | |
| Temperature resolution: | 0,1°C | 0,1°C |
| Timed operation: | Hours, minutes | Hours, minutes |
| Pump flow: | 10L/min | 10L/min |
| Case: | in epoxy painted steel | in stainless steel with side carrying handles |

Fig 1. Equipment piping

4. Information about safety in the use of the equipment

**Hazard
risk**

**Electric
risk**

**Risk of high
temperature contact**

4.1 Risks the operator is subjected to

- Risk of burns
- Electrical risk

4.2 Staff qualification

This equipment may only be used by properly trained personnel who have read and understood these instructions or have been properly trained in this equipment operation.

4.3 Use with hazardous substances

Do not treat flammable or explosive materials or those which chemically react with violence.

ATTENTION!

Burning risk

Take proper precautions
(use of thermal gloves, clamps,...)
to extract samples from inside the tank.

It can be at a high temperature.

4.4 Maintenance

Unplug the equipment from the mains before you handle inside it.

**The manipulation of the equipment internal electronic circuits by unauthorized personnel can cause irreparable damage.
Take it to one of J.P.SELECTA's authorized technical services.**

It is recommended to clean the accessible parts of the equipment with a cloth and a cleaner, or a similar solution. Do not use aggressive solvents.

4.5 Warranty

This product is guaranteed for one year. It does not cover damages caused by an incorrect use or causes beyond the control of J. P. SELECTA, s.a.u.

Any manipulation of the equipment by personnel not authorized by J.P. SELECTA, s.a.u. automatically cancels the guarantee.

Level switch inhibits the thermostat operation if no liquid level is detected.

5. Equipment Description

Immersion thermostats are widely used on laboratories for thermal treatment of samples, or heating/cooling other equipments.

The equipment includes a circulation pump that can be used to improve the homogeneity within a bath or to perform a closed liquid circulation circuit for cooling or heating.

The unit could work in a continuous or timed operation.

6. Installation

Unpack the equipment and visually check that the contents are in good conditions. Keep the package for some days.

To work in a comfortable way, please leave a free space of 50cm around the equipment and a minimum of 20cm at the rear part.

6.1 Piping connection

Connect the water inlet and outlet pipes (see Fig.1) from the thermoregulation circuit to the equipment.

Make sure the drain valve is closed.

6.2 Connection to the mains

Plug the equipment to a near socket suitable for the equipment power. Use the connection cable supplied or a similar one.

For your safety, the equipment must be earthed.

Before plugging the machine, check that the voltage supply and the frequency corresponds to the one indicated in the equipment characteristics plate.

7. Operation

All the equipment functions are set with the touch screen. Gently press with your finger in the center of each icon.

When starting the computer, the starting screen will be shown (Fig 7.1).

- **SETTINGS:**

- Time and date: menu to set the date and time.
- Temperature correction.
- Temperature alarm limits.

- **PROGRAMMING**

- Allows you to program the working temperature & time.

- **TASK**

- Allows you to select and launch the operating temperature, the working time and the pump flow.




Fig 7.1 Initial screen

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







1 Press over «TASK».

2 Press over temperature display (Frame changes to yellow color)

3 Set the temperature with (40°C on the picture)

4 Press over .

Notes:

In a similar way, hours and minutes of working time could be programmed.

The temperature set from this screen does not remain recorded once equipments are switched off. (See paragraph 7.2)


Fig 7.2 Setting screen

7.2 Setting screen

7.2.1 Set time and date

Press on Date/time button. Press successively on the box to select the day, month, year, hour, minute and second settings.

Use Up and Down icons to entry a value.

7.2.2 Set temperature, language and screen contrast

In a similar way, set the temperature, language and screen contrast by pressing the icons.

7.2.3 Working parameters protected under password

There are several parameters under password "0000" which affect the equipment operation.

. Press "Insert password" and follow the instructions:

password: 0000

. Press CAL to access these parameters:

log interval: 1 to 60min
 external sensor: yes/no
 calibration: -100 to +100
 lower temp. limit: -10
 upper temp. limit: -10

a) Log interval

It is the frequency of data collection for the logging.

b) External sensor

Select the value **YES** if you want to work with an external probe connected to the connector on the back side of the equipment

c) Calibration (OFFSET)

Use only this setting if you notice a deviation between the temperature indicated by the equipment and the one indicated by a standard thermometer.

For example: the equipment indicates 40.0°C, the pattern indicates 40.5°C. Enter a value of -0.5°C to correct the deviation.

d) The lower and the upper limit of temperature

These are the maximum and the minimum limits of temperature adjustments. This allows you to adjust the working temperature range of the equipment.

e) To update the firmware of the device, click on "Info" and the UPGRD button will appear at the bottom of the screen.

Contact the technical service to obtain the updated file.

7.3 Programming screen

Up to 10 programmes of temperature and time can be stored.

Each program has the following parameters:

Temperature: from ambient +5°C ... 200°C
 Time: HH:MM
 Undertemp: Alarm value for temperature drop.
 Overtemp: Alarm value for overtemperature.

From main menu, press over «PROGRAMMING».

Selected program. Program selector.


Fig 7.3. Programming screen


Fig 7.4. Selecting a program


Fig 7.5. Operation screen


Fig 7.6.1 INTERNAL circulation.


Fig 7.6.2 EXTERNAL circulation.

Fig 7.3 screen is displayed.

Select a program (from 00 to 09) with keys:

Set working temperature and the other parameters:

Press:

Keys will be shown to select temperature.

In a similar way press:

Sets de working time.
A time setting value of 00.00 means «continuous operation».

Sets the undertemperature alarm value.

Sets the overtemperature alarm value.

Press to save .

7.4 Selecting a program for operation

To select a stored program, press over SEL (Fig 7.4). (SEL indication changes its color to yellow).

Once selected, this program in placed on the TASK screen.

7.5 Operation screen

On the operation screen «TASK», the program selected (PROG 04) on the programming screen is shown (Fig 7.5) and it starts when pressing over .

This screen shows:

- Actual temperature
- Temperature set
- Working time set (SET)
- Elapsed time

7.6 Circulation pump

The liquid circulation pump could operate in two configuration modes:

- INTERNAL circulation:

Water (or other) will be pumped from horizontal pump output. (See fig 7.6.1).

Use this configuration to make water circulation into the bath vessel in order to homogenize the temperature from sample load.

- EXTERNAL circulation:

Water (or other) will be pumped from vertical pump output (See fig 7.6.2).

Use this configuration to made water circulation out of the bath vessel in order to temperate an external equipment. This is the factory setting.

7.7 External probe

The process can be controlled by means of an external probe.

Connect the probe to the rear connector and place it in the media you want to control the temperature. To adjust it, follow the steps in section 7.2.3.

7.8 Communications

7.8.1 USB

Access the SETTINGS > LOGS DOWNLOAD screen to download the logs of the processes carried out. Then follow the instructions that appear on the screen.

7.8.2 RS-232

RS232 port configuration. RS232 communication protocol based on technical specification SPT-0004-B.

| | |
|--------------------|-----------|
| Transmission speed | 56.000bps |
| Data | 8 bits |
| Stop | 1 bit |
| Parity | No parity |

7.9 Alarms

| Error | Screen | Description | Type |
|-------|-----------------------------|-------------------------------|------|
| 01 | Level | Safety or level thermostat | A |
| 02 | Overtemp | High temperature level | A |
| 03 | Undertemp | Low temperature level | A |
| 04 | External RTD opened | Open external probe | A |
| 05 | External RTD shortcircuited | Shortcircuited external probe | A |
| 06 | Internal RTD opened | Open internal probe | A |
| 07 | Internal RTD shortcircuited | Shortcircuited internal probe | A |
| 08 | Power Fail | Without net | B |

Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.

J.P. SELECTA, s.a.u. 2020