DIGIT-COOL

3001373

**ULTRATERMOSTATOS DE CIRCULACIÓN
CIRCULATION ULTRATHERMOSTAT**

INDICE GENERAL

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GENERAL INFORMATION

- 1) *Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributor immediately.*
- 2) *Do not install or use the equipment without reading this handbook before.*
- 3) *This handbook must always be attached to the equipment and it must be available for all users.*
- 4) *If you have any doubts or enquiries, please contact your supplier or J.P. Selecta's technical service.*
- 5) **IMPORTANT! J.P. SELECTA WILL NOT ACCEPT ANY EQUIPMENT TO BE REPAIRED IF IT IS NOT DULY CLEANED.**
- 6) *If any modification, elimination or lacking in maintenance of any device of the equipment by the user transgress the directive 89/655/CEE, the manufacturer is not responsible for the damage that can occur.*
- 7) *Do not use the equipment with liquids which can give off vapours capable of making explosive mixtures.*

PACKING LIST

The standard equipment consists of the following components:

Descripción / Description	Codigo / Code	Cant. / Qty.
DIGIT-COOL	3001373	1
Manual de instrucciones / Instruction Manual	80128	1

ACCESSORIES

Stainless steel tanks:

Internal dimensions cm:

Code	Litres	Depth	Length	Width
6003921	20	15	48	30
6003931	27	20	48	30

Overall dimensions cm:

Code	Depth	Length	Width
6003921	24	61	41
6003931	29	61	41

Ring set lid: 105, 80, 60 and 37mm diameter reduction rings and a thermometer location hole in the lid.

Code: 1001374

Gablet lid. Code: 3001295

Tube rack support. Capacity for 3 racks.

Code: 1001296

Tube racks:

Code	Width	Tubes Ø mm	No. tubes	Depth/Length/		
1001202		13	36	8	8.2	23.6
1001203		16	24	8	8.2	23.6
1001204		20	24	8	8.2	23.6

TECHNICAL FEATURES

Voltage supply 115-230V 50/60 Hz according to the characteristics plate indications.

Código / Code	3001373		
Rango de regulación °C / Regulation range °C	De 0 a 100 From 0 up to 100		
Estabilidad °C / Stability °C	± 0,1		
Homogeneidad / Homogeneity	± 0,1		
Error de consigna /	± 0,25	± 0,25	± 1
Resolución / Resolution	± 0,1		
Medidas exteriores en cm. Overall dimensions cm.	Alto / Height	45	
	Ancho / Width	21	
	Fondo / Depth	51	
Bomba Pump	Presión / Pressure	150 mbar	
Potencia HP / Compressor HP	Caudal H ₂ O / H ₂ O flow	12 l/min	
Consumo W. / Consumption W.		3/8	
Peso Kg. / Weight Kg.		1460	
		22	

EQUIPMENT DESCRIPTION

DIGIT-COOL immersion thermostat allow regulation with a stability of 0.1°C of the liquid temperature from 0°C up to 100°C.

The thermostats are equipped with:

- *Temperature sensor by Pt100 probe.*
- *Circulation pump with adjustable flow for agitation and external temper.*
- *Cooling coil for regulation at ambient temperatures through water tap connection.*
- *Safety thermostat and fluid level float.*
- *Adjustable upper and lower alarm temperature .*

- *Adjustable maximum and minimum working temperature.*
- *Correction of adjustable reading.*
- *RS232 Connection.*
- *Reading in Fahrenheit or Celsius degrees..*


FILLING THE BATHS

For the correct use of immersion thermostats, it is necessary first to fill the tank to be used with liquid

up to the adequate level, i.e. for covering at least the heating element and the circulation pump. Otherwise the thermostat does not work by the floating effect.

The liquid or mixture to be used must be adequate to the model of thermostat and the temperature range.

As a guidance, here are some recommended liquids, however, other similar liquids can be used:

***Silicone oil 50cps**

INSTALLATION

Place the DIGIT-COOL on the right part of the tank; and over a flat horizontal levelled surface, trying to leave a free space at the back of the equipment of 10 cms.

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Be sure that the voltage supply is the same as the one indicated on the characteristic plate of the equipment.

Do not use the equipment if it is NOT earthed.

If you have to change the plug, bear in mind the following:

Blue cable: Neutral.

Brown cable: Phase.

Yellow/green cable: Earth.

It is recommended the use of decalcified water. DO NOT USE DEIONIZED OR DEMINERALIZED WATER as it has a great tendency to rust. To prolong the life of the bath and avoid rusting of the parts, dissolve 1 gr of CO_3Na_2 (Sodium Carbonate) per litre of water used in the bath.

CIRCULATION

Fig 7.6.1 **INTERNAL** circulation.

Circulation pump

The liquid circulation pump could operate in two configuration modes:

- INTERNAL circulation:

Water (or other) will be pumped from horizontal pump output.
(See fig 7.6.1).

Use this configuration to made water circulation into the bath vessel in order to homogenize the temperature from sample load.

Fig 7.6.2 **EXTERNAL** circulation.

- EXTERNAL circulation:

Water (or other) will be pumped from vertical pump output.
(See fig 7.6.2).


Use this configuration to made water circulation out of the bath vessel in order to temperate an external equipment.

CAUTION!!! IMPORTANT FOR YOUR SAFETY
*When working at higher temperatures, once finished it is recommended to switch off the heater leaving the stirring on.
Thus, we avoid the system filtration on the machine.*

OPERATION

CONTROL PANEL:

1. Display.
2. Main switch.
3. "Enter" push button.
4. Parameters push button.
5. Increase set value push button.
6. Decrease set value push button.
7. Start /Stop push button.

DEFAULT CONFIGURATION PARAMETERS:

- | | |
|--|----------------------|
| - Overtemperature alarm | OVERTEMP = OFF. |
| - Undertemperature alarm | UNDERTEMP = OFF. |
| - Minimum working temperature | LOWER LIM = 0°C |
| - Maximum working temperature
(depending on model). | UPPER LIM = 100°C |
| - Scale selection ^o | SELECCIÓN = CELSIUS. |
| - Calibration | CALIBRATE = 0°C. |


STARTING UP:

1. Fix the thermostat at one side of the tank, on the right part.
2. Fill the tank with liquid until it pushes the level float.
3. Connect the thermostat into the mains.
4. Turn on the main switch(2):
 - Heating position only if one will work with temperatures above ambient temperature;
 - Hot / Cold position if one will work with temperatures lower or equal to ambient temperature.
5. Select the temperature (SV) with buttons (5) and (6).
6. Press (7) to start up. While working, PV: (Point value) shows the real fluid temperature and SV: (Set point) shows the desired temperature.

MECHANICAL SAFETY THERMOSTAT

The safety thermostat is adjusted at factory.

If it starts, stop the machine, let the bath cool and reset the thermostat by pressing the grey PVC washer placed at the rear part.

REAR PART

- 9. Safety thermostat control screw.
- 10. Safety thermostat manual reset button (grey).
- 11. Fuse carrier.

ACCESS TO CONFIGURATION PARAMETERS:

OVERTEMPERATURE

- Press (7) to stop the thermostat (STATUS: STOP).
- Press (3) once, and on the display you can see: **"OVERTEMP"**, by pressing (5) or (6) the overtemperature degrees are selected. The alarm is activated when the read temperature (PV) exceeds the overtemperature degrees of the selected temperature (SV). The heating is cut off.
- The following screen appears:

OVERTEMP	STATUS
PV	ALARM

where PV is the present temperature read.

- To remove the screen indications, press (7), let the fluid cool and press (7) again to start.

UNDERTEMPERATURE

- Press (7) to stop the thermostat (STATUS: STOP).
- Press (3) once, and on the display you can see: **"UNDERTEMP"**, by pressing (5) or (6) the undertemperature degrees are selected. The alarm is activated when the read temperature (PV) is below the undertemperature degrees of the selected temperature (SV). The heating is cut off.

- *The following screen appears:*

UNDERTEMP STATUS
PV ALARM

where PV is the present temperature read.

- *To remove the screen indications press (7), let the fluid cold down and press (7) again to start.*

MINIMUM WORKING TEMPERATURE

- *Press (4) until the display shows **LOWER LIM**.*
- *By pressing (5) and (6), select the desired value.*
- *Press (3) several times until the display shows the initial screen.*

MAXIMUM WORKING TEMPERATURE

- *Press (4) until the display shows **LOWER LIM**.*
- *Press (3) until the display shows **UPPER LIM**.*
- *By pressing (5) and (6), select the desired value.*
- *Press (3) several times until the display shows the initial screen.*


WORK WITH FARENHEIT OR CELSIUS SCALE

- *Press (4) until the display shows **LOWER LIM**.*
- *Press (3) until the display shows **SELECT**.*
- *By pressing (5) and (6), select the desired scale.*
- *Press (3) several times until the display shows the initial screen.*

CALIBRATION

- Press (4) until the display shows **LOWER LIM**.
- Press (3) until the display shows **CALIBRATE**.
- By pressing (5) and (6), select the desired value.
- Press (3) until the display shows the initial screen.

Access to the configuration parameters:

ALARMS THAT CAN APPEAR ON THE DISPLAY:

OVERTEMPERATURE: OVERTEMP (See page 7).

UNDERTEMPERATURE: UNDERTEMP (See pag 7).

LEVEL:

The screen shows:

LEVEL	STATUS
FAIL	ALARM

- 1- Cause: Low fluid level.
Effect: Cut heating off.
Action to solve:
 - 1) Stop the equipment by pressing (7).
 - 2) Add fluid.
 - 3) Start the equipment by pressing (7).

- 2- Cause: The mechanical thermostat starts
Effect: Cut heating off.
Action to solve:
 - 1) Stop the equipment by pressing (7).

- 2) *Reset the safety thermostat placed at the rear part of the machine. Review the temperature adjustment of the thermostat.*
- 3) *Start the equipment by pressing (7).*

PROBE:

1- The screen shows:

RTD FAIL	STATUS
INT OPEN	ALARM

Cause: Temperature probe open.

Effect: Cut heating off.

Action to solve:

- 1) *Stop the equipment by pressing (7).*
- 2) *Disconnect it from the mains.*
- 3) *Take it to the Technical Service.*

2- The screen shows:

RTD FAIL	STATUS
INT SHORT	ALARM

Cause: Shorted external temperature probe.

Effect: Cut heating off.

Action to solve:

- 1) *Stop the equipment by pressing (7).*
- 2) *Disconnect it from the mains.*
- 3) *Take it to the Technical Service.*

3- The screen shows:

RTD FAIL	STATUS
EXT OPEN	ALARM

Cause: External temperature probe cut or disconnected.

Effect: Cut heating off.

Action to solve:

- 1) *Stop by pressing (7).*
- 2) *Check that the external probe is correctly connected.*
- 3) *Start by pressing (7).*

If the failure persists:

- 1) *Stop the equipment by pressing (7).*
- 2) *Disconnect it from the mains.*
- 3) *Take it to the Technical Service.*

4- The screen shows:

RTD FAIL	STATUS
EXT SHORT	ALARM

Cause: Shorted external temperature probe.

Effect: Cut heating off.

Action to solve:

- 1) *Stop the equipment by pressing (7).*
- 2) *Check that the external probe is correctly connected.*
- 3) *Start the equipment by pressing (7).*

If the failure persists:

- 1) *Stop the equipment by pressing (7).*
- 2) *Disconnect it from the mains.*
- 3) *Take it to the Technical Service.*

MAINS FAILURE:

The screen shows:

PV...
SV...**POWER
FAIL***Cause: Mains failure.**Effect: The message **POWER FAIL** appears on the screen and the thermostat is still working.**Action to solve:*

- 1) Pressing button (7) the message **POWER FAIL** disappears.
- 2) Press button (7) again to start up the equipment.

SPARE PARTS

To guarantee the safety of the equipment, the spare parts must be purchased from J.P.SELECTA, s.a.u.

Descripción / Description	Código / Code
Resistencia calefactora / Heating element	39019
Circuito regulador / Control circuit	
Circuito display / Display circuit	
Circuito pulsadores / Button circuit	
Boya / Float	36020
Sonda Pt100 / Pt100 probe	43059
Motor / Motor	24017

GUARANTEE

This product is guaranteed for one year. The guarantee does not cover damages caused by an incorrect use or causes beyond the control of J.P. SELECTA, s.a.u.

Any manipulation of the equipment by unauthorized personnel by J.P. SELECTA, s.a.u. cancels the guarantee automatically.

COMUNICACION RS-232

RS-232 COMUNICACION OVERVIEW

- 9600 bauds, 7 data bits, 1 bit stop, without parity no handshake.
- Up to 15m of cable.
- Cable connection DTE to DTE 9 pins male.
- Protocol ASCII with LRC against transmission errors
- MASTER-SLAVE Communication (PC computer as a master).
- The device never sends anything if the computer does not ask for it.

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- Messages starts with character #.
- Device address must be included.
This field goes from 00 to 99, size always is 2 bytes in decimal-ASCII format.
- Data field is the message we want to send.
- LRC is a calculated value to check good transmission.
- At the end of the frame always a CR (0D) followed by LF (0A).
- Slave device answers always sending its address

LRC ERROR CHEKING FIELD

The frame includes a transmission error cheking field called LRC. This field sizes 2 bytes. This field is calculated by the sender device and includes it on the frame. The receiver device also calculates it and if they are equal, the transmission was good. If not, a communication error has occurred.

LRC calculation procedure:

1. Convert each character on its hexadecimal value.
2. Add them.
3. Take the smallest byte (the first two characters starting from the righthand side).
4. Subtract the previous result from FF.
5. Add 1.

LRC CALCULATION EXAMPLE

#01SVT<BL>+098,037<CR><LF>

ASCII	HEX		
#	23		
0	30		
1	31		
S	53		
V	56		
T	54		
<BL>	20		
+	2B		
0	30		
9	39		FF
8	38	-	C9
,	2C		36
0	+ 30	+	1
	2C9	LRC =	37

- «Get» commands always ends by «?»
- Characters always in CAPITAL.
- Use spaces only when indicated by <BL>.
- Decimal separator is «,».

CONNECTION CABLE

*It could be found is any specilized purchaser:
Cable Null - 9 pins modem*

PC	Device
1 _____	7, 8
2 _____	3
3 _____	2
4 _____	6
5 _____	5
6 _____	4
7 _____	1
8 ┌───┐	

Connectors in both cable sides are 9 pins female.

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MAINTENANCE

Before removing the side cover, disconnect the equipment from the mains.

If tapwater is used, there will be calcareous depositions on the bottom of the tank and on the heating element. Use of the heating element must be avoided when covered with these depositions because it loses calorific power and it is destroyed quicker.

It is advisable to clean the inner tank and the heating element periodically with a solution of 10% of hydrochloric acid at 28% or with a solution of water and vinegar (acetic acid). Then rinse it with clean water.