**UNITRONIC-ORBITAL  
UNITRONIC-ORBITAL C**

**6032011  
6001173**

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**GENERAL INFORMATION**

- 1) *Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributor immediately.*
- 2) *Do not install or use the equipment without reading this handbook before.*
- 3) *This handbook must always be attached to the equipment and it must be available for all users.*
- 4) *If you have any doubts or enquiries, please contact with your supplier or J.P. Selecta's technical service.*
- 5) **IMPORTANT! J.P. SELECTA WILL NOT ACCEPT ANY APPARATUS TO BE REPAIRED IF IT IS NOT DULY CLEANED.**
- 6) *If any modification, elimination or lacking in maintenance of any device of the equipment by the user transgress the directive 89/655/CEE, the manufacturer is not responsible for the damage that can occur.*
- 7) *Do not use the apparatus with liquids which can give off vapours capable of making explosive mixtures.*

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The standard equipment consist of the following components:

Descripción / Description	Código / Code	Cant. / Qty
UNITRONIC-ORBITAL	6032011	1
UNITRONIC-ORBITAL C	6001173	1
Muelle largo / Large spring	25024	4
Muelle corto / Short spring	25025	4
Manual de instrucciones / Instruction manual	80063	1

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**ACCESSORIES**

( They must be installed in the factory )

Denominación	Código	Denomination	Code
TAPA FORMA TEJADILLO En acero inox. AISI 304.	6000713	GABLE COVER In AISI 304 stainless steel.	6000713
ADAPTADORES PARA MATRACES ERLLENMEYER Construidos en acero inox. al temple.		ADAPTERS FOR ERLLENMEYER FLASKS In tempered stainless steel.	
* Adaptador para 50ml.	1001003	* Adapter for 50ml	1001003
* Adaptador para 100ml.	1000978	* Adapter for 100ml	1000978
* Adaptador para 250ml.	1000979	* Adapter for 250ml	1000979
* Adaptador para 500ml.	1000980	* Adapter for 500ml	1000980
* Adaptador para 1000ml.	1000981	* Adapter for 1000ml	1000981
PLATAFORMA PARA MATRACES ERLLENMEYER 3001004		PLATFORM FOR ERLLENMEYER FLASKS 3001004	

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In AISI 304 s.s. with multiadapter holes to hold. Erlenmeyer flasks of 50, 100, 250, 500 and 1,000 ml capacities.

Código Code	Medidas plataforma cm Platform dimensions cm		Capacidad de adaptadores en plataforma para: Adapter capacity in platform for:					Peso Weight Kg
	Ancho / Length	Fondo / Width	50ml	100ml	250ml	500ml	1,000ml	
3001004	23.5	34.5	16	12	6	5	3	1

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*Voltage supply 115-230V 50/60 Hz according to the characteristics plate indications.*

Código Code	Capacidad litros Capacity litres	Medidas utiles cm. Tank dimensions cm			Medidas de la cubeta móvil cm. Moving tank dimensions cm			Medidas exteriores cm. Overall dimensions cm			Consumo Consump. W	Peso Weight Kg
		Alto / Depth	Ancho / Length	Fondo / Width	Alto / Depth	Ancho / Length	Fondo / Width	Alto / Depth	Ancho / Length	Fondo / Width		
6032011	27	18	48	29	17	35	24	30	74	40.5	1500	25
6001173	27	18	48	29	17	35	24	30	74	63.5	1750	43

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*The thermostatic baths UNITRONIC-ORBITAL with swinging action of the mobile tank, are equipped with digital electronic regulator of temperature, time and speed. It includes a safety device which disconnects the heater when the circuit fails or low fluid level.*

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
*It has an RS-232 output for a full control with a computer (see communications protocol at the end of this instructions).*

*The baths are equipped with Pt100 thermoresistance as a temperature sensor, a drain tap AISI 304 s.s. in the mobile tank, with tensor springs permitting the adaptation of all types of flasks, racks and bottles. The bath UNITRONIC-ORBITAL C (Code 6001173) has a built-in compressor group fitted on to the antivibrating coupling with ventilated evaporator.*

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*To use the baths, fill the tank with liquid up to the required level, which must at least cover the heating element. Suit the liquid or mixture to the type of bath and range of working temperature. We herewith show details of some of these liquids as a guide:*

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Alarma: OVERTEMP 0 = OFF, +10  
 Alarma: UNDERTEMP 0 = OFF, -10

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**Temperature:** 0.0 °C +100.0 °C  
**Refrigeration:** Grupo compresor

**BATH 320 OR code 6032011**

**Temperature:** 10.0 °C +100.0 °C  
 (Regulation under ambient temperature requires a refrigerator)

**Resolution:** ± 0,1  
**Power:** 1.600 W  
**Voltage supply:** 230/50Hz – 60Hz  
 115/50Hz – 60Hz  
**Sensor:** Pt100  
**Time:** Timer from 00h01m to 99h59m.  
 Chronometer up to 99h59m.

**Regulation:** PID

**Output:** Serial port

**Alarms:** Overtemp off - +10  
 Undertemp off - -10  
 Tto-safety Safety thermostat  
 Rtd fail Shortcuted probe  
 Power fail Mains error  
 Motor Motor error

**Configuration:** Maximum temperature: 100,0 °C  
 Minimum temperature:  
     Code 6032011 = 10, 0  
     Code 6001173 = 0, 0  
 Maximum shaking: 100 opm.  
 Minimum shaking: 10 opm.  
 Temperature: °F °C  
 Calibrat: Correction of °C °F ±10.0°  
 Terminal ID: Number of periferic  
 serial port

**Alarm programing:**

Alarm: OVERTEMP 0 = OFF, +10  
 Alarm: UNDERTEMP 0 = OFF, -10

**INSTALLATION**

Place the bath on a flat, horizontal and level surface. Bearing in mind the weight of liquid when choosing the support. Provide for bath drainage.

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*Do not use the autoclave if it is not earthed.*

*Make sure that the voltage indicated on the description plate coincides with the mains power.*

*If you have to change the plug bear in mind the following:*

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*It is recommended the use of decalcified water. DO NOT USE DEIONIZED OR DEMINERALIZED WATER as it has a great tendency to rust. To prolong the life of the bath and avoid rusting of the parts, dissolve 1 gr of  $\text{CO}_3\text{Na}_2$  (Sodium Carbonate) per litre of water used in the bath.*

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**WARNING!!! WHEN WORKING AT TEMPERATURES ABOVE 60°C TAKE THE NECESSARY PRECAUTIONS TO AVOID POSSIBLE BURNS**

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**WARNING!!! TAKE PRECAUTIONS TO AVOID POSSIBLE LIQUID SPLASHES WHEN WORKING WITH STIRRERS**

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
*When pressing **IG**, the bath will start, the display shows:*

*Real temperature    °C / °F.  
 Real speed            o.p.m.  
 Time*

*Pressing **P3** the display shows the programmed **values**.  
 Pressing **P1** (3 sec.) we can start programming.  
 Pressing **P5** start or stop the control of the bath.  
 If a mains failure occurs, when it came back the bath starts with its old settings.*

**PROGRAMMING:**

*If we are in mode **stop** and press **P1** (3 sec.) we go to **temperature** programming.  
 Pressing **P3** or **P4** increase or decrease the temperature.*

*Pressing **P2** we go to **shaking** programming.  
 Pressing **P3** or **P4** increase or decrease the shaking.*

*Pressing **P2** we go to working **time** programming.  
 Pressing **P3** or **P4** increase or decrease the time.*

*Pressing **P2** we go start again the cycle temperature, shaking, time.*

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*Pressing **P1** (3 sec.) we go to **temperature alarm** programming.  
If no key is pressed for 5 seconds, the bath will go automatically to **stop** mode.*

#### OVERTEMP AND UNDERTEMP ALARM PROGRAMATION:

*If we are in mode **programming** and press **P1** (3 sec.) we go to **overtemp** temperature alarm programming.*

*Pressing **P3** or **P4** increase or decrease the overtemp alarm.*

*Pressing **P2** we go to **undertemp** temperature alarm programming.  
Pressing **P3** or **P4** increase or decrease the overtemp alarm.*

*Pressing **P2** we go back **overtemp** temperature alarm programming.  
If no key is pressed for 5 seconds, the bath will go automatically to **stop** mode.*

#### CONFIGURATION:

*If we are in mode alarm programming and press **P1** (3 sec.) we go to **configuration of maximum temperature**.  
Pressing **P3** or **P4** increase or decrease the maximum temperature.  
Pressing **P2** we go to configuration of **minimum temperature**.  
Pressing **P3** or **P4** increase or decrease the minimum temperature.  
Pressing **P2** we go to configuration of **maximum shaking**.  
Pressing **P3** or **P4** increase or decrease the maximum shaking.  
Pressing **P2** we go to configuration of **minimum shaking**.  
Pressing **P3** or **P4** increase or decrease the maximum shaking.  
Pressing **P2** we go to configuration of temperature (**°C** or **°F**).  
Pressing **P3** or **P4** selects **°C** or **°F**.  
Pressing **P2** we go to **calibration** of temperature.  
  
Pressing **P3** or **P4** increase or decrease the calibration.  
  
Pressing **P2** we go to **terminal ID**.  
Pressing **P3** or **P4** increase or decrease terminal ID.  
  
Pressing **P1** (3 sec.) we go back to **stop** mode.  
  
If no key is pressed for 5 seconds, the bath will go automatically to **stop** mode.*

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**RUN MODE:**

*Pressing **P3** we see the programmed values.*

**PAUSE:**

*Pressing **P4** the bath is in pause mode, which stops the shaking, the time and keep the temperature.*

*If we press **P4** again, go back to **start** mode.*

*Pressing **P5** stops the control of the bath.*

*Pressing **P1** (3 sec.) we go to temperature **programming**.*

**PROGRAMMING IN RUN MODE:**

*If we are in mode **start** and press **P1** (3 sec.) we go to **temperature programming**.*

*Pressing **P3** or **P4** increase or decrease the temperature.*

*Pressing **P2** we go to **shaking programming**.*

*Pressing **P3** or **P4** increase or decrease the shaking.*

*Pressing **P2** we go to working **time programming**.*

*Pressing **P3** or **P4** increase or decrease the time.*

*Pressing **P1** (3 sec.) we go back to **run** mode.*

*If no key is pressed for 5 seconds, the bath will go automatically to **run** mode.*

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**STOP AND START SCREENS:****HEATER:**

*Temperature programming. From 100.0 to 10.0 and OFF = stop the heating.*

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**SHAKING:**

*Shaking programming. From 100 to 10 rpm, and  
OFF = stop the shaking.*

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\* *If temperature = OFF and shaking = OFF, when we press **P5 (RUN-STOP)** the bath will not start and will give us a message **SETTINGS OFF**. This message will go out going back to programming.*

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**CALIBRATE:***Calibration adjustment of the bath.***TERMINAL ID:***Number of the serial port terminal.***ALARMS:****OVERTEMP:***Temperature too high.***Function:**

- Deactivate run mode.*
- Start the alarm.*
- Stop the heating.*
- Stop the motor.*

**Solution:**

- To stop the alarm press **P5** (START/ STOP) and go to STOP mode.*

**UNDERTEMP:***Temperature too low.***Function:**

- Deactivate run mode.*
- Start the alarm.*
- Stop the heating.*
- Stop the motor.*

**Solution:**

- To stop the alarm press **P5** (START/ STOP) and go to STOP mode.*

**LEVEL:***Safety or level thermostat.***Function:**

- Deactivate run mode.*
- Start the alarm.*
- Stop the heating.*
- Stop the motor.*

**Solution:**

- To stop the alarm press **P5** (START/ STOP) and go to STOP mode.*

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**RTD FAIL:***Open or crossed probe.***Function:**

- Deactivate run mode.*
- Start the alarm.*
- Stop the heating.*
- Stop the motor.*

**Solution:**

*To stop the alarm press **P5** (START/STOP) and go to STOP mode.*

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**POWER FAIL:***Mains failure.***Function:**

*The bath have a mains failure, and if there are not any other alarm, it will keep controlling the last programmed parameters.*  
*If another alarm appears, its own parameters will appear.*  
*The buttons **P1** ENTER, **P3** SETTINGS and **P4** PAUSE does not work.*

**Solution:**

*To stop the alarm press **P5** (START/STOP), the message POWER FAIL disappears and continue with the **RUN** mode.*  
*If another alarm appears, pressing **P5** (START/STOP) and go to STOP mode.*

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**EEPROM FAIL:***Eeprom error.***Function:**

*An EEPROM reading error with the lost of the last programmed parameters.*  
*The bath set its default working parameters automatically.*  
*The bath go into **STOP** mode.*

**Solution:**

*To stop the alarm press **P5** (START/STOP) and go to STOP mode.*

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**Function:**

*There is an problem with the motor or tachometer.  
Start the alarm.  
The bath go into **STOP** mode.*

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
**Solution:**

*To stop the alarm press **P5** (START/STOP) and go to **STOP** mode.*

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1. *Safety thermostat adjustment knob.*
2. *Safety thermostat manual RESET knob.*

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*The safety thermostat is at the rear of the bath casing.*

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1. *With a screwdriver turn the safety thermostat clockwise as far as possible by means of knob.*
2. *Select the desired temperature.*
3. *Start the machine up and leave the bath to settle.*
4. *Once the bath is stabilized, slowly turn the safety thermostat knob anticlockwise, until it activates. The reset knob will come out simultaneously.*
5. *Stop the bath and leave it to cool.*
6. *When the temperature is at 50% of the working temperature, reset the safety thermostat by pressing knob.*
7. *The bath is ready to work with.*

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**NOTE:**

*In the 6001173 bath, the main heating switch (1) has two positions, one is for the normal function of the bath, and the other is for the bath to work with compressor.*

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### RS-232 COMMUNICATION OVERVIEW

- 9600 bauds, 7 dat bits, 1 bit stop, without parity no handshake.
- Up to 15m of cable.
- Cable connection DTE to DTE 9 pins male.
- Protocol ASCII with LRC against transmission errors
- MASTER-SLAVE Communication (PC computer as a master).
- The device never sends nothing if computer does not ask for it.

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- *Messages starts with character #.*
- *Device address must be included. This field goes from 00 to 99, size allways 2 bytes in decimal-ASCII format.*
- *Data field is the message*
- *LRC is a calculated value to check good transmission.*
- *At the end of the frame allways a CR (0D) followed by LF (0A).*
- *Slave device answers allways sending its address*

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*Frame includes a transmission error cheking field called LRC. This field size 2 bytes. This field is calculated by the sender device and includes it on the frame. The reciver device also calculates and if are equal the transmission was good.*

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**LRC calculation procedure:**

1. Convert each character on it hexadecimal value.
2. Add
3. Take the smallest byte (two first characters starting from right hand side)
4. Subtract from FF on previous result
5. Add 1

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
It could be found is any specilized purchaser:  
Cable Null - modem de 9 pins

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Conectors in both cable sides are 9 pins female

## INSTRUCCIONES RS- 232 EQUIPO FRIGITERM - ULTRATERM / RS-232 INSTRUCTIONS FOR FRIGITERM - ULTRATERM

**PVT?** Solicitar temperatura de la sonda PV (Point Value) *Get probe temperature PV (Point value):*  
 Comando / *Command:* #01PVT?KK<CR><LF>  
 Respuesta / *Answer:* #01+000,0KK<CR><LF>

**SVT?** Solicitar temperatura de consigna SP (Set point) *Get Set point temperature:*  
 Comando / *Command:* #01SVT?KK<CR><LF>  
 Respuesta / *Answer:* #01+000,0KK<CR><LF>

**SVT** Establece temperatura de consigna SP (Set point) *Ser point temperature*  
 Comando / *Command:* #01SVT<BL>+000,0KK<CR><LF>  
 Respuesta / *Answer:* #01OKKK<CR><LF>

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## MAINTENANCE

Before removing the covers, disconnect the apparatus from the mains.

The manipulation of the internal electronic circuits of the plate by unauthorized personnel can cause irreparable damage. Make sure to take it to one of J.P.Selecta,S.A. authorized technical services.

### CLEANING:

If normal tapwater is used calcareous deposits will appear in the bath and in the heating element. A heater covered with lime should not be used as it can lose heating power and deteriorate more quickly. Periodical cleaning of the inside of the tank and the heating element with a solution of 10% of 28% Chlorhydric acid or a solution of water and vinegar (Acetic acid) is advised. Rinse with clean water afterwards.

## SPARE PARTS

To guarantee the safety of the equipment, the spare parts must be purchased from J.P.SELECTA,S.A.

Código / Code	6032011	6001173
SopORTE circuito / Circuit support	5260	5260
Alargo Schuko / Schuko cable	7001	7001
Circuito / Circuit	11023	11023
Cojinetes / Bearings	15447	15447
Motor / Motor	24053	24053
Retenedor cable / Cable retainer	28005	28005
Disco cuentavuelTas / Tachometer disk	30053	30053
Grifo baño / Bath tap	34032	34032
Racord termostato / Thermostat connector	34141	34141
Resistencia / Heating element	39012	39012
Termostato de seguridad / Safety thermostat	43031	43031
Sonda Pt100 / Pt100 probe	43065	43065
Transformador / Transformer	45020	45020

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## GUARANTEE

This product is guaranteed for one year. The guarantee does not cover damage caused by incorrect use or causes beyond the control of J.P. SELECTA,S.A.

Any manipulation of the apparatus by personnel not authorized by J.P. SELECTA,S.A. automatically cancels the guarantee.

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