PRECISBIG 6001240
PRECISBIG 6001241

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1. GENERAL INFORMATION

- 1) Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributor immediately.
- 2) Do not install or use the equipment without reading this handbook before.
- 3) This handbook must always be attached to the equipment and it must be available for all users.
- 4) If you have any doubts or enquiries, please contact with your supplier or J.P. Selecta's technical service.
- 5) **IMPORTANT! J.P. SELECTA WILL NOT ACCEPT ANY EQUIPMENT TO BE REPAIRED IF IT IS NOT DULY CLEANED.**
- 6) If any modification, elimination or lacking in maintenance of any device of the equipment by the user transgress the directive 2009/104/CE, the manufacturer is not responsible for the damage that can occur.
- 7) Do not use the equipment with liquids which can give off vapours capable of making explosive mixtures.

2. PACKING LIST

The standard equipment consist of the following components:

Piece	Code	Qty.
PRECISBIG	6001120	1
PRECISBIG	6001121	1
Instructions manual	80248	1

3. ACCESSORIES

Description	Code
GABLE LID For baths 6001240.	6000715
GABLE LID For baths 6001241.	6000716

4. TECHNICAL FEATURES

Voltage supply 115-230V 50/60 Hz according to the characteristics plate indications.

Code		6001240	6001241
Maximum temperature °C		100	
Capacity litres L.		160	220
Tank dimensions cm	Depth	40	
	Length	80	110
	Width	50	
Overall dimensions cm	Depth	96	
	Length	100	140
	Width	70	
Consumption W.		9000	9000
Weight Kg.		130	160

Temperature up to 100°C


Stability ±0,25°C

Homogeneity ±1%

Set error ±1%

Resolution 0.1°C

- * Remember that these are static baths.
- * Depending on the kind of fluid, the homogeneity can be altered.

Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.


5. EQUIPMENT DESCRIPTION

PRECISBIG baths with adjustable temperature from ambient +5°C up to 100°C. To use the baths fill the tank with liquid until the level covers at least the heater element. The liquid or mixture to be employed must be adequate to the model of bath and the working temperature.

Baths incorporate a drain tap for an easy pouring.

As a guidance, here are some liquids (to use these liquids, follow the manufacturer's instructions):

6. INSTALLATION

Place the bath on a flat, horizontal and level surface. Bear in mind the weight of the liquid when choosing the support. Provide for bath drainage.

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Be sure that the voltage supply is the same as the one indicated on the characteristic plate of the equipment.

Do not use the equipment if it is NOT earthed.

Power supply at 400V III + N 10A x phase.

It is recommended the use of decalcified water. DO NOT USE DEIONIZED OR DEMINERALIZED WATER as it has a great tendency to rust. To prolong the life of the bath and avoid rusting of the parts, dissolve 1 gr of CO₃Na₂ (Sodium Carbonate) per litre of water used in the bath.


7. OPERATION

7.1 CONTROL PANEL

General switch with signal lamp.

1. Heating indicators; open thermostat; failure
2. Reading display (temperature or time).
3. Setpoint value display (temperature or time).
4. Visualization mode change.
5. Increase value.
6. Decrease value.
7. Start / stop button.
8. Safety thermostat regulation screw (grey).
9. Safety thermostat reset button (white).
10. Control panel
11. Magnetothermic
12. Safety thermostat

7.2 WORKING MODE OF BATH

STARTING UP

1. Activate the main magnetothermic.
2. Press the start switch (1) and the screen (2) will light up at the same time.
3. Press keys (5 or 6) until you select the desired temperature value in the green display (the upper one).
4. Press key (4) if you want to select an operating time. Select time with keys (5 or 6).
5. Press key (7) to start the cycle.
6. The heating symbol will light and blink.
7. Press the key (7) again to stop the process. Note: during operation, the temperature or time adjustment value can be changed.

COMMENTS:

- * Do not program the equipment at lower temperatures than ambient +5°C, because the thermostat might have difficulties in regulating correctly.
- * When using the time function, the green display (the upper one) will indicate the remaining time. If you want to see the target temperature value, press key (4).
- * Place the load in the equipment trying to leave spaces for the air circulation. Do not load the incubator more than 70% of its useful chamber volume.
- * The device detects if the power fails. If, when the equipment's use has ended, the process is not stopped by pressing key (7), **P.FAIL** message can be avoided when the device is restarted.

ADJUSTING THE SAFETY THERMOSTAT:

From the factory, the thermostat is adjusted to its maximum temperature. To adjust the safety thermostat temperature, do as follows:

1. With the equipment connected to the mains, turn the safety thermostat screw (8) clockwise to the end of the scale, with the help of a flat screwdriver.
2. Select the desired safety temperature using the control keys. This temperature has to be minimum at 3°C above the working temperature.
3. Turn on the equipment with key (7).
4. When the equipment is stabilized at the above selected temperature, turn the safety thermostat screw (8) slowly to the left until it clicks, and the red indicator lamp lights up.

5. Stop the equipment and leave it cool with the lid open for 15 minutes.
6. Select the working temperature by using the temperature selection keys (4) & (5) (minimum at 3°C below the desired safety temperature).
7. Turn the equipment on again with key (7). Once the temperature is stabilized, the machine is ready to work.

NOTE: Safety thermostat is especially useful when the user works habitually at the same temperature.

8. MAINTENANCE

Before removing the cover disconnect the equipment from the mains.

The control elements and heating elements must only be manipulated by authorized personnel.

The distilled or demineralized water has a great tendency to rust. To prolong the life of the bath and avoid rusting of the parts, dissolve 1 gr of CO₃Na₂ (Sodium Carbonate) per litre of water used in the bath.

CLEANING:

If normal tapwater is used calcareous deposits will appear in the bath and in the heating element. A heater covered with lime should not be used as it can lose heating power and deteriorate more quickly. Periodical cleaning of the inside of the tank and the heating element with a solution of 10% of 28% Chlorhydric acid or a solution of water and vinegar (Acetic acid) is advised. Rinse with clean water afterwards.

For the cleaning of the different parts use the following products:

Cleaning of stainless steel: Alcohol

Cleaning of plastic: Alcohol with cotton duster.

9. ERROR MESSAGES

P.FAIL Power failure. The bath continues with the process.

Err will appear in green on the top line of the display.

SAFE /

TSEG The safety thermostat is open.
(press button 9 to reset it)

SENS Open temperature probe. The control does not detect the temperature probe. Probe connection fails.
(check probe connections or replace it)

HEAT The resistance does not work. The equipment detects that the temperature does not increase.
(check the resistance function. If necessary, replace it)

10. ADJUSTMENTS

The offset (A) and the °C/°F unit (B) can be adjusted.

A. Adjust the probe reading correction (offset)

Press keys (5) and (6) continuously and at the same time (about 5 seconds). The **OFFs** message will appear.

Use keys (5) / (6) to set the actual temperature reading value. Press key (7) to confirm.

B. Presentation unit: °C / °F

Press key (4). The message **Unit** will appear. Using keys (5) / (6), select the unit: °C / °F. Confirm with key (7).

11. SPARE PARTS

To guarantee the safety of the equipment, the spare parts must be purchased from J.P. SELECTA, s.a.u.

Control circuit	38500
Temperature probe	43044
Level circuit	29525
Heating element	39060
Safetu thermostat	43031
SSR	36028

12. GUARANTEE

This product is guaranteed for one year. The guarantee does not cover damages caused by an incorrect use or causes beyond the control of J.P. SELECTA, s.a.u.

Any manipulation of the equipment by unauthorized personnel by J.P. SELECTA, s.a.u. cancels the guarantee automatically.