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Circulación
Externa

1. General information

- Handle the package carefully. Unpack and check that the contents match what is indicated in the "Packing List" section. If a damaged component or the absence of any is observed, please advise the distributor immediately.
- Do not install nor use the equipment without previously read this manual. These instructions are an inseparable part of the equipment and should be available to all users.
- Check any doubt with J.P. SELECTA's technical service.
- Do not use the equipment in explosive, flammable or corrosive atmospheres.

2. Packing list

Recirculation unit:	6001491
SILICONE TUBE 14 X 8 (METRE)	0046041

3. Technical specification

Inner tank volume:	20L
External dimensions: height/width/depth (mm)	70 x 43 x 56
Power voltage:	230 VAC 50Hz
Electrical power (W):	2100 + 2100
External fuse:	20A Replaceable 10x5mm
Temperature range:	From... Ambient T ^a + 5°C ... up to... 90°C
Display resolution:	1°C
Temperature stability:	±0.5°C
Programmable operation time:	From 0h 1m to 99h
Programmable waiting time:	From 0h to 99h

4. Equipment description

Equipment for water recirculation, for feeding condensation circuits, distillation, reactor layers, viscosimeters, etc...

Water forced circulation gets a great temperature homogeneity and fast temperature rise.

Temperature controller includes the following functions:

- Temperature selection up to 90°C in increments of 1°C.
- Programming temperature profiles for starting increasing ramps in case of reactors with very viscous material.
- Autotuning function to optimize the temperature control.
- Two manual power levels to facilitate operation.

5. Installation

Once unpacked the equipment, visually check if it's in good condition. Keep the packing for a few days.

5.1 Location

Place the equipment on a flat, levelled and stable surface suitable for the equipment weight.

Preferably on the ground near the equipment to be thermostated.

For comfort, leave a free space of 50cm around the equipment.

5.2 Connection to the mains

Choose a plug near the equipment and adequate to the equipment power. Use the supplied connection cable or a similar one.

Make an installation with 2.5mm section cable if you are going to work at full power in a continuous way. Bear in mind that are 4000W, about 18 A.

For your safety, the power plug must be earthed.

Before plugging in, check that the voltage and frequency of the mains corresponds to the one indicated in the equipment.

Important:
The fluid must enter by the reactor's lower entrance.


Fig 6.2. Installation example to thermostat a reactor.



5.3 Fluid line installation

Install a similar circuit than the one in figure 6.2.

Use the proper hose for diameters and temperature. A suitable hose is:

1001540: Roll of 25 meters. Temperature up to 200°C.

6. Control panel

7. Operation

7.1 Starting

- MAIN SWITCH OFF.
- Check that the tank drain valve is closed.
- Fill the tank with the suitable fluid:

For instance:

- Water up to 90°C
- Silicone oil (50cPs) up to 150°C

- View the fluid level in the tank. (See 6)
- Connect the hoses to DELIVERY & RETURN. (See 6)
- Turn the CIRCULATION VALVE knob to INNER POSITION.
- Place the power selector to 50%.
- Check and review it is correct.

- MAIN SWITCH ON

The equipment starts the circulation pump.

Fluid must circulate only internally.

- Select the working temperature in the Temperature controller.
- Slowly move the CIRCULATION VALVE to EXTERNAL POSITION.
- View that the fluid begins to flow as expected and not leaking.
- In case of leakage, move the valve to INTERNAL position.

Internal
Circulation

External
Circulation