INTERHEATER – 180°C

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INDEX

1. General Info	18
2. Packaging Contents	18
3. Technical Specification	19
3.1 Electrical connection	19
3.2 Accesories	19
4. Description	20
4.1 Placement	20
4.2 Components	21
4.3 Schematics	22
5. Operation	23
5.1 Initial preparation	23
5.2 Touchscreen	23
5.3 Graphics	24
5.4 Reports export	24
5.5 Configuration	24
6. Maintenance tasks	25
6.1 General considerations	25
7. Safety thermostat:	25
7.1 Rearming the safety thermostat	25
7.2 Adjusting the safety thermostat	26
8. Recycling & equipment distmantling:	26
9. Spare parts:	27
10 . Error List & alarms:	29
11. Warranty:	30

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ENGLISH

1. General Info

Handle the package carefully. Unpack it and check that the contents match those indicated in the "Packaging Contents" section. If you notice any damaged or missing components, notify the distributor immediately.

Do not install or use the equipment without first reading this instruction manual.

These instructions are an integral part of the equipment and must be available to all users.

Any questions can be clarified by contacting the technical service of J.P. SELECTA, s.a.u.

ATTENTION! NO MACHINE WILL BE ACCEPTED FOR REPAIR THAT IS NOT PROPERLY CLEANED AND DISINFECTED.

Any modification, removal, or lack of maintenance of any device on the machine violates the User Directive 2009/104/EC, and the manufacturer is not liable for any resulting damage.

DO NOT USE THE EQUIPMENT WITH FLUIDS THAT MAY RELEASE VAPORS THAT ARE CAPABLE OF FORMING EXPLOSIVE MIXTURES.

2. Packaging Contents

Item	
Fluid Recirculator " Interheater-180°C"	(x1) cod. 6001497
Reinforced Glass Fiber Silicone Hose \varnothing 12x19 (2,5 mts)	(x2) cod. 46153
Stainless steel clamp \varnothing 12/22	(x4) cod.44288
Instructions Manual	M.80463
EU Declaration of Conformity	

3. Technical Specification

Item	
Required Voltage (V):	(400V / 50Hz) other values on demand
Installed Power (W):	9700
Max. Intensity (A):	20 A
Thermal Circuit Breaker (A):	4 poles 32 A
Power Cord:	3ph+N+T/T (6mm ²) x 2 mts. + CETAC plug 32A (other connections on demand)
Net Weight (Kg):	25
Thermal fluid tank capacity:	25L
External Dimensions: (cm depth. x width x height):	70x45x67
Temperature range (°C)	20..180°C (setpoint value must be 5°C higher than room temperature)
Tank pressure	atmospherical
Hose connections	Hose nipple 1/2" (ø12mm)
Thermal stability	+/-0.5°C at 100°C
Setpoint error	+/-1°C
Room conditions	Max. 40°C HR 70%
USB connection	yes

3.1 Electrical connection

Check the electrical installation in the room where the fluid recirculator will be installed.

For your safety, we recommend protecting the equipment with a residual current device.

Model	Voltage	Power (W)	Thermal Circuit Breaker
6001497	400V (3f+N+TT)	9700	32A

3.2 Accesories

To optimize the thermal performance of a reactor system, it is recommended to operate it in conjunction with a fluid cooler (chiller) and connect it to the recirculator inlets as described later in this manual.

JP Selecta recommends the following equipment:

Fluid Recirculator cod. JP Selecta: 6001492 / 6001493 / 6001494 / 6001496

4. Description

Thermal fluid recirculation equipment for feeding condensation and distillation circuits, reactor jackets, viscometers, etc.

Thermal fluid circulation achieves high temperature uniformity and rapid temperature rise.

The temperature controller includes the following functions:

- Temperature selection from 20°C to 180°C.
- Real-time monitoring of fluid temperature via a TFT touch screen.
- Possibility of downloading records by means of USB port.

4.1 Placement

Place the equipment on a flat, horizontal and levelled surface, ensuring that there is at least 50 cm of free space at the back and sides of the equipment for maintenance or adjustments.

Depending on the equipment's operating temperatures, water can be used as a thermal fluid below 90°C. For higher temperatures, a suitable fluid capable of reaching the desired temperature without degrading should be used.

When using silicone oils, it is normal for vapors to be generated due to heating above 150–160°C. It is recommended to install the equipment in a well-ventilated area or, alternatively, to have a duct to extract these vapors, connected to the device aeration connection.

Review the fluid's safety data sheet for more information.

Connect the recirculator preferably to the floor, near the equipment to be thermostated.

Please note that the recirculation fluid can reach temperatures of 180°C, so it is recommended to cover the hoses with an insulating cover to avoid burns from direct contact if the hoses are accessible to the operator.

The equipment is equipped with a 32 A, 5-pole CETAC type plug (3ph+N+T/T). The installation of a socket base appropriate for this connection is required, or if preferred, connect the cables directly to a terminal block located in a box conveniently prepared for this purpose.




Although the recirculation equipment is fitted with a thermal circuit breaker as an overload protection element, the installation where the equipment is to be located must also have a residual current device to protect against electric shocks.

Before plugging in, check that the voltage and frequency of the electrical network correspond to those indicated on the equipment's nameplate


4.2 Components.

Description	
01	Thermal Fluid filling cap
02	Deaeration registry tank
03	Main Switch
04	Aux Cooling Unit (Chiller) electrical socket enabling switch
05	Aux. Cooling Unit (Chiller) electrical socket. Max 16 A.
06	Recirculator Touchscreen
07	Thermal Circuit Breaker 32A
08	Usb connector to download registers
09	Safety Thermostat
10	Thermal Fluid tank drain valve
11	Outlet to aux. cooling unit (Chiller) Note: To connect to the cooling unit inlet (Chiller)
12	Inlet from aux. cooling unit (Chiller) Note: To connect to the cooling unit outlet (Chiller)
13	Outlet to reactor // Note: To connect to the reactor inlet
14	Bypass valve to the reactor inlet
15	Inlet from reactor // Note: To connect to the reactor outlet

4.3 Schematics.


5. Operation

5.1. Initial preparation.

- Connect the recirculator to the mains.
- Set the main switch (3) to the OFF position.
- Check that the tank drain valve (10) is closed.
- Fill the tank with the appropriate fluid through the filler cap (1).
 - Water up to 90°C
 - Silicone oil up to 190°C
- Check the fluid level in the tank.
- Connect the supplied DELIVERY and RETURN hoses and secure them with their cable ties, as indicated in chapter 4.
- Check that the vent (2) has a duct to the outside of the room or verify that the room is sufficiently ventilated to prevent vapor buildup.
- Set the BYPASS valve knob (14) to the RECIRCULATION POSITION (horizontal position of the valve knob).

5.2. Touchscreen.


Once the main switch is turned on and the screen is started (wait about 20-30 seconds), the appearance of the start screen is as shown in the attached image.

- Select the set temperature by clicking on the circle and moving it over the circular sector between 20°C and 180°C, until the desired value is displayed.
- Check that the bypass valve (14) is in the reactor RECIRCULATION position (bypass valve control in horizontal position).
- Press the key to start the equipment.
- The equipment will start the circulation pump.
- The fluid must circulate only internally (internal recirculation).
- Slowly move the bypass valve (14) to the reactor feed position (knob in vertical position).
- Visualize that the fluid begins to circulate as expected and without leaks.

- If you notice small leaks in any of the connections, move the valve to the RECIRCULATION position and check for leaks (tighten the clamps).

5.3. Graphics.

It is possible to monitor the temperature evolution in real time, using the key enabled on the screen.

5.4. Reports export.

It is possible to download the generated logs using the enabled key. To do this, insert a USB memory stick (pendrive), formatted FAT32, into the connector (8).

5.5. Configuration.

Using the config. key, you can modify:

- a) Date and time.
- b) Language.
- c) Calibration and offset application of the equipment's temperature control probes (inlet and outlet to the reactor).

6. Maintenance tasks

6.1 General considerations

Check the condition of the hoses and their connections monthly, replacing them if damaged.

Check the hose insulation covers monthly. Keep in mind that the equipment operates at high temperatures (around 180°C) and under pressure, so it is very important that these insulation covers are in perfect condition.


Keep the equipment clean of thermal fluid residue or general dirt by using a non-aggressive household cleaning product. Use a microfiber cloth to clean painted surfaces.

The equipment's touchscreen is a delicate element; weekly cleaning is recommended using a household cleaning product suitable for electrical and electronic equipment and a microfiber cloth.

Check the condition of the wiring and electrical connectors monthly, and repair or replace any damaged insulation.

7. Safety thermostat:

This safety device is located at the rear right hand of the vacuum oven, and limits the temperature inside the autoclave chamber.

8. Safety thermostat regulation screw (grey).
9. Safety thermostat reset button (white).

7.1 Rearming the safety thermostat

Safety thermostat is an element to prevent the temperature from exceeding a value. If by a fortuitous event the thermostat is triggered, that is, it opens the power circuit, the heating element will stop working. Usually there is an indicator light, or an error message, that indicates that the thermostat is in this position.

If the thermostat is triggered, check that the equipment looks good and it is not overheating at any point.

To change its status, simply press the white button (9), and you will hear a “click” that indicates that the thermostat has been reset. If the equipment has not cooled down sufficiently, the thermostat cannot be reset, and you must wait for the temperature to drop before it can be reset.

7.2 Adjusting the safety thermostat

In some equipment such as ovens and baths and with certain processes it may be useful to adjust the thermostat trigger value, but this operation should NEVER be performed to the rest of equipments.

From factory, the thermostat is adjusted to its trigger temperature, set over 200°C.

Stop the equipment and leave it to cool with the lid open for 15 minutes.

8. Recycling & equipment dismantling:

Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.

9. Spare parts:

		
<p>13005 Contactor CWB25</p>	<p>13019 Th. Circuit Breaker 4P 32A</p>	<p>16144 Valv.Solenoid 220V NC</p>
		
<p>16145 Valv. Solenoid 220V NO</p>	<p>16220 Pressure Switch</p>	<p>20134 Level switch</p>
		
<p>24353 Pump</p>	<p>43017 Thermostat</p>	<p>34190 Check valve 3/4"</p>
		
<p>16143 Power Supply</p>	<p>34412 Check valve 1/2"</p>	<p>34356 Ball valve 1/2"</p>

		
<p>85284 PLC Touchscreen</p>	<p>20017 Main Switch</p>	<p>36062 Solid State relay 3ph</p>
		
<p>85280 PLC</p>	<p>85281 Converter A/D</p>	<p>85287 Temperature transmitter</p>
		
<p>85266 Fan 24vdc</p>	<p>43170 Temperature Probe</p>	<p>39376 exchanger 90kw</p>
		
<p>39375 Heater 4,5 kW 3ph</p>	<p>15500 Schuko socket 50x50</p>	

10 . Error List & alarms:

Code	Name	Description
00	POWER FAILURE <u>Reason:</u> The equipment has suffered a power outage. <u>Solution:</u> revise cable & electrical installation	(POWER FAIL)
01	E2PROM FAIL <u>Reason:</u> Controller failure <u>Solution:</u> Contact aftersales dpt.	(E2PROM_FAIL)
02	SAFETY THERMOSTAT <u>Reason:</u> The safety thermostat is switched on. <u>Solution:</u> Check the set temperature and compare it with a standard probe to detect a possible fault in the heating system.	(SAFETY THERM)
03	FLUID LEVEL AT MINIMUM <u>Reason:</u> Fluid Level is too low. <u>Solution:</u> Check the fluid level or sensor signal	(LOW LEVEL ALARM)
04 & 06	TEMPERATURE PROBE FAIL <u>Reason:</u> Temperature probe is generating an abnormal reading <u>Solution:</u> Revise the temp. probe status (shorted or cut)	(TEMP FAIL)
05 & 07	OVERTEMPERATURE <u>Reason:</u> The temperature probe is detecting an abnormal reading, too high <u>Solution:</u> Revise the temp. probe status (shorted or cut)	(OVERTEMP)
08	LOW PRESSURE OF FLUID <u>Reason:</u> The pressure sensor is generating an abnormal reading. <u>Solution:</u> Check hoses, leakages, pump fail	(PRES FAIL)

11. Warranty:

This recirculation equipment has a one-year warranty.

The warranty does not cover damage caused by improper use or causes beyond the control of JP Selecta.

Any tampering with components or the equipment circuit by personnel not authorized by JP Selecta automatically voids the warranty.