PRECISDIG 6001120
PRECISDIG 6001121
PRECISDIG 6001122
PRECISDIG 6001123

PRECISDIG 6001220
PRECISDIG 6001221
PRECISDIG 6001222
PRECISDIG 6001223

CON TAPA MAKROLON / WITH MAKROLON FOLDING LID:

PRECISDIG 6001131
PRECISDIG 6001132


PRECISDIG 6001133

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
***WATER AND OIL BATHS WITH OR WITHOUT FOLDING LID
WITH DIGITAL REGULATION AND TEMPERATURE AND TIME READING***

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**Notice to customers:**

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.

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INDEX

1. GENERAL INFORMATION	12
2. PACKING LIST	12
3. ACCESSORIES	12
4. TECHNICAL FEATURES	13
5. EQUIPMENT DESCRIPTION	14
6. INSTALLATION	15
7. OPERATION	16
8. MAINTENANCE	18
9. ERROR MESSAGES	18
10. ADJUSTMENTS	19
11. SPARE PARTS	19
12. GUARANTEE	19

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1. GENERAL INFORMATION

- 1) Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributor immediately.
- 2) Do not install or use the equipment without reading this handbook before.
- 3) This handbook must always be attached to the equipment and it must be available for all users.
- 4) If you have any doubts or enquiries, please contact with your supplier or J.P. Selecta's technical service.
- 5) **IMPORTANT! J.P. SELECTA WILL NOT ACCEPT ANY EQUIPMENT TO BE REPAIRED IF IT IS NOT DULY CLEANED.**
- 6) If any modification, elimination or lacking in maintenance of any device of the equipment by the user transgress the directive 2009/104/CE, the manufacturer is not responsible for the damage that can occur.
- 7) Do not use the equipment with liquids which can give off vapours capable of making explosive mixtures.

2. PACKING LIST

The standard equipment consist of the following components:

Piece	Code	Qty.
Precisdig	6001120	1
Precisdig	6001121	1
Precisdig	6001122	1
Precisdig	6001123	1
Precisdig	6001220	1
Precisdig	6001221	1
Precisdig	6001222	1
Precisdig	6001223	1
Precisdig	6001131	1
Precisdig	6001132	1
Precisdig	6001133	1
Agibat	6032015	1
Agibat	6032016	1
Instruction manual	80379	1

3. ACCESSORIES

Description	Code
REMOVABLE CONSTANT LEVEL For baths to be used without perforated cover.	6001400
FIXED CONSTANT LEVEL It must be fitted in factory.	6000123

TUBE RACKS ADAPTABLE TO THE SUPPORT FRAME

For bath code	Cover with holes and reducing discs 105, 80, 60 and 37mm.Ø in stainless steel	Gable cover	Perforated shelf	Support-frame for holding racks
6001120 - 6001220	6003232 (2 holes)	6000710	6000022	1001199
6001121 - 6001221	6003234 (4 holes)	6000711	6000044	1001200
6001122 - 6001222	6003236 (6 holes)	6000712	6000066	1001201
6001123 - 6001223	2x6003236 (6 holes = 12 holes)	6000714	6000088	2 x 1001201
6032015 - 6032016*	6003236*	6000712*	6000066	1001201

*Only for model 6032016 without Makrolon lid

Code	For tubes Ømm.	Nº of tubes up to	Overall dimensions cm			Nº of racks for baths of capacity litres			
			Height	Width	Length	5	12	20	45
1001202	13	36	8	8,2	23,5	1	3	5	2x5
1001203	16	24	8	8,2	23,5	1	3	5	2x5
1001204	20	24	8	8,2	23,5	1	3	5	2x5

4. TECHNICAL FEATURES

Voltage supply 115-230V 50/60 Hz according to the characteristics plate indications.

Code		6001120	6001121 6001131	6001122 6001132	6001123 6001133	6001220	6001221	6001222	6001223
Maximum temperature °C		100				200			
Capacity litres L.		5	12	20	45	5	12	20	45
Tank dimensions cm	Depth	15							
	Length	30	30	48	62	30	30	48	62
	Width	14	30	30	50	14	30	30	50
Overall dimensions cm	Depth	20 / 40 (modelo con tapa MAKROLON)							
	Length	37	40	58	69	37	40	58	69
	Width	27	42	42	62	27	42	42	62
Consumption W.		600	1100	1500	2400	1100	1400	2100	3800
Weight Kg.		6	8 / 9	10 / 11	18 / 19	6	8	10	18

Code		6032015	6032016
Maximum temperature °C		80	
Capacity litres L.		20	
Tank dimensions cm.	Depth	15	
	Length	48	
	Width	30	
Overall dimensions cm	Depth	20	
	Length	58	
	Width	42	
Consumption W.		2100	
Weight Kg.		11	12

Temperature up to 200°C

Stability $\pm 0,5^{\circ}\text{C}$

Homogeneity $\pm 2,5\%$

Set error $\pm 1,5\%$

Resolution 1°C

- * Remember that these are static baths.
- * Depending on the kind of fluid, the homogeneity can be altered.

5. EQUIPMENT DESCRIPTION


PRECISDIG baths with adjustable temperature from ambient $+5^{\circ}\text{C}$ up to 100°C or up to 200°C , depending on the model. To use the baths fill the tank with liquid until the level covers at least the heater element. The liquid or mixture to be employed must be adequate to the model of bath and the working temperature.

AGIBAT baths allow liquid temperature regulation from ambient $+5^{\circ}\text{C}$ to 80°C and also inner stirring with adjustable flow. The equipment also has a temperature probe and water level, and a temperature and time control by digital electronic circuit. These baths come complete with drain tap for an easy drainage.

The baths incorporate a drain tap for an easy pouring.

As a guidance, here are some liquids (to use these liquids, follow the manufacturer's instructions):

6. INSTALLATION

Place the bath on a flat, horizontal and level surface. Bear in mind the weight of the liquid when choosing the support. Provide for bath drainage.

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Be sure that the voltage supply is the same as the one indicated on the characteristic plate of the equipment.

Do not use the equipment if it is NOT earthed.

If you have to change the plug bear in mind the following:

Blue cable: Neutral.

Brown cable: Phase.

Yellow/green cable: Earth.

It is recommended the use of decalcified water. **DO NOT USE DEIONIZED OR DEMINERALIZED WATER** as it has a great tendency to rust. To prolong the life of the bath and avoid rusting of the parts, dissolve 1 gr of CO_3Na_2 (Sodium Carbonate) per litre of water used in the bath.

7. OPERATION


7.1 CONTROL PANEL

General switch with signal lamp.

1. Heating indicators; open thermostat; failure
2. Reading display (temperature or time).
3. Setpoint value display (temperature or time).
4. Visualization mode change.
5. Increase value.
6. Decrease value.
7. Start / stop button.


REAR PART

8. Safety thermostat regulation screw (grey).
9. Safety thermostat reset button (white).

Note: only for **AGIBAT** models:

10. Stirring flow adjusting nut.
11. Temperature probe / water level sensor.

Adjusting stirring flow

The stirring flow can be adjusted by turning the nut (10). The nut's setup does not allow being totally off as, if this happened, the circulation pump would get damaged.

7.2 WORKING MODE OF THE BATH

STARTING UP

1. Press the start switch (1) and the screen (2) will light up at the same time.
2. Press keys (5 or 6) until you select the desired temperature value in the green display (the upper one).
3. Press key (4) if you want to select an operating time. Select time with keys (5 or 6).
4. Press key (7) to start the cycle.
5. The heating symbol will light and blink.
6. Press the key (7) again to stop the process. Note: during operation, the temperature or time adjustment value can be changed.

COMMENTS:

- * Do not program the equipment at lower temperatures than ambient +5°C, because the thermostat might have difficulties in regulating correctly.
- * When using the time function, the green display (the upper one) will indicate the remaining time. If you want to see the target temperature value, press key (4).
- * Place the load in the equipment trying to leave spaces for the air circulation. Do not load the incubator more than 70% of its useful chamber volume.
- * The device detects if the power fails. If, when the equipment's use has ended, the process is not stopped by pressing key (7), **P.FAIL** message can be avoided when the device is restarted.

ADJUSTING THE SAFETY THERMOSTAT:

From the factory, the thermostat is adjusted to its maximum temperature. To adjust the safety thermostat temperature, do as follows:

1. With the equipment connected to the mains, turn the safety thermostat screw (8) clockwise to the end of the scale, with the help of a flat screwdriver.
2. Select the desired safety temperature using the control keys. This temperature has to be minimum at 3°C above the working temperature.
3. Turn on the equipment with key (7).
4. When the equipment is stabilized at the above selected temperature, turn the safety thermostat screw (8) slowly to the left until it clicks, and the red indicator lamp lights up.
5. Stop the equipment and leave it to cool with the lid open during 15 minutes.

6. Select the working temperature by using the temperature selection keys (4) & (5) (minimum at 3°C below the desired safety temperature).
7. Turn the equipment on again with key (7). Once the temperature is stabilized, the machine is ready to work.

NOTE: Safety thermostat is especially useful when the user works habitually at the same temperature.

8. MAINTENANCE

Before removing the cover disconnect the equipment from the mains.

The control elements and heating elements must only be manipulated by authorized personnel.

The distilled or demineralized water has a great tendency to rust. To prolong the life of the bath and avoid rusting of the parts, dissolve 1 gr of CO_3Na_2 (Sodium Carbonate) per litre of water used in the bath.

CLEANING:

If normal tapwater is used calcareous deposits will appear in the bath and in the heating element. A heater covered with lime should not be used as it can lose heating power and deteriorate more quickly. Periodical cleaning of the inside of the tank and the heating element with a solution of 10% of 28% Chlorhydric acid or a solution of water and vinegar (Acetic acid) is advised. Rinse with clean water afterwards.

For the cleaning of the different parts use the following products:

Cleaning of stainless steel: Alcohol

Cleaning of plastic: Alcohol with cotton duster.

9. ERROR MESSAGES

P.FAIL Power failure. The bath continues with the process.

Err will appear in green on the top line of the display.

SAFE /

TSEG The safety thermostat is open.
(press button 9 to reset it)

SENS Open temperature probe. The control does not detect the temperature probe. Probe connection fails.
(check probe connections or replace it)

HEAT The resistance does not work. The equipment detects that the temperature does not increase.
(check the resistance function. If necessary, replace it)

10. ADJUSTMENTS

The offset (A) and the °C/°F unit (B) can be adjusted.

A. Adjust the probe reading correction (offset)

Press keys (5) and (6) continuously and at the same time (about 5 seconds). The **OFFs** message will appear.

Use keys (5) / (6) to set the actual temperature reading value. Press key (7) to confirm.

B. Presentation unit: °C / °F

Press key (4). The message **Unit** will appear. Using keys (5) / (6), select the unit: °C / °F. Confirm with key (7).

11. SPARE PARTS

To guarantee the safety of the equipment, the spare parts must be purchased from J.P. SELECTA, s.a.u.

DESCRIPTION	PART NUMBER									
Model	6001120	6001121 6001131	6001122 6001132	6001123 6001133	6001220	6001221	6001222	6001223	6032015	6032016
Circuit	38500								29416	
Heating element	39061	39051	39012	39123	39050	39060	39059	39124	39059	39059
Thermostat	43031				43017				43031	
Pt100 probe	43044								43044	

12. GUARANTEE

This product is guaranteed for one year. The guarantee does not cover damages caused by an incorrect use or causes beyond the control of J.P. SELECTA, s.a.u.

Any manipulation of the equipment by unauthorized personnel by J.P. SELECTA, s.a.u. cancels the guarantee automatically.