TERMOSTATOS DE BLOQUE METÁLICO HEATING BLOCS

**TEMBLOC
MULTIPLACES**

**7462200
7471200**

Indice

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1 Safety

The equipment comes complete with safety features.

This manual indicates areas of possible risk.

1.1 Safety Icons

It identifies risky situations and safety measures that should be adhered to. The icons make reference within the paragraph, marked as a lined grey.

Risk danger

Risk danger

Beware, follow instructions as indicated.

Electrical risk

Risk of electrical shock if accessing zones are shown with this icon.

Beware, follow instructions as indicated.

Risk of burns in high temperature areas

The temperature in the indicated zone with this icon can exceed 60°C. Use thermal protective gloves to complete described functions.

Beware, follow instructions as indicated

Important information

- Important information in obtaining the best results or optimum performance of the equipment.
- Important information to extend the life of the equipment and maintain its optimum performance.

1.2 Risks that the operator can be subjected to

- Possibility of touching surfaces that are over 60°C.
- Electrical risk.

1.3 Qualified users

This equipment should only be used by personnel that are qualified users.

This equipment should only be used by personnel that have understood this manual or have previously read instructions about the use of steam sterilisers.

2 General information

- 2.1 Handle the equipment with care. Unpack and check that all items coincide with the packing list. If you see any discrepancy or damage of any kind, notify your distributor as soon as possible.
- 2.2 Read this manual before using the equipment.
- 2.3 These instructions are an important part of the equipment and should be made available to all users.
- 2.4 If in any doubt or clarify on the use of this instrument, contact your local distributor or the technical service of J.P. Selecta, s.a.u.

- 2.5 ATTENTION, EQUIPMENTS NOT CLEAN AND DISINFECTED SHALL NOT BE REPAIRED.**
- 2.6 Any unauthorised modification, lack of maintenance or deliberate damage to the machine, transgresses the directive 89/655/CEE where the manufacturer will not be responsible for damages that could occur due to undue care and attention.
- 2.7 Do not use the equipment for fluids that can produce vapours that can form inflammable or explosive mixtures.

3 Technical specifications

Voltage supply 115-230V 50/60Hz according to the characteristics plate indications.

| Code | Temperature °C | Stability °C | Resolution °C | Overall dimensions cm | | | Consumption W | Weight Kg |
|---------|-------------------|-----------------|------------------|-----------------------|-------|-------|------------------|--------------|
| | | | | Height | Width | Depth | | |
| 7462200 | 30-200 | ± 0,75 | 1 | 11 | 19 | 29 | 350 | 3,7 |
| 7471200 | | | | 11 | 31 | 31 | 700 | 7,4 |

4 Contents list

4.1 Standard contents

The standard equipment consists of the following items:

- | | |
|--------------------------|---------|
| - Tembloc | 7462200 |
| - Multiplaces | 7471200 |
| - Removing blocks handle | 5014 |
| - Reading thermometer | 9742027 |
| - Instruction manual | 80076 |

5 Introduction

The metallic thermostatic blocks **Tembloc** and **Multiplaces** have been designed to work in a temperature range from +30°C to 200°C with a stability of $\pm 0.75^\circ\text{C}$. They have accessories such as blocks to be used with 1.5 ml Eppendorf tubes and with 12 mm, 16 mm and 20 mm diameter tubes as well as a block without holes ready to be tooled with the drills needed by the user. **Tembloc** can fit one round block meanwhile **Multiplaces** has a capacity for three rectangular blocks.

6 Equipment description

6.1 Equipment

The thermostatic blocks **Tembloc** and **Multiplaces** have the following features:

- Pt100 temperature regulation probe
- Thermocouple type K temperature limitation probe
- Backing plated heating element of 350W/700W, 800 depending on the model

6.2 Features

- Temperature regulation from +30°C to 200°C
- Stability $\pm 0.75^\circ\text{C}$
- Resolution 1°C
- Working cycles from 1 minute to 99h:59minutes
- Electronic temperature limiter

6.3 Accessories

TEMBLOC:

Metal blocks, size: 105 \varnothing x 55 mm high.


- 7000346 for 34 x 6 mm \varnothing tubes
- 7000714 for 18 x 1.5 ml. Eppendorf tubes
- 7001224 for 24 x 12 mm \varnothing tubes.
- 7001618 for 18 x 16 mm \varnothing tubes.
- 7000208 for 8 x 20 mm \varnothing tubes.
- 7000256 for 6 x 25 mm \varnothing tubes.
- 7000715 Blank block with no holes for customized demand.

MULTIPLACES:

Metal Blocks: 55 high x 95 wide x 75 mm deep.


- 7000306 For 30 x 6 mm \varnothing tubes (total 3 blocks = 90 tubes).
- 7000716 For 20 x 1.5ml, Eppendorf tubes. (total 3 blocks = 60 tubes.)
- 7031220 For 20 x 12 mm \varnothing tubes (total 3 blocks = 60 tubes).
- 7031612 For 12 x 16 mm \varnothing tubes (total 3 blocks = 36 tubes).
- 7003208 For 8 x 20 mm \varnothing tubes (total 3 blocks = 24 tubes).
- 7001256 For 6 x 25 mm \varnothing tubes (total 3 blocks = 18 tubes).
- 7000717 Blank block without holes, can be customised to customer requirements.
- 7001474 Block for 15 cuvettes 10 x 10. Capacity: 3 blocks.
- 7001475 Block for 96 well microplates of 1.2 ml. Capacity: 1 block.

6.3 Control panel

1. Temperature and time display
2. Modification temperatura/time key
3. Increasing value key
4. Decreasing value key
5. Start/Stop key

6.4 Rear part

6. Main switch
7. Fuse carrier
8. Power supply cord

7 Installation

7.1 Placement

Locate the equipment near a socket keeping in mind its amperage and on a leveled surface enough resistant to support the weight indicated on the table of section 3 (Technical features). Leave a space of 10 cm by the rear side and also both laterals.

DO NOT INSTALL THE EQUIPMENT IN EXPLOSIVE ATMOSPHERES OR IN PLACES WHERE COULD BE SOLVENTS OR OTHER PRODUCTS WHICH MAY PRODUCE EXPLOSIVE OR FLAMMABLE GAS.

7.2 Power supply connection


Make sure that the voltage given by the socket is the indicated on the feature label attached on the rear part.

Do not use the equipment unless it is correctly earthed.

If you replace the plug keep in mind the following:

| | |
|--------------------|---------|
| Blue wire: | Neutral |
| Brown wire: | Phase |
| Yellow/green wire: | Ground |

8 Operation

8.1 Control panel

1. Adjustment (change between temperature / time) button
2. Increase Value key
3. Reduce Value key
4. Stop / Start cycle button

8.2 Starting up

1. Press the start switch and the screen will light up at the same time.
2. Press keys (2)/(3) until you select the desired temperature value in the green display (the upper one).
3. Press key (1) if you want to select an operating time. If time is 00:00, regulation will be continuous until the STOP key is pressed.
4. Press START key (4) to start the cycle.
5. The heating symbol  will light and blink.
6. Press the key (4) again to stop the incubation process. Note: during operation, the temperature or time adjustment value can be changed.

COMMENTS:

- * Do not use the equipment at a lower temperature than ambient +5°C, because the thermostat might have difficulties in regulating correctly.
- * When using the time function, the green display (the upper one) will indicate the remaining time. If you want to see the target temperature value, press key (1).

9 Error messages

P.FAIL Power failure. The equipment continues with the process.

Err will appear in green on the top line of the display.

TSEG Safety thermostat is open.

SENS Open temperature probe. Control does not detect the temperature probe.

HEAT The heating element does not work. The equipment detects that the temperature does not increase.

10 Adjustments

The offset (A) and the °C/°F unit (B) can be adjusted.

A. Adjust the probe reading correction (offset)

Press keys (2) and (3) continuously and at the same time (about 5 seconds). The **OFFs** message will appear.

Use keys (2) / (3) to set the actual temperature reading value. Press key (4) to confirm.

B. Presentation unit: °C / °F

Press key (1). The message **Unit** will appear. Using keys (2) / (3), select the unit: °C / °F. Confirm with key (4).

Do not touch the blocs or its housing: There is risk of burnings. To place or remove the blocks into its housing user always the suplied handle.

11 Maintenance

Before of remove the cover to handle inside, unplug the the equipment from the socket.

The handle of the electronic circuit boards by unskilled personel may cause damages. Make sure to give the equipment to an authorised technical assistance.

11.1 Cleaning

To clean the different parts of the equipment recommend to use a cotton cloth wet by a water and alcohol solution.

11.2 Spare parts

The most frequent spares are the following:

| Spare description | 7462200 | 7471200 |
|-------------------------------|---------|---------|
| Main switch | 20017 | 20017 |
| Electronic temperature switch | 29156 | 29156 |
| Control circuit board | 38500 | 38500 |
| Heating element | 39147 | 39148 |
| Safety circuit probe | 43037 | 43037 |
| Control circuit probe | 43028 | 43028 |
| Mains filter | 45015 | 45015 |

12 Warranty

This product is guaranteed for one year. The guarantee does not cover damages caused by an incorrect use or causes beyond the control of J.P. SELECTA, s.a.u.

Any manipulation of the equipment by personnel not authorized by J.P. SELECTA, s.a.u. automatically cancels the guarantee.

Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.