**PERCOM N-M II**

**3001301**

**BOMBA PERISTÁLTICA**  
**PERISTALTIC PUMP**

## Indice

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## 1. General information

- Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributors immediately.

- Do not install or use the equipment without reading this handbook first. This instruction manual must be attached to the equipment and must be available for all users.

- If you have any doubts or enquiries, contact with J.P. SELECTA's technical service.
- Do not use the equipment in explosive, flammable or corrosive atmospheres.

## 2. Packing list

PERCOM N-M II Pump (with head)	<b>3001301</b>
Instructions manual	80357

## 3. Accessories

Remote Control Pedal	1001045
----------------------	---------

## 4. Technical specifications

Model:	<b>3001301</b>
Speed range:	0.1-350rpm, reversible
Speed resolution:	0.1rpm
Speed control mode:	Digital control
External control:	TTL level control start/stop & change direction
	0-5V/0-10V/4-20mA control speed
Communication:	RS232/RS485 interface. Supports standard MODBUS communication protocol
Power supply:	AC 220 ±10%, 50Hz/60Hz
Consumption:	22W
Working conditions:	Temperature 0-40°C
	Relative humidity < 80%
	No condensation
Drive dimensions (LxWxH):	208 x 131 x 194
Pump weight:	3.98Kg
IP rate:	IP31

### 4.1 Main Functions and Features

- Control speed through a keypad. Manually control the rotating speed or automatically by an external control interface. It is easy to operate.
- Match with a dispensing controller to achieve the filling function.
- Fill and empty key that activate the motor at maximum speed.
- High torque, low power loss and a good adaptability.
- Drive circuit has an excellent performance, a good radiating condition, a low working noise and a stable running. It has a memory function to automatically store the running parameters in time.
- It is not suitable for outdoor application.

## 5. Equipment Description

With an outstanding silent performance, this standard peristaltic pump has become the first choice of laboratory application. It usually supports various analytical instruments, and it is also highly appreciated by their costs performance.

With an easy-load pump head to easy load and fix the tubing, it can connect multichannel pump heads in series to use according to customers' requests. The pump housing has a perfect rigidity and structure, with a high temperature and corrosion resistance, especially in the strong chemical corrosion resistance, such as organic solvents, etc. This feature decreases effectively the damage ratio of the pump head, and also reduces maintenance costs for customers.

Rotating speed can be manually controlled or automatically through an external interface. The memory function, after connecting power again, works according to the last state of load. Timing test function is very convenient for users to test and calibrate flow rates. It provides 7 kinds of external control modes for an easy and flexible use.

These pumps are normally used for university laboratories and research institutes; for ion chromatography and titrators; for industry production.

## 6. Installation

Place the pump on a flat, horizontal and level surface, trying to leave a space of about 10 cm at the back and on each side.

**CAUTION!!! IMPORTANT FOR YOUR SAFETY**

**Do NOT use the pump to dry or for thermic treatment of products which could form explosive gases.**

**Be sure that the voltage supply is the same as the one indicated on the characteristics plate of the equipment.**

**Do not use the equipment if it is NOT EARTHED.**

**If you have to change the plug bear in mind the following:**

**Blue cable: Neutral.**


**Brown cable: Phase.**

**Yellow/green cable: Earth.**

## 7. Operation

### 7.1 Operation Panel

- **Direction Key**

Drive control running direction. By pressing left / right reversal key once, the drive running direction will change once. When the drive rotates left, left rotation display lamp lights; when the drive is rotating right, right rotation display lamp lights.

- **Start/ Stop Key**

To run and stop the drive. By pressing the key start / stop once, the drive running state will change once. When the drive is running, the start / stop display lamp lights; also, when the drive stops working, the start / stop display lamp turns off. If you press the key for a long time, you could set the RS485 address. If you press the key for more than 5 secs, a number will blink, press the acceleration / deceleration key to set the required address, and then press the start / stop key to save the address.

- **Speed + / - Key**

Use when increasing or decreasing the motor running speed. By pressing the acceleration / deceleration key once, the motor speed will increase or decrease 0.1r/min. If pressing the key in a continuous way, the drive rotating speed will quickly increase or decrease for shortening the adjustment time.

- **Full speed Key**

Under the internal control mode, if pressing the full speed key, the pump will run at full speed. The, "Full" state will be shown, and it will finish emptying quicker, during the filling or the washing operation; if pressing the full speed button again, the pump will return to the normal status. Under the full mode, the start / stop key will also stop the drive from running.

- **Menu Key**

Enter the menu mode with the unit stopped. With the + key you can select the corresponding options, press the menu key again to accept.

- **External control mode Key**

The external control mode key only has effect when the drive stops working. It is used for cut over the different kinds of external control modes.

There are 7 different kinds of external modes:

- 0) Internal control mode (the display lamp of external control start / stop, the external control direction & the external control speed are off).
- 1) External control start / stop mode (the display lamp of external control start / stop is on; the external control direction & the external control speed indicator are off).
- 2) External control commutation mode (the display lamp of the external direction is on; the external control start / stop and the external control speed light are off).
- 3) External control speed mode (the display lamp of the external speed control is on; the external control start / stop and the external control direction are off).


- 4) External control start/stop and direction mode ( the display lamps of the external control start/stop and the direction are on; the display lamp of the external control speed is off).
- 5) External control start/stop and speed mode (the display lamps of the external control start/stop and the speed are on; the display lamp of the external control direction is off).
- 6) External control direction and speed mode (the display lamps of the external control direction and the speed are on; the display lamp of the external control start/stop is off).
- 7) External control start/stop, direction and speed mode (the display lamps of the external control start/stop, the direction and the speed are on).

**Note: When the display lamp of the external control is on, the relevant key of the internal function will be locked.**

### 7.2 Remote Control

The 15 needle interface on the back of the machine is the interface which control the machine through the external signal. We will give customer standard external control circuit when leave factory, the pin definition of external control circuit is as below:

- 1 pin WK-SPEED (blue wire): The signal of external connecting control motor rotating speed. Standard 0-5V.
- 2 pin WK-RUN/STOP (green wire): The signal input part of external control start/stop (high level drive stop running, low level drive run).
- 3 pin WK-DIRECT (red wire): The signal input part of external direction control (high level and low level control separately the left-hand and right-hand).
- 15 pin GND (black wire): Ground wire.


PIN	E-CONTROL INPUT INTERFACE DEFINITION	INPUT WIRE COLOR
1	E-c Ground	Black
2	485 Interface A	Brown
3	E-c Ground	Black
4	Start/stop 2	Grey
5	Direction	Green
6	E-c Ground	Black
7	0-10khz Input	Blue
8	0-5v Input	Blue
9	485 Interface B	White
10	E-c Ground	Black
11	Start/stop 1	Red
12	E-c Enable	Yellow
13	E-c Ground	Black
14	4-20m A Input	Blue
15	0-10 V Input	Blue

Note:

- 1.External control receive 0-5V voltage signal for factory default, if other signal is needed, please contact distributors. To avoid damage the machine, remould voluntarily is forbidden.
- 2.External control start/stop,reversing signal, both passive signal
- 3.To avoid burnout the external control interface, the plug must be inserted and plucked in off-position.

DB15 external control port regular line sequence as below:

**RS232 Communication Pins Definition(standard):**

6 pin 232-T: upper computer send data, joint the 3 pin of DB-9 connector.  
 7 pin 232-R: upper computer receive data, joint the 2 pin of DB-9 connector.  
 15 pin GND: ground line, joint the 5 pin of DB-9 connector.


**RS485 Communication Pins Definition(optional):**

4 pin 485-B-(red wire): 485 communication signal end, joint 485 transverter B-.  
 5 pin 485-A+(yellow wire): 485 communication signal end, joint 485 transverter A+.  
 15 pin GND(black wire): ground line, joint 485 transverter GND.


**Segregated Type External Control Pins Definition:**

12 pin yellow wire: external control start/stop  
 Segregated external control signal input port, contact 24V to start by default.  
 13 pin red wire: external control reversing  
 Segregated external control signal input port, contact 24V to reverse by default.  
 14 pin black wire: Segregated interface signal input ground line.

**Segregated External Control Connection Diagram:**

Note: Percom N-M II pumps can match both foot pedal and dispenser.


**Notice to customers:**


The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.


## 8. Operating menus

When A0-0 (speed display) mode is selected, turning on A1 option menu (external control, system setting). When the A0-1 (flow display) mode is selected, turning on A1 option menu (external control, system settings), A2 (pump head, pump tubing setting), E0 (calibration function).

**Application one: speed mode transmission.**

Transmission fluid to rotate speed of 80rpm/min

1. Press button to enter A0-- (display of choice)
2. Press the knob to enter A0-0 (and rotate speed display mode)
3. Press the knob to determine the current model for the speed display mode.
4. To rotate ESC and press the knob (return to the main screen)
5. Rotate the knob to 80.0
6. Confirm operation of 80rpm/min: 80.0

**Application two: flow mode transmission**

Transmission fluid to rotate seed of 80ml/min

1. Press button to enter A0-- (display of choice)
2. Press the knob and rotate to enter A0-0 (flow)
3. Press the knob to determine the current model for the flow pattern.
4. To rotate and press the knob ESC (return to the main screen)
5. Rotate the knob to 80.0
6. Confirm operation of 80rpm/min: 80.0

**Application three: calibration fluid**

With 8ml/min speed transmission fluid

1. Turn the peristaltic pump to flowrate moderate
2. Need to choose pump head and tube
3. Return to main screen, adjust the flowrate as 8ml/min. If the flowrate is not precise, one can use calibration function to adjust.
4. Enter to EO-- displaying 8ml
5. Confirm that peristaltic pump will transfer and auto stop when reach the set time, and metering and store the liquid volume.
6. Pull to stop, displaying 8000. Use the digital knob to input the actual calibration volume.
7. Then press the knob to the main screen.

## 9. Three Roller Pump Head

### 9.1 Installation of Head

Take two special screws to insert into the open holes of the pump head separately. Insert the pivot bracket of the pump into the groove of the drive shaft joint, and put the two special screws to aim at the drive connecting holes at the same time. Wind up the screws separately. Use the long screw under series to wound many pump heads.

Note: The user try to use the same power to wind up the two screws and also not overwind for avoiding appearing the running noise because of under bracing changing the shape.

### 9.2 Installation of Tubing

Pull the lever 180° counterclockwise, and open the top board. Lift the cartridges of two sides separately, and put the tubing between the rollers and the top board. Pull the lever 180° clockwise and push the top board to the correct position. The cartridge jams the tubing stably for preventing the tubing movement in the process of work.

### 9.3 Tubing Technical Parameters

<b>Cross sections (1:1)</b>							
<b>Wall thickness</b>	1.6						
<b>ID (mm)</b>	0.8	1.6	2.4	3.1	4.8	6.4	7.9
<b>Max pressure (Mpa)</b>	continuous	0.17		0.14	0.1	0.07	
	intermittent	0.27		0.24	0.14	0.1	

## 10. Maintenance

Unplug the equipment from the mains before removing the oven's cover to handle inside.

**The manipulation of the equipment internal electronic circuits by unauthorized personnel can cause irreparable damage. Take it to one of J.P.SELECTA's authorized technical services.**

The pump is life time lubricated.

Check condition of power cord (and foot switch tube) annually for cracks and wear. When pump is not working, please loose the cartridges of pressing the tubing for avoiding changing the shape of tubing because of longtime extrusion.

Pump head cannot resist super corrosive liquid. Please pay attention to it when it is using.

Keep the rollers of pump head clean and dry. Do not use aggressive solvents.

If the surface of rollers is not clean, they can quicken the tubing wearing, and reduce the useful life of tubing.

If there is liquid on the rollers, please dry them. Longtime moisture can damage the rollers.

## 11. Troubleshooting

Trouble	Cause	Solution
Fan and display normal work but the pump does not wok	Motherboard blown fuse	Change it
	The unit is broken	Change it
Fan, motor and display normal work, but the pump does not work	The joint is broken	Change it
The pump is working but display cannot be seen	The connected line is lost or damaged	Reinstall the line or change it

## 12. Warranty

This product is guaranteed for one year. It does not cover damages caused by an incorrect use or causes beyond the control of J. P. SELECTA, s.a.u.

Any manipulation of the equipment by personnel not authorized by J.P. SELECTA, s.a.u. automatically cancels the guarantee.