SCRUBBER

4001611

BR-1612

4001612

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**Scrubber unit and
Water circulating vacuum pump**

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Table of content

- 1 Safety10**
- 1.1 Hazards10
- 1.2 Risks which the operator is put under10
- 1.3 Personnel training.....10
- 1.4 Proper uses of the equipment.....10
- 1.5 Modifications.....10
- 2. Packing list10**
- 3. Description11**
- 4. Technical specification:12**
- 5. Installation12**
- 6. Operation.....13**
- 6.1 Water circulating vacuum pump:13
- 6.2 Gas neutralizer: «Scrubber»14
- 6.2.1 Neutralization reagent addition14
- 6.3 Starting the complete digestion system:14
- 7. Maintenance16**
- 7.1 Scrubber maintenance:16
- 7.2 Vacuum pump maintenance.....16
- 7.3 Spare parts16
- 8. Warranty16**
- Notice to customers16*

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1 Safety

The **SCRUBBER** unit and the water circulating vacuum pump have the suitable measures for their use in laboratory. However, the fat analysis method can be hazardous and the user must be informed of some of those risks.

1.1 Hazards

It identifies the risk and the safety actions to be applied.

Hazard of contact with corrosive reagents

Hazard of contact with solvents. Follow the instructions for the described operation.

Hazard of electrical shock

Electrical hazard for entering in areas marked with this signal. Follow the instructions for the described operation.

Hazard of burns to be in contact with hot surfaces

Temperature in the heater surface exceeds 200°C. Use the describer tools to perform the operation on the unit.

1.2 Risks which the operator is exposed to

The process of determination of fat with the described equipment exposes to the operator to the following risky situations:

- Manipulation of strong acids
- Manipulation of glass parts.
- Possibility of touching parts to a temperature above 60°C.
- Electrical risk.

1.3 Personnel training

This equipment must only be operated by personnel trained and suitable qualified to work in a chemical laboratory and its common hazards and who have read and understood this manual or have been qualified on its manipulation and operation.

1.4 Proper uses of the equipment

This equipment is for laboratory use, for fat extraction sample. Mainly in soxhlet analysis.

The improper use of this unit (not following manual instructions) can endanger the protection prepared by the equipment.

This equipment is NOT designed and not manufactured to work in an explosive ambient area except the ones described on this manual.

1.5 Modifications

The operation modification or the manipulation of the safety systems unauthorized by the manufacturer, could expose the operator to hazards not indicated in this manual.

2. Packing list

Scrubber:	4001611
2x	Bottles with gas trap glassware.
1x	Inox holder.
1x	1m ISOVERSINIC Ø 10x14 hose
1x	Flask of 3Kg Neutralizing solution Na ₂ CO ₃ Cod.: 31242
Circulating pump:	4001612
1x	Pump
1x	Power cord

3. Description


Fume suction and neutralize system specially generated in the Kjeldahl digestion process.

The system uses a scrubber unit than retains the acid condensations under a Na_2CO_3 solution, and a water circulation pump that provides a big flow rate of vacuum to fume suction.


The vacuum pump 4001612 uses water circulation to provide high level of suction on the front two taps. This scrubbing application only uses one tap.

The 4001612 pump could be used in other laboratory application to liquid handling, cooling... as it could be configured for water circulation.

Components for a complete Kjeldhal digestion system

4. Technical specification

• Maximum vacuum level:	0.098MPa (0.98bar)
• Sucking air flow rate:	10 L/min (single tap)
• Suction:	Two aspiration taps.
• Tank capacity:	10 litres.
• Ports for external water circulation	
• Consumption:	280W
• Power supply:	230V 5A fused
• Noise level:	Lower than 65dB
• Made from:	ABS and stainless steel.

5. Installation

The digestion complete set will have to include the following equipment:

- Digestion Block
- Temperature controller.
- Fume collector.
- ISOVERSINIC Ø12x17 hose
- Scrubber.
- ISOVERSINIC Ø10x14 hose
- Water circulation vacuum pump.

Choosing the right working place

On the digestion process, strong acids and reagents are used and must handle with caution, also irritant fumes are generated. For this reason the working place must be chosen with care to avoid inconvenience on other laboratory activities.

The location must have the following conditions:


- Must be wide enough for safety work.

- Ventilation available. It is not necessary to install the whole chain into an extraction cabinet, but the room must have windows or air removal system to ventilate if some time any gas leakage is produced (The Kjeldahl digestion gases are not toxic but irritant).
- At least 2 electricity plugs must be available.
- Water tap and drainage are not necessary but suggested.
- The surface must be flat, stable and strong.

Hoses Connecting:

Connect the hoses according the following diagram:

Only VITON, TEFLON or suitable material must be used for hoses.

Fill the bottle indicated on the scrubber with the reagents according (see 6.2.1).

Fill the pump blue plastic tank with around 10-12 litres of low calcium content water.

Hoses for connection:

If the supplied hose set is too short to use, you can contact J.P. SELECTA, s.a.u. dealer with the following specification:

Fume collector - scrubber connection 12x17 Hose :


- Code **46085** ISOVERSINIC black 12x17

Scrubber - Vacuum pump connection 10x14 Hose :

- Code **46084** ISOVERSINIC black 10x14

Fig. 6.1 Configuration Water Pump

- A- Circulating Water: «OFF»
 B- Emptying tube: «Connected»
 C- Water IN /OUT: «N/Connection»

6. Operation

6.1 Water circulating vacuum pump

It is formed with a motor-pump that produces a water circulation through two water jet pumps which produce the vacuum.

The TAP1 and TAP2 ports aspirates an air flow rate around 10L/min each one. The suctioned air is melt on the tank water.

In order to increase the vacuum pump life, it is strongly suggested to use the scrubber unit before the pump to neutralize the gas that comes from the digester.

Apply this configuration for water pump 40001612:

- Switch «CIRCULATING WATER»: OFF position.
- Emptying tube: «CLOSED»
- Do not connect anything on IN and OUT water ports.

**WARNING:**

Use the scrubber in order to neutralize the sucked air before it arrives to the vacuum pump.

Do not directly connect the vacuum pump to the digestion unit. This situation damages the vacuum pump.

Pump setup:

- Check that the drain tube is correctly inserted in the tank upper connector.
- Fill the tank with about 10-12 litres of low calcium content water.
- Set the external «Circulating water switch» in the OFF position.

The pump is now ready to be used.

Pump start:

When switching «ON» the main switch, the pump starts suction by means of «TAP1» and «TAP2» ports.

To check if it is correctly working, please close with the finger «TAP1» and «TAP2» ports and check the vacuum-meter indicates around 0.098.

6.2 Gas neutralizer: «Scrubber»

The scrubber bottles and the vacuum pump tank must be filled with a sodium carbonate dissolution to neutralize the digestion acids gases.

(Part number: **31242**)

6.2.1 Neutralization reagent addition

Bottle 1: 600ml water + 150gr Na_2CO_3

Bottle 2: 600ml water + 150gr Na_2CO_3

Vacuum pump tank: 10gr Na_2CO_3 + 10-12L H_2O

It is important to check that there's no leak in the hose connections. To do this, once the hose is connected to the vacuum pump, check that the vacuum level in the pump is about 0.098 by closing with the finger the gas entry on the collector fume's side.


Fig 6.2.1 Addition of neutralization reagents.

6.3 Starting the complete digestion system

The suggested sequence for preparing and starting the complete digestion system is:

- Prepare the reagents and the samples in the digestion tubes. See the digester's instructions.
- Program the temperature controller for a three steps digestion cycle. According to the digester instructions manual.
- Insert the sample tubes in the block-digester and cover them with the fume collector.
- Check that the vacuum pump is ready. That is, the water tank and the external circulating are OFF.
- Check that the scrubber bottles have the anhydride sodium carbonate solution. (See 6.2.1)

WARNING!

Digestion process needs strong chemical products. Pay attention to its security sheets about how to use it appropriately.

Only staff suitable trained on chemical products operation should work on digestion process of samples.

- Assembly the hoses:
 - Ø 10 hose from the Scrubber to the vacuum pump.
 - Ø 12 hose from the collector to the Scrubber.
- Start the digestion cycle from the temperature controller.
- Start the vacuum pump.
- Visually check the digestion process every 30 or 40 minutes.
- When finishing the digestion process, the temperature process beeps. Then carry out the following operations:
 - Remove the tubes rack from the digester without disconnecting any hose and leave it on the holder.
 - Keep the vacuum pump on.
 - Carefully carry out this operation.

WARNING: Possible sulphuric acid spillage.

Use protective glasses and gloves for any manipulation of the rack with samples or if you just want to check the digester for a visual inspection.

- Let it cool for about 15-30 minutes with the fume collector system connected.
- Stop the vacuum pump.
- Carefully disconnect the hose from the fume collector if you want.

WARNING: Use gloves.

CAUTION: The drops that can fall from the collector or from the hose are sulphuric acid drops. Manipulate with care to avoid chemical burns on the skin.

- Remove the fume collector from the tube rack and leave it on the holder.

• **Fume collector cleaning:**

You can simply clean the fume collector with water.

• **Water tank emptying:**

To empty the tank, just remove the «emptying tube» from the upper connector.

• **Water circulating function:**

The vacuum pump can also work as a water closed circuit.

- Connect the hose to the «water circulating out port».
- Place the «water circulating switch» in the ON position.
- When starting the pump, this will move the tank water to the «water circulating out port».

7. Maintenance

7.1 Scrubber maintenance

Along the time, the Na_2CO_3 solution could damage the glass bottle.

Clean the Scrubber bottles with abundant water after use.

Monthly visually check the hoses status and, if losses are found, change them.


7.2 Vacuum pump maintenance

Never left the vacuum pump filled of used water after the digestion process.

On a three month basis, visually check the status of the internal hosing to substitute damaged ones.

To do it, remove the upper cover. Before removing the cover, unplug the pump from the power supply.


7.3 Spare parts

8. Warranty

This product is guaranteed for one year from date of purchase against faulty workmanship. The guarantee does not cover damages caused by an incorrect use or causes beyond the control of J.P. SELECTA, s.a.u.

Any mistreatment of the equipment by unauthorised personnel not approved by J.P. SELECTA, s.a.u. cancels the guarantee automatically.

Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.