**BOMBA PERISTÁLTICA PR-2004**  
**PERISTALTIC PUMP PR-2004**

**4129051**

## Indice

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| PIN | E-C INPUT INTERFACE DEFINITION | E-C OUTPUT INTERFACE DEFINITION |
|-----|--------------------------------|---------------------------------|
| 1   | E-c Ground                     | Suspended                       |
| 2   | 485 Interface A                | E-c output enable output B      |
| 3   | E-c Ground                     | Start/stop output B             |
| 4   | Start/stop 2                   | Direction Output B              |
| 5   | Direction                      | E-c Ground                      |
| 6   | E-c Ground                     | 4 ~ 20mA output negative        |
| 7   | 0-10kHz Input                  | E-c Ground                      |
| 8   | 0-5V Input                     | E-c Ground                      |
| 9   | 485 Interface B                | E-c output enable output A      |
| 10  | E-c Ground                     | Start/stop A                    |
| 11  | Start/stop 1                   | Direction Output A              |
| 12  | E-c Enable                     | 0~10kHz output                  |
| 13  | E-c Ground                     | 4 ~ 20mA output positive        |
| 14  | 4-20mA Input                   | 0~10V output                    |
| 15  | 0-10V Input                    | 0~5V output                     |

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## 1. General information

- Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributors immediately.
- Do not install or use the equipment without reading this handbook first. This instruction manual must be attached to the equipment and must be available for all users.
- If you have any doubts or enquiries, contact with J.P. SELECTA's technical service.
- Do not use the equipment in explosive, flammable or corrosive atmospheres.


## 2. Packing list

|                          |                |
|--------------------------|----------------|
| Peristaltic pump PR-2004 | <b>4129051</b> |
| Instructions manual      | 80418          |
| Remote Control Pedal     | 1001046        |

## 3. Technical specifications

|                           |   |
|---------------------------|---|
| Model:                    | <b>4129051</b>  |
| Speed range:              | 0.1 to 600 rpm  |
| Speed resolution:         | 0.1rpm  |
| Control method:           | TFT screen  |
| External control:         | Start/Stop, direction, speed control (optional), RS485 serial communication |
| Power supply:             | AC 220V±10%, 50Hz / 60Hz  |
| Consumption:              | <50W  |
| Working conditions:       | Temperature 0-40°C<br>Relative humidity <80%                                |
| Drive dimensions (LxWxH): | 254x184x145mm   |
| Pump weight:              | 4.8Kg   |
| IP rate:                  | IP31  |
| Rango del caudal:         | 0.0015 to 2280 ml/min   |

### 3.1 Main Functions and Features

- A - Main switch
- B - External control interface (RS485)
- C - Fuse
- D - Power connector

## 4. Equipment Description

- Equipment with 3.5 inches color LCD screen. Running states real-time displaying function, showing the flow rate and motor revolving speed at the same time, touch screen controlling.
- With intelligent calibration function, could calibrate the flow rates and infusion volume; suitable for high precision transferring fluids.
- Customer data saving function: can save 9 sets of data for customers.
- Accurate angle controlling function ensures to high precision distribution and metering.
- Memory function available in case of power off, reliable and safety.

## 5. Installation

Place the pump on a flat, horizontal and level surface, trying to leave a space of about 10 cm at the back and on each side.

**CAUTION!!! IMPORTANT FOR YOUR SAFETY**

**Do NOT use the pump to dry or for thermic treatment of products which could form explosive gases.**

**Be sure that the voltage supply is the same as the one indicated on the characteristics plate of the equipment.**

**Do not use the equipment if it is NOT EARTHED.**

**If you have to change the plug bear in mind the following:**

**Blue cable: Neutral**  
**Brown cable: Phase**  
**Yellow/green cable: Earth**

## 6. Operación

### 6.1 Function buttons

1. Home key
2. Start/Pause key
3. Direction key
4. Full speed key
5. Display

**1. Home key**

After user's entry the interface, you can press this button to quickly return to the main interface without saving the data.

**2. Start/Pause key**

Control drive's run and stop, first time you press the button the drive will run, press again the drive will stop, and then press it, the drive continues to run.


**3. Direction key**

Each time you press the button, the revolving direction of the drive will be changed.

**4. Full Speed key**

Also known as emptying and cleaning key, press the button the drive will run at full speed, press again to return to the original state.



**6.2 Areas in control panel**


1. Graphical menu area: Settings & Data Operation / Lock & Unlock / Sound on/off / Remote Control
2. Filling parameter input area: Filling Volume / Filling Time / Pause Time / Copy Number
3. State parameter display area: These parameters are displayed: Animation Instructions Operation / Direction / Pump Head Model / Pump Tubing Model / Rotate Speed / Flow rate / Filling Volume/Filling Number

1. System settings
2. Data storage
3. Lock /unlock panel to prevent operate
4. Sound on / off
5. Remote control mode

### 6.3 Graphical menu

## 7. Parameters settings

Instructions for filling parameter area:

### 7.1 Set the Filling Volume

1. Click the filling quantity value, will show the small keyboard.
2. Click the Volume Range 0.001ml-9999L (show the max. and min. value for the using pump head and tubing, can only input the value in this range).
3. Unit: ml / L
4. When the filling stop can click on the input value, but if the drive running, you can't click on the input value.
5. Set the keyboard:

| Filling volume: |   |              |     |
|-----------------|---|--------------|-----|
| Max: 9999L      |   | Min: 0.001ml |     |
| 124.5 ml        |   |              |     |
| 1               | 2 | 3            | C   |
| 4               | 5 | 6            | U   |
| 7               | 8 | 9            | Esc |
| 0               | . | OK           |     |

Enter the filling volume you want

Delete

Unit setup

Exit

Confirm

### 7.2 Set the Filling Time

1. Click the Filling time quantity value, will show a small keyboard, then input your filling time
2. Range: 0.1sec - 9999min; Unit: sec, min
3. When the filling stop you can click on the input value, but if the drive running, you can't click on the input value
4. With countdown function
5. Set the keyboard:

| Filling time: |   |           |     |
|---------------|---|-----------|-----|
| Max: 9999min  |   | Min: 0.1s |     |
| 234.5 s       |   |           |     |
| 1             | 2 | 3         | C   |
| 4             | 5 | 6         | U   |
| 7             | 8 | 9         | Esc |
| 0             | . | OK        |     |

Filling time range

Delete


Unit setup

Exit

Confirm


### 7.3 Set the Stop Time

1. Click the Stop Time quantity value, will show a small keyboard, then input your interval time
2. Time Range: 0.1sec-9999min
3. Units:sec(second), min (minute)
4. Stop time value is to set the interval time between two fillings; it can be only set when the filling is stopped.
5. During the filling period, the stop time cannot be set
6. With countdown display function
7. Set the keyboard:

### 7.4 Set the Filling Number

1. Click on the Filling Number quantity value, will show a small keyboard, and then input your filling number
2. Range: 0-9999 times, 0 on behalf of the infinite time
3. When the filling stop you can click on the input value, but if the drive running, you can't click on the input value
4. Set the keyboard:

### PUMP TUBE SELECTION

Select the current use of the pump tubing models

Turn on the power, enter the initial interface, Click Settings, enter the setup menu; Click System Set, enter the system setup menu; Click Tube, enter the pump tube selection menu; then choose your current pump Tube. Click Enter, save and exit; Click Cancel, exit without saving. (if you just want to change the tubing which the same pump head use, you can just enter the Pump tube selection menu, to choose the suitable model).

## 8. Configuration

### External control input

Select External control input parameters.

Turn on the power, enter the initial interface, Click Settings, enter the setup menu; Click System Set, enter the system setup menu; Click External control Input, enter the external control Input menu; choose the suitable control items; Click Enter, saved and exit; Click Cancel, exit without saving.

### Communication Settings

Select the pump communication parameters.

From initial interface press the System icon to enter the setup menu; click System Settings to enter the system setup menu; click Communication Settings to enter the communication settings menu; choose the appropriate control items (including: broadcast address, local address, baud rate); click Enter to save and exit; click Exit to exit without saving.

The pump also has 485 serial communication bus interfaces to be connected to the host computer (computer, PLC, SCM). PC can simultaneously connect up to 30 devices with 485 serial communication bus interface functions; when using with multiple devices to communicate with the host computer, it must know the machine number for each device, this device ID is the only number that should be connected together, and each machine number of the devices is not the same. The factory default value for each pump is "1"; it could be reset by remote controlling through the host computer or by manually changing from the setting menu. Change the following steps: Press the Enter key to enter the menu options, move the cursor to adjust the parameter settings, press the Enter key to enter, find the machine number option and press the Enter key to enter, you can adjust the machine number.

(see chapter **10.2 Settings > external control**)

### Suck-back Set

Select Suck-back parameters

Turn on the power, enter the initial interface, Click Settings, enter the setup menu; Click System Set, enter the system setup menu; Click Suck-back, enter the suck-back set menu; Choose the suitable control items (angel and speed); Click Enter, saved and exit; Click Cancel, exit without saving.

### The gear icon

Show the running state of the pump.

When the pump is running, the gear will run; When the pump stop, the gear stop.

### Direction & Full speed icon

By pressing the arrow keys to change the running direction of the motor, so different direction arrows will change.

Press the Full Speed button will display double arrows, indicating motor full speed running.


### Other status values shown

- Pump Head show
- Pump Tube show
- Speed display, the unit is rpm, the system automatically displays based on flow rate
- Filling times showed a cumulative function
- Real conveying fluid volume display with countdown function
- Real filling time display with countdown function

## 9. Calibration

Pump accuracy calibration, in order to reach the high precision dispensing, when the pump replace a new hose, external conditions changed or other conditions causing flow rate deviation, the flow rate can be corrected to meet users' requirements of flow accuracy. Turn on the power, enter the initial interface, Click Settings, enter the setup menu, Click Calibration, enter the calibration guide menu:


Calibration Wizard:

Pipe:14#, Flowrate:224.8ml/min

Amount of flowrate: 320ml

(pls enter the actual value you need)

Next Exit

Calibration Data:

Liquid volume 1: 321.3ml


Liquid volume 2: 321.3ml

Liquid volume 3: 321.3ml

(Tested repeatedly input the results)

Back Next Exit

The interface shows the current using pump tube, conveying fluid volume and advice calibration volume. If the suggested calibration volume is big, you can reduce the value according to the filling condition; input the calibration value and then click Next, enter the calibration menu, or click Exit Back to previous menu.

Running:

Time reaminging:85. 41s

Cancel

The pump will automatically stop after finishing setting, or click Cancel to cancel the calibration, and return to the calibration interface. After 3 times calibration and enter the measured values, click Next to save the data, click Exit, to finish the calibration.

Mode Selection

The system default allocation mode

Data

Preserve the existing filling data; you can save at most nine filling data

It has Delete Save Use Exit Options

## 10. Settings

### 10.1 3 Roller Pump Head

#### Micro flow tube

| Tube                 | 0.5X0.8    | 1X1 | 2X1 | 3X1 | 2.4X0.8 | 3.2X0.8 |
|----------------------|------------|-----|-----|-----|---------|---------|
| Hose section         |            |     |     |     |         |         |
| Wall thickness (mm)  | 0.8        | 1   | 1   | 1   | 0.8     | 0.8     |
| Inside diameter (mm) | 0.5        | 1   | 2   | 3   | 2.4     | 3.2     |
| Pressure (Mpa)       | Continuous |     |     | 0.1 |         |         |
|                      | Interval   |     |     | 0.1 |         |         |

#### The basic flow tube

| Tube                 | 13#        | 14# | 19#  | 16# | 25#  | 17# | 18#  | 15# | 24#  | 35# | 36#  |  |
|----------------------|------------|-----|------|-----|------|-----|------|-----|------|-----|------|--|
| Hose section         |            |     |      |     |      |     |      |     |      |     |      |  |
| Wall thickness (mm)  | 1.6        |     |      |     |      |     | 2.4  |     |      |     |      |  |
| Inside diameter (mm) | 0.8        | 1.6 | 2.4  | 3.1 | 4.8  | 6.4 | 7.9  | 4.8 | 6.4  | 7.9 | 9.6  |  |
| Pressure (Mpa)       | Continuous |     | 0.17 |     | 0.14 |     | 0.1  |     | 0.07 |     | 0.17 |  |
|                      | Interval   |     | 0.27 |     | 0.24 |     | 0.14 |     | 0.1  |     | 0.27 |  |

#### Industrial hose

| Tube                 | 73#        | 82# | 86#  | 90# | 88#  | 92#  |
|----------------------|------------|-----|------|-----|------|------|
| Hose section         |            |     |      |     |      |      |
| Wall thickness (mm)  | 3.3        |     | 6.3  |     | 4.8  |      |
| Inside diameter (mm) | 9.6        | 13  | 9.5  | 19  | 12.7 | 25.4 |
| Pressure (Mpa)       | Continuous |     | 0.17 |     | 0.1  |      |
|                      | Interval   |     | 0.27 |     | 0.1  |      |

#### 10.1.1 Head installation

Take two special screws to insert them into the open holes of the pump head separately. Insert the pump pivot bracket into the groove of the drive shaft seal and, at the same time, fit the two special screws into the connection holes of the disc. Screw them separately.

Use the long screw at the bottom to attach multiple pump heads.

Note: The user must try to use the same power to tighten both screws.

#### 10.1.2 Installing the tubes

Pull the lever 180° to the left and open the top plate.

Lift the cartridges from both sides separately and place the tube between the rollers and the top plate.

Pull the lever 180° clockwise and push it into the correct position. The cartridge presses the tube stably to prevent movement of the tube during the working process.

#### 10.1 External control

Drive external control interface (DB-15) introduction:

The 15-pin interface on the back of the machine is the controlling interface to operate the machine through the external signals.

Firstly prepare a DB15 connector (with holes) and various colors signal wires, then open the 15-pin interface, weld on the signal wires according to the below sheet, and fix the wires together using the clip on the interface, you may inject some glue to reinforce these wires; and at last install the shell and screws.


### External control interface technical parameters

| PIN | E-C INPUT INTERFACE DEFINITION | E-C OUTPUT INTERFACE DEFINITION |
|-----|--------------------------------|---------------------------------|
| 1   | E-c Ground                     | Suspended                       |
| 2   | 485 Interface A                | E-c output enable output B      |
| 3   | E-c Ground                     | Start/stop output B             |
| 4   | Start/stop 2                   | Direction Output B              |
| 5   | Direction                      | E-c Ground                      |
| 6   | E-c Ground                     | 4 ~ 20mA output negative        |
| 7   | 0-10kHz Input                  | E-c Ground                      |
| 8   | 0-5V Input                     | E-c Ground                      |
| 9   | 485 Interface B                | E-c output enable output A      |
| 10  | E-c Ground                     | Start/stop A                    |
| 11  | Start/stop 1                   | Direction Output A              |
| 12  | E-c Enable                     | 0~10kHz output                  |
| 13  | E-c Ground                     | 4 ~ 20mA output positive        |
| 14  | 4-20m A Input                  | 0~10V output                    |
| 15  | 0-10V Input                    | 0~5V output                     |

#### Notes:

- 1.The external control interface is a selective accessory.
2. The using method of the external interface: (A) Enabled wire and Ground wire connect or shut, control the entry of the external control. (B) Start/ Stop 1 wire and Ground wire connect or shut, control the start and stop of the pump; Start/Stop 2 wire and Ground wire control the start and stop through the voltage signal (less than 30V). (C) Direction wire and Ground wire connect or shut, control the running direction of the pump. (D) Between Speed wire and Ground wire, join up 0-5V, 0-10V, 4-20mA, 0-10kHz, etc. controlling wire signal.

If you need other functions, please contact the manufacturer for external control.

## 11. Maintenance

Before removing the covers of the equipment to manipulate inside, disconnect it from the mains.

**Manipulation of electronic circuits by unauthorized personnel can cause damages difficult to repair. Check with the authorized technical service.**

The pump must be permanently lubricated.

Annually check the condition of the power cord (and foot switch tube) for cracks and wear.

When the pump is not running, release the pressure cartridges from the tube to prevent the tube from changing shape due to a prolonged extrusion.

The pump head cannot withstand super corrosive liquids. Please pay attention when you are using it.

Keep the rollers on the pump head clean and dry. Do not use aggressive solvents.


If the rollers surface is not clean, it can accelerate the wear of the tube and reduce its useful life.

If there is liquid on the rollers, wipe them dry. Prolonged humidity can damage them.

## 12. Warranty

This product has a warranty of one-year. The warranty does not cover damages caused by a misuse or causes beyond the control of J.P. Selecta, s.a.u.

Any manipulation of the equipment by unauthorized personnel cancels the benefits of the guarantee.

### Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.