<b>FIBROMAN HT-W 100ml</b>	<b>3031470</b>
<b>FIBROMAN HT-W 250ml</b>	<b>3031471</b>
<b>FIBROMAN HT-W 500ml</b>	<b>3011472</b>
<b>FIBROMAN HT-W 1000ml</b>	<b>3031473</b>
<b>FIBROMAN HT-W 2000ml</b>	<b>3031474</b>

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## General Information

Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributor immediately.

Do not install or use the equipment without reading this handbook before.

This instructions manual must always be attached to the equipment and it must be available for all users.

If you have any doubts or enquiries, please contact your supplier or J.P. Selecta.

If any modification, elimination or lacking in maintenance of any device of the equipment by the user could be dangerous, the manufacturer is not responsible for the damage that can occur.

Do not use the equipment with liquids which can give off vapours capable of making explosive mixtures.

## Packing List

The standard equipment consists of the following components:

Heating mantle	1
Power supply cord	1
External temperature supply	1
Instruction Manual	1

## Accessories

This equipment has the following accessories:

Part	Code
SUPPORT BAR 12mm diameter x 750mm long	6000270
DOUBLE BOSSHEAD	1000069
SPRING clamp to hold thermometers	7000025
CLAMP to support and hold condensers	7000161

## Technical Specifications

Code	3031470	3031471	3031472	3031473	3031474
Flask volume (ml)	100 250	500	1000	2000	
Ø Flask (aprox) (mm)	65/67 83/87	102/105	132/139	164/170	
Height / Width / Depth (cm)	11/14/18	12/18/24	12/18/24	13/21/27	15/25/31
Power (W)	160 160	300	420	760	
Ø External probe	4 4	4	4	4	
Range de T <sup>a</sup> (°C)	From ambient temperature +5 to 400°				
Temperature variation (°C)	From 0.5 to 1 sample depending				
Weight (Kg)	1.5 3	3.5	4	5	

## Equipment Description

**FIBROMAN HT-W** are a family of heating mantles for round bottom flasks which have an electronic circuit for temperature control and measure, with the possibility of being controlled from an external probe directly placed into the flask. The heating mantles are mounted inside a stainless steel case.

## Installation

The equipment must be installed over a flat surface, slightly ventilated and free of dust and moisture. Please leave a free space of about 20cm around the equipment for your safety.

Check that the voltage supply and frequency coincides with the heating characteristic plate.

Foresee a mains plug at least at 1m away.

Be sure that the electric supply grid is properly connected to the earth.

If you have to change the plug, bear in mind that the yellow/green cable is the earth, and the brown and the blue ones are indistinctly phase and neutral.

Please avoid direct light from sun and air flows.

The equipment must be located under a fume hood if the sample releases irritating or corrosive fumes.

## Operation

### Starting up

1. Connect the equipment into the mains by means of the main switch.
2. The temperature red display (Current temperature) shows 400°C for a while and then indicates the probe temperature. Then, the LED is lit and the mantle begins heating till it reaches the programmed temperature, which is the one indicated in the green display (Selected temperature).
3. The LED will light up and off indicating that the heater element is activated or deactivated.

### IMPORTANT: Mantle «Initial burning»

When the mantle is started-up for the first time, the heating mantle produces smokes and «burning» smell. This effect is completely usual because of the mantle components burning. After a few minutes it will disappear.

Its also posible that this step was run during the factory testing. In this case the mantle will show a burnt color as the «Initial burning» is already done. In any case, never, used mantles are delivered.

### Selecting working temperature

Program the working temperature by pressing the increase-decrease keys (see Figure 1).

Any pulse will increase the temperature in steps of 1°C, and if the key is kept pressed, variation will accelerate. The programmed temperature is kept in memory, regardless whether the equipment has stopped.

### Please use flaks that fit the heating mantle size.

The equipment can operate in two different ways:


Fig. 1  
Front panel.

Fig. 2  
Temperature probe connector location.

If you don't use the external probe, the «bridge connector» must be placed.

Fig. 3  
External probe operation.

### Operating without external probe:

The temperature controller shows and controls the mantle temperature. This temperature is slightly higher than the sample temperature.

This operating mode is appropriate when there's no need of a very accurate control of the sample temperature.

1. Place the «bridge connector» (see Fig. 2). It normally comes placed by default.
2. Select the working temperature.

### Operating with external probe:

The temperature controller shows and controls the probe temperature immersed in the sample. This allows an accurate control of the temperature at which the sample is being treated.

#### ATTENTION!

- The sample liquid must be compatible with the probe stainless steel.
- Do not submerge the probe in hard acids.

This operating mode is appropriate when you need a very accurate control of the sample temperature.

1. Extract the «bridge connector» (see Fig. 2) and connect the external probe (see Fig. 3).
2. Fix the probe with the help of a bar and a clamp. Then submerge it about 2 cm in the sample.
3. Select the working temperature.

#### DANGER

- Be sure that the probe is immersed in the flask liquid before switching the mantle on.**
- If not, the sample temperature will rise in an uncontrolled way.**

## Maintenance

### CLEANING:

It is recommended the use of the following products for the cleaning of the equipment different parts:  
Cleaning of casing, plastic and painted metals: Alcohol with cotton or a duster.

## Guarantee

This product is guaranteed for one year. The guarantee does not cover damage caused by incorrect use of causes beyond the control of J.P. SELECTA, S.A.

Any manipulation of the equipment by unauthorized personnel by J.P. SELECTA, S.A. cancels the guarantee automatically.