**PLACA CIRCULAR DOBLE DUPLAC  
DOUBLE CIRCULAR HOTPLATE DUPLAC**

**4000634**

**PLACAS CALEFACTORAS  
HOTPLATES**

## INDICE

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# ENGLISH VERSION

## GENERAL INFORMATION

- 1) Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributor immediately.
- 2) Do not install or use the equipment without reading this handbook before.
- 3) This handbook must always be attached to the equipment and it must be available for all users.
- 4) If you have any doubts or enquiries, please contact with your supplier or J.P. Selecta's technical service.
- 5) **IMPORTANT! J.P. SELECTA WILL NOT ACCEPT ANY EQUIPMENT TO BE REPAIRED IF IT IS NOT DULY CLEANED.**
- 6) If any modification, elimination or lacking in maintenance of any device of the equipment by the user transgress the operative directive, the manufacturer is not responsible for the damage that can occur.
- 7) Do not use the equipment with liquids which can give off vapours capable of making explosive mixtures.

### Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.

## PACKING LIST

The standard equipment consist of the following components:

Piece	Code	Quantity
Double circular hotplate	4000634	1
Instructions manual	80111	1

## EQUIPMENT DESCRIPTION

Double circular hotplate are designed to work with a big range of temperatures up to 400°C.

The temperature regulation is done by an independent energy regulator in each plate.

## TECHNICAL FEATURES

Voltage supply 115-230V 50/60 Hz according to the characteristics plate indications.

Code	N° of places	Ø Plate cm	Overall dimensions cm			Consumption per plate W	Maximum temperature °C	Weight Kg
			Height	Width	Depth			
4000634	2	19	15	47	24	800	400	8.5

## INSTALLATION

Place the equipment on a flat, horizontal, level surface, trying to leave a space of about 10 cm at the back and on each side.

For a longer life of the equipment provided with the hotplate, it is recommended to cover its surface as much as possible by the container to be heated. This will avoid overheating and deformation of the plate by absorbing more heat.

**CAUTION!!! IMPORTANT FOR YOUR SAFETY**

**Be sure that the voltage supply is the same as the one indicated on the characteristic plate of the equipment.**

**Do not use the equipment if it is NOT earthed.**

**If you have to change the plug, bear in mind the following:**


**Blue cable: Neutral.**

**Brown cable: Phase.**

**Yellow/green cable: Earth.**

## OPERATION

### CONTROL PANEL:

### STARTING UP:

1. Move the control knob of power regulator (1) from the "0" position (disconnected), to the "12" position (maximum power).
2. It will be noted that according to the control position, the ON-OFF cycle switch of the indicator lamp for the working of the heater (2) will be longer or shorter, which will permit adjusting to the desired working power.
3. To stop, move the control button to "0" position.

### **WARNING!!!**

**Should any liquid spill onto the plate, simply clean it with a damp cloth.**

### **WARNING!!!**

**Caution with the hotplate surface, serious burns can be produced, even after switching off the equipment.**

## MAINTENANCE

Before removing the cover, disconnect the equipment from the mains.

The manipulation of the internal electronic circuits of the equipment by unauthorized personnel can cause irreparable damages. Take it to one of J.P.SELECTA's authorized technical services.

### CLEANING:

For the cleaning of the different parts, use the following products:

Cleaning of stainless steel: Alcohol

Cleaning of plastic: Alcohol with a cotton duster.

## SPARE PARTS

To guarantee the safety of the equipment, the spare parts must be purchased from J.P.SELECTA, s.a.u.

Description	Code
Control knob	6001
Simerstat regulator	43021
Hotplate	32005

## WARRANTY

This product is guaranteed for one year. The guarantee does not cover damages caused by an incorrect use or causes beyond the control of J.P. SELECTA, s.a.u.

Any manipulation of the equipment by unauthorized personnel by J.P. SELECTA, s.a.u. cancels the guarantee automatically.