**AGITADOR ORBITAL Y DE VAIVÉN
con plataforma calefactora
ORBITAL MIXER AND SHAKER
with heating platform**

ROTATERM

3000435

Indice

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ENGLISH VERSION

1 Safety

The equipment comes complete with safety features.

This manual indicates areas of possible risk.

1.1 Safety Icons

Identifies risky situations and safety measures that should be adhered to.

The icons make reference within the paragraph, marked as a grey line.

Risk danger

Risk danger

Beware, follow instructions as indicated.

Electrical risk

Risk of electrical shock if accessing zones are shown with this icon.

Beware, follow instructions as indicated.

Risk of burns in high temperature areas

The temperature in the zones indicated with this icon can exceed 60°C. Use thermal protective gloves to complete described functions.

Beware, follow instructions as indicated.

Important information

- Important information in obtaining the best results or optimum performance of the equipment.
- Important information to extend the life of the equipment and maintaining its optimum performance.

1.2 Risks that the operator can be subjected to

- Possibility of trapping or strokes.
- Electrical risks.
- Possibility of touching surfaces that are over 60°C.

1.3 Qualified users

This equipment should only be used by personnel who are qualified users.

This equipment should only be used by personnel that have understood this manual or have previous adequate instruction of the use of steam sterilizers.

2 General information

- 2.1 Handle the equipment with care. Unpack and check that all items coincide with the delivery note. If you see any discrepancy or damage of any kind, notify your distributor as soon as possible.
- 2.2 Read this manual before using the equipment.
- 2.3 These instructions are an important part of the equipment and should be made available to all users.
- 2.4 If in any doubt or clarity on the use of this instrument, contact your local distributor or the service department of J.P. Selecta, s.a.
- 2.5 **ATTENTION, ANY EQUIPMENT NOT CLEAN AND DISINFECTED SHALL NOT BE REPAIRED.**
- 2.6 Any unauthorized modification, lack of maintenance or deliberate damage to the machine, transgresses the directive 89/655/CEE where the manufacturer will not be responsible for damages that could occur due to undue care and attention.
- 2.7 Do not use the equipment for fluids that can produce vapours that can form inflammable or explosive mixtures.

3 Technical specifications

Voltage supply 115-230V 50/60Hz according to the characteristics plate indications.

Code:	3000435
Oscillation:	20-230 u/min
External dimensions:	
Height	30 cm
Width	46 cm
Depth	38 cm
Platform dimensions:	
Width	41 cm
Depth	31 cm
Power:	2300 W
Weight:	32 kg

4 Contents list

4.1 Standard contents

The standard equipment consists of the following items:

- | | |
|--|---------|
| - ROTATERM (including tensile springs) | 3000435 |
| - Instruction manual | 80033 |
| - EU Declaration of Conformity | |
| - Hex socket key n°5 | |
| - Spring clips | 5900 |

5 Introduction

The **ROTATERM** shakers are designed to give a multiuses, both mixer and shaking working. Allows electrical control of the speed to 20 or 230 u/min and the oscillation amplitude from 15 to 20 mm. The equipment is able to work with 35 kg.

The heating plate allows regulate from room temperature +5°C to 250°C on its surface. Electronic temperature regulator and digital reading (3 figures) with Pt100 probe and resolution 1°C.

Top platform made from aluminium with heating elements uniformly distributed, provided by a device with tensile spring clips which allows look flasks, bottles, etc.

6 Equipment description


6.1 Equipment

The **ROTATERM** are equipped with the following features:

- Digital circuit board which controls the speed and the working time.
- Digital circuit which controls the temperature.
- Mechanical device to select orbital or shaking movement.
- Mechanical device to select oscillation movement amplitude.
- Legs with suckers which have great holding on a even surface.


6.2 Features

- Speed regulation from 20 to 230 u/min.
- Temperature regulation from room temperature +5°C to 250°C.
- Working time from 1 min to 999 min (or for no defined time).
- Mixer or shaker movement.
- Movement amplitude of 15 or 20 mm.
- Maximal load 35 kg.

6.4 Control panel

1. Main switch
2. Speed and time display
3. Speed selection indicator led
4. Time selection indicator led
5. Change set key
6. Increasing speed/time set key
7. Decreasing speed/time set key
8. Shaking start/stop key
9. Increasing temperature set key
10. Decreasing temperature set key
11. Heating start/stop key
12. Overtemperature indicator pilot lamp
13. Temperature display

6.5 Rear part

14. Pt100 probe (control)
15. Thermocouple probe (safety)
16. Fuse carrier
17. Control circuits/platform connector
18. Power supply cord

7 Installation

7.1 PLACEMENT

Place the piece of equipment near a socket, keeping in mind its amperage on a leveled surface enough resistant to support the weight indicated on the table of the section 3 (Technical features) and the oscillations produced by the working of the equipment. Do not use wheeled tables. Leave at least 30 cm at each side of the equipment.

THERE ARE RISK OF BLOWS PRODUCED BY THE PLATFORM MOVEMENT AND FINGER TRAPPINGS UNDER IT.

DO NOT INSTALL THIS EQUIPMENT IN EXPLOSIVE ATMOSPHERES OR IN PLACES WHERE THERE ARE SOLVENTS OR OTHER PRODUCTS, WHICH MAY PRODUCE EXPLOSIVE OR FLAMMABLE GAS.

7.2 POWER SUPPLY CONECTION

Make sure that the voltage given by the socket is the indicated on the feature label attached on the rear part .


Do not use the equipment unless it is correctly earthed.

If you replace the plug, keep in mind the following:

Blue wire:	Neutral
Brown wire:	Phase
Yellow/green wire:	Earth

8 Operation

8.1 STARTING UP

- Connect the equipment to the power supply.
- Place and lock with the tensile spring clips the material to be shaken on the platform.
- Turn on the main switch (1).
The display (13) shows «888» and then «-A-» while the display (2) shows «-A-» for some seconds.
Keep in mind that if the equipment stops by the switch before ending the work (unless leaving finished the set time or pressing the start/stop key), the platform will start to work at the last speed, time and temperature set.
- Normally, when the equipment is switched on the display shows «000» and the display (13) remains switched off (the plate does not heat).
- To start pressing the start/stop key (8). The work stops on pressing the start/stop key (8) or the set time finishes.
- To start the heating for a not defined time, press the start/stop key (11).
- If a working time was set, the heating also stops when it finishes.

8.1.1 SETTING THE SPEED


After doing what is explained on the section 8.1:

- Select the led (3) lighted up with the key (6).
- Press the key (5) and, keeping it pressed, with the keys (6) and (7) select the set speed which will be showed on the display (2).

8.1.2 SETTING THE WORKING TIME

After doing what is explained on the section 8.1:

- Select the led (4) lighted up with the key (7).
- Press the key (5) and, keeping it pressed, with the keys (6) and (7) select the set time which will be showed on the display (2).

8.1.3 SETTING A WORKING TEMPERATURE

After doing what is explained on the section 8.1:


- Activate the heating by pressing the key (11). The display (2) lights up.
- Select the set temperature with the key (10) and (11). The selection will be shown on the display (13) for a few seconds. Then the display will show the real temperature read by the control probe.
- The heating can work independently of the shaking except in temporized mode in which both stop at the same time.

8.2 MOVEMENT CHANGES

To access the mechanism which allows movement changes is necessary to remove the platform. To do this, do the following:


- Disconnect the equipment from the power supply.
- Unplug the connector of the rear part (see picture).
- Loosen the 4 screws of the rear part and the frontal part with an allen key size 5 mm (see picture).
- Remove the platform.

Once finished the change in the mechanism, assemble the platform again and tighten the allen screws.


Placement of the different shafts involved in the movement changes.

8.3 HOW TO CHANGE FROM ORBITAL INTO SHAKER MOVEMENT

- * Remove the platform by unscrewing the 4 Allen screws (and 4 washers) with a 5 mm Allen key by lifting it vertically.
- * Loosen and remove the nut and washer of the eccentric shaft (4) with a 10 mm key.
- * Loosen and remove the nut and dowel screw collar (5).
- * Remove the dowel O-ring gasket (1).
- * Remove P1 piece from its housing (2-4) and place it between 1 and 5, trying to keep its initial position, (swinging collar in 2 and the cavity produced because of the thickness difference between the piece and the radial bearing at the top).
- * Place a nut and washer on the dowel screw (5) and tighten slightly with a 10 mm key.
- * Place the collar and the remaining nut in the eccentric shaft (4) and tighten slightly with a 10mm key.
- * Place the O-ring gasket in the dowel (2).
- * Replace the platform and tighten the 4 Allen screws with their washers.


3.4 HOW TO CHANGE FROM SHAKER INTO ORBITAL MOVEMENT

- * Remove the platform by unscrewing the 4 Allen screws (and 4 washers) with a 5 mm Allen key by lifting it vertically.
- * Loosen and remove the nut and the dowel screw washer (5).
- * Loosen and remove the nut and the top eccentric shaft collar (4) (the lower collar must remain in the shaft 4).
- * Remove the dower O-ring gasket (2).
- * Remove P1 piece from its housing (1-5) and place it between 2-4, trying to keep its initial position (the swinging collar in 1 and the cavity produced because of the thickness difference between the piece and the radial bearing at the top).
- * Place the nut and washer on the dowel screw (5) and tighten slightly with a 10 mm key.
- * Place the O-ring gasket at the dowel (1).
- * Replace the platform and tighten the 4 Allen screws with their washers.

8.5 HOW TO CHANGE FROM LONG TO SHORT STROKE IN ORBITAL MOVEMENT

- * Remove the platform by unscrewing the 4 Allen screws (and 4 washers) with a 5 mm Allen key by lifting it vertically.
- * Loosen and remove the nut and the eccentric shaft washer with a 10 mm key.
- * Remove P1 piece from its housing 2 and 4.
- * Remove the collar from shaft 4.
- * Remove P2 piece from its housing 3 and 4.
- * With a 14 mm key loosen and remove the hexagonal extension of the shaft 4 and place it in the nearest thread hole to the centre of the shaft.
- * Place P2 piece between 3 and 4, trying to keep its initial position (the swinging collar in 3 and the cavity produced by the thickness difference between the piece and the radial bearing at the top).
- * Place the collar in 4.
- * Place P1 piece between 2 and 4, trying to keep its initial position (swinging collar in 2 and the cavity produced by the thickness difference between the piece and the radial bearing at the top).
- * Place the washer and the nut in 4 and tighten slightly with a 10 mm key.
- * Replace the platform and tighten the 4 Allen screws with their respective washers.

The cavity produced because of the thickness difference is on the top part.

8.6 HOW TO CHANGE FROM LONG TO SHORT STROKE IN SHAKER MOVEMENT:

- * Remove the platform by unscrewing the 4 Allen screws (and 4 washers) with a 5 mm Allen Key.
- * Loosen and remove the nut and collars of the eccentric shaft 4 with a 10 mm key.
- * Remove P2 piece from its housing 3 and 4.
- * With a 14 mm key loosen and remove the hexagonal extension of shaft 4 and place it in the nearest thread hole to the centre of the shaft.
- * Place P2 piece between 3 and 4, trying to keep its initial position (swinging collar in 3 and the cavity produced by the thickness difference between the piece and the radial bearing at the top).
- * Place the short collar in 4.
- * Place the long collar in 4.
- * Screw the nut on shaft 4 and tighten slightly with a 10 mm key.
- * Replace the platform and tighten the 4 Allen screws with their washers.

9 Maintenance

Before removing the cover to handle inside, unplug the equipment from the socket.

The handle of the electronic circuit boards by unskilled personnel may cause damages. Make sure to give the equipment to an authorized technical assistance.

9.1 CLEANING

To clean the different parts of the equipment, we recommend to use a cotton cloth wet in a water and alcohol solution.

9.2 SPARE PARTS

The most common spares are the following:


Spare description	Code
Bar W H6 10x295	5019
Bar W H6 10x270	5020
Bar holder	5022
Holder with bearing	5023
Bearing GE 6PB	5037
Bearing 6201.2ZR.J22C	15217
Bearing 626.2ZR	15302
Power supply cord	7001
Main switch	20100
Motor control circuit board	28328
Heating control circuit board	29359
Tachometer sensor circuit board	29013
Fuse carrier	15318
Motor reducer	24053
Mains filter	45022
Transformer	45020
Overtemperature circuit board	29161
Connecting circuit boards cord	7318
Pt100 probe	43073
Thermocouple probe	43038
Heater element	39200
Alarm pilot lamp	22004

10 Guarentee

This product is guaranteed for one year. The guarantee does not cover damages caused by an incorrect use or causes beyond the control of J.P. SELECTA, S.A.

Any manipulation of the equipment by personnel not authorized by J.P. SELECTA,S.A. automatically cancels the guarantee.

11 Electric diagram

Pos	Description	Code
1	Power supply cord	7001
2	Mains filter	45022
3	Fuse carrier	15318
4	Main switch	20100
5	Transformer	45020
6	Tachometer sensor circuit board	29013
7	Motor control circuit board	29359
8	Ferrites	45027
9	Motor reducer	24053
10	Heating control circuit board	29328
11	Overtemperature circuit board	29161
12	Circuit boards connecting cord	7318
13	Pt100 probe	43073
14	Themocouple	43038
15	Heater element	39200
16	Pilot lamp	22004

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