**AGITADOR DE NORIA /
ROTARY DISK MIXER**

3000446 - ORBIT-II

INDICE

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1. PRODUCT DESCRIPTION

1.1. Introduction

This rotary disk mixer is equipped with a long-life BLDC motor and compact size. The rotational motion provides smooth oscillation and offers a pulse function for thorough mixing. It allows you to configure up to 9 user-defined programs via the digital display, ensuring precise and easy adjustment. The removable disc plate can be adjusted to fit tubes of different diameters. The disc platter can be easily removed and replaced for cleaning.

1.2. Intended Use

This rotary disk mixer is used in laboratories to mix blood samples, viscous samples and liquid-solid suspensions.

NOTE: Before using the instrument, please read this user manual carefully. This manual is intended to assist with the operation and care of the unit only and not with its repair. For repairs, please contact the supplier.

1.3. Symbology

1.4. Packing List

- Rotary Disk Mixer
- Power Extension Cable
- Instruction Manual
- EU Declaration of Conformity

1.5. Features

- Smooth rotation motion for optimal mixing Adjustable speed from 10 to 80 RPM
- Adjustable timer from 1 minute to 99 minutes and infinite mode
- Brushless DC motor, maintenance-free, for quiet operation and long service life
- Detachable disc plate for proper cleaning and better protection against contamination
- Optional disc platters available to accommodate tubes of different diameters Pulse mode function to alternate clockwise and counterclockwise movement with 2-second pauses for better, thorough mixing
- Automatic ignition option in case of power failure; When reset, the computer automatically restarts and runs for the remaining time
- Digital display for easy setting of values Compact design optimizes the use of space in the laboratory
- Suitable for mixing blood samples and viscous substances
- Closed-loop operation for precise speed control under different load conditions
- 24 VDC Low Voltage Instrument
- Operates in up to 80% humidity and 40°C

1.6. Technical Specification

Motor Type	Brushless DC
Movement Type	Rotation
Speed Range	10..80 RPM
Speed Accuracy	+/- 1 RPM
Speed Regulation	+/- 1 RPM
Timer Range	0..99 minutes and infinite mode
Disc Angle	Flexible, 0-90°
Maximum Capacity	2kg
Removable Disc Plate	Yes
Dimensions	395x336x326mm
Protection Class	IP21
Permissible Room Temperature	5°C.. 40°C
Permissible Relative Humidity	<80%
Weight	4.5 kg
Altitude	Up to 2000 meters above sea level.
Pollution Degree	2
Outdoor / Indoor Use	Indoor use only
Power adapter	110/230Vac 50/60Hz – 24V 1.5A
Power Consumption	4.5W

2. SAFETY CONSIDERATIONS

Please read all safety and usage information provided in this manual carefully before using this device.

- Never use the instrument in any manner not specified in this manual.
- Use of the equipment in a manner not specified in this manual or by the manufacturer may result in loss of warranty.
- Do not attempt to stop the platter by locking or holding it. Do not attempt to remove the disc platter during operation.
- Never move the instrument while the equipment is in operation.
- Never use the instrument in any manner not specified in this manual.
- Repairs should only be performed by an authorized service technician.

- Use only recommended genuine spare parts for best results and product safety.
- If liquids are spilled on the disc platter, the instrument should be cleaned carefully and properly before being used again.
- Prior to mixing, the tubes should be visually inspected for damage to the material. Damaged tubes should not be used, as this can lead to sample loss and product contamination.
- The maximum capacity of 2 kg must not be exceeded.
- This product should be used only for the specified applications. It must not be operated in hazardous or flammable environments or used to mix explosive or highly reactive substances.
- Do not place potentially hazardous material within the safety area or enclosure.
- Use your personal protective equipment according to the risk category of the medium to be processed.
- There is a risk of liquid splashes.
- Place the instrument on a flat, level, solid surface. Make sure the feet of the instrument are clean and in good condition.
- Load the tubes safely and evenly.
- Always distribute samples evenly across the device.
- Do not move the instrument while it is in operation or connected to the mains.
- Make sure the knob is firmly tightened before performing any operation.
- Make sure only closed tubes are used for mixing.
- Do not lean on the equipment, you could damage the equipment or cause injury to the operator. Do not fill the tubes near the device. Spilling liquids can damage the device.
- In the event of contamination caused by aggressive agents, the device should be cleaned immediately using a neutral cleaning liquid.
- If any damage is observed, contact the service technician.

3. INSTALLATION

Gently remove the top packaging and pull out the disc rotor by holding it from the bottom. When this equipment is used for the first time, make sure that all packaging accessories have been removed from the product. Keep all packaging in a safe place for warranty purposes.


3.1. Location & Mounting

Place the device on a flat, level surface and make sure all four legs rest firmly on the surface. Avoid installing on slippery or vibration-prone surfaces.

- The ideal ambient temperature is $20^{\circ}\text{C} \pm 5^{\circ}\text{C}$; Avoid placing the device in direct sunlight.
- Keep a clearance of at least 30 cm on all sides to ensure efficient cooling.
- Keep the device away from heat sources or water to avoid problems with the sample temperature.
- Do not place the device in a place where it will be difficult to operate.


3.2. Plug in the power adaptor

- Connect one end of the power adaptor to the back of the device and the other to the power source, as shown in the figure below.
- Make sure the power switch is in the OFF position while connecting the power adaptor.


4. PART LISTING


Front Side

Back Side

5. USER INTERFACE AND DISPLAY


No.	Description
1. Display	The display shows the set speed as well as the time value. It also displays selected Program when in Prog Mode
2. Speed	Press 'Speed' to select speed value and use increment or decrement button to set require speed
3. Time	Press 'Time ' to select time value and use increment or decrement button to set require run time
4 / 5. Increment / Decrement	Press ' Increment / Decrement ' to set require value of speed, time or program (Prog Mode)
6. Start/Stop	Press " Start/Stop " to start the operation or to stop the ongoing operation.
7. Pulse & LED	Press 'Pulse' once to start pulse operation and press it again to stop operation. Pulse LED will glow to indicate Pulse is activated.
8. Prog	Hold down 'Prog" to enter in program mode and use other require buttons to set parameter for selected program.

6. REMOVING AND INSTALLING THE DISC PLATE


6.1. Removing the Disc Plate

The disc plate can be removed or changed by turning the knob counterclockwise. Turn the knob counterclockwise, and then remove the disc plate as shown in the representative image below.


6.2. Installing the Disc

Follow the reverse process exactly to insert back the disc. The following image shows the step-by-step procedure for insert the disc plate back into the unit.

7. OPERATION

7.1. Powering the device ON

After connecting the power adaptor to the back side of the device, turn on the main power supply by setting it to "ON" and then turn on the device from the rear. Make sure the tubes are properly loaded before use.

7.2. Speed adjustment

Press "Speed" to select the speed value, and then press "Increase/Decrease" to set the required speed. Operating speed ranges from 10 RPM to 80 RPM, with increments/decrements of 1 RPM. Speed values are saved after blinking 5 times.

7.3. Time adjustment

Press "Time" to select the time value and use the increment or decrease buttons to set the run time. Run time can be adjusted between 1 minute and 99 minutes, with increments/decrements of 1 minute. It can also be set to continuous or infinite mode, indicated by the " ∞ " symbol on the time display. The value will be saved after 5 flashes. The countdown timer is activated as soon as the operation begins. The time cannot be changed when the unit is in operation.

7.4. START/STOP Operation

Once the required speed and time have been set, press "Start/Stop" to START the trade. The speed gradually increases from a low RPM to the set RPM, and the time starts counting as soon as the operation starts. Press "Start/Stop" again to STOP the trade at any time. The speed gradually decreases from the set RPM to zero (0) and stops in a few seconds.

NOTE: Speed and time cannot be changed or modified during operation.

7.5. Default Pulse Mode


Pulse mode is used for vigorous mixing of samples using a rotational motion with momentary clockwise and

counterclockwise pauses every 30 seconds. The default cycle time for pulse mode is 30 seconds. Set the required speed and total operating time before activating pulse mode, so that the device starts operation in pulse mode with that specific speed and total time. Once the speed and total operating time have been set, press the "Pulse Button" to START operation in pulse mode with the selected values. In pulse mode, the pulse LED will illuminate indicating that pulse mode is activated.


Example: Speed set to 50 RPM and total run time set to 2 minutes. There will be a clockwise rotation movement for 30 seconds (cycle time). Once the 30 seconds are complete, the disc platter will have 2 cycles of momentary pause clockwise and counterclockwise. After these 2 cycles, the disc platter will rotate for 30 seconds counterclockwise. This entire process constitutes 1 pulse cycle. Press "Pulse" again to stop the pulse mode operation in progress. The user can also press "Start/Stop" to STOP the ongoing pulse mode operation.

7.6. Program Mode


The program mode operation function allows the user to save 9 programs, each with different speed and time parameters. Users can use any of these programmed parameters by simply activating the program mode. To activate program mode, press and hold "Prog" until "Prog on" appears on the screen, indicating that program mode is activated.

After a few seconds, the display will show "P1" and the speed and time values alternately. This indicates that the current program is P1, with a certain RPM and time value.


Press the "Program" button to select any program from P1 to P9. Programs can also be selected by pressing the increase or decrease buttons. Once you have selected the desired program, press the "Speed" button to select and adjust the speed value. Press "Time" to select and adjust the time value.

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
Empieza a parpadear

Empieza a parpadear

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

The speed and time value for that particular program is now saved.


Press the "Start" button to start operation in program mode. Once the configured RPM is reached, the program "Px" (where x = 1 to 9) and its values of "Speed" and "Time" will be displayed alternately.

To stop the operation, press "STOP"

To exit program mode, press and hold the "Program" button again to DISABLE program mode. "Prog off" will appear on the screen, indicating that the program mode is disabled.

NOTE: On each restart, the device automatically exits program mode. Displays the last speed and time defaults used.

7.7. Power Failure Mode

The function of this feature is that, in the event of a power outage, the equipment saves the SPEED, TIME and other details configured by the user. When power is restored, the equipment continues to operate for the remaining time, with a cycle time variation of ± 5 minutes. The equipment saves the operating time every 5 minutes of operation.

To select this feature, press and hold the "TIME" button. While this function is not activated, the disc rotor will not resume operation after a power failure.

Example: The user has set 45 RPM and a time of 35 minutes. After 23 minutes of operation, a power outage occurs and operation stops. Just before the cut, the remaining time displayed on the screen was 23 minutes. When power is restored, the equipment resumes operation after 25 minutes, with the RPM set and the other operating values.

8. REMOVING DISC TUBES FROM DISK PLATE

To remove the tubes from the disc platter, do not attempt to turn the platter with your bare hands. The procedure for removing the tubes attached to the bottom of the disc platter is as follows:

Press the PULSE button once to rotate the platter so that the tubes at the bottom of the platter are at the top.

Next, remove the tubes from the disc platter clamps one by one.

Similarly, follow the process described above to remove the tubes, in order to avoid damage to the disc platter and tubes.

9. MAINTENANCE AND CLEANING

Warning: Make sure the unit is disconnected from the power source before attempting any maintenance or cleaning.

Clean the unit with a soft, damp cloth or mild non-corrosive detergent ($\text{pH} \leq 8$) and do not directly spray the unit.

Wear clean gloves during cleaning and make sure the unit is completely dry before putting it into operation.

Cleanup is made easier if spills are taken care of immediately.

In any case, acid and alkali spills should be removed immediately, as these chemicals can attack and damage the finish of the case. Make sure proper safety precautions are observed.

10. WARRANTY STATEMENT

This product is warranted to be free from defects in materials and workmanship. Your product will be duly repaired upon notification, in accordance with the conditions of the warranty.

Conditions:

This warranty is valid only if the product is used for its intended purpose and within the guidelines specified in this instruction manual. This warranty does not cover damage caused by accident, neglect, misuse, improper service, acts of God, or other causes not resulting from defects in materials or original workmanship. This warranty does not cover incidental or consequential damages, business losses, or any other damages arising from the use of this product.

The warranty will be void for any modification not made at the factory, which will immediately terminate all liability on our part with respect to the products or damage caused by their use. The customer will be responsible for the product or the use of the product, as well as for any supervision necessary for safety. If requested, the products must be returned to the distributor properly packed and insured, and all shipping costs must be paid.

Products received without proper authorization will not be processed for warranty or technical service. All items returned for service must be shipped postage-prepaid, in the original packaging or other suitable packaging, padded to prevent damage. We are not responsible for damage caused by improper packaging.

11. PRODUCT DISPOSAL

Notice to customers about recycling:

The product is composed of various components and various materials that must be recycled or, failing that, deposited in the corresponding disposal sites. debris when the life of the product has been completed or when, in the event of a otherwise, it must be discarded. To do this, the fine user who acquires the product must know the current regulations of each municipality and/or locality in function of electrical and electronic waste. The user who acquires this product Pipeline must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if it is a prior dismantling and knowledge of the components incorporated. If you do not know the procedure to follow, consult with the city council for more information.