**3000980 AGITADOR DE PLAQUETAS  
PLATELET SHAKER**

**APL-48  
APL-48**

**3000981 AGITADOR DE PLAQUETAS  
PLATELET STIRRER**

**APL-96  
APL-96**

# Indice

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**1. General Information**

- Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributor immediately.
- Do not install or use the equipment without reading this handbook first. This instruction manual must be attached to the equipment and must be available for all users.
- If you have any doubts or enquiries, contact with J.P. SELECTA's technical service.
- Do not use the equipment in explosive, flammable or corrosive atmospheres.


**2. Packing List**

Model	<b>3000980</b>	3000981
Stirrer with support trays	1	1
Trays	8	16
Instructions manual	80270	80270

**3. Technical Specifications**

Useful measures: height/width/depth (cm):	385/445/390	385/840/390
Displacement length (mm):	40	40
Voltage:	See characteristic plate	
Electric power (W):	50 W	50 W
Weight (Kg):	22	39
Sound level (dB):	20	21
Loading and unloading programmable time:	99s	
Loading and unloading break counter:	99	
Fail counter by supply failure:	99	
Date and time indication:	hh/mm/ss dd/mm/yy	
Parameters protection with password:	Fixed password 1111	

## 4. Installation and Start-up

### 4.1 Packaging

Once the equipment is unpacked, visually check that it is in good conditions. Keep the package for some time just in case you will need it.

### 4.2 Levelling

Level the equipment if necessary, by using the four adjustable level rubber stops on the stirrer base.

Please bear in mind the equipment weight already loaded, to determine the weight the table must support.

### 4.3 Connections to the mains

Plug the equipment to a near socket suitable for the equipment power. Use the connection cable supplied or a similar one.

For your safety, the equipment must be earthed.

Before plugging the equipment, check that the voltage supply and the frequency corresponds to the one indicated in the equipment characteristics plate.

If you change the plug, bear in mind the following:

Blue cable:	Neutral.
Brown cable:	Phase.
Yellow/Green cable:	Earth.

The audible alarm of power failure is activated during the first connection to the mains.

For a correct operation, the equipment connection voltage must be  $\pm 10\%$  of nominal voltage. To avoid problems and failures, it is important that the equipment works within these margins.

### 4.4 Parameters selection rotatory button

This equipment has an only button in the control panel, which has a couple of movements: circular and rotatory.

## 5. Maintenance

- Before cleaning, disconnect the computer from the mains.
- Do not spill liquids on the equipment.
- The equipment external parts can be cleaned with a little alcohol and water solution.

## 6. Warranty

This product is guaranteed for one year. The guarantee does not cover damages caused by an incorrect use or causes beyond the control of J. P. SELECTA, s.a.u

Any manipulation of the equipment by personnel not authorized by J.P. SELECTA, s.a.u. automatically cancels the warranty.

## **7. Operation**

The platelet stirrer has been designed for a quiet and continuous blood bag operation.

To enable a continuous operation monitoring, the equipment records all break times, both due to loading and unloading, as unforeseen ones, due to power failure.

### **7.1 Starting-up**

Start the equipment by pressing the main switch.

By a (short) press of the button, it stops for loading and unloading.

The stop lasts for a determined «T» time. After this time, the stirrer starts automatically.

With a long press you can access the record menu. See chapter 8.

### **7.2 Break counter**

One can access the break counter by the single button (chapter 8).

This counter can be reset to zero.

The break lasts for a determined «T» time. After this time, the stirrer starts automatically.

### **7.3 Power supply fails**

The computer registers the power failures.

For every stop, the computer registers date, time and duration of the stop.

Registration of stops can be reset.

### **7.4 Password**

Access to the above parameters is password protected.

Password is 1111.


Password can be momentarily unlocked, to facilitate access to counters and registers before initiating treatment.

### **7.5 Settings**

The agitation speed is factory set. If it is necessary to modify this speed, remove the blue rear plug and adjust the speed by using a flat screwdriver.

## 8. Menu Diagram

**Notice to customers:**

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.

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