7000385

**Agitador de varilla SE-100
SE-100 rod stirrer**

Indice

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1 Safety

The equipment includes appropriate security measures.

Risky situations that must be respected are indicated along this manual.

1.1 Safety icons

The safety icons identify the risky situations and the safety measures that must be taken.

Risk of danger

Electrical risk

Respect the indicated instructions to perform the operation.

Electrical risks

Risk of electrical accident when accessing the areas indicated by this sign or when performing the operations indicated by this icon in this manual.

Respect the indicated instructions to perform the operation.

Important information

- Important information for good results or for an optimum equipment operation.
- Important information to extend life of the equipment or to prevent degradation of any of the components.

1.2 Personnel qualification

This equipment can only be used by personnel who have been properly qualified.

This equipment can only be used by personnel who have read and understood these instructions or who has been properly qualified in the operation of this equipment.

2 General information

Handle the package with care. Unpack and check that the contents coincides with the "Packing List". If any part is damaged or missing, please advise the distributor immediately.

Do not install or use the equipment without previously reading this instructions manual. They are an inseparable part of the equipment and should be available to all users of the machine. Any questions can be clarified by contacting J.P. SELECTA's technical service.

Any modification, deletion or lack of maintenance of the device could cause damages. The manufacturer is not liable for the damages that may arise.

Do not use the equipment with fluids which can give off vapours or form explosive or flammable mixtures.

Maintenance

Manipulation of electronic circuits by unauthorized persons may cause damages difficult to repair. Check with J.P. SELECTA's technical service.

To clean the accessible parts, use a cloth stained with alcohol cleaner or a similar product. Do not use aggressive solvents.

Warranty

This product is guaranteed for one year. The warranty does not cover damages caused by an improper use or causes beyond J. P. SELECTA, s.a.u.

Any manipulation of the equipment by an unauthorized personnel cancels the warranty benefits.

3 Technical specifications

Voltage supply 230V 50 Hz as indicated in the machine's features plate.

Code:	7000385
External dimensions (H x W x D) mm	175 x 75 x 135
Power	100W
Speed range	60 ... 2000 rpm
Stirring maximum viscosity (mPas)	100.000 cPs
Weight	2.8 Kg
Rod diameter (mm Ø)	Up to 9.5 (through axis)
Maximum torque	450Ncm
Rear clamping bar	Ø16mm x 185mm

4 Packing list

4.1 Standard package content

The standard equipment consists of the following components:

- Stirrer	700385
- Instructions manual	80320

5 Introduction


Stirrer 7000385 is a desktop rod stirrer for an accurate and reproducible stirring.

6 Installation and operation

- Place the equipment near a power socket and plug it.
- Place the equipment over a flat level surface, strong enough to support the equipment weight.
- Select the appropriate accessories (see section 7).
- Activate the main switch ON/OFF.
- Select the appropriate speed agitation.

CAUTION:

It is advisable to start the agitation with the speed knob at minimum to avoid splashing.

Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.

7 Accessories

7.1 Base supports

DOUBLE BASE SUPPORT T "MAXI"
Rod of AISI 304 stainless steel of 20 Ø x 800 mm high. Designed to support weights up to 20 Kg at 300 mm height from the base. Useful for stirrers.
Weight: 7 Kg.
Part No. **6000410**

DOUBLE BASE SUPPORT T TELESCOPE
Stainless steel base AISI 304 of 32 mm Ø. Telescopic elevation from 725 to 1025 mm. Maximum support weight 20 Kg. Actual weight: 7.7 Kg. Compatible with Part No. 5057022 (See page 166).
Part No. **5057011**

7.2 Bossheads


Bosshead "MAXIM"
For large diameters. Chrome plated with one opening to hold rods up to 30 mm Ø. Second opening up to 20 mm Ø.
Part No. **7000408**

BOSSHEAD "SUNTECH"
For retort bars of 5 to 22 mm Ø. With screw clamp joints. Ideal for supporting overhead stirrers.
Part No. **5900600**

7.3 Rotors

STIRRING SHAFTS WITH ROTORS
Shaft with AISI 304 stainless steel rotor. Shaft: 7 mm Ø.

Part No.	Rotor type	Ø shaft mm	Total length mm	Features
6002891	1	90	400	Double fly vane
6002892	2	60	400	Fly vane
6002893	3	50	400	Fixed blade
6002894	4	50	400	Double cross vane
6002895	5	58	400	Helix
6002896	6	70	450	Six hole paddle
6002897	7	70	450	Three hole paddle
5050921	8	46	400	Radial turbine
5050932	9	45	400	Turbo helix 45 mm Ø
5050933	10	66	400	Turbo helix 66 mm Ø

SHAFTS WITH PTFE PADDLES
Shaft of 8 mm Ø.

Part No.	Rotor type	Ø shaft mm	Total Length mm	Features	Speed range r.p.m.
5050934	11	50	500	Helix PTFE	500
5050935	12	80	500	Anchor PTFE	500

ADJUSTABLE HEIGHT SHAFT STIRRERS
Shafts made of AISI 304 stainless steel coated in PTFE and PTFE stirrer paddles. Shaft: 8 mm Ø. Adjustable height paddles using a locking thread. Can be used with differing paddles on the same shaft to create a multiple mixing action.

Part No.	Rotor type	Ø shaft mm	Total Length mm	Features
5050936	13	70	400	Shaft PTFE
5050937	14	40	400	Helix PTFE
5050938	Shaft PTFE Spare or additional paddle for multiple stir effect			
5050939	Helix PTFE Spare or additional paddle for multiple stir effect			