DOSI-FIBER 4 PLAZAS / 4 UNITS
DOSI-FIBER 6 PLAZAS / 6 UNITS

4000599
4000623

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GENERAL INFORMATION

- 1) *Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributor immediately.*
- 2) *Do not install or use the equipment without reading this handbook before.*
- 3) *This handbook must always be attached to the equipment and it must be available for all users.*
- 4) *If you have any doubts or enquiries, please contact with your supplier or J.P. Selecta's technical service.*
- 5) **IMPORTANT! J.P. SELECTA WILL NOT ACCEPT ANY EQUIPMENT TO BE REPAIRED IF IT IS NOT DULY CLEANED.**
- 6) *If any modification, elimination or lacking in maintenance of any device of the equipment by the user transgress the directive 89/655/CEE, the manufacturer is not responsible for the damage that can occur.*
- 7) *Do not use the equipment with liquids which can give off vapours capable of making explosive mixtures.*

PACKING LIST

The standard equipment consists of the following components:

Piece	Code	
Dosi-Fiber 4 units / Dosi-Fiber 6 units	4000599	4000623
Main unit (4 or 6 units)	1	1
Crucible rack (4 or 6 units)	1	1
Crucible handle (4 or 6 units)	1	1
P-2 porositu crucible	4	6
PVC 3m hose of 8mm inner Ø	1	1
PVC 2m hose of 5mm inner Ø	1	1
3/8' racord	2	2
1/4' racord	2	2
Connection cable	1	1
Water Jet Pump code 7000293	1	1
Instruction manual (80055)	1	1

ACCESSORIES

Description	Code
Double hotplate for reagents	4000634
Glass jar for reagents	1000563
Rack for 4 crucibles	4000600
Rack for 6 crucibles	4000624
P-2 porosity crucible	4000601

TECHNICAL FEATURES

Voltage supply 115-230V 50/60 Hz according to the characteristics plate indications.

Code	4000599	4000623	
Consumption	700	1000	
Nr. units	4	6	
Weight Kg	18.5	25	
Overall dimensions cm	Height	56	56
	Length	43	57
	Width	32	32
Cooling water consumption	1l./minute		

EQUIPMENT DESCRIPTION

Extraction system for the determination of crude fibre.


It is possible to do extractions using: Nitric, Acetic or Trichlore-acetic acids.

INSTALLATION

Place the equipment on a flat, horizontal and level surface, taking care of leaving a space of about 10 cm. at the rear part and both sides.

- 1- Connect the cooler water inlet «1» (racord 3/8') with the 8 mm inner dia. hose to a water tap with a minimum flow of 1 litre/minute.
- 2- Connect the cooler water outlet «2» (racord 3/8') with the 8 mm inner dia. hose to a drainage.
- 3- Connect the vacuum inlet «3» (racord 1/4') to a water jet pump with the 5 mm inner dia. tube.
- 4- Connect the drainage residual outlet of the tray «4» (racord 1/4') to a drainage.
- 5- Connect it into the mains.

CAUTION!!! IMPORTANT FOR YOUR SAFETY

Be sure that the voltage supply is the same that the on indicated on the characteristics plate.

Do not use the equipment if it is not earthed.

If you change the plug bear the following in mind:

Blue cable: Neutral
Brown cable: Phase
Yellow/green cable: earth

IMPORTANT NOTE FOR THE CRUCIBLE HANDLING

The porous crucibles are made of borosilicate glass, with a standard filter of porosity no.2 (40-70 micras).

A cold crucible must not be put directly into a muffle furnace at 500°C since the crucible would be immediately broken. To avoid it, please follow one of the two options:

- Place the crucible in the cold furnace and then switch it on.

- Place the crucible on an thermal insulating material and put altogether in the hot furnace.

Crucibles must NOT be put on a cold surface after taking them from the oven. They must be cooled slowly.

It is recommended to handle the crucible with the help of the handle for crucible.

OPERATION

Sample and reagents preparation:

1. Grind the sample to a thick size (normally 1 mm. sieve).
2. Heat the reagent on a hotplate up to 95/100 °C (accessory 4000634 or similar).
- 3- Fill the crucible with grinded samples and put them on the crucible rack (4). The rack can be fixed to the front of the main unit. Use the handle for crucibles to take them and insert them in the main unit facing the heating elements (6). Lower the fixing lever (2) and reflector.
- 4- Put all valves (5) in «OFF» position.
- 5- Open the cooler water inlet tap. Flow to 1 or 2 litres/minute.
- 6- Turn on the main switch (POWER) (9), the amber pilot lamp comes on. Put the potentiometer (7) in «OFF» position.

Hot extraction procedure:

- 7- Lift the top lid (1) and add the reagent in each column (10). To the necessary quantity, (observe the quantity by the graded column).
- 8- Turn the potentiometer (7) clockwise up to a position of 80 or 90%. The heating element starts up.
- 9- Add foam-suppressor in each column.
- 10- When the reagent begins to boil, lower the heater by turning the potentiometer (7) unclockwise up to 20-30%.
- 11- During the extraction, the second reagent or distilled water can be heated.
- 12- At the end of the extraction, turn off the heater switch (9).
- 13- Open the tap of the water jet pump (if this system has been used to produce the vacuum pressure). Place the valves (5) in «Aspirar» position. When the filtration is finished, close the valves.
- 14- If during the filtration, it is necessary to dissolve residue, turn on the air pump switch (8) (PRESSURE) and place the valves in «soplar» position, returning it to «aspirar» position afterwards. The power of the blow pump is adjustable inside. See «MAINTENANCE» section.

DO NOT STOP THE PUMP (PRESSURE) WITH THE VALVE IN «SOPLAR» POSITION.

- 15- Wash the samples with warm distilled water using the spray. The water is introduced in each column inlet. Place the valves in «aspirar» position to dry the sample. Close the valves again. If the method requires more extractions repeat the whole cycle.
- 16- To remove the crucibles from the extraction unit use the crucible handle. Insert it in the crucibles and immediately free them using the lever on the left.
- 17- Transfer them to the rack «4».

MAINTENANCE

Before removing the cover disconnect the equipment from the mains.

The manipulation of the internal electronic circuits by unauthorized personnel can cause irreparable damage. Take it to one of the J.P.SELECTA s.a.u. authorized technical services.

- 1- Remove the rear lid of the main unit and find the «soplar» pump (See page 9, spare part 2)
- 2- Adjust the pump power throughout the grooved button:
 - Clockwise: more power.
 - Unclockwise: less power.
- 3- This pump is adjusted at maximum power in factory.

GUARANTEE

This product is guaranteed for one year. The guarantee does not cover damage caused by incorrect use or causes beyond the control of J.P. SELECTA, s.a.u.

Any manipulation of the equipment by personnel not authorized by J.P. SELECTA, s.a.u. automatically cancels the guarantee.

APPLICATION NOTES

1- CRUDE FIBER DETERMINATION BY WEENDE METHOD.

2- CRUDE FIBER DETERMINATION BY WIJKSTRON METHOD.

3- FIBER COMPOSITION DETERMINATION BY VAN SOEST METHOD.

- NEUTRAL DETERGENT FIBER (N.D.F)

- ACID DETERGENT FIBER (A.D.F)

- FRACTIONAL DETERMINATION:

- HEMICELLULOSE

- CELLULOSE

- LIGNINE

- SILICE

1- CRUDE FIBER DETERMINATION BY WEENDE METHOD.

Required reagents:

*Sulphuric acid H_2SO_4 0.128M
(7.1 ml 96% to 1 l with distilled water)*

*Potassium hydroxide KOH 0.223M
(12.5 g to 1 l with distilled water)*

Foam-suppressor such as octanol

Acetone

Required equipment:

Dosi-Fiber equipment

Scale ± 0.1 mg

Vacuum pump

«Kitasatos» flask

Porous crucibles

Muffle furnace 500°C

Oven at 150°C

Vacuum dryer

Procedure:

- 1- *Weight and transfer (precision $\pm 0.1\text{mg}$) of 1 to 1.5 g sample into a porous crucible. The quantity of sample is W_0 .*
- 2- *Place the crucible into the Dosi-fiber unit.*

Hot acid hydrolysis

- 3- *Ensure that the valve are in «OFF» position.*
- 4- *Add to each column 100-150 ml of preheated sulphuric acid and some drops of foam-suppressor.*
- 5- *Open the cooling circuit and turn on the heating elements (power at 90%).*
- 6- *Wait until it boils, reduce the power at 30% leave boiling during the extraction time (30 min. to 1 hour depending on the material). For more effective hydrolysis, put air pump in «soplar» position.*
7. *To stop heat. Open the vacuum circuit and place the valves in «ABSORPTION» . Fill with distilled water with the spray in each column. Repeat this process 3 times.*

Hot basic hydrolysis:

- 8- *Repeat the points 3 to 7 by using KOH instead of H_2SO_4 .*

Cold extraction with Acetone.

Do not make the cold extraction with acetone in the Dosi-fiber equipment.











- 9- *Prepare the «Kitasatos» flask with the vacuum pump. Place the crucible in the Kitasatos inlet and add acetone at the same time that the vacuum circuit is absorbing towards the flask. Repeat it 3 times.*
- 10- *Put the samples to dry in an oven at 150°C during 1 hour.*
- 11 - *Leave them to cool inside the vacuum dryer.*
- 12 - *Weight the samples with a precision of $\pm 0.1\text{ mg}$. The weighed quantity is W_1 .*
- 13- *Ash crucible samples in a muffle furnace at 500°C during a minimum of 3 hours.*
- 14- *Leave them inside the vacuum dryer to cool. Take note of the recommendations given for the crucible handling.*
- 15- *Weigh the crucibles with a precision of $\pm 0.1\text{ mg}$. The weighed quantity is W_2 .*
- 16- *Carry out the following calculation:*

CRUDE FIBRE DETERMINATION. Weende method.

J.P. SELECTA, S.a.

Reactive A: H_2SO_4 0.128M
Reactive B: KOH 0.223M

Repeat with reagent «B»


$$\% \text{ Crude fiber} = \frac{W_1 - W_2}{W_0} \times 100$$

David Pecanins

$$\%Fiber = \frac{W_1 - W_2}{W_0} \times 100$$

2- CRUDE FIBER DETERMINATION BY WIJKSTRON METHOD.

Required reagents:

Sulphuric acid H_2SO_4 0.32M
(17.76 ml 96% to 1 liter with distilled water)

Potassium hydroxide KOH 0.556M
(31.02 g to 1 liter with distilled water)

Foam-suppressor as octanol

Acetone

Required equipment:

Dosi-Fiber equipment

Scale ± 0.1 mg

Vacuum pump

«Kitasatos» flask

Porous crucibles

Muffle furnace 500°C

Oven at 150°C

Vacuum dryer

Procedure:

- 1- Weigh and transfer (precision $\pm 0.1\text{mg}$) of 1 to 1.5 g sample into a porous crucible. The quantity of sample is W_0 .
- 2- Place the crucible into the Dosi-fiber unit.

Hot acid hydrolysis:

- 3- Ensure that the valves are in «OFF» position.
- 4- Add to each column 100-150 ml of preheated sulphuric acid and some drops of foam-suppressor.
- 5- Open the cooling circuit and turn on the heating elements (power at 90%)
- 6- Wait until it boils, reduce the power at 30% leave boiling during the extraction time (10 min.). For more effective hydrolysis, put air pump in «soplar» position.
7. To stop heat. Open the vacuum circuit and place the valves in «ABSORPTION» . Fill with distilled water with the spray in each column. Repeat this process 3 times.

Hot basic hydrolysis:

- 8- Repeat the points 3 to 7 by using KOH instead of H_2SO_4 .

Cold extraction with Acetone**Do not make the cold extraction with acetone in the Dosi-fiber equipment**


- 9- Prepare the «Kitasatos» flask with the vacuum pump. Place the crucible in the Kitasatos inlet and add acetone at the same time that the vacuum circuit is absorbing towards the flask. Repeat it 3 times.
- 10- Put the samples to dry in an oven at 150°C during 1 hour.
- 11 - Leave them to cool inside the vacuum dryer.
- 12 - Weight the samples with a precision of $\pm 0.1\text{ mg}$. The weighed quantity is W_1 .
- 13- Ash crucible samples in a muffle furnace at 500°C during a minimum of 3 hours.
- 14- Leave them inside the vacuum dryer to cool. Take note of the recommendations given for the crucible handling.
- 15- Weight the crucibles with a precision of $\pm 0.1\text{ mg}$. The weighed quantity is W_2 .
- 16- Carry out the following calculation:


CRUDE FIBRE DETERMINATION. Wijkstrom method.

J.P. SELECTA, S.a.


Reagent A: H₂SO₄ 0.32M
Reagent B: KOH 0.556M


Triture and homogenize sample
Weigh from 1 to 1.5g of sample (precisión ±1mg)
W₀


Put crucibles on the grid


Put crucibles on Dosi-Fiber

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
Dry the samples in an oven at 100°C for 1 h.
Allow to cool in a desiccator.


Place a crucible in the kitasato. Fill with acetone and filter it.
Repeat 3 times.

Weigh the crucibles
W₁


Ignite in muffle furnace at 500°C for 3 h.
Allow to cool in desiccator.


Weigh the crucibles
W₂


Repeat with reagent <>

Start heating.
Add 120ml of reagent on each column.
Wait for boiling.
Add some drops of antifoam.

Keep boiling for 10m.


Stop heating.
Start the vacuum tube.
Spray distilled water to drag the particles adhered to the glass.
Repeat 3 times.

$$\% \text{ Crude fibre} = \frac{W_1 - W_2}{W_0} \times 100$$

$$\% \text{Fiber} = \frac{W_1 - W_2}{W_0} \times 100$$

3- FIBER COMPOSITION DETERMINATION BY VAN SOEST METHOD.

- NEUTRAL DETERGENT FIBER (N.D.F):

The neutral detergent fiber is the insoluble residue in neutral detergent (sodium sulphate lauril 3% pH 7). It includes Hemicellulose, Lignine, Cellulose, Protein and Silice. Required reagents: -Neutral detergent solution (N.D.S): Make a solution containing 30 g of sodium sulphate lauril in 500 ml of distilled water and add 10 ml of 2-Etoxiethanol P.A.. Dissolve 18.61 g of EDTA $\text{Na}_2 \cdot \text{H}_2\text{O}$ and 6.81 g of cahidrate sodium tetraborate P.A. in 250 ml of hot distilled water. Add this solution to the above. Dissolve 4.56 g of disodium phosphate anhydre P.A. in 250 ml of hot distilled water. Add this solution to the above wich contains all the other reagents. The final pH must be between 6.9 and 7.1, if it is not, adjust with NaOH or HCl solution.

- Foam-suppressor such octanol

- Acetone

Required equipment:

Dosi-Fiber equipment
 Scale $\pm 0.1 \text{ mg}$
 Vacuum pump
 «Kitasatos» flask
 Porous crucibles
 Muffle furnace 500°C
 Oven at 150°C
 Vacuum dryer

Procedure:

- 1- *Weight and transfer (precision ± 0.1 mg) of 1 to 1.5 g sample into a porous crucible. The quantity of sample is W_0 .*
- 2- *Place the crucible into the Dosi-fiber unit.*
- 3- *Ensure that the valves are in «OFF» position.*
- 4- *Add 100 ml of neutral detergent solution (N.D.S.) to each column.*
- 5- *Open the cooling circuit and turn on the heating elements (power at 90%)*
- 6- *Wait until it boils, reduce the power by 30% and leave boiling during the extraction time (60 min.).*
7. *To stop heat. Open the vacuum circuit and place the valves in «ABSORPTION» . Fill with distilled water with the spray in each column. Repeat this process 3 times.*

Cold extraction with Acetone.

Do not make the cold extraction with acetone in the Dosi-fiber equipment

- 8- *Prepare the «Kitasatos» flask with the vacuum pump. Place the crucible in the Kitasatos inlet and add acetone at the same time that the vacuum circuit is absorbing toward the flask (like in the 1 and 2 methods). Repeat it 2 times.*
- 9- *Put the samples to dry in an oven at 150°C during 1 hour.*
- 10- *Leave them to cool inside the vacuum dryer.*
- 11- *Weight the samples with a precision of ± 0.1 mg. The weighted quantity is W_3 .*
- 12- *Ash crucible samples in a muffle furnace at 500°C during a minimum of 3 hours.*
- 13- *Leave them inside the vacuum dryer to cool. Take note of the recommendations given for the crucible handling.*
- 14- *Weight the crucibles with a precision of ± 0.1 mg. The quantity weighted is W_4 .*
- 15- *Carry out the following calculation:*

$$\% \text{ F.N.D.} = \frac{W_3 - W_4}{W_0} \times 100$$

**- ACID DETERGENT FIBER (A.D.F):**

The acid detergent fiber is the insoluble residue in acid detergent (cethyl trimethyl ammonium bromure). It includes lignine, cellulose, and silice.

Required reagents:

- Acid detergent solution (A.D.S): To add 20 g of cethyl trimethyl ammonium bromure to 1 l of sulphuric acid 1N.
- Acetone.
- Foam-suppressor such as octanol.

Required equipment:

Dosi-Fiber equipment.
Scale ± 0.1 mg.
Vacuum pump.
«Kitasatos» flask.
Porous crucibles.
Muffle furnace 500°C.
Oven at 150°C.
Vacuum dryer.

Procedure:

- 1- Weight and transfer (precision ± 0.1 mg) of 1 to 1.5 g sample into a porous crucible. The quantity of sample is W_0 .
- 2- Place the crucible into the Dosi-fiber unit.
- 3- Ensure that the valve are in «OFF» position.
- 4- Add to each column 100 ml of acid detergent solution (A.D.S).
- 5- Open the cooling circuit and turn on the heating elements (power at 90%)
- 6- Wait until it boils, reduce the power by 30% and leave boiling during the extraction time (60 min.).
7. To stop heat. Open the vacuum circuit and place the valves in «ABSORPTION». Fill with distilled water with the spray in each column. Repeat this process 3 times.

Cold extraction with Acetone

Do not make the cold extraction with acetone in the Dosi-fiber equipment

- 8- Prepare the «Kitasatos» flask with the vacuum pump. Place the crucible in the Kitasatos inlet and add acetone at the same time that the vacuum circuit is absorbing toward the flask (like in the 1 and 2 methods). Repeat it 2 times.

- 9- Put the samples to dry in an oven at 150°C during 1 hour.
- 10- Leave them to cool inside the vacuum dryer.
- 11- Weight the samples with a precision of ± 0.1 mg. The weighed quantity is W_5 .
- 12- Ash crucible samples in a muffle furnace at 500 °C during a minimum of 3 hours.
- 13- Leave them inside the vacuum dryer to cool. Take note of the recommendations given for the crucible handling.
- 14- Weight the crucibles with a precision of ± 0.1 mg. The weighed quantity is W_6 .
- 15- Carry out the following calculation:

$$\% \text{ F.A.D.} = \frac{W_5 - W_6}{W_0} \times 100$$

- FRACTIONALLY DETERMINATION:

Required reagents:

- Acid detergent solution (A.D.S): To add 20 g of cethyl trimethyl ammonium bromure to 1 l of sulphuric acid 1N.
- HBr 48% solution.
- Sulphuric acid solution 72% in weight.
- Foam-suppressor such as octanol.
- Acetone.

Required equipment:

Dosi-Fiber equipment.
Scale ± 0.1 mg.
Vacuum pump.
«Kitasatos» flask.
Porous crucibles.
Muffle furnace 500°C.
Oven at 150°C.
Vacuum dryer.

Procedure:

- 1- Determinate the A.D.F until the residue W_5 is obtained, as in the previous procedure.
- 2- Place the crucible with the residue W_5 into the Dosi-fiber unit.

- 3- *Ensure that the valves are in «OFF» position.*
 - 4- *Add 25 ml of sulphuric acid solution 72%.*
 - 5- *Put switch of air pump in “Soplar” position. Leave the extraction during 3 hours.*
 - 6- *Open the vacuum circuit and place the valves in «ABSORPTION». Fill with distilled water with the spray in each column. Repeat this process 3 times.*
 - 7- *Put the samples to dry in an oven at 150°C during 1 hour.*
 - 8- *Leave them to cool inside the vacuum dryer.*
 - 9- *Weight the samples with a precision of ± 0.1 mg. The weighed quantity is W_7 .*
 - 10- *Ash crucible samples in a muffle furnace at 500°C during a minimum of 3 hours.*
 - 11- *Leave them inside the vacuum dryer to cool. Take note of the recommendations given for the crucible handling.*
 - 12- *Weight the crucibles with a precision of ± 0.1 mg. The weighed quantity is W_8 .*
- W_7 = Lignine and minerals
 W_8 = Minerals= Silice and acid detergent ashes.
- 13- *Place the crucible with the residue W_8 into the Dosi-fiber unit.*
 - 14- *Ensure that the valves are in «OFF» position.*
 - 15- *Add 4 ml of HBr 48% solution.*
 - 16- *Put switch of air pump in “Soplar” position. Leave the extraction during 1.5 hours.*
 - 17- *Open the vacuum circuit and place the valves in «ABSORPTION». Fill with distilled water with the spray in each column. Repeat this process 3 times.*

Cold extraction with Acetone.

Do not make the cold extraction with acetone in the Dosi-fiber equipment.

- 18- *Prepare the «Kitasatos» flask with the vacuum pump. Place the crucible in the Kitasatos inlet and add acetone at the same time that the vacuum circuit is absorbing toward the flask (as in the 1 and 2 methods). Repeat it 2 times.*

19- Ash crucible samples in a muffle furnace at 500 °C during minimum 3 hours.

20- Leave them inside the vacuum dryer to cool. Take note of the recommendations given for the crucible handling.

21- Weight the crucibles with a precision of ± 0.1 mg. The weighted quantity is W_9 .

$$W_9 = \text{Silice}$$


22- Carry out the following calculation:

$$\% \text{Hemicellulose} = \frac{W_3 - W_5}{W_0} \times 100$$

$$\% \text{Cellulose} = \frac{W_5 - W_7}{W_0} \times 100$$

$$\% \text{Lignine} = \frac{W_7 - W_8}{W_0} \times 100$$

$$\% \text{Silice} = \frac{W_9}{W_0} \times 100$$

Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.