EQUIPOS DE DIGESTIÓN

Digestion systems

DQO 6		4000638
DQO 12		4000639
DQO 20		4000640
Bloc-Digest Macro	6	4000629
Bloc-Digest Macro	12	4000630
Bloc-Digest Macro	20	4000631
Bloc-Digest micro	12	4001047
Bloc-Digest micro	24	4001048
Bloc-Digest micro	40	4001049
RAT-2		4001538

Indice

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1. Information and warranty

- Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributors immediately.
- Do not install or use the equipment without reading this handbook first. This instruction manual must be attached to the equipment and must be available for all users.
- If you have any doubts or enquiries, contact with J.P. SELECTA's technical service.
- Do not use the equipment in explosive, flammable or corrosive atmospheres.

1.1 Warranty

This product is guaranteed for one year. Warranty does not cover damages by an improper use or causes beyond the control of J. P. SELECTA, s.a.u.

Any manipulation of the equipment by personnel not authorized by J.P. SELECTA, s.a.u. automatically cancels the guarantee.

2. Packing List

Packing list contents for each product:

4000629 BLOC DIGEST 6x tube Ø42

Heater block 6p	4000507
Temperature profiler:	4001538
Tube rack	4005071
Digestion macro tubes	4042300 6x
Fume collector	4005072
Control wire	0403507
Power supply cord.	0007001

4000630 BLOC DIGEST 12x tube Ø42

Heater block 12p	4000508
Temperature profiler:	4001538
Tube rack	4005081
Digestion macro tubes	4042300 12x
Fume collector	4005082
Control wire	0403507
Power supply cord.	0007001

4000631 BLOC DIGEST 20x tube Ø42

Heater block 20p	4000509
Temperature profiler:	4001538
Tube rack	4005091
Digestion macro tubes	4042300 20x
Fume collector	4005092
Control wire	0403507
Power supply cord.	0007001

4001047 BLOC DIGEST 12x micro tube Ø26.

Heater block 12p	4001050
Temperature profiler	4001538
Tube rack	4001045
Macro digestion tubes	4042300 6x
Fume collector	4001056
Control wire	0403507
Power supply cord.	0007001

4001048 BLOC DIGEST 24x micro tube Ø26.

Heater block 20p	4001051
Temperature profiler	4001538
Tube rack	4001045
Macro digestion tubes	4042300 6x
Fume collector	4001057
Control wire	0403507
Power supply cord.	0007001

4001049 BLOC DIGEST 40x micro tube Ø26.

Heater block 6p	4001052
Temperature profiler	4001538
Tube rack	4001045
Macro digestion tubes	4042300 6x
Fume collector	4001058
Control wire	0403507
Power supply cord.	0007001

4000638 Thermoreactor DQO 6x tube Ø42.

Heater block 6p	4000507
Temperature profiler	4001538
Tube rack	4005071
Tube D.Q.O.	1000641 6x
Reflux condenser	1000642 6x
Condenser tube stand 6x	4000643
Control wire	0403507
Power supply cord.	0007001

4000639 Thermoreactor DQO 12x tube Ø42.

Heater block 12p	4000508
Temperature profiler	4001538
Tube rack	4005081
Tube D.Q.O.	1000641 12x
Reflux condenser	1000642 12x
Condenser tube stand 12x	4000644
Control wire	0403507
Power supply cord.	0007001

4000640 Thermoreactor DQO 20x tubo Ø42.

Heater block 20p	4000509
Temperature profiler	4001538
Tube rack	4005091
Tubes D.Q.O.	1000641 20x
Reflux condenser	1000642 20x
Condenser tube stand 20x	4000645
Control wire	0403507
Power supply cord.	0007001

3. Equipment description

Digestion unit BLOC-DIGEST is used for digestion in Kjeldahl nitrogen determination and measurement of the COD. Other possible application is hydrolysis of samples in determining fat.

The digestion equipment usually uses a heating block, a control unit, a rack for sample tubes manipulation and a fumes removal system.

For smokes neutralization some other equipments are needed:

4001611 Scrubbing unit

4001610 Water circulation pump

The heater block houses the digestion tubes, in order to get an optimal heating transfer to samples. It only uses tubes with suitable diameter according to the heater.

The heater is protected against overtemperature with a device that cuts the power supply if temperature exceeds 550°C.

The temperature of the heating block is regulated by the control unit that is separated from the heating block to prevent spoilage splashes, spills, fumes, etc. This unit (RAT 2) allows up to 20 programs of 4 steps for temperature and time.

In the Kjeldahl digestion or hydrolysis our for heating in several steps allows optimal control of the foams produced, to be able to evaporate water from the sample before digestion at 400 ° C

For extracting the fumes produced by the digestion, fume collector may be used with a jet pump, a sucker or scrubber.

Binding of steam hood with the sample tube is not airtight. This allows absorbing small explosions without damaging the glassware. The extractor fan leads to the drain due to the vacuum or suction produced at the end of the collector outlet tube. Therefore it is essential to the proper functioning of the jet pump must produce a negative pressure of 0.6 bar.

For D.Q.O. digestion fumes that are not removed no fumes are extracted but few glasses condense the steam produced is used and return the liquid to the tube sample. (Reflux condenser).

The Kjeldhal digestion systems evolve irritating acid fumes.

They need a proper installation for treating released fumes.

3.1 Main parts and assembly

Assemble the digester's system according to this guide.

Note: The fumes have, for storing and transportation, a small metal part to protect the white gasket.




Remove this parts if they are from metal.

Fumes removal

Allocated over the sample tubes, they drive the fumes to scrubbing and neutralizer system.

Digestion tubes

Allocates the samples.


Sample tubes rack

Accommodates the sample tubes to allow its transportation and provides a thermal isolation when they are on the heater block.




Temperature controller RAT

Controls the temperature profiles applied to samples.

Heater

Heats the samples.

Grid holder

Acomodates the sample tube rack when it's outside of the heater.

Fumes removal holder

Allocates the fumes collector caps when its outside the heater system.

3.2 Technical specification (Digester blocks)

Block:	6x MACRO	12x MACRO	20x MACRO	12x MICRO	24x MICRO	40x MICRO
Tubes:	6x Ø42	12x Ø42	20x Ø42	12x Ø42	24x Ø42	40x Ø42
Electric power: (W)	1600	2100	2500	1600	2350	2500
Weight (kg)	18	25	31	16	22	27

Process control:

Temperature range:	45..450°C
Load maximum power:	10A
Temperature stability:	±1°C
Memory:	20 programs of 4 steps
Maximum time per step:	600 minutes
RS-232 serial communication:	Optional


4. Installation


Installing a Kjeldahl digestion system requires a suitable place for it specifically, due to acid fumes generated during digestion.

The installation within cabins, hoods, ... should consider the compatibility of the materials used with acidic fumes.


It is best to neutralize the acid fumes through a "scrubber" and a suction pump. The scrubber that contains a concentrated solution of calcium carbonate, the effect of the vacuum produced by the "Pump" digestion gases pass through the calcium carbonate solution that neutralizes the gases released by the evaporation of sulfuric acid digestion.

"Scrubber" unit 4001611

Assembly with fume extraction by vacuum tube (Its effectiveness depends on the vacuum tube).


"Scrubber" unit 4001611

- Place all the equipment on an appropriate construction (table, bench, ...) to the weight of it.
- Remote or isolated from the presence of personnel.
- Need several outlets.

Details for fume extraction:

1. Kjeldahl digestion generates highly corrosive acid fumes and irritating to humans and animals. They are not toxic.
2. Although means for extracting fumes are installed, it should be set in an isolated area of the laboratory.
3. Avoid unnecessary exposure to fumes.
4. If installed in a fume hood, this must be constructed of suitable materials, plastics, glass, ...
5. The draw-fumes with a jet vacuum pump may be sufficient, but within an enclosed with external ventilation. Its effectiveness depends on the actions of the vacuum tube. A water tap (with a minimum pressure of approximately 2.5 bar) and a drainage are needed.
6. Fumes extraction with a scrubber unit or two and a vacuum pump water circulation is recommended. It can be done without hoods or booths.
7. For MACRO 20 places systems, it is recommended installing two Scrubber units both connected to the two vacuum inputs.

Vacuum pump for the absorption of fume.

ATTENTION:

Possible splashes of
sulfuric acid.

Use

gloves and

protective goggles

in the manipulation of the rack with
samples

or to approach the digester for visual
inspection.

**WARNING:**

Use gloves.

**The drops that falls
from the
collector or the hoses
are from
sulphuric acid.**

**Manipulate with caution
to avoid chemical burns on the
skin.**

4.1 Digestion process operation

The suggested sequence for preparing and starting the complete digestion system is:

- Prepare the reagents and the samples in the digestion tubes. See the digester's instructions.
- Program the temperature controller for a three steps digestion cycle. According to the digester instructions manual.
- Insert the sample tubes in the block-digester and cover them with the fume collector.
- Check that the vacuum pump is ready. That is, the water tank and the external circulating in OFF.
- Check that the scrubber bottles have the sodium carbonate solution.
- Assembly the hoses:
 - Ø 10 hose from the Scrubber to the vacuum pump.
 - Ø 12 hose from the collector to the Scrubber.
- Start the digestion cycle from the temperature controller.
- Start the vacuum pump.
- **Visually check every 30 or 40 minutes, the digestion process.**
- When finishing the digestion process, the temperature process beeps.

Then carry out the following operations, with caution:

- Remove the tubes rack from the digester without disconnecting any hose and leave it on the holder.
- Keep the vacuum pump on.
- Remove the fume collector from the tube rack and leave it on the holder.

- **Fume collector cleaning:**

You can simply clean the fume collector with water.


- **Water tank emptying:**

To empty the tank, just remove the «emptying tube» from the upper connector.

- **Water circulating function:**

The vacuum pump can also work as a water closed circuit.


- Connect the hose to the «water circulating out port».
- Place the «water circulating switch» in the ON position.
- When starting the pump, this will move the tank water to the «water circulating out port».

The RAT temperature profiler perform "profiles" as a sequence of "steps" defined with time and temperature.

5. Programming temperature profiler

5.1. Information shown on display

Note: Pay attention where the cursor is located.


5.2. Profiler states

There are four possible states of the RAT controller.


- STP: "Stop" which is the default state, waiting for the user to select a temperature profile to execute or modify.
- RUN: The controller executes a temperature profile.
- PRG: Menu for programming a temperature profile.
- CFG: Configuration options, units and calibration.

Note: State: PRG and CFG, the controller returns to the state STP if no key is pressed for more than 10s.

Note: Once the program running is stopped by pressing the R / S key.

5.3. Main menu


Use ▲ ▼ to set Celsius o Fahrenheit temperature units.

Use ▲ ▼ to set an offset value for temperature calibration if necessary.

Use ▲ ▼ for interface language set.

5.4. Programming menu

Select a profile to configure (default PR 1).

Select a step (step) to configure it (default Step 1) by using keys ▲ ▼.

Temperature selection keys using ▲ ▼. The temperature should always be higher than the previous step (by default, the computer displays the temperature of the previous step).

Select duration (minutes) step using the keys ▲ ▼.

Continue to the next step.

⋮

The temperature profile that ends in step lasts "OFF".

5.5. Review a temperature profile

In stop state (STP) and using the arrows, program steps are displayed one to query parameters.

6. Profiles for digestion process

The possibility of performing a process of digestion temperature profile has several advantages, especially in the control of foam produced by the control samples and smoke generation.

Control foam production.

The foam production during digestion has to be controlled, otherwise the foam may overflow tube height of digestion and go outside, with two consequences:

- Damaging the block heater. (The foam contains sulfuric acid)
- Causing loss of sample.

The foam is caused by the moisture (water) sample. So the first step temperature profile is dedicated to the removal of moisture (water) sample.

Control of production Smokes.

Fuming during digestion has to be controlled, otherwise these fumes, highly acidic, are sent to the atmosphere.

The second step of the temperature profile is dedicated to smoke control.


The Kjeldahl digestion has an area of "white smoke" to 300 ° C.

Profile typical Kjeldahl digestion.

Step 1: 125° 30 'Extract moisture.

Step 2: 300° 30 'Checking white fumes.

Step 3: 400° 60 'ammonium mineralization.

RAT controller executes temperature "profiles", as a sequence of "steps" defined with time and temperature.

- Dried Samples: Example: cereals: Step 1 may be 15 '
- Normal Samples: Example: Meat: Step 1 may be 30 '
- Samples "water": Example: residual waters: Step 1 may be 90 '

Examples of other "steps that may be necessary according to the sample:

- Low recovery of ammonium: Increase step 3
- Foam: Add a step to 150
- Excess smoke: add a step to 280.

Overall for better control, use more steps in the profile.

Profile for COD digestion.

Step 1: 150° 15 '

Step 2: 170° 60 'Digestion.

Profile Objective: providing a gentle rise in temperature.

Profile for hydrolysis (in determining fat SOXHLET).

Step 1: 150° 15 '

Step 2: 170° 60 ' Digestion.

Profile Objective: providing a gentle rise in temperature.


7. Spare parts

To ensure safety of the equipment spare parts must be purchased from, J.P. SELECTA, s.a.u.

Descripción / Description	Código/Code
Resistencia/Heater element 1600W (Macro 6)	39070
Resistencia/Heater element 1600W (Micro 12)	
Resistencia/Heater element 2100W (Macro 12)	39069
Resistencia/Heater element 2100W (Micro 24)	
Resistencia/Heater element 2500 W (Macro 20)	39068
Resistencia/Heater element 2500 W (Micro 40)	
Sonda tipo K / K termocuple probe	43026
Junta vitón blanco(Macro) /White viton gasket (Macro)	21157
Junta tórica vitón verde (Macro)/Green viton o-ring gasket (Macro)	21158
Junta vitón negro (Micro)/Black viton gasket	21167
Tubo viton / Viton tube	46084
Rampa vidrio (Macro 6)/Ramp glassware (Macro 6)	47003
Colector vidrio (Macro 6)/Colector glassware (Macro 6)	47007
Rampa vidrio (Macro 12)/Ramp glassware (Macro 12)	47002
Colector vidrio (Macro 12)/Colector glassware (Macro 12)	47006
Rampa vidrio (Macro 20)/Ramp glassware (Macro 20)	47001
Colector vidrio (Macro 12)/Colector glassware (Macro 20)	47005
Rampa vidrio (Micro 12)/Ramp glassware (Micro 12)	47013
Colector vidrio (Micro 12)/Colector glassware (Micro 12)	47014
Rampa vidrio (Micro 24)/Ramp glassware (Micro 24)	47015
Colector vidrio (Micro 24)/Colector glassware (Micro 24)	47016
Rampa vidrio (Micro 40)/Rampa glassware (Micro 40)	47011
Colector vidrio (Micro 40)/Colector glassware (Micro 40)	47012

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9. Kjeldahl analysis step by step

9.1 Sample Preparation

- Grind, mix thoroughly and mix the sample.
- Weigh between 1 and 2 grams of sample.
- In samples with very small nitrogen content (sewage, etc) take enough samples containing at least 5 mg of nitrogen.

9.2 Digestion

- Add between 10 and 15 ml (tube macro) 96-98% H₂SO₄ and 15-20 g of catalyst (for the micro tube is maximum 5ml H₂SO₄).
 - Set up a system for the extraction of fumes or scrubber with NaOH.
 - Perform digestion in three steps:
 1. Depending on the water content of the sample digestion start evaporating water at 150 for 20 to 60 minutes.
 2. Perform a second step at 280°C for 30 minutes to reduce the production of white fumes.
 3. Continue digestion at 400 for 60-90 minutes.
- Visual Control: The result is a clear liquid with clear, green or blue coloration depending on the catalyst used. There should be no black debris attached to the tube wall.

Note: During digestion foam production in the samples must be controlled. If this is excessive, should be extended step # 1.

9.3 Dilution

- Remove the tubes shows the block digester and allow to cool to room temperature (It may be forced immersing the tubes cautiously in a little water).
- Add about 25ml of water to each tube (10 Tube MICRO).
- Add the water slowly and moving the tube if allowed to solidify the sample. If necessary, gently heat the tube (for example introducing it into the heating block still hot).
- Allow to cool back to room temperature.
- Avoid loss of nitrogen and violent reactions to insert the distillers still hot tube.

9.4 Distillation

- Check the level of NaOH and Boric Acid in the tank. Check that the metering pumps are primed and meter the correct volumes.
- Place a 250ml Erlenmeyer flask to the condenser with 50 ml of boric acid and a few drops of indicator.
- Schedule a dosage of 50-75 ml of NaOH (25ml for MICRO).
- Insert the sample tube into the distiller.
- Distill to collect in 250ml Erlenmeyer (50ml distillate Boric 200ml).

Visual Control: Once the NaOH was added, the sample must have a bluish color, if not, add more NaOH.

9.5 Tritation and calculation

- Titrate the distillate with HCl or H₂SO₄ 0.1 or 0.25 N to change colour (Endpoint: pH 4.65).
- Perform the calculation: Nitrogen mg = 14 * Volumen HCl * Normality.
- To move to protein corrected by the appropriate factor depending on the nature of the sample (6.25 by default).
- Periodically perform a blank test and subtract the result.

10 Hydrolysis for fat determination

Material: (for 6 samples)

- 6 glass funnels Ø100mm
- 6 filter papers Ø100 normal (not slow one)
- 6 digestion tubes (SELECTA code 4042300)
- Complete digestion system.
- 6 Erlenmeyer 500ml.

Reagents needed:

- 1L of HCl 4N (350ml de HCl 35% and add until 1L. with distilled water)
- Granules of pumice stone.

Procedure:

- Number the digestion tubes.
- Weight around 3gr. of sample and put it into the digestion tube.
- Add 100 ml of HCl 4N in each tube.
- Add a spoonful of granules of pumice stone.
- Place this into the block and digest at 100°C - 130°C during 1 hour (see page 27 table). It must boil slowly.
- Retire the digestion block and let it cool down.
- Prepare the Erlenmeyer with the funnels and the filter paper. See the photo below.
- Number the Erlenmeyer.
- Empty the content of the sample tube into the funnel with filter paper.
- Clean the tube with water and empty into the funnel. Repeat it until the water from the tube is transparent.
- Put water on the funnel.
- Empty the Erlenmeyers.
- Repeat the operation until the water from Erlenmeyer is transparent (about 5 times).
- Let completely dry the paper filter (1 hour).
- You will observe that the filter keeps the sample fat.
- Proceed to the extraction with the det-gras equipment. Let dry for 5h, and to extract in two phases: Boiling and rinsing and let dry for 30 minutes.

11 Digestion on C.O.D. determination

- Add 5.0ml of potasic dichromat and some glass balloons to regulate ebullition to the sample and it is homogenized carefully.
- Add slowly 15ml. of silver sulphate - sulphuric acid, shaking the tube carefully with a circular movement, cooling down with cool water of the tap or an ice bath, in order to avoid any loss of organically volatile substances.
- Place the tubes dully prepared into the rack tubes holder on their support. With the help of a rack, insert it into the heating block.
- Join on each tube a refrigerant in order to fit correctly the polished gaskets.
- Program the operating temperature (170-200°C) and the length time (around 2 hours) of the back flowing into the temperature regulator and RAT time.
- Press on «START» and at this moment the metallic block will heat where the tubes are placed. This will be done in a progressive way and when the preselectionated temperature is reached, the time of back flowing will be counted.
- Once the time is finished, the heating in the metallic block stops automatically and an acoustic sound is heard. At this moment, you should press on «STOP» in order to finish the process.
- Let the tubes cool down for a while and clean the projections that have been produced in the inner wall of the refrigerants into the sample recipient with a little volume of distilled water.
- Separate the refrigerants from the tubes, and the tubes are extracted with the help of the rack, that is placed on the support. The sample is diluted with 75 ml. of distilled water and it is cooled down at ambient temperature.

Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.