DESTILADOR ENOLÓGICO
ALCOHOL DISTILLATION UNIT

DE-1626

4001626

Indice

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1 Packing list and accessories

1.1 Standard package content

- Distillation unit
- Tube 1.5m water input. 10 bar hose, connector 3/8". (1)
- Tube 1.5m water outlet. Silicone 8x14, connector 3/8". (2)
- Tube 1.5m remain steam. Silicone 8x14, connector 3/8". (3)
- Tube 1.5m boiler empty. Silicone 8x14. (4)
- Tube 1.5m water tank. Silicone 8x14, conenctor de 1/4". (5)
- Mains cord (6)
- Distillation Tube Ø 42 4042300 (7)
- Distillation Tube Ø 52 1001422 (8)
- Distillation Tube Ø 80 1000646 (9)
- 200ml graduated flask 1001636 (10)
- Drop collector tray (11)
- Tank 2l for H₂O (12)
- Instruction manual (12)

1.2 Accessories

Tubes, racks & tongs:

- Distillation tube Ø 42 4042300
- Distillation tube Ø 52 1001422
- Distillation tube Ø 80 1000646
- Graduated flask 200 ml 1001636

Calibrated alcoholimeters



- 0 - 10 1001629
- 10 - 20 1001630
- 20 - 30 1001631
- 30 - 40 1001632
- 40 - 50 1001633
- 50 - 60 1001634
- 90 - 100 1001635
- Alcoholimeter holder 1000015

Rack :

- Rack for 6 tubes Ø42. 4000648
- Rack for 6 tubes Ø52. 4001473
- Rack for 6 tubes Ø80. 4001613

Safety Goods:

- Tongs for tubes 7000532
- Heat protection gloves 1001255
- Safety glasses 5028192

2 Safety

The alcohol distillation unit DE-1626 has the suitable measures for its use in laboratory. However, working with it can be hazardous and the user must be informed of some of those risks. In order to be respected, the risk of danger situations are remarked across this manual.

2.1 Safety symbols

It identifies the risky and safety actions to be applied. The symbols make reference to the paragraph marked with a grey line.

Warning

Warning of a dangerous operation.

Indicated instruction must be followed.

Hazard of electrical shock

Electric hazard for entering in areas marked with this signal, or for following the manual instruction marked by this icon. Follow the instructions for the described operation.

Hazard of burns to be in contact with hot surfaces

Temperature in this area can exceed 60°C. Use heat protected gloves for the described action. Follow the instructions for the described operation.

Important information

- Important information to obtain good results or for an optimal operation of the equipment.
- Important information to extend the life of the equipment or to avoid the degradation of some of its components.

2.2 Risks which the operator is put under

The process of determination of Kjeldahl Nitrogen with the described equipment exposes to the operator to the following situations of risk:

- Manipulation of glass parts.
- Possibility of touching parts to a temperature above 60°C.
- Possible exhibition to irritating acid steam.
- Electrical risk.

2.3 Personnel training

This equipment must only be operated by personnel trained and suitable qualified to work in a chemical laboratory and its common hazards.

This equipment must only be operated by personnel who have read and understood this manual or have been qualified on its manipulation and operation.

2.4 Proper uses of the equipment

This equipment is for laboratory use, for sample distillation by steam pulling. Mainly in Kjeldahl nitrogen analysis.

The improper use of this unit (not following manual instructions) can endanger the protection prepared by the equipment.

Do not use to distillate explosive or flammable substances.

This equipment its not designed and not manufactured to work in an explosive ambient area.

2.5 Modifications

The operation modification or the manipulation of the safety systems unauthorized by the manufacturer could expose the operator to hazards not indicated in this manual.

3 Overview

The DE-1626 is a steam distillation, mainly designed to obtain the distilled from sample for later measuring the %Vol (alcohol contents). The unit is validated according to CEE nº 2676/90 (wine analysis) and CEE Nº 2870/2000 (high degree drinks analysis) for alcohol degree determination.

The DE-1626 yields a fast steam distillation and an efficient recovery of volatile substances, as alcohol and acetic acid.

As the DE-1626 will essentially be used for alcohol degree determination, this procedure is given in high detail on the chapter 11. Chapters 5 and 6 describe the general operation of the unit and how to obtain the distilled from samples.

The unit has the following systems:

The steam generator

It generates steam from the water reservoir placed outside the unit.

The steam generator includes a **Safety thermostat** that cuts the electrical supply to the steam generator when overtemperature is detected (probably due to lack of water) and a **Safety pressure switch** that interrupts the electrical supply to the steam generator when there is an overpressure.

It also includes a conductivity probe to maintain a constant water level in the steam generator.

The distillation circuit

It is mainly made by glassware:

- The rectifier: with a special design to concentrate the volatile substances.
- The two stage coolant: it condensates the steam leaked from the sample. This glassware is cooled with a external water circulation that is automatically open/closed.


Automatic stop and addition for leveling

Made by a level sensor device and a valve that closes the distilled output and allows collector flask leveling.

Control panel

Addition valve to adjust level.

Level sensor for automatic stop.

Collector flask

Distillation tube

Plate, fuses and connector to power supply.

4 Installation

Once the unit is unpacked, visually check the packing list (see 1.1) and that there are no parts damaged during transportation. Save the packing for some days for later sending.



4.1 Suitable location

Place the unit over a flat level and stable surface, suitable for the unit weight.

No longer that 1.5m from a water tap, waste output and power supply plug.

To work comfortably, let a free space of 50cm on each side and not least than 20cm on the rear part.

Prepare to refill the water reservoir from the top of the unit.

4.2 Connection to power supply

Locate a power supply socket (rated to 16A) near the unit.

Use the power cord supplied or another of similar specification.

For safety reason, the socket must have ground terminal, since the unit must work grounded.

Before plugging on the socket, check that its voltage match the voltage marked on the unit plate.

4.3 Connection to tap water and waste

Connect the water entry hose to a water tap (1). Connect the water output hose from (2) to waste.


The hoses «remaining steam» (3) and «Steam generator empty» (4) must go to waste.

Do not forget the rubber white nuts on the hose's connections.

Press plastic connections by hand. Do not use any tool.

The tap water must have, at least, 1 bar pressure.

Keep the tap opened. The DE-1626 will open / close the cooling water automatically.

Note: It is critical in order that the unit works properly, the hose «Remaining steam» (3) MUST NOT BE IMMersed under water, or blended or blocked by any reason.

4.4 Water for steam generator

The unit includes a 2 litres water reservoir, at the top, to store the water for the steam generator.

As the system used to control the level of water on the steam generator works by conductivity, water of at least 20 microSiemens must be used. Water, also, must have a low content of calcium to increase the life of the steam generator.

Water with this conductivity could be obtained from a demineralised unit. Or from distilled water and add 5 mg/l Sodium Chloride.

5 Operation

The control panel allows function selection and shows the distiller alarms status

Once the unit is plugged in, put it on, by pressing the main switch. The unit goes to stand by status (ON).

Note: If main switch is pushed consecutively with an interval of less than 3 seconds, the unit could start in abnormal mode. So wait for a few seconds to restart the unit.

On starting the unit for the first time, the water pump starts to work, automatically, to fill the steam generator. This operation usually takes near 40 seconds and it will stop when the steam generator is filled. Restart the unit if the alarm blocks it.

If, after this 40 seconds, this operation is approx not finished, please check that the water tank (visual check on the right hand side level), and also check the water conductivity. (See 4.4)

5.1 Insert/ Remove the distillation tube

- Insert the distillation tube starting from the PTFE tube.
- Push down the fixing lever.
- Check the upper black rubber closing and the tube fixing.

Proceed in a similar way to remove the distillation tube.

ATTENTION

Use thermal protection gloves.
The distillation tube is above 60°C
once distillation is finished.

Check properly
close

Check
fixing

5.2 Place the collector flask.


The collector flask must be placed carefully to assure that the automatic stop system will work properly.

The flask wall must be in contact with the level sensor face.


Note:

To assure that the automatic stop system will work properly:

The front side adjustable legs must be adjusted to left the unit slightly up from the front side.


Placing 200ml graduated flask

Placing the Erlenmeyer

START

5.3 Distilling

On pushing the «START» key, the unit starts sample distillation.

This opens coolant water circuit and initiates the steam generation.

After 30s approx the sample starts to boil and its steams are driven to the cooler glassware for its condensation. The distilled is collected on the graduated flask or on the Erlenmeyer.

When the level on the collected distilled will be near the sensor, the distillation will be stopped automatically and the outlet valve closed to avoid that the level exceeds the mark of the graduated flask. Approximately 180ml will be collected.

Push the «STOP» key if you want to stop the sample distillation process.

5.4 Levelling

The level sensor will stop the distillation approximately at 180 ml.

Remove a little bit (see picture), the collector flask to add some drops of distilled to adjust up to near 200ml (an accurate adjust to 200ml must be done later when solution will be at 20°C).


On pushing «ADD» key, the outlet valve is opened to allow passing of some drops of distilled from the cooler glassware. Push it some times to adjust the level near the flask mark.

Proceed similar with the Erlenmeyer flask.

Placement for the collector flask to add drops to near the graduated mark.

ADD

5.5 Alarms

To avoid dangerous situations for the operator or to protect the unit life, it includes some safety systems and some alarms:

(1) Alarm: «No water for steam generator»

- If the lamp (1) lights on sporadically :This means that water level control has detected a low level of water on the steam generator and automatically starts the water pump to get the suitable level of water.
- If the lamp (1) lights on permanently: This means that water level control has detected a low level of water on the steam generator for a long time and it automatically refills with water. This situation stops distillation. Probably the water reservoir is empty.

If the alarm blocks the unit, reset the unit by ON/OFF button.

(2) Heater status lamp.

This lamp lights when the heater is ON.

(3) Alarm: «Overtemperature» on the steam generator

This red lamp lights on when an overtemperature on the steam generator has switched the safety thermostat.

It inhibits all distiller functions. And cuts the power supply to the heating element.

Once it is switched, it must be let to get cool for some minutes. Solved the problem that causes the switching and rearmmed the safety thermostat by pushing the central button. The safety thermostat is located on the electrical box, right hand side cover must be removed to get it.

The safety thermostat switching is due to low level of water on the steam generator, probably because the water reservoir is empty or due to water pump malfunction.

6 Sample distillation

To obtain the distilled with the maximum content of alcohol or volatile substances and to get a good performance, the following working routine must be followed. Once the distilled is obtained, the alcohol degree must be measured (see chapters 11 and 12).

6.1 Check the level of external water tank

Refill the water tank with distilled water stabilized with 5 mg/l Sodium Chloride to get a minimum conductivity of 20 microSiemens.cm (see 4.4).

6.2 Unit warm-up

To assure that the unit gives maximum efficiency on the distillation, a warm-up step must be run previous sample distillation batch. To do it, just insert a distillation tube (any) with 25ml (aprox) of water and collect the distiller during 5 minutes (use an Erlenmeyer flask as collector).

Run this warm-up step after 1h of inactivity.

6.3 Choose the properly distillation tube

The distillation tube must be properly chosen to avoid the rectifier contamination and to get an optimal energy transfer to sample.

Use the following guide to choose the distillation tube:

- Tube Ø42 for samples of volatile acidity.
- Tube Ø52 for samples of volatile acidity.
- Tube Ø52 for samples of 100ml for alcohol degree.
- Tube Ø80 for samples of 200ml for alcohol degree.

Collector flasks.

- 200ml graduated flask to collect distilled for alcohol contend.
- 250ml Erlenmeyer flask to collect distilled for volatile acidity.

6.4 Distillation

Please follow this steps: (on the same order) (see also chapter 11 for accurated details)

- Visually check that the solenoid cooler is completely empty from rests of previous sample. If necessary, remove it on a separate flask by pushing «ADD» key.
- Add near 10ml of distilled water on the collector flask. This is to avoid the loss of alcohol on the first condensed drops.
- Place the collector flask on the flask holder. The flask wall must be in contact with the sensor face. The end of the tube from the cooler must be immersed under the level of the 10 ml of water.
- Place the sample on the distillation tube and insert it on the DE-1626.
- Push «START» key.
- When the level of distilled obtained reaches the sensor face, the unit will stops distillation automatically.
- First of all, remove the distillation tube.
- Add distilled drops by pushing «ADD» key until levelling the collector flask.

7 Unit verification (OQ)

Once the unit is installed and the operator trained with its operation, the unit could be qualified.

This process could be used as a Operation qualification (OQ).

There are two basic qualification operations:

7.1 Functional checking

Run this operation during the installation steps and daily for unit warm-up.

Follow chapter 5 instructions to distillate until the unit stops automatically (use Erlenmeyer or graduated flask).

When the unit is on the steady state, the distillation rate is between 30 - 40 ml per minute.

7.2 Checking the maximum etanol loss

This checking is widely used to certify that the unit works properly. Run some of this tests for training purposes before analyzing valuable samples.

It consists on preparing some 10% dilution of distilled water and etanol.

Prepare this solutions by using flask, and solutions at 20°C and room temperature at 20±3°C.

Obtain the distiller (see chapter 6 and chapter 11) 5 consecutive times and measure the alcohol degree of the distilled. The distilled must have an alcohol degree at least of 9.9%Vol. In other words, alcohol losses must be at east of 0.02% per distillation.

After the distillation, compare the original sample with the distilled one.

Note:

To get these results, you must work with high precision:

- Keep part of the hidroalcoholic initial solution to compare with the final distilled.
- Store all flasks stopped.
- Use only certified volumetric material.
- Work at room temperature 20±3°C.
- Do meassure only at 20°C 20±0.5°C

See chapter 11 for more details.

8 Maintenance

The distillation unit DE-1626 needs a minimal maintenance operations. Run visual inspections inside to locate possible losses of liquid or steam on the joints, connections, or hoses.

To avoid metal parts corrosion, dry any occasional leakage of water inside the unit.

8.1 Cleaning the unit

8.1.1 Cleaning the distillation glassware

Every day clean the distillation glassware after samples distillation. Just distillate 50ml of distilled water for 5 or 6 minutes.

8.1.2 Chamber cleaning

The chamber and the drops collect tray must be cleaned every day just with a humid cloth and dried, to avoid deterioration.

8.2 Inspection of hoses (weekly)

Visually check the hose condition and its connections to different components. To avoid premature deterioration of the unit due to liquid leakage, replace as soon as possible any hose that has any leakage or damage.

Take these following cautions for hose replacing:

- Unplug the unit from the electrical power supply.
- Try, as far as possible, (see 8.3) that the hose to remove is empty.
- Use only hoses of material indicated on 11.2.
- Use clamps as if factory assembled.
- Once the hose is replaced, it must be checked for leakages.

8.3 Preparing for transport or storage

Once the unit has been used, proper actions should be taken for transport or storage. Initially clean the unit according to the former instructions.

8.3.1 Remove completely the water from the unit

- **Emptying the water reservoir and the steam generator**

1. Open the tap (4).
2. Start the unit, push «START» and then wait until the water pump empties the water reservoir.

- **Emptying the water from the cooling circuit**

1. Disconnect the hose for water **entry** from the tap water.
2. Switch ON the unit and push «START».

Note: From this point ,actuate in a fast way to avoid overtemperature on the steam generator.

3. Blow for the end of the **output** coolant water. The water will **come out** on the end cap of the **water in** hose.

CAUTION

Do not ship off the unit if water has not been completely removed.

During transportation, water leakages could fall over important components and cause severe damages.

During transportation, water leakages could fall outside of the package.

9 Technical specification

Mains voltage:	230V 50/60Hz (See plate)
	Other power supply voltages could be available (Ask to your dealer)
Electrical power:	2450W (See plate)
Steam generator power:	2400W
Maximum consumption:	11 A
Fuses:	15 A fast, dimensions: 20x5 mm
Power supply cord:	Phase, Neutral, Earth (section 1mm ²)
Weight (Net):	35 Kg
Dimensions: W x L x H	300 x 300 x 900 mm
Water reservoir volume: (external)	2 liters
Steam generator water consumption:	1 ... 1.25 ml per ml of distilled obtained.
Cooling water flow:	6 l/min at a temperature lower than 20°C
Ethanol losses:	<0.02% Vol per distillation
Distillation rate:	approx. 40 /minute
Short distillation time:	6 minutes, approx 200ml.
Long distillation time:	8 minutes, approx 250ml
Cooling water entry connector:	3/4"
Cooling water output connector:	3/4"
Remaining steam output connector:	3/4"
Sound level:	Lower than 60 dBA
Ambient conditions:	Inner use. Altitude up to 2000m Temperature from 5°C to 40°C. Relative maximum damp 80% for above 31°C temperatures, decreasing linealy to 50% for 40°C.
Overvoltage category:	Category II
Pollution degree:	2
Max pressure of boiler:	4 bar
Pressure switch cut:	0.1 bar
Glass material:	Borosilicate 3.3
Carriage material:	Iron with chemicals resistant coating
Chamber material:	Aluminium with RILSAN coating chemical resistant
Steam generator material:	Stainless steel.
Hosepipe material:	Silicone

32050

28129

10 Technical service

For repair and substitution of deteriorated parts, please contact with J.P. SELECTA, s.a.u. or your distributor. They will inform you about the technical services authorized in your area.

Find a spare parts list for easy buying and identifying.


10.1 Spare parts

- 16213** Electrovalve N. A.
- 16067** Water electrovalve.
- 16025** Distilled outlet clamp valve.
- 24278** Water bomb support.
- 28129** Lockgum.
- 28061** Drop tray.
- 32050** Steam generator 2400W.
- 47043** Cooler glass DE-1626
- 47044** Rectifier glass DE-1626
- 47045** Cooled union tube glass DE-1626
- 47259** Anti feedback valve.
- 29520** Relay electronic PCB.
- 29511** Control panel electronic PCB

Hoses:

- 46030** **TEFLON TUBE 8X10MM** Tube immersed on the sample.
- 46041** **Silicone tube 8x14** Water and Steam.
- 46042** **Silicone tube 6x10** Water and Steam.
- 46143** **3/4 WATER INLET HOSE** Connection to the water network
- 28264** **LFL Tygon tube** Distilled outlet.

Adjusting screw

10.1.1 Adjusting the level sensor

It could be possible that along the time the collector level sensor changes its sensibility and never detects the distilled level or it actuates just by sensing the glass flask.

Turn the adjusting screw (see picture) to adjust it.

Place a flask with water on the collector and turn the screw until the yellow LED lights.

WARNING
Before any work inside the equipment
disconnect it from the mains.

11 Measuring the alcohol degree

11.1 Alcoholometer


The volumetric alcohol degree (written as %Vol) is defined as the number of litres of Ethanol and his homologous (Metanol, etc...) content in 100 litres of wine (or alcoholic drink in general), all measured at room temperature of 20°C.

The determination basis is easy: the alkalinized sample is distilled and the alcohol % is measured by means of special graduated aerometers named alcoholometers.

There are other methods to measure the alcohol content on the distilled (hydrostatic scale, refractometer, ...) but the alcoholometer is the more accurate and repetitive method.

To obtain accurate results, you must work with a method that considers all factors that could have influence on the volumetrics and on the alcohol losses.

A complete method is presented here. It could be simplified according to laboratory requirements.

Alcoholimeter reading

11.1.1 Material

- Cylinder 250 ml (30cm high)
- Alcoholometers (see «accessories» on Chapter 1).
- Thermometer 0-35°C resolution 0,5°C.
- 200ml Calibrated flask (certified) with stopper.
- Bath or oven to work at 20°C.

11.1.2 Reagents

- Calcium Hydroxide 2 mol/l: PANREAC 625409
- Anti-foam liquid silicone PANREAC 211628
- Distilled water.

11.1.3 Production of distilled

In order to not have a negative effect on the results, room temperature must be $20 \pm 3^\circ\text{C}$.

Material preparation


Store all the material that will be used at $20 \pm 0,5^\circ\text{C}$ (cylinder, alcoholometers, calibrated flask, wine, samples..)

Removing CO₂

If the sample is a young or sparkling wine, proceed to remove the CO₂ by stirring the sample or with a ultrasonic bath.

Distilling

- Fill the 200ml calibrated flask up to his mark with the sample. Remove bubbles (all at 20°C).
- Pass the content to the Ø80 distillation tube.
- To be sure that all the sample has been passed to the distillation tube, clean 2 or 3 times the calibrated flask with around 10ml of distilled water and add it to the distillation tube.
- Add around 10ml of calcium hydroxide solution.
- Add some drops of anti-foam silicone solution.
- Use the same calibrated flask to collect the distilled.
- Add around 10ml of distilled water on the calibrated flask and place it on the DE-1626. Check that the output of the tube remains immersed under the water level (see figure).
- Start distillation.

Check that the output of the tube remains immersed under the water level.

levelling

- Once the distillation is finished, add some drops of distilled (see 5.4) up to near (without arrive) the 200ml mark.
- Put the top on the flask and let it on the 20°C bath or cabinet for at least 3h.
- Add distilled water at 20°C until the 200ml mark.

Measuring with the alcoholometer

Previously to start this operation, all material that will be used (flask, cylinder, alcoholometers,...) must be stored at $20 \pm 0.5^\circ\text{C}$ for 4 h.

- Put distilled into 250ml cylinder.
- Place the alcoholometer (with a little spin movement).
- Wait 1 minute. Read %Vol from the bottom part of the mechanism.
- Repeat the reading of %Vol, 3 times at least and average. Write the result by using 2 decimals figures. Round last figure at 0 or 5.

Notes:

- If it is necessary, use a magnifier to read better the alcoholometer.
- To assure that the liquid temperature do not change during %Vol measurement, the room temperature must be $20 \pm 3^\circ\text{C}$.
- If the material used has not been stored at 20°C , measure the temperature for 1 minute before introducing the alcoholometer. The measured temperature must be round. (Example: 20.4°C will be round at as 20°C ; 20.6°C as 21°C) and then modify the %Vol according to table 1. (This will produce a leak of accuracy on the result)

	Alcoholic strength by volume,% vol.												
	8	9	10	11	12	13	14	15	16	17	18	19	20
14°C	+0.85	+0.91	+0.97	+1.04	+1.12	+1.20	+1.29	+1.39	+1.49	+1.58	+1.68	+1.78	+1.88
15°C	+0.73	+0.77	+0.83	+1.89	+0.95	+1.02	+1.09	+1.16	+1.24	+1.32	+1.40	+1.48	+1.56
16°C	+0.60	+0.63	+0.67	+0.72	+0.77	+0.82	+0.88	+0.94	+1.00	+1.06	+1.12	+1.19	+1.25
17°C	+0.46	+0.48	+0.51	+0.55	+0.59	+0.62	+0.67	+0.71	+0.75	+0.80	+0.84	+0.89	+0.94
18°C	+0.31	+0.33	+0.35	+0.37	+0.40	+0.42	+0.45	+0.48	+0.51	+0.53	+0.56	+0.59	+0.62
19°C	+0.16	+0.17	+0.18	+0.19	+0.20	+0.21	+0.23	+0.24	+0.25	+0.27	+0.28	+0.30	+0.31
20°C	+0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00
21°C	-0.17	-0.18	-0.19	-0.19	-0.20	-0.22	-0.23	-0.25	-0.26	-0.28	-0.29	-0.30	-0.31
22°C	-0.34	-0.36	-0.37	-0.39	-0.41	-0.44	-0.47	-0.49	-0.52	-0.55	-0.57	-0.60	-0.62
23°C	-0.51	-0.54	-0.57	-0.60	-0.63	-0.66	-0.70	-0.74	-0.78	-0.82	-0.86	-0.90	-0.93
24°C	-0.70	-0.73	-0.77	-0.81	-0.85	-0.89	-0.94	-0.99	-1.04	-1.10	-1.15	-1.20	-1.25
25°C	-0.89	-0.93	-0.97	-1.02	-1.07	-1.13	-1.19	-1.25	-1.31	-1.37	-1.43	-1.49	-1.56
26°C	-1.08	-1.13	-1.18	-1.24	-1.30	-1.36	-1.43	-1.50	-1.57	-1.65	-1.73	1.80	-1.87

Table 1.: Alcohol degree correction according to temperature.

11.2 Scale + Immerser pending


11.2.1 Material

- 100 ml Flask
- Immerser with its conversion table.
- Thermometer 0-35°C resolution 0.5°C
- Electronic balance

11.2.2 Procedure

Before start, check that ambient temperature is stable inside $20 \pm 3^\circ\text{C}$ and the distilled liquid temperature is stable inside $20 \pm 0.5^\circ\text{C}$.

- Put the distilled in the 100ml flask.
- Insert the thermometer.
- Wait 1 minute and read the temperature.
- Remove the thermometer. Insert completely the immerser.
- Wait 1 minute. Take a read from electronic balance and calculate the %Vol from the conversion table supplied with the immerser.
- Repeat the reading, at least 3 times and average. Write the result, only with 2 tenth. Round second tenth at 0 or 5.
- If the measured temperature, once round it is not 20°C (p.e. 20.4°C is round at 20°C ; 20.6°C is round at 21°C) the %Vol must be corrected according table 1.

11.3 ABBE Refractometer

11.3.1 Material

- Abbe Refractometer range 1.330 to 1.346.

11.3.2 Procedure

There are 2 types of refractometers:

- Without temperature measurement.
- With temperature measurement.

On the first case, before using the refractometer and measuring the %Vol, remain refractometer and the sample in a chamber at $20 \pm 0.5^\circ\text{C}$ for a 4h period. Then measure it in a room at $20 \pm 3^\circ\text{C}$.

On the second case probably, the refractometer has the capacity to climate the measuring area and the sample at 20°C . Before reading the refraction index, check that temperature is in $20 \pm 0.5^\circ\text{C}$.

Follow the refractometer instructions sheet to read the refraction index.

By means of table 2, look for the %Vol by introducing the refraction index.

Remember to check the refractometer calibration.

The presence of methanol on the distilled produces a lower refraction index in front of a pure hydroalcoholic solution (Ethanol - Water). The hydroalcoholic solutions are used for calibration standards.

Table 2. Correspondence between the refractive index at 20°C and the GAV

Index of refraction at 20°C	Alcohol levels at 20°C		Index of refraction at 20°C	Alcohol levels at 20°C	
	Hydroalcoholic mixtures	Distillates		Hydroalcoholic mixtures	Distillates
1,33628	6,54	6,48	1,34222	16,76	16,65
1,33642	6,79	6,74	1,34236	16,99	16,88
1,33656	7,05	7	1,3425	17,22	17,12
1,3367	7,3	7,27	1,34264	17,44	17,34
1,33685	7,58	7,54	1,34278	17,68	17,56
1,33699	7,83	7,79	1,34291	17,89	17,78
1,33713	8,09	8,05	1,34305	18,12	18,01
1,33727	8,34	8,3	1,34319	18,36	18,23
1,33742	8,62	8,56	1,34333	18,59	18,46
1,33756	8,87	8,81	1,34347	18,82	18,7
1,3377	9,12	9,06	1,34361	19,05	18,92
1,33784	9,36	9,3	1,34375	19,28	19,17
1,33799	9,63	9,55	1,34389	19,51	19,4
1,33813	9,87	9,81	1,34403	19,75	19,62
1,33827	10,12	10,05	1,34417	19,98	19,86
1,33841	10,35	10,29	1,34431	20,22	20,09
1,33856	10,61	10,54	1,34445	20,44	20,33
1,3387	10,86	10,78	1,34458	20,65	20,54
1,33884	11,1	11,02	1,34472	20,89	20,76
1,33898	11,33	11,26	1,34486	21,11	20,99
1,33912	11,47	11,5	1,345	21,34	21,21
1,33926	11,81	11,74	1,34513	21,55	21,44
1,3394	12,05	11,98	1,34527	21,78	21,65
1,33955	12,3	12,22	1,34541	22	21,87
1,33969	12,53	12,46	1,34555	22,23	22,1
1,33983	12,76	12,69	1,34568	22,44	22,31
1,33997	13	12,92	1,34582	22,67	22,54
1,34011	13,23	13,15	1,34596	22,9	22,75
1,34025	13,47	13,4	1,3461	23,13	22,96
1,34039	13,7	13,62	1,34623	23,33	23,17
1,34053	13,93	13,86	1,34637	23,57	23,4
1,34067	14,16	14,09	1,34651	23,81	23,61
1,34081	14,41	14,32	1,34665	24,04	23,85
1,34096	14,66	14,57	1,34678	24,26	24,09
1,3411	14,89	14,81	1,34692	24,48	24,31
1,34124	15,13	15,06	1,34706	24,72	24,56
1,34138	15,36	15,28	1,3472	24,95	24,78
1,34152	15,59	15,5	1,34733	25,16	25
1,34166	15,83	15,74	1,34747	25,4	25,23
1,3418	16,06	15,96	1,3476	25,62	25,45
1,34194	16,29	16,19	1,34774	25,86	25,7
1,34208	16,52	16,41	1,34788	26,1	25,93

12 Volatile acidity determination

12.1 Equipment

Distilling unit DE-1626
Tritation Burette
Kitasatos flask or Ultrasonic bath
Vacuum pump

12.2 Reagents

Saturated solution of borax
Potassium iodide crystallized
Hydrochloric acid (d = 1,8 ó 1,9 g/l)
Tartaric acid crystallized
Sodium hidroxide 0,05N
Fenolftaleine 1%
Iodyne solution 0,005 M (0,01N)
Solution Starch paste 5g/l

12.3 Sample preparation

Remove the carbonic gas bubbles from sample by using ultrasonic bath or kitasatos flask and vacuum pump.

12.4 Distillation

Use the distillation tube 250ml (small one) and add 20 ml sample on it.

Add about 0,5 g of tartaric acid and start the steam generator.

Distill up to 250 ml on the erlenmeyer flask.

12.5 Tritation

Tritate distilled with sodium hidroxide 0.1M with two drops of fenolftaleine. Consumption will be «V».

12.6 Free SO₂ Tritation

Add on the 12.5 tritated solution, 4 drops of hidrochloric acid, and 2 ml de starch paste and some glasses of Potassium iodide.

Tritate the free SO₂ with 0.005 M Iodyne. Consumption will be «V'».

12.7 Combined SO₂ tritation

Add the Saturated solution of borax (Around 20ml) up to a pink solution appears. Tritate with 0.005 M Potasium Iodyne. Consumption will be «V''».

12.8 Calculation

V = NaOH consumption (ml) on tritation paragraph 12.5

V' = Iodyne 0.005M consumption on tritation paragraph 12.6

V'' = Iodyne 0.005M consumption on tritation paragraph 12.7

$$\text{Volatile acidity (meq/l)} = 5 \times (V - V'/10 - V''/20)$$

$$\text{Volatile acidity (g/l sulphuric acid)} = 0.49/2 \times (V - V'/10 - V''/20)$$

$$\text{Volatile acidity (g/l acetic acid)} = 0,3 \times (V - V'/10 - V''/20)$$

13 Determination of sorbic acid (SO₂)

(Method of determination by spectrophotometry of ultraviolet absorption)

Principle:

Sorbic acid separated by steam dragged distilling is determined in the distillate by spectrophotometry of ultra violet absorption. The substances which interfere in the measure of the absorption in the ultra violet are eliminated by evaporation to dryness of the distillate of the sample, slightly alkalinized by a solution of calcium hydroxide.

Necessary reagents:

- Crystallized tartaric acid.
- Solution of hydroxide of calcium Ca(OH)₂, approximately 0.02M.
- Reference solution of sorbic acid C₆H₈O₂ of 20 mg per litre: Dissolve 20 mg of sorbic acid in approximately 2 ml of solution 0.1M of sodium hydroxide. Pour into a calibrated flask of 1000 ml, level with water. 26,8 mg of potassium sorbate C₆H₇KO₂ can be dissolved in water and filled with water up to 1000 ml.

Necessary material:

- DE-1626 distilling unit by steam dragging.
- Water bath at 100°C.
- Zirconium or porcelain capsule of 55 mm in diameter.
- Spectrophotometer permitting the measuring a wavelength of 256 nanometers with 1 cm quartz tanks.

Procedure:

DISTILLATION

1- Place 10 ml of wine into the small bubble tube, add 1 or 2 g of tartaric acid. Put it in the DE-1626 distilling unit and collect 250 ml of distillate in the gauged flask.

PATTERN CURVE

2- By dilutions, prepare with water - from the reference solution - four pattern solutions which contain 0.5, 1, 2.5 and 5 mg of sorbic acid per litre respectively.

3- With a spectrophotometer, measure the respective absorbance at 256 nm of pattern solutions with regards to the distilled water.

4- Trace the absorbance curve regarding the concentration of solutions. The relation is lineal.

DETERMINATION

5- Introduce 5 ml of distillate in a capsule of 55 mm diameter, add 1 ml of calcium hydroxide solution and a drop of copper sulphate solution. Evaporate to dryness in a bath of boiling water.

6- Collect the residuum with some millilitres of distilled water, drag quantitative in a 20 ml gauged flask and drag with the rinsing water, measure the absorbance at 256 nm in comparison to witness solution obtained with 1 ml of solution of hydroxide of calcium and diluted at 20 ml with water.

7- Take the value of the measured absorbance and subtract the **C** concentration of the solution in sorbic acid.

Note: In current practice, this evaporation or dryness may not be necessary. Measure the absorbance directly in the diluted distillated at 1/4 against the distilled water.

EXPRESSION OF THE RESULTS

The concentration in sorbic acid of the wine expressed in milligrams per litre equals: $100 \times C$

C = concentration in sorbic acid of the analysed solution by spectrophotometry, expressed in milligrams per litre.

14 Guidelines for IQ and OQ

Some laboratories need the installer makes an IQ document to certificate the good installation of this equipment. As well as an OQ document to certificate the good working of this equipment.

The installer can elaborate those documents following the examples:

Checking list for IQ (Installation Qualification):

- Installer name.
- Manufacturer name.
- Unit identification: Code, model and serial number.
- Unit location: Laboratory, Person in charge.
- Checking packing list:
- Unit location match the requirements?
Water and electricity requirements are suitable?
- Electrical supply.
Socket:
 - Ground terminal?
 - Suitable to power consumption. (16A)
 - Voltage according plate $\pm 10\%$
- Water connection.
 - Water pressure bigger than 1bar.
 - Hose to waste are properly fixed and not blocked.
- Water for steam generator: according to chapter 4.4
- Personnel training: Operator personnel understand the instruction manual?.
- Date and name of who has made the Installation (Signed).

Checking list for OQ (Operation Qualification):

- Installer name.
- Manufacturer name.
- Unit identification: Code, model and serial number.
- Unit location: Laboratory, Person in charge.
- Check the equipment distillate correctly. See: 7.1
- Check alcohol recovery by means of distillate hidroalcoholic solutions prepared with volumetrics starting from Etanol and distilled water, or distilling certified standards.

15 Warranty

The warranty covers the defective components or faulty assembling of the unit for 1 year since the purchasing date.

Keep purchasing bill as the only valid document for the purchasing date.

The warranty does not cover damages caused for a unit malfunction. The unit must be suitably qualified before being used for valuable samples processing.

Return goods are not accepted without previous J.P. SELECTA, s.a.u. approval.

16 Troubleshooting table

- **The green indicator «ON» does not light.**

Check fuses and wiring connection cables.

Check main switch.

- **The equipment still adding water till it breaks out off the sample.**

The equipment do not detect the water level for low conductivity. Add about 250ml of tap water in the water tank.

Check the cable (white) of the level probe of the steam generator.

- **At the end of distillation, the sample tube content is absorbed to the steam generator.**

Check the steam hose is not obstructed.

Locate the aerated electrovalve:

Check it actives, check the voltage with a tester (230V).

Check (detaching it) it is not obstructed.

- **The refrigerating water does not enter.**

Check the tap is open.

Locate the refrigerating water input valve, check the voltage with a tester (230V). If there is voltage and the valve does not open, hit it. That is due to a small lime rest blocking the water pass. If hitting is not enough, detach the valve to clean the water pass.

- **Security thermostat Alarm warned and the equipment is blocked.**

This alarm warns when there is an excessive temperature in the steam generator. In this case, the heating element stays disconnected by the thermostat and the control electronics blocks all the unit functions.

To rearm the thermostat, wait some minutes to let the temperature of the generator fall down to a more normal level, and pulse the rearm switch.

The cause of the thermostat disconnection is that the steam generator is lack of water. Check the hydraulic wiring to generate water, tank, pump, hoses, level probe.

- **Distillation does not stops automatically.**

Check that collector flask is properly placed on his holder (flask wall must be in contact with sensor face).

Check, from inside the unit, that the sensor detects the water level, this situation is indicated by a yellow light.

Adjust the sensibility.