DESTILADOR KJELDAHL MANUAL

MANUAL KJELDAHL DISTILLER

PRO - NITRO M


4002627

Indice

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Index

1 Packing list & accessories	27
1.1 Standard package content.....	27
1.2 Accessories	27
2 Safety	28
2.1 Safety symbols	28
2.2 Risks which the operator is subjected to.....	29
2.3 Personnel training.....	29
2.4 Proper uses of the equipment.....	29
2.5 Modifications.....	29
3 Working principle overview	30
4 Installation	31
4.1 Suitable location	31
4.2 Connection to power supply	31
4.3 Connection to tap water and waste.....	31
4.4 Water for steam generator	31
4.5 Connection to automatic tritator	31
5 Operation.....	32
5.1 Insert/ Remove the sample tube.....	32
5.2 Adding NaOH	32
5.3 Distillation.....	33
5.4 Alarms.....	33
6 Sample distillation.....	34
6.1 Check levels of tanks.....	34
6.2 Unit warm-up	34
6.3 Distillation	35
6.4 Titration.....	35
6.5 Protein contents (%) calculation.....	35
7 Reagents	36
7.1 Prepared reagents	36
7.2 Proper reagents set up.....	36
8 Checking the distiller (OQ).....	37
8.1 Functional checking	37
8.2 Checking Nitrogen Recovery.....	37
8.3 Blank test.....	37
8.4 Checking Nitrogen recovery with Acetanilide.....	38
9 Maintenance	39
9.1 Cleaning the unit:	39
9.2 Hose inspection (weekly)	39
9.3 Preparing the unit for shipping or storing.....	40
9.4 Level sensor cleaning.....	41
10 Technical specification	42
11 Technical service	43
11.1 Spare parts	43
11.2 Wiring diagram.....	44
11.3 Hydraulic wiring diagram.....	45
12 Kjeldahl analysis step by step	46
12.1 Sample preparation.....	46
12.2 Digestion.....	46
12.3 Dilution	46
12.4 Distillation	46
12.5 Titration and calculations.....	46
13 Guidelines for IQ y OQ.....	47
14 Warranty	47
15 Troubleshooting table	48

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
Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.

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
1 Packing list & accessories

1.1 Standard package content

- Distilling unit
- MACRO size sample tube Ø 42 (6)
- Tray
- Tube 1.5m incoming water. High pressure, coupling 3/8" (1)
- Tube 1.5m water outlet. Silicone 8x14, coupling 3/8". (2)
- Tube 1.5m remaining steam. Silicone 8x14, coupling 3/8". (3)
- Tube 1.5m generator emptying. Silicone 8x14. (4)
- Power supply wire (5)
- 2 litre tank for NaOH
- Instructions Manual
- Flask 100g of Ammonia sulphate. (7)

1.2 Accessories

Tubes, racks and tongs:

- MACRO sample tube Ø 42 4042300
- MICRO sample tube Ø 26 4001045

Racks:

- Rack MACRO 6 tubes. 4005071
- Rack MACRO 12 tubes. 4005081
- Rack MACRO 20 tubes. 4005091
- Rack MICRO 12 tubes. 4001053
- Rack MICRO 24 tubes. 4001054
- Rack MICRO 40 tubes. 4001055

Safety goods:

- Tongs for tubes 7000532
- Protection gloves 1001255
- Anti-acid gloves 5000044
- Safety glasses 5028192

Complete digestion sets

- Digestion unit MACRO 6 tubes. 4000629
- Digestion unit MACRO 12 tubes. 4000630
- Digestion unit MACRO 20 tubes. 4000631
- Digestion unit MICRO 12 tubes. 4001047
- Digestion unit MICRO 24 tubes. 4001048
- Digestion unit MICRO 40 tubes. 4001049

Smokes removal aids:

- Scrubber. 4001611
- Water circulation vacuum pump. 4001612

2 Safety

PRO-NITRO "M" unit has the suitable measures for its use in laboratory. However, the nitrogen analysis by the Kjeldahl method can be hazardous and the user must be informed of some of those risks.

In order to be respected, across this manual, the risk of danger situations are remarked.

2.1 Safety symbols

It identify the risk and the safety actions to be applied.

The symbols make reference to the paragraph marked with a grey line.

Warning

Warning of a dangerous operation.

Indicated instruction must be followed.

Hazard of contact with corrosive substances

Hazard of contact with sodium hydroxide and strong acid, substances that attack skin, clothes and other materials easily.

Follow the instructions for the describe operation.

Hazard of electrical shock

Electric hazard for entering in areas marked with this signal, or for following the manual instruction signaled by this icon.

Follow the instructions for the describe operation.

Hazard of burns to be in contact with hot surfaces

Temperature in this area can exceed 60°C. Use heat protected gloves for the described action.

Follow the instructions for the describe operation.

Wear gloves mandatory

The use of gloves is obligatory to run the operation.

In the section «Accessories» an advised model is indicated.

Mandatory use of safety glasses

The use of safety glasses or complete face protection is obligatory to conduct the described operation.

In the section «Accessories» an advised model is indicated.

Important information

- Information important to obtain good results or for an optimal operation of the equipment.
- Information important to extend the life of the equipment or to avoid the degradation of some of its components.

2.2 Risks which the operator is subjected to

The process of determination of Kjeldahl Nitrogen with the described equipment exposes to the operator to the following situations of risk:

- Manipulation of concentrated sulphuric acids and Hydrochloric.
- Sodium manipulation Hydroxide 40%.
- Manipulation of glass parts.
- Possibility of touching parts to a temperature above to 60°C
- Possible exhibition to irritating acid steam.
- Electrical risk.

2.3 Personnel training

This equipment must only operated by personnel trained a suitable qualified to work in a chemical laboratory and its common hazards.

This equipment must only operated by personnel who have read and understood this manual or have been qualified on its manipulation and operation.

2.4 Proper uses of the equipment

This equipment is for laboratory use, for sample distillation by steam pulling. Mainly in Kjeldahl nitrogen analysis.

The improper use of this unit (not following manual instructions) can endanger the protection prepared by the equipment.

Do not use to distillate explosive or flammable substances.

This equipment is not designed and not manufactured to work in an explosive ambient area.

2.5 Modifications

The operation modification or the manipulation of the safety systems unauthorized by the manufacturer, could expose the operator to hazards not indicated in this manual.

3 Working principle overview

The PRO-NITRO M is a steam distillation unit, designed, mainly for, protein determination according to direct Kjeldahl Nitrogen method.

It provides a sample steam distillate and an alkali adding in safety way.

The unit has the following systems:

The steam generator

Generates steam from the water reservoir placed at the top of the unit.

The steam generator includes two safety devices:

- Safety thermostat: Cuts the electrical power to the heater when overtemperature is detected on the steam generator. (Probably for poor level of water)
- Pressure switch: Cuts the electrical power to the heater when there is an overpressure on the steam generator.

Also includes a conductivity sensor to keep a constant level of water in the steam generator.


The distillation circuit



It is made, mainly, by two glassware parts:

- The splash protector: Its special design avoids that NaOH splash reach the cooler glassware.
- The cooler: It condensates the steams that are leaking from the sample. This glassware is cooled by a water circuit that comes from an external tap. An internal valve opens and closes this circuit automatically.

NaOH addition system

Composed by a properly specified pump to dispense NaOH directly to sample tube when distillation is running.

4 Installation

Once the unit is unpacked, visually check, the packing list (See 1.1) and that there are no parts damaged during transportation. Save the packing for some days for later sending.

4.1 Suitable location

Place the unit over a flat level and stable surface, suitable for the unit weight.

No longer that 1.5m from a water tap, waste output and power supply plug.

To work comfortably, let a free space of 50cm on each side and not least than 20cm on the rear part.

Prepare to refill the water reservoir from the top of the unit.

4.2 Connection to power supply

Locate a power supply socket (Rated to 16A) near to the unit.

Use the power cord supplied or another of similar specification.

For safety reason the socket must have ground terminal, since the unit must work grounded.

Before to plug on the socket check that its voltage match the voltage marked on the unit plate.

4.3 Connection to tap water and waste


Connect the water entry hose to a water tap (1). Connect the water output hose from (2) to waste.

The hoses «remaining steam» (3) and «Steam generator empty» (4) must go to waste.

Do not forget the rubber white nuts on the hose's connections.

Press plastic connections by hand. Do not use any tool.

The tap water must have, at least, 1 bar pressure.

Note: It is critical in order that the unit works properly, the hose «Remaining steam» (3) MUST NOT BE IMMERSed under water, or blended or blocked by any reason.

4.4 Water for steam generator

The unit includes a 6 litres water reservoir, at the top, to store the water for the steam generator.

As the system used to control the level of water on the steam generator works by conductivity, water of at least 20 microSiemens must be used. Water also must have a low content of calcium to increase the life of the steam generator.

Water with this conductivity could be obtained from a demineralised unit. If conductivity is low, add near 50mg of salt (NaCl o Na₂CO₃) for 10L. of water.

4.5 Connection to automatic tritator

The unit includes a TTL level output to connect it to an automatic tritator as the DL22 from METTLER TOLEDO.

Once the distillation step finishes the PRO-NITRO «M» sends a signal to the tritator in order to initiate the tritation. See chapter 10.

Moment indicator : «ON»

5 Operation

The control panel allows function selection and shows the distiller alarms status.

Once the unit is plugged in, put it on, by pressing the main switch. The unit goes to stand by status (ON).

Note: If main switch is pushed consecutively with an interval of less than 3 seconds, the unit could start in abnormal mode. So wait few seconds to restart the unit.

5.1 Insert/ Remove the sample tube

- Insert the tube starting to introduce the PTFE tube.
- Push down the fixing lever.
- Close the safety door.

Proceed similar to remove the sample tube.

CAUTION

Wear Anti-acid, Anti-thermal gloves.

**The white PTFE tube could have rest of Sodium hydroxide or Acid.
The sample tube's temperature is over 60°C after distillation.**

NaOH preselection

5.2 Adding NaOH

There are two modes for the NaOH adding function:

Note: The NaOH adding function is inhibited when there is no sample tube.

Automatic addition

On this mode, the NaOH volume is pre-selected before to start distillation.

The pre-selected volume will be dispensed, automatically, few seconds after distillation starts.

With the unit on stand-by status (Not distilling) keep pressed the «NaOH» key until it makes a different number of beeps that indicates the NaOH volume to be added. Release the key once you hear the number of beeps that you want.

Approximately each beep means 25ml.

Note: The NaOH volume set, is stored, only, until the unit is ON. Volume must be set again if the unit is switched off.

Note: The NaOH volume is set by "working time" of dosing pump. This means that the volume is only approximated, and could change with temperature, NaOH concentration and pump life time.


Manual NaOH
adding during distillation

Start / Stop distillation

To change the distillation time, keep pressed the "STEAM" key during 5 seconds:

- 1 beep: time is 6 minutes.
- 2 beeps: time is 8 minutes.

Alarm signals

Manual addition

During distillation NaOH could be added on the sample. For it, just push the «NaOH» key. A volume of around 25ml will be dispensed.

Note: Wait until the dispensing is finished before to push «NaOH» for another dispensing.

5.3 Distillation

On pushing the «STEAM» key the unit starts the sample distillation.

When an automatic addition of a NaOH volume has been set, it starts after a few seconds.

Approximately after 30s the sample starts to boil and the steams are driven to the cooler for its condensation. The distilled is collected on the Erlenmeyer.

Generally speaking it is accepted that once 150ml of distilled are collected distillation will be finished.

(50ml of Boric + 150ml of distilled; **Total:** 200ml)

Push «STEAM» to stop the distillation.

Note: Distillation stops automatically after 8 minutes (approx. 200ml of distilled). But it is suggested not to leave the unit unattended.

Note: By Keep pushed the key "STEAM" for 5 seconds, the distillation time is modified from 8 to 6 minutes. (It must be selected again when the unit is disconnected)

5.4 Alarms

To avoid dangerous situations for the operator or to protect the unit life, it includes some safety systems and some alarms:

(1) Alarm: «No water for steam generator»

- If the lamp (1) lights on sporadically : This means that water level control have detected a low level of water on the steam generator and automatically starts the water pump to get the suitable level of water.
- If the lamp (1) lights on permanently: This means that water level control have detected a low level of water on the steam generator for a long time and it automatically refills with water. This situation stops distillation. Probably the water reservoir is empty

(2) Alarm: «Safety door open or tube not present»

This red lamp lights on when the safety door is open or if there is no sample tube inserted on the unit.

This alarm inhibits the distillation and the NaOH adding.

(3) Alarm: «Overtemperature» on the steam generator

This red lamp lights on when an overtemperature on the steam generator has switched the safety thermostat.

It inhibits all distiller functions. And cuts the power supply to the heating element.

Once it is switched, it must be left to get cool for some minutes. Solve the problem that causes the switching and rearm the safety thermostat by pushing the central button. The safety thermostat is located on the electrical box, right hand side cover must be removed to get it.

The safety thermostat switching is due to low level of water on the steam generator, probably because the water reservoir is empty or due to water pump malfunction.

6 Sample distillation

To get good analysis' results and a good performance of the distiller it is suggested to follow a working routine.

6.1 Check levels of tanks

Before starting a batch analysis, check the following levels:

Steam generator water tank level

Refill it with demineralised water.

If the conductivity is low, add about 50mg of Salt (NaCl or Na₂CO₃) per 10 L. of water.

The water reservoir is 6 litres capacity, enough for analysing more than 20 samples.

Level of NaOH

Refill with NaOH 40% from the 2 litres flask. It is placed inside the unit.

Get the NaOH flask from the door placed on the left hand side cover.

The NaOH solution must be free of ammonia salts, to avoid its interference on the result.

Boric acid preparation

Prepare enough quantity of Boric acid 4% with mixing of indicators. Remember that a volume of 50 ml will be used on each distillation.

Tritation reagent preparation

Prepare enough quantity of tritration reagent, usually HCl or H₂SO₄ 0.25 Normality.

Note: See chapter 7 about the reagents.

Note: The precision of the results are affected directly by the precision of the acid normality used for tritration. **This is a common cause of error**, because, along the time, opened bottles tends to increase its normality.


Also, often there are errors in concentration calculation to prepare the acid for tritration.

6.2 Unit warm-up.

In order to get the best results with your Kjeldahl distiller, you have to run a warm-up step before starting samples' analysis. For it just insert a sample tube with 25ml of clean water and let distillate around 5 minutes. (Put an Erlenmeyer flask on the cooler output)

Run also this warm up, if the unit has been 2 or 3 hour without working.

Check the tube output is immersed under Boric's level.

6.3 Distillation.

Please follow these steps (on the same order as indicated):

- Open the safety door.
- Place an Erlenmeyer flask with around 25ml with some drops of indicator mix under the cooler output. **Check the tube output is immersed under the Boric's level.**
- Insert the sample tube.
- Pre-select the NaOH volume.
- Push «START»
- If necessary, more NaOH could be added to the sample by pressing «NaOH» button.

Visual Checking: When the NaOH volume is sufficient, the sample becomes blue. (If not, add more NaOH)

- Once 200ml of distillate are obtained on the Erlenmeyer (50ml Boric + 150ml of distilled press «START» button to stop the distillation.
- Open the safety door.
- Remove, first, the sample tube.
- Remove the Erlenmeyer. (wait drip)

6.4 Titration

Prepare the titration burette with the reagent usually HCl 0.25N.

Slowly dispense the HCl 0.25N on the Erlenmeyer stirring it at the same time.

Titrate until the colour changes (end point at pH 4.65).

6.5 Protein contents (%) calculation.

From the HCl volume (ml) consumption on the titration the Kjeldahl Nitrogen (mg) is calculated:

$$\text{Nitrogen (mg)} = 14 \times \text{Volume}_{\text{HCl}} \text{ (ml)} \times \text{Normality}_{\text{HCl}}$$

And, from the nitrogen we calculate the sample protein content:

Main protein factors

Almonds	5,18
Nuts	5,30
Nut-peanuts, nuts Brazil	5,41
Jello	5,55
Soya	5,71
Cebada, avena, centeno	5,83
wheat, whole flour	5,83
Flours (except the whole)	5,70
Rice	5,95
Corn	6,25
All the other foods	6,25
Saved	6,31
Milk and milky products	6,38

$$\% \text{ Protein} = \frac{P_2}{P_0} \times 100 \times F$$

P2: Nitrogen (mg).

P0: Sample weight (mg).

F: Protein factor.

7 Reagents

The reagents to be used must be free of ammonia.

7.1 Prepared reagents

The best thing is to use reagents factory prepared, mainly the HCl for titration. Any error on its preparation have a direct influence on the result.

As an example the following reagents could be used, or the equivalent one in other brands:

Boric acid 4% (For ammonia titration)	RV	PANREAC 282222
--	----	----------------

Mixed indicator 5 (Red of metile + Green Bromocresol)	RV	PANREAC 283303
--	----	----------------

Mixed indicator 4,4 (Red of metile + Blue of metilene)	RV	PANREAC 282430
---	----	----------------

HCl 0.1N		PANREAC 171023
----------	--	----------------

HCl 0.25N	SV	PANREAC 182318
-----------	----	----------------

H ₂ SO ₄ 0.1N	SV	PANREAC 181061
-------------------------------------	----	----------------

H ₂ SO ₄ 0.2N	SV	PANREAC 182011
-------------------------------------	----	----------------

Acetanilide 99% (Validation standard)		PANREAC 151005
---------------------------------------	--	----------------

Ammonia sulphate (Validation standard)		PANREAC 131140
--	--	----------------

Kjeldahl Catalyst 6.25% Cu tabl. 8g		PANREAC 174428
-------------------------------------	--	----------------

7.2 Proper reagents set up

As it is usual on most laboratories, the reagents to work with PRO-NITRO «M» could be prepared starting from concentrated ones.

The NaOH 40%, acid Boric 4% are not critical on the concentration and it is not necessary to work with precision.

Reagents for titration (HCl or H₂SO₄) must be high precision and calculation must be accurate. Any error on its preparation has a direct influence on the result.

8 Checking the distiller (OQ)

Once the unit is properly installed and the operator has acquired the knowledge to work with it. The unit must be checked.

This process is known as Operation Qualification.

There are four basics checking to do:

8.1 Functional checking

Run this checking after installation and also for some maintenance.

Follow the instructions given on the chapter 5 to distillate for 5 minutes.

Verify the distillation rate. When the unit is at steady state distillation rate is between 30 - 40 ml per minute.

8.2 Checking Nitrogen Recovery.

This checking is widely used to certify that the unit works properly. Run some of these tests for training purposes.

It consists to prepare some samples that have a known Nitrogen content. It is distilled, titrated and the nitrogen detected is calculated.

The rate (%) of nitrogen detected versus the sample Nitrogen content gives us the **Nitrogen Recovery**.

To allow this checking, a 100gr ammonia sulphate flask is delivered together with the distiller.

- Weight around 100mg (50mg reagent 0.1N) of ammonia sulphate. The exact weight will be P_0
- The exact quantity of Nitrogen is: $P_1 \text{ (mg)} = P_0 \times 0,212$
- Distillate adding 25ml of NaOH
- Once distilling and titration, see 6.5, we obtain the detected nitrogen (mg) : P_2

$$P_2 \text{ (mg)} = 14 \times \text{Volume}_{\text{HCl}} \text{ (ml)} \times \text{Normality}_{\text{HCl}}$$

- The recovery is calculated: $R(\%) = P_2 / P_1 * 100$
- Acceptable values are from 99.0 to 101.5 %.

Note: Some literature sets that the acceptable values for recovery must be above 99.5%. To get this result you must work with high precision and use reagents with high degree of purity.

8.3 Blank test.

This checking must be run on each batch in order to quantify the absence of ammonia on the reagents.

It consists on distilling and titrating a sample of 25ml of distilled water.

Add 15ml of NaOH on the distillation.

Run some blank tests in order to get a constant value of blank. If this value is repeated it could be subtracted from the detected Nitrogen.

8.4 Checking Nitrogen recovery with Acetanilide.

This test is widely used to certify the whole Kjeldahl analysis process, including the digestion, distillation and titration steps.

It consists in preparing some samples that have a known Nitrogen content. The samples are digested, distilled, titrated and the nitrogen detected is calculated.

The rate (%) of nitrogen detected versus the sample Nitrogen content gives us the **Nitrogen Recovery**.

To run this test Acetanilide is used.

Preparing the sample:

Nitrogen content of an acetanilide sample:

Formulae: C_8H_9NO

Molecular weight: 135.17

Factor: $14 / 135.17 = 0.1035$ mg of Nitrogen = $0.1035 * \text{mg Acetanilide}$.

- Weight around 250mg (100mg) of acetanilide. The exact weight P_0
- The Nitrogen content is: $P_1 \text{ (mg)} = P_0 \times 0,1035$

Sample digestion:

- Place the acetanilide on a MACRO tube, add 10ml of Sulphuric acid 98% (Example PANREAC 173163) and one tablet of Kjeldahl catalyser.

- Digest at 400°C for 1h. (The result must be a transparent blue liquid)

- Let cool at ambient temperature. Add 25ml of distilled water on each tube. Operate with caution to avoid splash. Water reaction over the sulphuric acid could be violent.

Distillation:

- Distillate by adding 75ml of NaOH. Use HCl 0.25N (0.1N) for titration.

- Once sample is distilled and titrated following chapter 6, the detected Nitrogen is calculated as: P_2

- Recovery is calculated as: $R(\%) = P_2 / P_1 * 100$

- Acceptable recovery is from 99.0 to 101.5 %.

Note:

If MICRO size tube is used, weight 0.3g of acetanilide, add 5ml of sulphuric acid, and one tablet of Kjeldahl catalyser. On the dilution add 10ml water.

Distillate for 6 minutes adding 25ml of NaOH.

Note:

Some literature sets that the acceptable values for recovery must be above 99.5%. To get these results you must work with high precision and use reagents with high degree of purity.

Clean regularly rest of NaOH between the closing joint and the splash protector.

Clean periodically the NaOH tank hose connector.

9 Maintenance

The distiller PRO-NITRO M uses NaOH and concentrated acids for the analysis process. These substances are very aggressive and will degrade some parts of the unit.

To avoid derating and increase distiller's life a suitable cleaning routine should be done.

9.1 Cleaning the unit:

9.1.1 Clean the distillation circuit:

After any session of analysis, distillate around 50ml of water for 5 minutes in order to remove any rest of NaOH on the distillation glassware.

9.1.2 Cleaning the chamber:

The chamber is protected with a chemical resistant coating. But it should be cleaned, daily, to increase its life.

Use a well damp cloth to clean the chamber and its components to remove any rest of NaOH.

Clean the tray with abundant water.

9.1.3 Cleaning of closing joint.

Near each three months clean, with water, the possible rest of NaOH observed on the closing joint.

Extract the closing pulling it and turning it downwards.

9.1.4 Cleaning the NaOH hose connection.

Weekly clean with a wet cloth the NaOH tank hose connection.

9.1.5 Cleaning the steam generator (weekly)

Switch off the equipment. Open the «Steam generator drainage valve» and wait till it will be completely empty. Close the valve and start the equipment. See how the steam generator is filled in chapter 4 of this manual.

9.2 Hose inspection (weekly)

Visually check the hose condition and its connections to components.

To avoid premature derating of the unit due to liquid leakage, replace as soon as possible any hose that have any leakage or damage.

9.2.1 Hose replacing:

Take these cautions to hose replacing:

- Unplug the unit from the electrical power supply.
- Try, as you can, (See 9.3) that the hose to remove is empty.
- Use only tubes material as indicated on 11.2.
- Use clamps as is factory assembled.
- Once the hose is replaced, it must be checked for leakages.

9.3 Preparing the unit for shipping or storing.

Once the unit has been used with NaOH, proper actions should be taken for its transportation or storing.

Initially clean the unit according to chapter 9.1

9.3.1 Remove completely the NaOH from the unit:

- **Clean the NaOH tank:**

Remove the NaOH 2 litres tank from the unit.
Clean it with abundant water for some minutes.

- **Cleaning the NaOH dosing circuit:**

Fill the NaOH 2 litres tank with water.
Connect to NaOH circuit.

Insert a MACRO tube on the unit. Push «STEAM» then, push «NaOH» several times. Remove the MACRO tube when filled. Insert an empty one up the 2 litres of water have circulated across the NaOH dosing circuit.

CAUTION

**The NaOH is very corrosive for the skin
Use chemical resistant gloves.
Take proper cautions procedures for NaOH handling.**

9.3.2 Remove completely the water from the unit:

- **Emptying the water reservoir and the steam generator:**

1. Open the tap (4).

2. Start the unit, push «STEAM» and then wait that the water pump empties the water reservoir.

- **Emptying the water from the cooling circuit.**

1. Disconnect the hose for water **entry** from the tap water.

2. Switch ON the unit and push «STEAM»

Note: From this point actuate in a fast way to avoid overtemperature on the steam generator.

3. Blow for the end of the **output** coolant water. The water will **come out** on the end cap of the **water in** hose.


CAUTION

Do not ship off the unit if NaOH and water have not been completely removed.

During transportation NaOH or water leakages could fall over important components and cause severe damages.

During transportation NaOH or water leakages could fall outside of the package.

Steam Generator

9.4 Level sensor cleaning


Inside the steam generator a constant water level is necessary.

To provide this water level there is a conductivity based level probe (a steel rod) that when it loses contact with the water, it shows a high resistance to the electronics and then the microprocessor orders the water pump to feed the steam generator with water from the reservoir up to the probe level.

After some time working with this unit (some months), some residues are settled on the probe causing a malfunction.

Effects of this problem:

The residues short-circuit the path between the probe and the steam generator walls (connected to earth) and so the sensor shows a LOW resistance to the electronics. It is the same signal as if the steam generator has the correct water level. This makes that the microprocessor does not supply water and this will probably cause the safety thermostat action.

Probe faston connector

Probe nut

This means that the probe must be removed and clean. You will need a cross screwdriver and a 12mm spanner and a 19mm one.

The cleaning frequency depend on two factors:


- The hardness of the water used.
- The usage of the unit.

As a general rule, please do it every **6 months**.

Probe 19mm nut

- Remove the left hand side cover (6 screws).
- Locate the steam generator.
- Unplug the probe connector (white cable with FASTON connector).
- Unscrew the sensor with the 19mm spanner.
- Unscrew the probe with a 12mm spanner.
- Remove the steel rod from the PTFE (white part).
- Clean the steel rod and the white PTFE part from the calcium deposits. Use diluted HCl or any home cleaner product for calcium deposits.
- Assemble it again.

Clean inside

Clean it

10 Technical specification

Mains voltage:	230V/50/60Hz (See plate) Other power supply voltages could be available (Ask to your dealer)
Electrical power:	1900W (See plate)
Maximum consumption:	8.5 A
Fuses:	10 A fast. dimensions: 20x5 mm
Power supply cord:	Phase, Neutral, Earth (section 1mm ²)
Weight (Net):	35 Kg
Dimensions: W x L x H	500 x 490 x 740 mm
Water reservoir volume:	6 litres
NaOH tank volume:	2 litres
Nitrogen recovery:	>99,5%
Detection range:	0.2 a 200 mg
Reproducibility: Relative Standard deviation)	< 1%
Distillation rate:	approx. 30 / minute
Short distillation time:	6 minutes, approx 150ml.
Long distillation time:	8 minutes, approx 200ml
Cooling water entry connector:	3/4"
Cooling water output connector:	3/4"
Remaining steam output connector:	3/4"
Dosification for NaOH pump:	from 0 to 100 steps of 20ml
Sound level:	Lower than 60 dBA
Ambient conditions:	Inner use. Altitude up to 2000m Temperature from 5°C to 40°C. Relative maximum damp 80%for above 31°C temperatures, decreasing lineally to 50% for 40°C.
Overvoltage category:	Category II
Pollution degree:	2
Glass material:	Borosilicate 3.3
Door and front material:	ABS injected
Carriage material:	Iron with overlap
Chamber material:	Aluminium with RILSAN coating chemical resistant.
Steam generator material:	Stainless steel.
Hosepipe material:	Silicone, synthetic fluorated elastomer.
Connection for an automatic titrating:	TTL outlet (from 1 to 10) in Pin 1 of 200ms to indicate the titrator the end of distillation. Ground in Pin5.

11 Technical service


For repair and substitution of deteriorated parts please contact J.P. SELECTA S.A. or its distributor. They will inform you about the technical services authorized in your area.

Find a spare parts list for easy buying and identifying.

11.1 Spare parts

32052

47046

47047

28129

24218

- 16087** Electrovalve N. A.
- 16067** Water electrovalve.
- 24278** KNF pump (NaOH)
- 24218** Water pump.
- 28129** Lockgum.
- 28403** Shelf.
- 32052** Steam generator 1800W.
- 47046** Splash protector.
- 47047** Refrigerant.
- 47259** Anti feedback valve.
- 29512** Relay electronic plate.
- 29519** Monitoring panel electronic plate.
- 29530** Level

Hoses:


- | | | |
|--------------|---------------------------------|------------------------------|
| 46030 | Teflon tube 8x10 | Tube immersed on the sample. |
| 46041 | Silicone tube 8x14 | Water and Steam. |
| 46042 | Silicone tube 6x10 | Water and Steam. |
| 46143 | Pressure hose | Tap water entry. |
| 46085 | ISO VERSINIC Black 12x17 | Distilling Glassware joint. |
| 46139 | ISO VERSINIC Black 6x10 | NaOH. |

WARNING

Before any work inside the equipment disconnect it from the mains.

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
12 Kjeldahl analysis step by step

12.1 Sample preparation

- Triture, homogenize and mix the sample.
- Weight between 1 and 2 grams of the sample.
- In low nitrogen content sample, (waste water, etc...) take enough samples to get 5mg nitrogen.

12.2 Digestion

- Add between 10-15ml (macrotube) of H_2SO_4 (for microtube maximum H_2SO_4 is 5ml).
- Rig up a system for smoke extraction or scrubber with NaOH.
- Make this digestion in 3 steps:

1. In function of the level of water in the sample. Begin digestion at 150°C for 20 or 60 minutes.
2. Make a second step by 280°C for 30 minutes to reduce white smoke production.
3. Continue the digestion with 400°C for 60 - 90 minutes.

Visual checking: the result is a neat transparent liquid with a clear blue coloration or green or yellow depending on the catalyser used. It must not still black rests stuck on the tube's wall.

Note: During digestion you must control foam production in the samples. If there is too much foam, length step 1.

12.3 Dilution

- Take out the sample tube off the digester block and let it cool to ambient temperature.
- Add about 25ml of water in each tube (10ml for microtube).
- Slow add water, moving the tube to not let solidify the sample. Heat the tube if necessary.(e.g. introducing the tube into the digester block)
- Let it cool again.
- To avoid nitrogen loss and violent reactions do not introduce the hot tube inside the distiller.

12.4 Distillation

- Check the tank level of NaOH and Boric acid. Check the dispensing bombs are primed and give the correct volume.
- Put a 250ml Erlenmeyer in the outlet cooler with 50ml Boric acid and some indicator drops.
- Program a 50 to 75 ml NaOH dosage (25ml for microtube).
- Insert the tube with the sample in the distiller.
- Distillate till recollecting 250ml in the Erlenmeyer (50ml Boric acid+200ml of distillate).

Visual checking: Once the NaOH has been added, the sample must get a blue coloration, if not add more NaOH.

12.5 Titration and calculations

- Titrate the distilled with HCl or H_2SO_4 0.1 or 0.25 N till colour changes. (end point: pH 4.65)
- Calculate: mg of Nitrogen = $14 * Volume_{HCl} * Normality$
- Calculate protein contents correct with the right factor based on the sample nature. (6.25 default value)
- Periodically make a blank test and subtract it from the result.

13 Guidelines for IQ y OQ

Some laboratories need the installer make an IQ document to certificate the good installation of this equipment. As well as an OQ document to certificate the good working of this equipment.

The installer can elaborate those documents following the examples:

Checking list for IQ (Instalation Qualification):

- Installer name.
- Manufacturer name:
- Unit identification: Code, model and serial number.
- Unit location: Laboratory, Person in charge.
- Checking packing list:
- Unit location matches the requirements?

Water and electricity requirements are suitable?

- Electrical supply.
Socket:
 - Ground terminal?
 - Suitable to power consumption. (16A)
 - Voltage according plate $\pm 10\%$
- Water connection.
 - Water pressure bigger than 1bar.
 - Hose to waste are properly fixed and not blocked.
- Water from steam generator: Low level of calcium, conductivity above 20 MicroSiemens.
- Personnel training: Operator personnel understand the instruction manual?.
- Date and name of who has made the Instalation (Signed).

Checking list for OQ (Operation Qualification):

- Installer name.
- Manufacturer name:
- Unit identification: Code, model and serial number.
- Unit location: Laboratory, Person in charge.
- Check the equipment distillate correctly. See: 8.1
- Check nitrogen recuperation with ammonia sulphate sample. See 8.2.
Note : ammonia sulphate weight, HCl normality, HCl consumption and obtained recovery.
- Check the white test See: 8.3
Note: Ammonia sulphate weight, HCl normality, HCl consumption and obtained recovery.
- You can check the nitrogen recovery with acetanilide sample. See : 8.4.
Note: Acetanilide weight, HCl normality, HCl consumption and obtained recovery.
- Date and name of who has made the Installation (Signed).


14 Warranty

The warranty covers the defective components or faulty assembling of the unit for 1 year since the purchasing date.

Keep purchasing bill as the only valid document for the purchasing date.

The warranty does not cover the damages caused for a unit malfunction. The unit must be suitably qualified before to be used for valuable samples processing.

Return good are not accepted without previous J.P. SELECTA, s.a. approval.

15 Troubleshooting table

- **The green indicator «ON» does not light.**

Check fuses and wiring connections cables.

Check main switch.

- **The equipment still adding water till it breaks out of the sample.**

The equipment do not detect the water level for low conductivity. Add about 250ml of tap water in the water tank.

Check the cable (white) of the level probe of the steam generator.

- **At the end of distillation, the sample tube content is absorbed to the steam generator.**

Check the steam hose is not obstructed.

Locate the aerated electrovalve:

Check it actives, check the voltage with a tester (230V).

Check (detaching it) it is not obstructed.

- **The alarm (2) for hose presence and lock door is always warned and block the unit running.**

Check the door is correctly closed.

Check the microswitches are o.k. (you will hear a click) when you close the door or lower the lever.

Check the cable (grey) continuity.

- **The refrigerating water does not enter:**

Check the tap is open.

Locate the refrigerating water input valve, check the voltage with a tester (230V). If there is voltage and the valve does not open knock it. That is due to a small lime rest blocking the water pass. If knocking is not enough, detach the valve to clean the water pass.

- **Security thermostat alarm warned and the equipment is blocked.**

This alarm warns when there is an excessive temperature in the steam generator. In this case the heating element stay disconnected by the thermostat and the control electronics block all the unit functions.

To rearm the thermostat, wait some minutes to let the temperature of the generator falls down to more normal level, and pulse the rearm switch.

The cause of the thermostat disconnection is that the steam generator lack of water. Check the hydraulic wiring to generate water, tank, pump, hoses and level probe.