
## **EXTRACTOR DE GRASAS**

### ***FAT EXTRACTOR***

**DET . GRAS N 2p            4002841**  
**DET . GRAS N 6p            4002842**

# Indice

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# Index

<b>1 Packing list and accesories</b> .....	<b>24</b>
1.1 Standard package contents 4002842 (6 samples) .....	24
1.2 Standard package contents 4002841 (2 samples) .....	24
1.3 Accesories.....	24
<b>2 Safety</b> .....	<b>25</b>
2.1 Hazards .....	25
2.2 Risks which the operator is put under .....	25
2.3 Personnel training.....	25
2.4 Proper uses of the equipment.....	25
2.5 Modifications.....	25
<b>3 Working principle</b> .....	<b>26</b>
<b>4 Installation</b> .....	<b>27</b>
4.1 Suitable location .....	27
4.2 Connection to power supply .....	27
4.3 Connection to tap water and waste.(Fig 4.3) .....	27
4.4 Selecting the locking gasket. ....	27
<b>5 Operation</b> .....	<b>28</b>
5.1 Control panel.....	28
5.2 Programming .....	28
5.3 Running a program.....	29
5.4 Sample and manual valves position.....	31
5.5 Up / Down the sample into the column. ....	31
5.6 Mounting the extraction thimbles. ....	31
5.7 Inserting the extraction thimbles into the unit. ....	32
5.8 Inserting the beakers.....	32
<b>6 The fat extraction process</b> .....	<b>33</b>
6.1 Sample preparation.....	33
6.2 Drying the sample .....	33
6.3 Weighing the beakers. ....	34
6.4 Placing the samples in the extractor unit. ....	34
6.5 Placement of the beakers and Solvent into the extractor. ....	34
6.6 Boiling step .....	35
6.7 Rising step .....	35
6.8 Solvent recovery phase .....	35
6.9 Solvent Evaporation step and weighting .....	35
6.10 Solvent adding.....	35
<b>8 Checking the unit. (OQ)</b> .....	<b>36</b>
8.1 Functional checking. ....	36
8.2 Solvent recovery checking. ....	36
<b>9 Maintenance</b> .....	<b>37</b>
9.1 Unit cleaning:.....	37
9.2 Replacing gaskets:.....	37
9.3 Hose replacing : .....	37
9.4 Fuse replacing:.....	37
9.5 Storing the unit: .....	37
9.6 Preparing the unit for transportation: .....	37
9.7 Spare parts .....	37
9.8 Diagnostics menu .....	38
9.9 Problem solving.....	39
9.10 Error meassage on the display:.....	39
<b>10 Technical specification</b> .....	<b>40</b>
<b>11 Warranty</b> .....	<b>40</b>
<b>12 Electrical wiring diagram</b> .....	<b>41</b>
<b>13 Process design</b> .....	<b>42</b>
<i>Notice to customers</i> .....	<b>43</b>

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
## 2 Safety

The DET GRAS N unit has suitable measures for its use in laboratory. However, the fat analysis method can be hazardous and the user must be informed of some of those risks.

### 2.1 Hazards

It identifies the risk and the safety actions to be applied.

The symbols make reference to the paragraph marked with a grey line.

#### Hazard of contact with corrosive and flammable solvents

Hazard of contact with solvents. Follow the instructions for the described operation.

#### Hazard of electrical shock

Electric hazard for entering in areas marked with this signal, or for following the manual instructions marked by this icon.

Follow the instructions for the described operation.

#### Hazard of burns to be in contact with hot surfaces

Temperature in the heater surface exceeds 60°C. Use the describer tools to perform the operation on the unit.

Do not touch the heater surface and the beakers during and after extraction process.

### 2.2 Risks to which the operator is subjected

The process of determination of fat with the described equipment exposes the operator to the following situations of risk:

- Manipulation of solvents
- Manipulation of glass parts
- Possibility of touching parts hotter than 60°C
- Electrical risk.

### 2.3 Personnel training

This equipment must only be operated by personnel trained and suitable qualified to work in a chemical laboratory and its common hazards.

This equipment must only be operated by personnel who have read and understood this manual or have been qualified on its manipulation and operation.

### 2.4 Proper use of the equipment

This equipment is for laboratory use, for sample fat extraction. Mainly in soxhlet analysis.

The improper use of this unit (not following manual instructions) can endanger the protection prepared by the equipment.

This equipment its NOT designed and not manufactured to work in an explosive ambient area except the ones described on this manual.

### 2.5 Modifications

The operation modification or the manipulation of the safety systems unauthorized by the manufacturer, could expose the operator to hazards not indicated in this manual.

In particular, do not operate the unit without the cover on the IP55 waterproof boxes. (Covers for electronics parts)

### 3 Working principle

The DET-GRAS N is a solvent fat extractor based on the classical «Soxhlet» method.

The units works by electrical heater. All the electrical and electronics are under IP55 protection.

The units includes the following main systems:

#### The heater

Built from an aluminium heat sink and 2 or 6 heating elements to get an homogenous heater system.

#### The extraction columns

Designed to yield an efficient fat extraction and to condensate a big rate of the solvent used.

#### Locking system

It moves the glass columns upwards / downwards and gives a reliable watertight lock between the glassware and the aluminium beakers.

#### Solvent cycle (See Fig 3.1)

The path that solvent runs into the glass column, yields the sample fat extraction.


At the beginning the solvent is on the aluminium beaker. Placed at the column bottom.

Due to the heating, the solvent evaporates once the boiling point is reached.

The solvent steam rises, through the sample, to the top of the column.

At the column top there is a double condenser, so the solvent becomes liquid again and returns back to the bottom, through the PTFE manual valve, and falls over the sample, carrying solved its fat.

**Fig 3.1 The solvent cycle**

When the solvent + fat falls into the aluminium beaker, the solvent is evaporated again, but the fat remains in the beaker.

After a few cycles, it is expected that all fat from the sample is moved into the beaker.

If the PTFE manual valve is closed, the solvent is retained at the top of the column. This step is called «Solvent recovery».

The sample can be placed in two positions:

- Top position. On this position the sample is placed at the solvent steam action and also at the solvent drops that come from the top of the column. This position is called «Rinsing».
- Down position: On this position the sample is immersed under the boiling solvent. This position is called «Boiling».


Fig 4.1.1: Free space on the rear part


Fig 4.1.2: Free space to let on both sides

Fig 4.3: Inlet and outlet water connections

## 4 Installation

Once the unit is unpacked, visually check the packing list (See 1.1) and that there are no parts damaged during transportation. Save the packing for some days for later sending.

### 4.1 Suitable location

Place the unit over a flat level and stable surface, suitable for the unit weight. No longer that 1.5m from a water tap, waste output and power supply plug.

To work comfortably, let a free space of 50cm on each side and not least than 20cm on the rear part.

### 4.2 Connection to power supply

Locate a power supply socket (rated to 16A) near the unit. Use the power cord supplied or another of similar specifications.

For safety reason the socket must have ground terminal, since the unit must be grounded.

Before plugging on the socket, check that its voltage matches the voltage marked on the unit plate.

### 4.3 Connection to tap water and waste (Fig 4.3)

Connect the water entry hose to a water tap (E). Connect the water output hose from (S) to waste.

Do not forget the rubber white nuts on the hose's connections.

Press plastic connections by hand. Do not use any tools.

Water pressure and temperature must be enough to cool the condensers. High cooling water temperatures lead to a poor condensation of the solvent.

### 4.4 Selecting the locking gasket

The unit is supplied with 2 sets of gaskets. One set is made from VITON and the other from BUTYL. The VITON ones are factory assembled.

**Note:** Take care to not mix the gaskets. If the gasket material and the solvent are no compatible, it will deteriorate in a short time.

Assembly the suitable gaskets according to this table:

#### Compatibility between solvent and gaskets material

Solvent	Ta(°C)	Viton	Butyl
Acetone	130 - 150		*
Acetonitrile	140 - 160		*
Carbon Tetrachloride	120 - 185	*	
Chloroform	100 - 140	*	
Ciclohexane	140 - 180	*	
Diclorometane	100 - 130	*	
Ether di-etilic	80 - 100	*	
Dioxane	180 - 200		*
Ethylate acetate	150 - 180		*
Hexane	140 - 155	*	
Petroleum Ether 40-60	100 - 135	*	
Petroleum Ether 60-80	120 - 160	*	
Toluene	200 - 220	*	
Trichloroethylene	150 - 180	*	

## 5 Operation

The control panel allows function selection and shows the unit status.

The DET-GRAS N does not have main switch. After plugging the unit, it starts to operate. By pressing «ESC» key from the main menu, the unit goes to "Sleep" status (Similar as when you stop the TV with the remote control).

Press any key to restart the unit.

**Note:** For a long time of no operation, please unplug the unit from mains.


### 5.1 Control panel

Figure 5.1 shows the DET GRAS N control panel that has an LCD display and 4 keys.

The control software has this main menu:

Fig 5.1 Control panel.

### 5.2 Programming

To perform a fat extraction with DET GRAS N, a program must be used. Up to 15 programs can be stored, in order to be used on a different extraction process.

A program has the following parameters:


- Program Name: 9 characters name. Choice a name that will be useful to remember, for example the type of sample.
- Solvent: On depending the solvent chosen, it automatically sets the heating temperature. This avoids the danger of using a temperature that could be dangerous according to the solvent used.
- «Boiling» time: Time (minutes) that the sample will remain immersed under the boiling solvent.
- «Rinsing» time: Time (minutes) that the sample will be under the action of the solvent steam and the condensed one.

**Parameters shortening:**  
 Bo: «Boiling» time  
 Ri: «Rinsing» time  
 Re: Solvent recovery time.  
 T: Temperature.

**Program number 0**  
 This program is only for checking purposes. It is not possible to change any parameter on it.

**Program 0 parameters:**

Bo: 10m  
 Ri: 10m  
 Re: 10m  
 T: 100°C

**Working with different temperatures**

In order to avoid dangerous situations, the DET GRAS N sets automatically the extraction temperature for any solvent, it can not be modified. (The most used solvents have been taken into account)

This avoids the dangerous to work at temperatures near to auto ignition solvent point.

However, if you need to change this temperature it can be changed.

To do it, when the solvent list is showed select the last one: "FREE"

**5.3 Running a program**

**CAUTION**  
 FIRE HAZARD  
 Be sure that you are using a program suitable to the solvent used on the extraction

**CAUTION**  
 FIRE HAZARD  
 Do not forget to open the coolant tap water before starting a program



To run a program:

- Open the water tap for coolant water.
- Select the menu «OPERATION» to run an stored program. See 5.2.
  - Bear in mind that program number "0" is only for checking purposes.
  - After «OPERATION» is selected, the first program shown on the display is the last program used or edited.
  - To start a program, press key ENTER/START and could be stopped at any time by pressing key ESC/STOP.
- The program initiates the «BOILING» step, so place the thimbles on the low position (See 5.5). Once the «BOILING» step finished, the unit makes some BEEP and the display flashes.

**Do not let the unit completely unattended.  
Pay attention to the end of each step to change, at time,  
the extraction thimbles position.**

- The unit jumps directly to «RINSING» step. Move the thimble to the up position. Once the «RINSING» time finishes, the unit makes some BEEP and the display flashes.
- The unit jumps directly to «RECOVERY» step. Close the PTFE manual valves for solvent recovery. Once the «RECOVERY» time finishes, the unit makes some BEEP and the display flashes.

#### 5.4 Sample and manual valves position

On each program step: , the sample position must be changed and also the PTFE manual valve. These changes are not automatic, they must be done manually.

At the end of each step a beep sounds and the display flashes.

**Note:** The device automatically changes phase. At the end of the «Boiling» phase, it automatically goes to «Rinsing» and the same happens with «Solvent recovery».

This table, summarizes the sample positions and manual valves on each step:

Step	Sample position	PTFE valves
Boiling	Down	Open
Rinsing	Up	Open
Recovery	Up	Closed

#### 5.5 Up / Down the sample into the column

Use the side knobs to move the sample up and down.

Place to the up position to:


- Insert the sample (thimble) to the unit
- Run «Rinsing» step.

Place to the down position to:

- Remove the sample (Thimble).
- Run «Boiling» step.

Mounting the thimble holder

Storage of extraction thimbles on the rack

#### 5.6 Mounting the extraction thimbles

Each extraction thimble must be mounted to a thimble-holder. It is made from a magnetic steel which allows hanging it on the column magnet.

- Clean the thimbles (if reused) with ether or some other degreasing solvent.
- Insert each thimble to a holder. Be sure to well tight.
- Use latex gloves to not add fat from your hands on the sample.
- Store it on the rack.

**Note:** Along the process. Use gloves, pincers or tongs to not add fat from your hands to the extraction thimbles.


### 5.7 Inserting the extraction thimbles into the unit

Whit the package comes a set of accessories that must be used for thimble manipulation and to insert all sample together into the extractor.

- Place the column magnets to the Up position (See 5.5).
1. Hang on the «Alignment tubes rack» on the front side of unit.
  2. Place the «alignment tubes» into the «Alignment tubes rack» holes.
  3. Place the «extraction thimbles» into the «alignment tubes».
  4. "Catch" the «alignment tubes» with the «insertion holder».
  5. Place the «insertion holder» under the columns.
  6. Take down the magnets until you capture the extraction thimbles.
  7. Raise the magnets (that carries the thimbles).
  8. Remove «Insertion holder».

The thimbles are ready to the extraction.

### 5.8 Inserting the beakers


The aluminium beakers collect the extracted fat from the sample. They are placed at the bottom of each column.

The locking system, between the beaker and the column is made by a gasket (Viton or Butyl) (See 5.3).

The operation to place the beakers must be done with care, because a wrong placement can damage the gaskets and produce solvent leakage.


Placing the beakers:

- Fit the beakers into the «Beakers rack».
- Used the side knob to raise the extraction thimbles.
- Place the «Beakers rack» under the glass columns.  
Press down the locking lever, and, at the same time, align the beakers and check that all fits under the white part of each column.
- Press the locking lever until it becomes engaged.
- To release the lever, pull the lever lock.


Placing beakers into the beaker rack.

Lever lock


Mounting the extraction sample.


Storing the extraction thimble on the rack

Handling the extraction thimbles.

Dry the sample around 105°C

Beakers handling

## 6 The fat extraction process

The fat extraction process includes the following steps:

1. Sample preparation.
2. Sample dry.
3. Beakers weighting.
4. Inserting samples to the extractor.
5. Placing aluminium beakers and solvent to the extractor.
6. «Boiling»
7. «Rising»
8. Solvent recovery step.
9. Evaporation solvent step

### 6.1 Sample preparation

The fat extraction by a solvent (Soxhlet) is efficient when the solvent can reach all the points of the sample. For this reason if the sample is triturated, the extraction will be more efficient. However, you must have in mind, that this operation can loose fat, that remains stucked on the mill walls. To minimize this error prepare at least 500gr of sample.

Weigh approx. 3 grams of sample. If it is necessary, crush or grind the sample (always after weighing it) clean off the crusher receptacle with some cotton wool soaked in solvent, to ensure there is no fat left in the bottom.

Place inside the extraction cartridge (thimble) the sample and the cotton wool used to clean off all the fat.

Clean the thimbles with ether or some other degreasing solvent and insert them in the extraction cartridges. Place them into the support rack.

To work accurately, during the manipulation of the cartridges it is advisable not to touch the cartridge with your fingers as grease from your hand could affect the result. It is recommended that rubber gloves and tongs be used.

Depending on the sample type perform a hydrolysis. (See application notes)

### 6.2 Drying the sample

It is frequently advised to dry the sample at 105°C for 1 to 24 hrs, according to the amount of water in the samples. Place the samples in the rack into the oven.

This stage can take a considerable amount of time. Depending on the sample type it can be omitted. To determine if drying is necessary, several assays with and without drying should be performed to determine any variation. If the results are similar, no drying process is required for similar samples.

If the samples are not being assayed immediately, it is recommended that they are stored in a desiccator.

Extraction thimbles ready to insert into the extractor.

Placing thimbles on the extraction unit.

Placing the beakers on the extraction unit

Right position for beakers

### 6.3 Weighing the beakers

The aluminium beakers should be clean, dry and free of fat. It is recommended to wash them with degreaser and volatile solvent to ensure they are dry.

Once clean, handle them with latex gloves or tongs.

Number the beakers with an indelible marker.

Accurately weigh (advisable to 0.0001 g) each beaker and record the result.

Place the beakers and the rack «beaker rack».

### 6.4 Placing the samples in the extractor unit

Insert the «insertion holder» with the «alignment tubes» and the extraction thimbles on the cut out profile in the front of the extraction unit.

Place in the rack «tubes support-cartridges»

With the help of the «magnetic clip» transfer the cartridges of the «extraction thimbles rack» to the «alignment tubes».

In the extraction unit locate the levers of each column into position «RISING».

With the «insertion holder» grab the «alignment tubes» and introduce them into the extraction unit taking care that each cartridge adheres to the magnet on each column.

### 6.5 Placement of the beakers and Solvent into the extractor


In the extraction unit, place the magnets of each column in position «RINSING» (Up)

Add into each aluminium beaker 50ml of Solvent.

By means of the «beaker rack», introduce the beakers under each column. Lower the lever and at the same time align the beakers so that they are inserted under the piece of white Teflon of each column (See 5.8).

The unit is now ready for the extraction.


From this point, the solvent cycle begins: boiling - evaporation - condensation, subjecting the sample to the action of the solvent (in a gaseous or liquid state) extracting the fat from the sample into the aluminium beaker.

On the «Boiling» step the sample is immersed under the solvent.

On the «Rinsing» step the solvent drags the fat from the sample.

PTFE manual valve position for solvent recovery

Adding extra solvent

## 6.6 Boiling phase

Place the magnets+thimbles of each extraction column at «BOILING» position.

In this position the samples are immersed under the boiling solvent. This phase accelerates the extraction process. However depending on the nature of the sample a longer «boiling» time may be required which also requires a longer «Rising» time.

A typical time for this phase is 30 - 45 minutes.

This time can be shortened if after several trials with less time provide the same results as the original longer period.

## 6.7 Rising phase

Place the magnets+thimbles of each extraction column into the «RINSING» position.

In this position the samples are subjected to the solvent vapours and to the liquid solvent condensed from the cooled reflux column. This phase is usually longer than the «BOILING» phase.

A typical time for this phase is 40 - 60 minutes.

This time can be shortened if the same results can be obtained with a reduced time.

## 6.8 Solvent recovery phase

Rotate to the horizontal position the Teflon key of each column.

In the extraction unit, make sure that the extraction cartridges are in the «RINSING» position.

In this phase, the solvent that comes from the condenser does not go back into the aluminium beaker but it is retained in the upper part of the collar.

A time of 10 to 30 minutes is enough to recover the 50ml of solvent.

The conclusion of this phase is that the aluminium beakers only contain the fat extracted from the sample.

## 6.9 Solvent evaporation phase and weighting

As this is a gravimetric method, any rests of solvent that remain into the beaker, after the extraction, will count as fat.

So it is important to guarantee the absolute absence of solvent on the beaker after extraction.

To do it first, remove the beakers with the beaker rack and dry in an oven at 105°C for 15 minutes (Raise the locking lever, to extract the aluminium beakers).

Weigh them. The increment difference in weigh is the weigh of the extracted fat.

## 6.10 Solvent adding

Extra volume of solvent can be added during the extraction process. Use a syringe to add solvent from the top of the glass column.

## 7 Checking the unit (OQ)

Once the user knows how the unit works, proceed to check it. This operation is known as Operation Qualification (OQ).

There are two levels of qualification:

### 7.1 Functional checking

This checking is without solvent. Do not place any beaker and thimble into the extractor.

Diagnostics utility (See chapter 9) can be used for this functional check.

Set a temperature of 100°C and verify that the temperature of the heating surface stabilizes at about 100°C ± 10°C, it remains stable and is uniform throughout the surface.

### 7.2 Solvent recovery checking

This OQ operation measures the solvent recovery rate.

- Prepare 2 or 6 aluminium beakers (on depending equipment model)
- Measure with a cylinder 50ml of solvent, for each beaker.
- Insert the beakers on the DET GRAS N.
- Run a P0 program. If done with another solvent, the recovery time may vary, becoming necessary in some cases up to 30 minutes:

Temperature:	According to solvent used.
Boiling:	15 minutes.
Rinsing:	15 minutes.
Recovery:	10 minutes.

- Let the unit get cold for 20 minutes.
  - Let PTFE manual valves close.
  - Let coolant water open.
- Remove the beakers set.
- Insert another cold set of aluminium beakers.
- Open PTFE manual valves. All solvent must fall into the beakers.
- Measure the volume of solvent recovered on any beaker (Use a cylinder). Work quickly to avoid solvent evaporation.
- Calculate the recovery rate:  $R(\%) = V_2 / V_1 * 100$
- Normal recovery values are above 50 %.


Fig 8.1 Gasket

## 8 Maintenance

Maintenance operation on the DET GRAS N (2 or 6 places) must be properly run in order to keep the safety degree, and specially to avoid any solvent leakage.

### 8.1 Unit cleaning

Use only a humid wash-cloth (small quantity of soap) to general clean the unit.

Dry quickly any leakage of solvent.

### 8.2 Replacing gaskets

Review the gasket each month or whenever a solvent leakage appears.

Change it if they are not flat and flexible.

### 8.3 Hose replacing

Use silicone hose Ø 14x8 to replace the coolant water hoses. Usually this hose must be replaced after changing a glass column.

### 8.4 Fuse replacing

Fuse are into the rear IP65 waterproof box. (See Fig 8.2).

Unplug the unit from the mains before open the IP65 box.

### 8.5 Storing the unit

Do not store the unit for a long time with solvent inside.

### 8.6 Preparing the unit for transportation

Do not transport the unit with solvent inside.


Fig 8.2 Remove the 4 screws to open the IP65 box and change the fuse.

### CAUTION

Before any work inside the equipment, disconnect it from the mains

28118


28120

28191


28105


28107

29513

39281

29514 / 21

### 8.7 Spare parts

For repair and replace deteriorated parts, please contact J.P. Selecta s.a.u. or your distributor. They will inform you about the technical services authorized in your area.

Find a spare parts list for easy buying and identifying.

- 28105** PTFE condensation disk.
- 28107** Magnet for extraction thimble.
- 28118** VITON gasket.
- 28120** BUTYL gasket.
- 28191** PTFE spout.
- 29513** Control board.
- 29514** 6 places power board.
- 29521** 2 places power board.
- 39281** Heating element 100W 230V.
- 46041** Silicone 8x14 hose
- 47049** Glass extraction column.

47049

### 8.8 Diagnostics menu

The DET GRAS N includes a menu to identify any problem with the measure and temperature controller.

This operation serves to homogenize the temperature of the different equipment internal probes.

- **Auto-calibration:** This operation allows heater surface temperature calibration against an external thermometer. (It must be provided with a suitable probe for surface measurement).

The process is as follows:


- On pressing ENTER, the unit starts for near 30 minutes to stabilize the temperature around 80°C.
- When this point is reached, a repeated beep sounds and the self-calibration process will be over.

- **Probes reading:** Use this operation to check the temperature probes and the heating elements.

The process is as follows:

- Set a working temperature using the Up/Down keys. Validate with ENTER.
- Toggle with the ESC key if the readings are shown with (C) or without (N/C) the correction.
- Start the unit by pressing ENTER (a "1" is shown on the up right display corner).
- Wait until the set temperature stabilizes. The reading of each probe appears on the display.
- Stop with ESC. To return to the main menu, keep ESC pressed for a few seconds.

From the main menu, keep ENTER key pressed for 10 s until the display changes. 

## 8.9 Problem solving

### 1. Solvent drain from between the seal and the aluminium beaker.

- **Damaged gasket.** Check that it looks flat not rippled, flexible and is well placed, and that it is solvent-related material.
- **Differences of height between the aluminium beakers.**
  - Use beakers of similar height.
  - Check that the top of the beaker is dirt free.
  - Check that the top edge of the beaker is not deformed.

### 2. Negative fat determination

- **Error in weighting of the aluminium beaker.**
  - Clean and dry the beakers before weighing them. Number the beakers before weighing them.
- **Uneven heat distribution through the surface heater.**
  - Use the diagnostics menu to check that all probes and heating elements work.
- **Faulty contact (heat transfer) between the beaker bottom and the heater.**
  - Check that the beaker does not have any dirt on the contact surface.
- **Unequal solvent content**
  - Check for solvent leakages.
  - Solvent contents should be 50ml at least, if necessary add solvent.

### 3. Solvent takes a long time in boiling.

- **Uneven heating on the heater surface.**
  - Use the diagnostics menu to check that all probes and heating elements work.

### 4. The solvent produces foam over the top of the extraction cartridge.

- **Sample with too much fat contents**
  - Reduce sample size
  - The aluminium beakers have fat residue from previous tests.
  - Sample with a high water content. Dry the sample for longer.

### 5. Insufficient fat quantity detected.

- **Lack on the process design.**
  - Check sample preparation.
  - Check times to drying, BOILING and RINSING.
  - Choice the solvent and working temperature.

## 8.10 Error message on the display

### **ERROR 1: AVERAGE TEMPERATURE DOES NOT GO UP**

Causes: Solid relay on power PCB, wiring, heater / temperature probes.  
Action: Error message on the display and disconnection to sleep status.

### **ERROR 2: OVERTEMPERATURE ON A TEMPERATURE PROBE (> 15% average T<sup>a</sup>)**

Causes: Solid relay on power PCB, heater/temperature probes, differences on solvent volume, beaker ausence.  
Action: Error message on the display.

### **ERROR 3: DANGER OVERTEMPERATURE (>25% T<sup>a</sup> set)**

Causes: Solid relay on power PCB, wiring, heater / temperature probes.  
Action: Error message on the display and disconnection to sleep status.

### **ERROR 4: PROBE READING DOES NOT GO UP (Reading < 50°C, with T<sup>a</sup> avg > 70°C)**

Causes: Heater / temperature probes.  
Action: Error message on the display and disconnection to sleep status.

### **ERROR 5: TEMPERATURE GROWS TO LOW (< 80% set temperature after a few minutes)**

Causes: Heater / temperature probes.  
Action: Error message on the display.

### **ERROR 6: IN STAND-BY STATUS, AVERAGE TEMPERATURE DOES NOT GO DOWN.**

Causes: Solid relay on power PCB, wiring  
Action: Error message on the display and disconnection to sleep status.

**ERROR 7: SOME PROBE IS ABOVE 250°C**

Causes: Probe damaged, triac of the shorted probe.

Action: The device warns and goes to rest immediately

**ERROR 8: IN OFF STATE AND WITH A MIDDLE TEMPERATURE LESS THAN 70°C, SOME PROBE IS ABOVE 80°C**

Causes: Probe damaged, triac of the shorted probe.

Action: The device warns and goes to sleep after a few seconds

**9 Technical specification**

Mains voltage:	230V 50/60Hz
(See plate)	Other power supply voltages could be available (Ask to your dealer)
Electrical power:	600W (6 places) 200W ( 2 places)
Maximum consumption:	3 A 1.5 A
Fuses:	4 A fast. dimensions: 20x5 mm
Power supply cord:	Phase, Neutral, Earth (section 1mm <sup>2</sup> )
Weight (Net):	25 Kg 19Kg
Dimensions: W x L x H	400x750x700 mm 400x450x700 mm
Thimbles size	External diameter 26mm /height 60mm
Max working temperature	80 to 240°C
Solvent recovery:	Better than 60% - 80%
Reproducibility:	Relative Standard deviation)< 1%
Sound level:	Lower than 40 dBA
Protection degree:	IP55
Ambient conditions:	Inner use. Altitude up to 2000m Temperature from 5°C to 40°C. Relative maximum damp 80% for above 31°C temperatures, decreasing linealy to 50% for 40°C.
Overvoltage category:	Category II
Pollution degree:	2
Glass material:	Borosilicate 3.3
Carriage material:	Iron with overlap with RILSAN coating chemical resistant.
Hosepipe material:	Silicone.
Number of programs:	15
"Boiling" extraction time:	from 0 to 999 min
"Rinsing" extraction time:	from 0 to 999 min
Solvent recovery time:	from 0 to 99 min

**10 Warranty**

The warranty covers the defective components or faulty assembling of the unit for 1 year since the purchasing date.


Keep purchasing bill as the only valid document for the purchasing date.

The warranty does not cover the damages caused for the unit malfunction. The unit must be suitably qualified before being used for valuable samples processing.

Return goods are not accepted without previous J.P. SELECTA, s.a.u. approval.

# 11 Electrical wiring diagram

## 12 Process design

The analysis process can vary depending on the nature of the sample. Before a process can be established, it is necessary to carry out a series of tests of known fat determinations. Once the process is designed, it is advisable to validate it with inter-laboratory exercises.

The following steps show the general process for the determination of fats:

1. Representative sampling of the samples
2. Homogenize, grind, grind ... the sample
3. Weigh the samples and tarnish the glasses. Number them.
4. Drying of the samples
5. Hydrolysis (optional)
6. Extraction
7. Solvent evaporation
8. Weigh and calculate the fat% on the sample weight

A first guidance on the process to follow can be obtained from the official tests or trials already designed. Although these are designed using the classic soxhlet method they are easily modifiable to execute them with the DET-GRAS extraction equipment. Initially, respect the operations described in the test outside the extraction.

Homogenization and grinding must be carried out without losing fat during handling. For this, it is advisable to prepare a large sample quantity (Ex: 1Kg) so that the effect of fat loss is negligible.

Another technique is to crush the sample in a container that can be completely cleaned with a solvent impregnated cotton so that the lost fat in the walls of the container is part of the sample.

In general, in the extraction, the "Rinsing" phase must be longer than "Boiling", since it is the proper stage of extraction. The "Boiling" phase can be understood as "preparation" to facilitate extraction.

Normally the meat tests recommend a hydrolysis prior to extraction. Hydrolysis complicates the test so you can study the steps of crushing, drying, extraction to avoid it.

Hydrolysis consists in the removal of organic matter by boiling the sample in HCl from 3 to 5N. Subsequently, each sample is washed by filtration. In the filter paper, after several washes, only fat remains. The paper is inserted into the extraction cartridge and then inserted into the extractor. Since fat is manifestly exposed to the action of the solvent extraction is very fast and complete.

Hydrolysis can be carried out with usual laboratory material such as heating plates, Erlenmeyer jars, filtration funnels, filter paper, ...


### **Choosing the solvent and the extraction temperature:**

Check in standards, specialized books, etc. which is the appropriate solvent for the extraction that is wanted to carry out. In case you do not have any information, oil ether is used 60-80°C.

It should also be taken into account that the working temperature does not spoil the extracted grease.

Normally, use the extractor's own temperatures for each solvent.

**Notice to customers:**

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.