DESTILADOR KJELDAHL SEMIAUTOMATICO
SEMIAUTOMATIC KJELDAHL DISTILLER


PRO - NITRO S 4002851

Indice

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Index

1 Packing list & accessories	31
1.1 Standard package content.....	31
1.2 Accessories	31
2 Safety	32
2.1 Safety symbols	32
2.2 Risks which the operator is put under	33
2.3 Personnel training.....	33
2.4 Proper uses of the equipment.....	33
2.5 Modifications.....	33
3 Overview	34
4 Installation	36
4.1 Suitable location	36
4.2 Connection to power supply	36
4.3 Connection to tap water and waste.....	36
4.4 Water for steam generator	36
4.5 Connection to reagents tanks	36
5 Operation.....	37
5.1 Insert/ Remove the sample tube.....	37
5.2 Operation modes	37
5.3 Parameters setting.....	37
5.4 MANUAL operation.....	39
5.5 AUTOMATIC operation	40
5.6 Alarms	40
6 Sample analyses.....	41
6.1 Check levels of tanks	41
6.2 Unit warm-up and cleaning.....	41
6.3 Distillation. (Sample analysis)	41
6.4 Analysis interruption	42
6.5 Titration and Protein contents (%) calculation.	42
7 Reagents	43
7.1 Prepared reagents	43
7.2 Proper reagents set up.....	43
8 Checking the distiller. (OQ).....	44
8.1 Functional checking	44
8.2 Checking Nitrogen Recovery.	45
8.3 Blank test.....	45
8.4 Checking Nitrogen recovery with Acetanilide.....	46
9 Maintenance	47
9.1 Cleaning the unit:	47
9.2 Hose inspection (monthly).....	47
9.3 Preparing the unit for shipping or storing.	48
9.4 Level sensor cleaning.....	49
10 Technical specifications	50
11 Technical service	51
11.1 Spare parts	51
11.2 Hydraulic tubing diagram.....	52
11.3 Wiring diagram. (1 of 2).....	53
12 Kjeldahl analysis step by step.....	54
12.1 Sample preparation.....	54
12.2 Digestion.....	54
12.3 Dilution	54
12.4 Distillation	54
12.5 Titration and calculations.....	54
13 Guidelines for IQ and OQ	55
14 Warranty	55
15 Troubleshooting table	56

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
1 Packing list & accessories

1.1 Standard package content

- Distilling unit (1)
- MACRO size sample tube \varnothing 42 (1)
- Tray (2)
- Tube 1.5m incoming water. High pressure, coupling 3/8" (3)
- Tube 1.5m water outlet. Silicone 8x14, coupling 3/8". (3)
- Tube 1.5m remaining steam. Silicone 8x14, coupling 3/8". (3)
- Tube 1.5m generator emptying. Silicone 8x14. (3)
- Tube 1.5m waste outlet. ISOVERSINIC 7x12 (3)
- Power supply wire (4)
- 2 litre tank for NaOH (5)
- 2 litre tank for Fixative solution (7)
- Instructions Manual (9)
- Flask 100g of Ammonia sulphate. (10)

1.2 Accessories

Tubes, racks and tongs:

- MACRO sample tube \varnothing 42 4042300
- MICRO sample tube \varnothing 26 4001045

Racks:

- Rack MACRO 6 tubes. 4005071
- Rack MACRO 12 tubes. 4005081
- Rack MACRO 20 tubes. 4005091
- Rack MICRO 12 tubes. 4001053
- Rack MICRO 24 tubes. 4001054
- Rack MICRO 40 tubes. 4001055

Safety goods:

- Tongs for tubes 7000532
- Protection gloves 1001255
- Anti-acid gloves 5000044
- Safety glasses 5028192


Complete digestion sets

- Digestion unit MACRO 6 tubes. 4000629
- Digestion unit MACRO 12 tubes. 4000630
- Digestion unit MACRO 20 tubes. 4000631
- Digestion unit MICRO 12 tubes. 4001047
- Digestion unit MICRO 24 tubes. 4001048
- Digestion unit MICRO 40 tubes. 4001049

Smokes removal aids:

- Scrubber. 4001611
- Water circulation vacuum pump. 4001612

Digester

Scrubber

Water circulation
vacuum pump

2 Safety

PRO-NITRO "S" unit has the suitable measures for its use in laboratory. However, the nitrogen analysis by the Kjeldahl method can be hazardous and the user must be informed of some of those risks.

In order to be respected, across this manual, the risk of danger situations are remarked.

2.1 Safety symbols

It identify the risk and the safety actions to be applied.

The symbols make reference to the paragraph marked with a grey line.

Warning

Warning of a dangerous operation.

Indicated instruction must be followed.

Hazard of contact with corrosive substances

Hazard of contact with sodium hydroxide and strong acid, substances that attack skin, clothes and other materials easily.

Follow the instructions for the describe operation.

Hazard of electrical shock

Electric hazard for entering in areas marked with this signal, or for following the manual instruction signalled by this icon.

Follow the instructions for the describe operation.

Hazard of burns to be in contact with hot surfaces.

Temperature in this area can exceed 60°C. Use heat protected gloves for the described action.

Follow the instructions for the describe operation.

Wear gloves mandatory

The use of gloves is obligatory to run the operation.

In the section «Accessories» an advised model is indicated.

Mandatory use of safety glasses

The use of safety glasses or complete face protection is obligatory to conduct the described operation.

In the section «Accessories» an advised model is indicated.

Important information

- Information important to obtain good results or for an optimal operation of the equipment.
- Information important to extend the life of the equipment or to avoid the degradation of some of its components.

2.2 Risks which the operator is exposed

The process of determination of Kjeldahl Nitrogen with the described equipment exposes to the operator to the following situations of risk:

- Manipulation of concentrated sulphuric acids and Hydrochloric.
- Sodium manipulation Hydroxide 40%.
- Manipulation of glass parts.
- Possibility of touching parts to a temperature above to 60°C
- Possible exhibition to irritating acid steam.
- Electrical risk.

2.3 Personnel training

This equipment must only be operated by personnel trained and suitable qualified to work in a chemical laboratory and its common hazards.

This equipment must only operated by personnel who have read and understood this manual or have been qualified on its manipulation and operation.

2.4 Proper uses of the equipment

This equipment is for laboratory use, for sample distillation by steam pulling. Mainly in Kjeldahl nitrogen analysis.

The improper use of this unit (not following manual instructions) can endanger the protection prepared by the equipment.

Do not use to distillate explosive or flammable substances.

This equipment is not designed and not manufactured to work in an explosive ambient area.

2.5 Modifications

The operation modification or the manipulation of the safety systems unauthorized by the manufacturer could expose the operator to hazards not indicated in this manual.

3 Overview

The PRO-NITRO "S" is a SEMI-Automatic steam distillation unit. That means that performs the sequence : Boric adding, NaOH adding, timed distillation and sample empty to waste.

This automation degree assures more reliable and robust results.

It is very easy to use but its operation must be perfectly known in order to get good results and to avoid losing valuable samples.

Before starting to analyse valuable samples do some training sessions running the operations related on the chapter 8, until the operator be familiar with the unit.

Please, start to identify the name of the analyser components:

The steam generator

Generates steam from the water reservoir placed at the top of the unit.

The steam generator includes two safety devices: **Safety thermostat**: Cuts the electrical power to the heater when overtemperature is detected on the steam generator. (Probably for poor level of water), and the: **Pressure switch**: Cuts the electrical power to the heater when there is an overpressure on the steam generator.

Also includes a conductivity sensor to keep a constant level of water in the steam generator.

The distillation circuit






It is made, mainly, by two glassware parts: the **splash protector**: Its special design avoids that NaOH splash reach the cooler glassware, and the: **cooler**: It condensates the steams that are leaking from the sample. This glassware is cooled by a water circuit that comes from an external tap. An internal valve opens and closes this circuit automatically.



NaOH and ammonia fixative solution addition system

Composed by properly specified pumps to dispense NaOH directly to sample tube when distillation is running and fixative solution (Boric acid with mix of indicators) directly to Erlenmeyer collector before the distillation.

Emptying systems

It consists in some pumps and valves to empty the sample tube and the collector once the analysis is finished.

4 Installation

Once the unit is unpacked, visually check, the packing list (See 1.1) and that there are no parts damaged during transportation. Save the packing for some days for later sending.

4.1 Suitable location

Place the unit over a flat level and stable surface, suitable for the unit weight.

No longer that 1.5m from a water tap, waste output and power supplied plug.

To work comfortably, let a free space of 50cm on each side and not least than 20cm on the rear part.

Prepare to refill the water reservoir from the top of the unit.

4.2 Connection to power supply

Locate a power supply socket (Rated to 16A) near to the unit.

Use the power cord supplied or another of similar specification.

For safety reason the socket must have ground terminal, since the unit must work grounded.

Before to plug on the socket check that its voltage match the voltage marked on the unit plate.

4.3 Connection to tap water and waste.


Connect the water entry hose to a water tap (1). Connect the water output hose from (2) to waste.

The hoses «remaining steam» (3) and «Steam generator empty» (4) and "waste" (5) must go to waste.

Do not forget the white rubber white nuts on the hose's connections.

Press plastic connections by hand. Do not use any tool.

The tap water must have, at least, 1 bar pressure.

Note: It is critical in order that the unit works properly, the hose «Remaining steam» (3) MUST NOT BE IMMersed under water, or blended or blocked by any reason.

4.4 Water for steam generator

The unit includes a 6 litres water reservoir, at the top, to store the water for the steam generator.

As the system used to control the level of water on the steam generator works by conductivity, water of at least 20 microSiemens must be used. Water, also, must have a low content of calcium to increase the life of the steam generator.

Water with this conductivity could be obtained from a demineralised unit. If the conductivity is lower, add (for 6l of water) near 50mg of salt or sodium carbonate.

Reagents tanks

4.5 Connection to reagents tanks

At the left hand side, there is a door to access reagent tanks, NaOH, fixative solution and HCl. Connect each one to the hose indicated with his name.

5 Operation

The control panel allows function selection and shows the distiller alarm status.

Once the unit is plugged in, put it on, by pressing the main switch. The unit goes to stand by status (ON).

Note: If main switch is press consecutively with an interval of less than 3 seconds, the unit could start in abnormal mode. So wait few seconds to restart the unit.

On starting the unit for the first time, the water pump starts to work, automatically, to fill the steam generator.

This operation, usually, takes near 40 seconds and it will stops, also automatically when the steam generator is filled. At 15s an alarm sounds and the unit becomes blocked. Stop the unit by Main switch, wait 4s and restart. Repeat this operation 3 or 4 times. until the steam generator is filled.

If, after this 1 minute approx., this operation is not finished, please check that the water tank has water (Visual check on the right hand side level), also be sure that the water has the suitable conductivity. (See 4.4)

5.1 Insert/ Remove the sample tube

- Open the safety door
- Insert the tube starting to introduce the PTFE tube.
- Push down the fixing lever.
- Close the safety door.

Proceed similar to remove the sample tube.

CAUTION
Wear Anti-acid,
Anti-thermal gloves.

The white PTFE tube could have rest of
Sodium hidroxide or Acid.

The sample tube's temperature is
above 60°C after the distillation.

5.2 Operation modes

A Kjeldahl distillation includes the following operations sequence:

- Boric adding,
- NaOH adding,
- Timed steam distillation
- Sample empty to waste

The PRO-NITRO «S» can do this as an automatic sequence («AUTOMATIC» mode) or run, separately, each operation (mode «MANUAL»).


5.3 Parameters settings

These settings must be done in both modes: «MANUAL» and «AUTOMATIC»

- **Setting the distilling time:**

- Press  .
- Use keys   to set the distilling time (minutes).

- **Setting the "Low water level alarm" time (low level of water into steam generator)**

- Press  .
- Keep pressed keys   for a 10s until display shows dashes.
- Press  to select steam.
- Use   to set the alarm steam time (factory setting 36).
- Use  to validate the setting.
- Press   to save the setting and exit from the settings mode.


Fig. 5.1 Control panel

- **Setting the heating power on stand-by status:**

- Press  .
- Keep pressed keys   for a 10s until display shows dashes.
- Press  to select % stand-by heating.
- Use   to set to set the % heating (factory setting 5).
- Use  to validate the setting.
- Press   to save the setting and exit from the settings mode.

- **Setting distillation power:**

- Press the keys   for 5 seconds until the display shows dashes.
- Press  and  until the display shows the current value.
- Modify the value by means of keys   .
- Press  and  to validate the setting.
- Press   to save and exit from the settings mode.

PARAMETERS PROGRAMMING

Factory settings:

Low level water alarm:	36
% heating:	5
Distillation power:	80%
Volume of Boric acid:	4
Volume of HaOH:	12


Fig. 5.1 Control panel

Status indication

1. Heater ON
2. Water feed fail.
3. Safety door open or sample tube not present.
4. Overtemp on steam generator.

Display and selection keys.

5. Display.
5. Up/down keys.

Function keys.

6. Operation function mode selection between «Manual» and «Automatic».
7. **H3BO3 / START**
«MANUAL» mode: Boric dose.
«AUTO» mode: cycle start.
8. **NaOH**
Modo «MANUAL»: NaOH dose.
Modo «AUTO»: No function.
9. **STEAM**
«MANUAL» mode: Distillation.
«AUTO» mode: No function.
10. **EMPTY/STOP**
«MANUAL» mode: Empty tube.
«AUTO» mode: STOP


Setting the dosing volume of NaOH:

- Press .
- Keep pressed keys   for a 10s until display shows dashes.
- Press  to select NaOH.
- Use   to set the NaOH pump working time (factory setting 12).
- Use  to validate the setting.
- Press  to save the setting and exit from settings mode.
- Approximately each second means 5ml.

Setting the dosing volume of Boric acid:

- Press .
- Keep pressed keys   for a 10s until display shows dashes.
- Press  to select H3BO3.
- Use   to set to set the Boric pump working time (factory setting 4).
- Use  to validate the setting.
- Press  to save the setting and exit from the settings mode.
- Approximately each second means 10ml.

5.4 MANUAL operation

- Select MANUAL mode. Key  toggles working mode. See Fig 5.1
- Set parameters according 5.3.
- Insert a sample tube (See 5.1) and an Erlenmeyer flask.
- **Adding Boric:** Press  Boric acid will be added on the Erlenmeyer flask.
- **Adding NaOH:** Press  NaOH will be added on the sample tube.

- **Distilling start:** Press . See Fig. 5.1 key «9»
Distillation starts and stops after the distilling time.
- Distillation can be stopped, on any time, by pressing  key.
- **Emptying sample tube:** Keep press  » for a few seconds.



Fig 5.2 Toggle key and indicator for mode set.

5.5 AUTOMATIC operation

- Select AUTOMATIC mode. Key «6» toggles working mode. See Fig 5.2
- Set parameters according 5.3.
- Insert a sample tube (See 5.1) and an Erlenmeyer flask.
- **Starting the sequence:** Press .
The following sequence will start:
 - **Adding H3BO3** on the Erlenmeyer.
 - **Adding NaOH** on the sample tube.
 - **Distilling until** the distilling time.

Distillation can be stopped, on any time, by pressing .

- **Emptying sample tube:** Keep press  for a few seconds. The sample tube will be emptied to the waste output.

Note: the parameters keep stored on the unit.

5.6 Alarms

In order to avoid dangerous situations or to protect the analyser life, some alarm systems are provided on the unit:

(2) Low level of water on the steam generator alarm

Indicates that a low level of water was detected on the steam generator for a long time. The microprocessor stops the distillation.

Probably the water reservoir is empty or the water does not have the suitable conductivity.

(3) «Tube/Door» alarm

Safety door not closed or sample tube not inserted.

This alarm blocks the distillation and the NaOH adding.


(4) «Overtemp» alarm

An excess temperature on the steam generator has switched the safety thermostat.

Block the distillation function.

Once is switched, it must be left to get cool for some minutes. Solve the root problem that causes the switching and reset the safety thermostat by pushing the central button. The safety thermostat is located on the electrical box, on the right hand side. A cover must be removed.

The safety thermostat switching is due to low level of water on the steam generator, probably because the water reservoir is empty or water pump malfunction.

Status indicator & alarms.

1. Heater ON
2. Water tank empty
3. Safety door open
4. Overtemperature on Steam generator.

6 Sample analyses

To get good analysis' results and a good performance of the distiller it is suggested to follow a working **routine**.

Checking water level

Before starting to analyse valuable samples do some training sessions running the operations related on the chapter 8, until the operator be familiar with the unit.

6.1 Check levels of tanks

Before starting a batch analysis, check the following levels:

Steam generator water tank level

Refill it with demineralised water.

Assure its conductivity adding around 0.1g per litre of sodium carbonate or NaCl.

The water reservoir is 6 litres capacity, enough for analysing more than 20 samples.

Level of NaOH tank

Refill with NaOH 40% from the 2 litres flask. It is placed inside the unit.

Get the NaOH flask from the door placed on the left hand side cover.

The NaOH solution must be free of ammonia salts, to avoid its interference on the result. In order to avoid NaOH pump obstruction, take care that NaOH becomes completely solved on water, no particles must be on the solution.

Level of Ammonia fixative solution tank

Refill with Fixative solution (Boric acid 1% + indicator mixed 5) 2 litres flask. It is placed inside the unit.

Remember that a volume of 60 ml will be used on each distillation.

Note: See chapter 7 for fixative solution preparation.

Priming the dispensers

The tubing for NaOH and fixative solution reagent must be properly primed, air free, to provide a good dispensing.

Carry on the refill operations to avoid air entry on the tubing. See 5.1.2 and 8.1 to know how to carry on the prime operations.

6.2 Unit warm-up and cleaning.

In order to get the best results with your Kjeldahl distiller, you have to run a warm-up step before starting samples analysis. This procedure is also used for cleaning the distillation circuit after analysis batch.

For it just insert a sample tube with 25ml of clean water. Use MANUAL mode to distillate for 5 minutes.

Check the reagents tanks levels before starting a batch analysis.

6.3 Distillation. (Sample analysis)

Please follow these steps: (On the same order as indicated)

- Open the safety door.
- Place an Erlenmeyer flask on the cooler output to collect the distilled.

Check, visually that the distilled output tube remains immersed under the level of Boric acid.

- Insert the sample tube. See 5.1
- Select mode «MANUAL» or «AUTOMATIC» and initiate the sequence.
- Wait for the end of distillation. Buzzer sound.
- Empty the sample tube.
- Open the safety door.
- Remove the sample tube. See 5.2
- The unit is completely prepared to analyse the following sample. Any operation of cleaning is needed, unless any malfunction has occurred during the analysis.

Visual Checking during the analysis:

Some visual checking can be applied from the colour that takes the solution on the collector or in the sample tube.

This checking is useful at the beginning to set the optimal parameters or to solve unit malfunctions.

Visual Checking on the sample tube: (Using CuSO₄ Catalyst)

When the NaOH volume is sufficient, the sample becomes blue. (If not, program more NaOH on the following a analysis)

Visual Checking on the collector:

The initial colour must be the fixative solution colour (Violet with mixed indicator or Red with mixed indicator 5). If it is grey or green, it is because of not neutralized rests of the previous sample remains on the collector. On this case stop the analysis and run a cleaning operation. See 6.2

6.4 Analysis interruption

Press «STOP» to abort the analysis.

Note: The analysis interruption, left on the unit rests not neutralized that will corrupts the following analysis. Run a cleaning operation (See 6.2) before to run new analysis.

6.5 Titration and Protein contents (%) calculation

Prepare the titration burette with the reagent usually HCl 0.25N.

Slowly dispense the HCl 0.25N on the Erlenmeyer stirring it at the same time.

Titrate until the colour changes (end point at pH 4.65).

From the HCl volume (ml) consumption on the titration, the Kjeldahl Nitrogen (mg) is calculated:

$$\text{Nitrogen (mg)} = 14 \times \text{Volume}_{\text{HCl}} \text{ (ml)} \times \text{Normality}_{\text{HCl}}$$

And, from the nitrogen we calculate the sample protein content:

$$\% \text{ Protein} = \frac{P_2}{P_0} \times 100 \times F$$

P2: Nitrogen (mg).

P0: Sample weight (mg).

F: Protein factor.

Common protein factors

Almonds	5,18
Nuts	5,30
Nut- peanuts	5,41
Jello	5,55
Soya	5,71
Barley, trenches, rye	5,83
wheat, whole flour	5,83
Flours (not whole one)	5,70
Rice	5,95
Corn	6,25
General	6,25
Saved	6,31
Milk and milky products	6,38

7 Reagents

The reagents to be used must be free of ammonia.

7.1 Prepared reagents

The best thing is to use reagents factory prepared, mainly the HCl for titration. Any error on its preparation have a direct influence on the result.

As an example the following reagents could be used, or the equivalent one in other brands:

Boric acid 4% (For ammonia titration)	RV	PANREAC 282222
--	----	----------------

Mixed indicator 5 RV (Red of metile + Green Bromocresol)		PANREAC 283303
---	--	----------------

Mixed indicator RV (Red of metile + Blue of metilene)		PANREAC 282430
--	--	----------------

HCl 0.1N SV		PANREAC 171023
-------------	--	----------------

HCl 0.25N	SV	PANREAC 182318
-----------	----	----------------

H ₂ SO ₄ 0.1N	SV	PANREAC 181061
-------------------------------------	----	----------------

H ₂ SO ₄ 0.2N	SV	PANREAC 182011
-------------------------------------	----	----------------

Acetanilide 99% (Validation standard)		PANREAC 151005
---------------------------------------	--	----------------

Ammonia sulphate (Validation standard)		PANREAC 131140
--	--	----------------

Kjeldahl Catalyst 6.25% Cu tabl. 8g		PANREAC 174428
-------------------------------------	--	----------------

7.2 Proper reagents set up

As it is usual on most laboratories, the reagents to work with PRO-NITRO «S» could be prepared starting from concentrated ones.

The NaOH 40%, acid Boric 4% are not critical on the concentration and it is not necessary to work with precision.

Reagents for titration (HCl or H₂SO₄) must be high precision and calculation must be accurate. Any error on its preparation has a direct influence on the result.


Fig. 5.1 Control Panel

8 Checking the distiller (OQ)

Once the unit is properly installed and the operator has acquired the knowledge to work with it. The unit must be checked.

This process is known as **Operation Qualification**.

There are four basics checking to do:

8.1 Functional checking

Run this checking after installation and also for some maintenance. Use MANUAL mode.

8.1.2 Checking the steam generator filling.

On starting the unit for the first time, the water pump starts to refill, automatically, the steam generator. This operation usually takes around 1 minute and stops automatically. At 15s an alarm sounds and the unit becomes blocked. Stop the unit by main switch, wait 4s and restart. Repeat this operation 3 or 4 times until the steam generator is filled.

If, after this 1 minute approx, the operation is not finished, switch the equipment off. Check that maybe the water tank has no water or that the water has not the suitable conductivity. (See 4.4)

If so happens, check, empty the steam generator by opening the drain and repeat the operation.

8.1.3 Checking the ammonia fixative solution

Insert a MACRO tube and an Erlenmeyer under the coolant glassware. Program a Boric dispense volume. See paragraph 5.

Press «H3BO3», this starts the dispensing pump, it takes liquid from the 2 litres tank to Erlenmeyer flask. Check, visually that the liquid gets the Erlenmeyer flask. And the tube outlet is immersed under the Boric surface.

Run this operation when the dispensing tube must be primed. If tubing is not primed the dispensed volume will be wrong.

8.1.4 Checking the NaOH pump and tubing

Insert a MACRO tube. Program a dispense NaOH volume. See paragraph 5.

Press «NaOH», this starts the NaOH dispensing pump, it takes liquid from the 2 litres tank to sample tube. Check, visually that the liquid gets the sample tube.

Run this operation when the dispensing tube must be primed. If tubing is not primed the dispensed volume will be wrong.

8.1.5 Checking distillation «Cleaning»

This operation is often used when working with the PRO-NITRO «S» because it performs a cleaning of the distillation glassware. It is also used for unit warm-up.

Place a MACRO sample tube, and the Erlenmeyer flask, with near 25ml of distilled water.

Select MANUAL mode. Press «STEAM», this must start the steam production and the coolant water circulation. This is the distillation.

After the warm-up, the distillation rate is between 30 - 40 ml / minute.

Left 5 minutes distilling and then press «STOP» to stop distillation.

8.1.6 Checking sample tube empty function

Place a MACRO sample tube with around 100ml distilled water.

Select MANUAL mode. Press «EMPTY/STOP». This starts the pump to empty the sample tube to the waste output.

8.2 Checking Nitrogen Recovery

This checking is widely used to certify that the unit works properly.

Run some of these test for training purposes before starting to analyse valuable samples, until the operator will be familiar with the unit.

It consists to prepare some samples that have a known Nitrogen content. It is distilled, titrated and the nitrogen detected is calculated.

The rate (%) of nitrogen detected versus the sample Nitrogen content gives us the **Nitrogen Recovery**.

To allow this checking, a 100gr ammonia sulphate flask is delivered together with the distiller.

- Weight near 100mg (50mg reagent 0.1N) of ammonia sulphate. The exact weight will be P_0
- The exact quantity of Nitrogen is: $P_1 \text{ (mg)} = P_0 \times 0,212$
- Add around 25ml of H_2O
- Distillate adding 15ml of NaOH
- Use HCl 0.25N (0.1N) as a titration reagent.
- Once distilling and titration (see 6.5) are finished, calculate the detected nitrogen (mg) : P_2
- The recovery is calculated: $R(\%) = P_2 / P_1 * 100$
- Acceptable values are from 99.0 to 101,5 %.

Note: Some literature sets that the acceptable values for recovery must be above 99.5%. To get this result you must work with high precision and use reagents with high purity degree.

8.3 Blank test

This checking must be run on each batch in order to quantify the absence of ammonia on the reagents.

It consists on distilling and titrating a sample of 25ml of distilled water.

Add 15ml of NaOH on the distillation.

Run some blank tests in order to get a constant value of blank. If this value is repeated it could be subtracted from the detected Nitrogen.

How many ammonia sulphate must I put on the test sample?

The HCl volume (Vol_{HCl}) for a quick analysis must be between 6 and 12 ml.

From the expression:

$$P_1 \text{ (mg)} = 14 \times Vol_{HCl} \text{ (ml)} \times N_{HCl}$$

We arrange it as:

$$P_0 \text{ (mg)} = \frac{14 \times Vol_{HCl} \text{ (ml)} \times N_{HCl}}{0.212}$$

And yields the table:

Normality	Sample amount Sufate ammon.
0,05	20 ... 40 mg
0,1	40 ... 90 mg
0,25	100 ... 200 mg
0,5	200 ... 400 mg

The volume consumption of HCl (Vol_{HCl}) on the titration can be forecasted from the expression:

$$P_1 \text{ (mg)} = 14 \times Vol_{HCl} \text{ (ml)} \times N_{HCl}$$

thus we obtain that HCl must be:

$$Vol_{HCl} \text{ (ml)} = 14 / [P_1 \text{ (mg)} \times N_{HCl}]$$

8.4 Checking Nitrogen recovery with Acetanilide

This test is widely used to certify the whole Kjeldahl analysis process, including the digestion, distillation and titration steps.

It consists in preparing some samples that have a known Nitrogen content. The samples are digested, distilled, titrated and the nitrogen detected is calculated.

Run this test some times to get familiar with the equipment and the kjeldahl analysis before analysing valuable samples.

The rate (%) of nitrogen detected versus the sample Nitrogen content gives us the **Nitrogen Recovery**.

To run this test acetanilide is used.

Preparing the sample:

Nitrogen content of an acetanilide sample:

Formulae: C_8H_9NO

Molecular weight: 135.17

Factor: $14 / 135.17 = 0.1035$ mg of Nitrogen = $0.1035 * \text{mg Acetanilide}$.

- Weight around 1g of acetanilide. The exact weight P_0
- The Nitrogen content is: $P_1 \text{ (mg)} = P_0 \times 0,1035$

Sample digestion:

- Place the acetanilide on a MACRO tube, add 10ml of sulphuric acid 98% (Example PANREAC 173163) and one tablet of Kjeldahl catalyser.

- Digest at 400°C for 1h. (The result must be a transparent blue liquid)

- Let cool at ambient temperature. Add 25ml of distilled water on each tube. Operate with caution to avoid splash. Water reaction over the sulphuric acid could be violent.

Distillation:

- Distillate by adding 75ml of NaOH. Use HCl 0.25N (0.1N) for titration.

- Once sample is distilled and titrated following chapter 6, the detected Nitrogen is calculated as: P_2

- Recovery is calculated as: $R(\%) = P_2 / P_1 * 100$

- Acceptable recovery is from 99.0 to 101.5 %.

Note:


If MICRO size tube is used, weight 0.3g of acetanilide, add 5ml of sulphuric acid, and one tablet of Kjeldahl catalyser. On the dilution add 10ml water.

Distillate for 6 minutes adding 25ml of NaOH.

Note:

Some literature sets that the acceptable values for recovery must be above 99.5%. To get these results you must work with high precision and use reagents with high degree of purity.

Clean regularly rest of NaOH between the closing joint and the splash protector.

Clean periodically the NaOH tank and fixative solution tank hose connector.

9 Maintenance

The distiller PRO-NITRO «S» uses NaOH and concentrated acids for the analysis process. These substances are very aggressive and will degrade some parts of the unit.

To avoid derating and increase distiller's life a suitable cleaning routine should be done.

9.1 Cleaning the unit

9.1.1 Clean the distillation circuit

After any session of analysis, distillate around 50ml of water for 5 minutes in order to remove any rest of NaOH on the distillation glassware. See 6.2

9.1.2 Cleaning the chamber

The chamber is protected with a chemical resistant coating. But it should be cleaned, daily, to increase its life.

Use a well damp cloth to clean the chamber and its components to remove any rest of NaOH.

Clean the tray with abundant water.

9.1.3 Cleaning of closing joint

Near each three months clean, with water, the possible rest of NaOH observed on the closing joint.

Extract the closing pulling it and turning it downwards.

9.1.4 Cleaning the NaOH and Boric pipping (monthly)

Clean with distilled water rests of NaOH and fixative solution the tank hose connectors

- Cleaning the 2 Litres tank: Clean it with abundant water for a few minutes.
- Cleaning the NaOH: Fill the NaOH tank with water, connect the hose connector, Insert a sample tube, and press the NaOH button until fill the sample tube.

Repeat with Boric pump.

9.1.5 Cleaning the Steam generator

Stop the unit by the main switch. Open the manual valve (Empty steam generator valve) and wait until it becomes empty. Close valve and restart the unit. The steam generator refills automatically as explained on paragraph 4.

9.2 Hose inspection (monthly)

Visually check the hose condition and its connections to components.

To avoid premature derating of the unit due to liquid leakage, replace as soon as possible any hose that have any leakage or damage.

9.2.1 Hose replacing

Take these cautions to hose replacing:

- Unplug the unit from the electrical power supply.
- Try, as you can, (see 9.3) that the hose to remove is empty.
- Use only tubes material as indicated on 11.2.
- Use clamps as is factory assembled.
- Once the hose is replaced, it must be checked for leakages.

9.3 Preparing the unit for shipping or storing

Once the unit has been used with NaOH, proper actions should be taken for its transportation or storing.

Initially clean the unit according to chapter 9.1

9.3.1 Remove completely the NaOH from the unit

- **Clean the NaOH tank:**

Remove the NaOH 2 litres tank from the unit.
Clean it with abundant water for some minutes.


- **Cleaning the NaOH dosing circuit:**

Fill the NaOH 2 litres tank with water.
Connect to NaOH circuit.
Insert a MACRO tube on the unit. Press «NaOH» several times. Remove the MACRO tube when filled. Insert an empty one, up the 2 litres of water has circulated across the NaOH dosing circuit.

Repeat the operation with the Fixative solution tank (on the menu press «H3BO3»).

CAUTION
The NaOH is very corrosive for the skin.
Use chemical resistant gloves.
Take proper cautions procedures for NaOH handling.

9.3.2 Remove completely the water from the unit

- **Emptying the water reservoir and the steam generator:**

1. Open the tap "Steam generator empty valve".
2. Start the unit, and wait that the pump empties all the water from the water reservoir. Alarm «STEAM» will be activated. Stop the unit by main switch, Wait 4s and restart.

- **Emptying the water from the cooling circuit.**

1. Open the "Coolant glassware valve" See picture.
2. Blow for the end of the **output** coolant water. The water will **come out** from the valve.

Coolant glassware valve

CAUTION
Do not ship of the unit if the chemicals (NaOH, H3BO3) and water have not been completely removed.

During transportation, chemicals or water leakages could fall over important components and cause severe damages.

During transportation, chemicals or water leakages could fall outside of the package.

9.4 Level sensor cleaning

Inside the steam generator a constant water level is necessary.

To provide this water level there is a conductivity based level probe (a steel rod) that when it loses contact with the water, it shows a high resistance to the electronics and then the microprocessor orders the water pump to feed the steam generator with water from the reservoir up to the probe level.

After some time working with this unit (some months), some residues are settled on the probe causing a malfunction.

Effects of this problem:




The residues short-circuit the path between the probe and the steam generator walls (connected to earth) and so the sensor shows a LOW resistance to the electronics. It is the same signal as if the steam generator has the correct water level. This makes that the microprocessor does not supply water and this will probably cause the safety thermostat action.

This means that the probe must be removed and clean. You will need a cross screwdriver and a 12mm spanner and a 19mm one.

The cleaning frequency depends on two factors:

- The hardness of the water used.
- The usage of the unit.

As a general rule, please do it every **6 months**.

- Remove the left hand side cover (6 screws).
- Locate the steam generator.
- Unplug the probe connector (white cable with FASTON connector).
- Unscrew the sensor with the 19mm spanner.
- Unscrew the probe with a 12mm spanner.
- Remove the steel rod from the PTFE (white part).
- Clean the steel rod and the white PTFE part from the calcium deposits. Use diluted HCl or any home cleaner product for calcium deposits.
- Assemble it again.

10 Technical specifications


Mains voltage:	230V 50/60Hz (See plate) Other power supply voltages could be available (Ask to your dealer)
Electrical power:	1900W (See plate)
Maximum consumption:	8.5 A
Fuses:	10 A fast. Dimensions: 20x5 mm
Power supply cord:	Phase, Neutral, Earth (section 1mm ²)
Weight (Net):	35 Kg
Dimensions: W x L x H	500 x 490 x 740 mm
Water reservoir volume:	6 litres
NaOH, H3BO3 tank volume:	2 litres
Nitrogen recovery:	>99,5%
Detection range:	0.2 to 200 mg
Reproducibility:	Relative Standard deviation) < 1%
Distillation rate:	approx. 30 / minute
Analysis time for 50mg of N:	6 minutes. approx
Analysis time for 100mg of N:	9 minutes, approx
Cooling water entry connector:	3/4"
Cooling water output connector:	3/4"
Remaining steam output connector:	3/4"
Dosification for NaOH pump:	0 to 250 ml ±5%
Dosification for H3BO3 pump:	50 to 65 ml ±5%
Sound level:	Lower than 60 dBA
Ambient conditions:	Inner use. Altitude up to 2000m Temperature from 5°C to 40°C. Relative maximum damp 80% for above 31°C temperatures, decreasing linealy to 50% for 40°C.
Overvoltage category:	Category II
Pollution degree:	2
Glass material:	Borosilicate 3.3
Door and front material:	ABS injected
Carriage material:	Iron with overlap
Chamber material:	Aluminium with RILSAN coating chemical resistant.
Steam generator material:	Stainless steel.
Hosepipe material:	Silicone, synthetic fluorated elastomer.
RS-232 output:	Connector 9 pins male. Pin 1 (Signal) and pin 5 (GND).

11 Technical service


For repairing and substitution of deteriorated parts please contact J.P. SELECTA s.a.u. or its distributor. They will inform you about the technical services authorized in your area.

Find a spare parts list for easy buying and identifying.

11.1 Spare parts

32052

47046


47047

28129

- 16087** Electrovalve N. A.
- 16067** Water electrovalve.
- 24218** Water pump. sample emptying.
- 24278** KNF pump (NaOH and fixative solution)
- 28129** Lock black rubber.
- 28403** Shelf.
- 32052** Steam generator 1800W.
- 47021** Collector flask.
- 47046** Splash protector.
- 47259** Anti feedback valve.
- 29336** Power plate
- 29337** Control / front plate
- 29530** Level

Hoses:

- | | | |
|--------------|---------------------------------|------------------------------|
| 46030 | Teflon tube 8x10 | Tube immersed on the sample. |
| 46041 | Silicone tube 8x14 | Water and Steam. |
| 46042 | Silicone tube 6x10 | Water and Steam. |
| 46143 | Pressure hose | Tap water entry. |
| 46084 | ISO VERSINIC Black 10x14 | Distilling Glassware joint. |
| 46139 | ISO VERSINIC Black 6x10 | NaOH, Boric. |


24218


WARNING

Before doing any works inside the equipment, please disconnect it from the mains.

11.2 Hydraulic tubing diagram

11.3 Wiring diagram

12 Kjeldahl analysis step by step

12.1 Sample preparation

- Triture, homogenize and mix the sample.
- Weight between 1 and 2 grams of the sample.
- In low nitrogen content sample, (waste water, etc...) take enough samples to get 5mg nitrogen.

12.2 Digestion

- Add between 10-15ml (MACRO tube) of H_2SO_4 (for MICRO tube maximum H_2SO_4 is 5ml).
- Rig up a system for smoke extraction or scrubber with NaOH.
- Make this digestion in 3 steps:

1. In function of the level of water in the sample. Begin digestion at 150°C for 20 or 60 minutes.
2. Make a second step by 280°C for 30 minutes to reduce white smoke production.
3. Continue the digestion with 400°C for 60 - 90 minutes.

Visual checking: the result is a neat transparent liquid with a clear blue coloration or green or yellow depending on the catalyser used. It must not still black rests stuck on the tube's wall.

Note: During digestion you must control foam production in the samples. If there is too much foam, length step 1.

12.3 Dilution

- Take out the sample tube off the digester block and let it cool to ambient temperature.
- Add about 25ml of water in each tube (10ml for MICRO tube).
- Slow add water, moving the tube to not let solidify the sample. Heat the tube if necessary. (e.g. introducing the tube into the digester block)
- Let it cool again.
- To avoid nitrogen loss and violent reactions do not introduce the hot tube inside the distiller.

12.4 Distillation

- Check the tank level of NaOH and Boric acid. Check the dispensing bombs are primed and give the correct volume.
- Put a 250ml Erlenmeyer in the outlet cooler and add 50ml Boric acid and some indicator drops.
- Program a 50 to 75 ml NaOH dosage (25ml for microtube).
- Insert the tube with the sample in the distiller.
- Distillate till recollecting 250ml in the Erlenmeyer (50ml Boric acid+200ml of distillate).

Visual checking: Once the NaOH has been added, the sample must get a blue coloration, if not add more NaOH.

12.5 Titration and calculations

- Titrate the distilled with HCl or H_2SO_4 0.1 or 0.25 N till colour changes (end point: pH 4.65).
- Calculate: mg of Nitrogen = $14 * Volume_{HCl} * Normality$
- Calculate protein contents correct with the right factor based on the sample nature. (6.25 default value)
- Periodically make a blank test and subtract it from the result.

13 Guidelines for IQ and OQ

Some laboratories need the installer make an IQ document to certificate the good installation of this equipment. As well as an OQ document to certificate the good working of this equipment.

The installer can elaborate those documents following the examples:

Checking list for IQ (Installation Qualification):

- Installer name.
- Manufacturer name:
- Unit identification: Code, model and serial number.
- Unit location: Laboratory, Person in charge.
- Checking packing list:
- Unit location matches the requirements?

Water and electricity requirements are suitable?

- Electrical supply.
Socket:
 - Ground terminal?
 - Suitable to power consumption. (16A)
 - Voltage according plate $\pm 10\%$
- Water connection.
 - Water pressure bigger than 1bar.
 - Hose to waste are properly fixed and not blocked.
- Water from steam generator: Low level of calcium, conductivity above 20 MicroSiemens.
- Personnel training: Operator personnel understand the instruction manual?.
- Date and name of who has made the Installation (Signed).

Checking list for OQ (Operation Qualification):

- Installer name.
- Manufacturer name:
- Unit identification: Code, model and serial number.
- Unit location: Laboratory, Person in charge.
- Check the equipment distillate correctly. See: 8.1
- Check nitrogen recuperation with ammonia sulphate sample. See 8.2.
Note : ammonia sulphate weight, HCl normality, HCl consumption and obtained recovery.
- Check the white test See: 8.3
Note: Ammonia sulphate weight, HCl normality, HCl consumption and obtained recovery.
- You can check the nitrogen recovery with acetanilide sample. See: 8.4.
Note: Acetanilide weight, HCl normality, HCl consumption and obtained recovery.
- Date and name of who has made the Installation (Signed).

14 Warranty

The warranty covers the defective components or faulty assembling of the unit for 1 year since the purchasing date.

Keep purchasing bill as the only valid document for the purchasing date.

The warranty does not cover the damages caused for a unit malfunction. The unit must be suitably qualified before to be used for valuable samples processing.

Return good are not accepted without previous J.P. SELECTA, s.a.u. approval.

15 Troubleshooting table

- **The green indicator «ON» does not light.**

Check fuses and wiring connections cables.

Check main switch.

- **The equipment still adding water till it breaks out of the sample.**

The equipment does not detect the water level for low conductivity. Add about 250ml of tap water in the water tank.

Check the cable (white) of the level probe of the steam generator.

- **At the end of distillation the sample tube content is absorbed to the steam generator.**

Check the steam hose is not obstructed.

Locate the AIR electrovalve:

Check it actives, check the voltage with a tester (230V).

Check (detaching it) it is not obstructed.

- **The alarm (2) for hose presence and lock door is always warned and block the unit running.**

Check the door is correctly closed.

Check the microswitches are o.k. (you will hear a click) when you close the door or lower the lever.

Check the cable (grey) continuity.

- **The refrigerating water does not enter**

Check the tap is open.

Locate the refrigerating water input valve, check the voltage with a tester(230V). If there is voltage and the valve does not open knock it. That is due to a small lime rest blocking the water pass. If knocking is not enough, detach the valve to clean the water pass.

- **Security thermostat Alarm warned and the equipment is blocked.**

This alarm warns when there is an excessive temperature in the steam generator. In this case the heating element stay disconnected by the thermostat and the control electronics block all the unit functions.

To rearm the thermostat, wait some minutes to let the temperature of the generator falls down to more normal level, and pulse the rearm switch.

The cause of the thermostat disconnection is that the steam generator lack of water. Check the hydraulic wiring to generate water, tank, pump, hoses, level probe.