DESTILADOR L-3
WATER DISTILLER L-3

4903000

INFORMACIÓN GENERAL

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GENERAL INFORMATION

- 1) Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributor immediately.
- 2) Do not install or use the equipment without reading this handbook before.
- 3) This handbook must always be attached to the equipment and it must be available for all users.
- 4) If you have any doubts or enquiries, please contact with your supplier or J.P. Selecta's technical service.
- 5) **IMPORTANT! J.P. SELECTA WILL NOT ACCEPT ANY EQUIPMENT TO BE REPAIRED IF IT IS NOT DULY CLEANED.**
- 6) If any modification, elimination or lacking in maintenance of any device of the equipment by the user transgresses the directive 89/655/CEE, the manufacturer is not responsible for the damage that can occur.
- 7) Do not use the equipment with liquids which can give off vapours capable of making explosive mixtures.

PACKING LIST

The standard equipment consists of the following components:

Piece	Qty.
Main unit	1
Heating element	1
Water inlet tube	1
Instruction manual	1

EQUIPMENT DESCRIPTION

L-3 Water Distiller characteristics:

- Continuous and completely automatic working.
- Pressure switch, set at 1 bar of entry pressure water. If the pressure is lower, the pressure switch cuts the heating to avoid the risk that an insufficient cooling may cause breakage of the glass distillation.
- Maximum water use due to the circuit design. The incoming water is used for cooling and maintaining the water level into the boiler, and the surplus exits through the outlet (19).

TECHNICAL DATA

- Distilling capacity: 3 litres/hour
- Conductivity: 3 μ S
- Water consumption: 60 litres/hour
- Consumption: 2.200 W.
- Voltage: 220 V / 50 Hz
- Weight: 4.5 Kg.
- Size (Height Length Width): 345x560x190 mm.

Supply voltage could be different. Check the characteristics plate.

ATTENTION

Since distilled water is in contact with the atmosphere before the conductivity measurement can be carried out, its value is actually much less (up to 10 times less) than what it is actually measured. The measurement of any ultrapure water in contact with atmosphere has a conductivity of 1-2 μ S/cm and a pH of 5, due to a small amount of CO₂ (0.05 ppm) absorbed from the atmosphere.

At a general level, features on distilled water are:

	Conduc- tivity (uS/ cm)	Resistiv- ity (M Ω m · cm)	Silicates (mg/l)	Heavy metals (mg/l)	Reduced Perman- ganate (min)	Sodium (mg/l)	Calcium hardness (mg/l)	Ammo- nium (mg/l)	Bacteria (u.f.c./ml)	pH (at 25 °C)
Once distilled	10 – 2	0,1 – 0,5	1 – 0,5	1 – 0,5	30	5 – 2	3 – 1	0,01	<10	5 – 6,5
Bidistilled	2 – 1	0,5 – 1	0,7 – 0,1	0,8 – 0,1	30	1 – 0,5	0,3 – 0,1	0,01	<10	5 – 6,5

ACCESORIES

Description

- Polyethylene reservoir up to 50L.
- Decalcifier

Code

5903050
703052

INSTALLATION

Assembly the heater element, put it into the glass part (10), thread the adjusting nut on the glass thread and connect to the base of the furniture (14).

The water supply, will be made by a 3/4" GAS thread tap (type washing-machine), where the filter and the thread (2) of the union tube will be placed. The zone (4) is connected to the thread of the water entry (5) located on the base (14).

Connect the mains wire (14) to the power supply.

The frequency of the cleaning is proportional to the salts, metals, lime or other elements that may be found in the used water. We recommend you to incorporate a decalcifier between point 1 and 2, in order to reduce the risk of breaking of the heating element and to lengthen the time of cleaning.

Never use the heating element connector to move the distiller. The heating element can be broken.

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Be sure that the voltage supply matches the power indicated on the characteristics plate.

Do not use the equipment if it is NOT EARTHED.

If you have to change the plug, bear in mind that the yellow/green cable is the protective earth and it is connected into the metallic envelope of the equipment.

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ATTENTION!: This equipment do not have fuse. A circuit breaker and a differential switch must be installed at mains.

	Tensión/Voltage	Magnetotérmico/Break circuit	Diferencial/Differential
4903000	230 V	16 A	16 A/30 mA
	115 V	32 A	32 A/30 mA

OPERATION


Before starting, check that the water level is over the heating element (15) and then pour the surplus from the tube (19) towards the nearest drain tube.

Once the distiller is plugged onto the mains, push the switch (12) on to start the unit. A green light indicates that the distiller is ready for distilling.

When the water tap (1) is opened, distillation starts. If the amber light is on, it means that there is not enough water pressure. Open the tap at maximum to get more pressure. If the tap is already at maximum, it means that the installation is not suitable for the distiller.

After ten minutes, the boiling and distillation starts.

Distilled water gets out from the tube (18).

- 1- Tap.
- 2- Upper tube thread.
- 4- Inferior tube thread.
- 5- Water entry.
- 6- Feed tube.
- 7- Refrigerant.
- 8- Level.
- 10- Glass part.
- 12- Switch on/off.
- 14- Furniture base.
- 15- Heater element.
- 16- Steam trap.
- 17- Steam way out.
- 18- Distillate tube.
- 19- Drain tube.
- 20- Cleaning orifice.
- 22- Level orifice.

MAINTENANCE

Before removing the covers, disconnect the equipment from the mains.

Manipulation of the plate internal electronic circuits by unauthorized personnel can cause damages of difficult solution. Make sure to take the equipment to a J.P. Selecta, s.a.u. authorized technical services.

- The impurities of the water can be built up in the boiler (10). To remove them, just we have to insert through the entrance (20), a 10% aqueous citric acid, acetic acid, hydrochloric acid mixture or another decaling product, up to the resistance (15). Let convenient time until the waste has been attacked.

- To empty, disconnect the distiller from the main supply, and overturn the boiler to empty through the entrance. Rinse with water 5 times and dispense the 3 first litres after the washing.

GUARANTEE

This product is guaranteed for one year. The guarantee does not cover damages caused by an incorrect use or causes beyond the control of J.P. SELECTA, s.a.u.

Any manipulation of the equipment by unauthorized personnel by J.P. SELECTA, s.a.u. automatically cancels the guarantee.

SPARE PARTS

To guarantee the safety of the equipment, all spare parts must be purchased from J.P. SELECTA, s.a.

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