**DESTILADOR R-4 RESER
WATER STILL R-4 RESER**

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1.1. Table of contents

1.1. Table of contents	2
1.2. Getting started	2
1.3. Product overview	3
1.4. Product specifications	4
1.5. Technical features	5
1.6. Operating	6
1.7. Warning & cleaning	7
1.8. Trouble shooting	8
1.9. Water flow diagram	8
1.10. Water purity	9

1.2. Getting starting

Your water still has been designed with function, reliability and safety in mind. It is your responsibility to install it in conformance with local electrical codes. For safe operation, please pay attention to the alert symbols through the manual.

This manual contains important operating and safety information. You must carefully read and understand the contents of this manual prior to the use of this equipment.

Warning

Warning alert you to a possibility of personal injury.

Caution

Caution alert you to a possibility of damage to the equipment.

Note

Notes alert you to pertinent facts and conditions.

Hot surface

Hot surface alert you possibility of burning injury by hot instrument surface.

Explosive

Explosive alerts you to possibility of explosion by high pressure.

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1.3. Product overview

R-4 RESER water still produces distilled water from tap water with space saving design and effective distilled water production.

Unique design and safety feature with two water float switches provide automatic operation. By lack of water, the heating is automatically switched off to protect heater from over heating.

Hydraulic thermostat that switched off the apparatus to avoid the over heating.

Exclusive system for limiting the flow of water from the tap water optimizing the flow down to a necessary minimum.

Water still is easy to clean-up by just removing upper part of the body (condenser unit). By controlling cooling water flow, you can save a lot of water wasted before.

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A -> ancho; F -> fondo; H -> altura

1.4. Product Specifications

Model	R-4 RESER
Capacity	4.0 l / hour
Tank capacity (l)	8.0 litres
Dimensions (mm)	620 (W) x 330 (D) x 460 (H)
Weight (Kg)	24 Kg
Heater	3,00 KW
Cooling Water	Approx. 55 l / hour
Material Inner	Stainless Steel
Material Outer	Galvanized Steel & Powder Coating
Electrical Requirements	220 V, 50 / 60 Hz
Safety	Low Water & Over Heating Cut-Out

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
1.5. Technical features

- This is an automatic water still with built-in tank for reserving distilled water.
- With boiling and reservoir tanks made of stainless steel.
- Heater of stainless steel.
- Housing made of galvanized steel and coated with epoxy paint.
- Cable Schuko with ground wire.
- Lighting switch ON/OFF.
- Lighting pilot of distillation function.
- Magnetic switch of `Distillation starting`
- Magnetic switch of `Full tank`
- Over heating hydraulic thermostat.
- Electrovalve that controls the inlet of water to the equipment.
- Built-in tap to supply distilled water.
- Built-in tap to empty the boiling tank for cleaning.
- Exclusive system for limiting the flow of water from the tap water optimizing the flow down to a necessary minimum.

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
1.6. Operation

- 1) Connect the water still to a plug, suitable for an intensity of 13 A.
- 2) Connect to a tap water with thread of 3/4" using the supplied tube.
- 3) Connect this tube to the drain outlet (6) and the another end to the drains, in a lower level than the water still.
- 4) The stopcock (5) must be closed during the usual operation of the water still. It must only opened for cleaning the boiling tank.
- 5) Switch on the switch (1), it will be lighting during all the process. Tap on the water tap. When the water is over the heaters, and it connects the magnetic switch, the lighting pilot will light and will start the distilling process.
- 6) The distilled water is stored in the reservoir tank of the water still. When it reaches the maximum level, the water still is disconnected automatically in order not to flood the reservoir tank. When the level of the distilled water in the reservoir tank falls, the water still will start automatically.

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- 7 -

1.7. Warning & cleaning

1. Make sure the electric supply is in accordance to the given voltage on the name-plate, 50 / 60 Hz, single phase.
2. Place water still on a flat surface.
3. Do not touch top of the instrument. During operation, surface of the top of the water still and steam exhausted from the hole is very hot.
4. Quit the upper lid of the water still and the individual lid of the boiler tank to access to the tank and the heater.
4. A mixture of 10% formic acid, 10% acetic acid and 80% distilled water is recommended for cleaning the apparatus. Fill the solvent into the evaporator to just above the highest lime deposit and switch on the water still. The evaporator is clean after approx 10 to 15 minutes. Pour out the lime and solvent mixture through the outlet (5) and thoroughly rise with water.
6. In the event of overheating, it activates the security device and you should reset the water still pushing RESET (5).

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
Trouble	Check first	Trouble Shooting
Power failure	Check Electric Supply	Plug firmly into the electric supply
Power on but water do not boil	Check water supply	Supply water enough and appropriate presure
	Check heater	Replace heater
Contact sales representative or customer service department		

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ATTENTION

Since distilled water is in contact with the atmosphere before the conductivity measurement can be carried out, its value is actually much less (up to 10 times less) than what it is actually measured. The measurement of any ultrapure water in contact with atmosphere has a conductivity of 1-2 $\mu\text{S}/\text{cm}$ and a pH of 5, due to a small amount of CO_2 (0.05 ppm) absorbed from the atmosphere.

At a general level, features on distilled water are:

	Conductivity ($\mu\text{S}/\text{cm}$)	Resistivity ($\text{M}\Omega \cdot \text{cm}$)	Silicates (mg/l)	Heavy metals (mg/l)	Reduced Permanganate (min)	Sodium (mg/l)	Calcium hardness (mg/l)	Ammonium (mg/l)	Bacteria (u.f.c./ml)	pH (at 25°C)
Once distilled	10 – 2	0,1 – 0,5	1 – 0,5	1 – 0,5	30	5 – 2	3 – 1	0,01	<10	5 – 6,5
Bidistilled	2 – 1	0,5 – 1	0,7 – 0,1	0,8 – 0,1	30	1 – 0,5	0,3 – 0,1	0,01	<10	5 – 6,5