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## 1.1. GENERAL INFORMATION

- The package must be handled carefully. When you unpack it, you should check that what it contains inside is what is indicated in the Packing List section. If any component is missing or otherwise damaged, contact the dealer immediately.
- Please, read the instructions before installing the equipment. The instruction manual will be an inseparable part of the product; for this reason, keep it in good condition for future reference.
- The modification of the components of the product transgresses the directive of use 2009/104/CE; and therefore, the manufacturer is not responsible for the damages that this brings as damage.
- The product must be installed on a horizontal surface and never in places where there are vapours of contaminating or flammable substances.
- Please, dispose of at least the first distillation before using the equipment.
- Keep the packaging during the warranty period.
- Warning and Observation are indicated as follows:

**Warning: This must be strictly followed to avoid physical damage.**

**Note: It contains suggestions and restrictions for the product.**


## 1.2. PACKING LIST

Description	Qty.
Metallic water still	1
Feed water tube	1
Filter	1
Silicone tube 12x18	2
Instruction manual	1

## 1.3. PRODUCT SPECIFICATONS

Specifications	AC-L4		AC-L8
	Distillation capacity (L/H)	4	4
Cooling water (L/MIN)	1	1	1,4
Conductivity 20 <sup>0</sup> C (μS/cm)	2	2	2
Wattage (W)	3000	3000	6000
Voltage (V)	230	110	230
Intensity (A)	13	27	26
Weight (Kg)	8		12
Dimensions, height x width x depth (mm)	435 x 370 x 220		635 x 370 x 260
Safety	Thermostat, level switch		

1.3.1. WATER DISTILLER. FRONT PANEL

1.4. INSTALLATION

1. Connect the equipment to a 230V or 110V power outlet using the power supply cable.
2. Connect the grey hose to the water supply socket. The other end of the hose, with the filter previously installed, will be connected to a 3/4" tap.
3. The silicone tube must be connected to the outlet for distilled water and the other end to a collecting container whose level is below the outlet by about 25 cm.

## 1.5. START-UP

1. Check that the installation has been carried out properly.
2. Open the water inlet tap previously connected to a 3/4" outlet, feed water with a flow of 60 L/H for AC-L4. The flow should be 80 L/H for the AC-L8 still.
3. Turn switch S to the ON position. The green light will come on and after about 10 minutes the distillation process will begin.

**Warning: Always disconnect the power supply before carrying out the cleaning operation**

## 1.6. MAINTENANCE

1. With the switch in the OFF position, and the power cable disconnected from the mains remove the refrigerant.
2. Pour in hydrochloric acid, strong water, or hot water until it comes out of the drain. Leave on for 15 minutes.
3. Pour all liquid and impurities into a drainpipe.
4. Rinse with plenty of water.
5. Repeat steps 2, 3 and 4 if necessary.
6. Carry out a distillation that must be discarded.

## 1.7. SPARE PARTS


Description	Symbol	Cat. Number	Qty.
Contacteur, 230 V	C1	0100196	1
Contacteur, 110 V	C1	YE.DCONT10001	1
Safety thermostat	T	9959056	1
Heater, 110 V	R	YE.DRDES00110	1
Heater, 230 V	R - R1, R2	9959055	1 - 2
Silicone gasket (AC-L8)	---	9959058	1
Silicone gasket (AC-L4)	---	9959077	1
Plastic switch level	N	9959057	1

## 1.8. TROUBLE SHOOTING

Problem	Action
The power cord is connected to the mains, the electrical switch S is ON and the green light does not turn on.	Measure and check that the voltage in the network is 230V. (or 110V)  Check the status of the switch.
Water supply is adequate and distillation does not happen	Following this order: Check level switch N. The circuit closes upwards and activates contactor C1. Replace if necessary. There is current supply to the input of contactor C1 and it does not work: replace it. There is current supply to the resistance: R, R1, R2 and it does not heat. Replace the resistance.
In the event of overheating that activates the safety mechanism.	Press the RESET button.


**Note: It will not be allowed to repair any water still that is not properly cleaned.**

### Notice to customers:


The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.


1.9. ELECTRIC DIAGRAM

- S > Interruptor con piloto neon / Neon light switch
- T > Termostato de seguridad / Security thermostat
- N > Interruptor de nivel / Level switch
- C1 > Contactor / Contactor
- R > Resistencia calefactora / Heating resistance

S > Interruptor con piloto neon / Neon light switch  
 T > Termostato de seguridad / Security thermostat  
 N > Interruptor de nivel / Level switch  
 C1 > Contactor / Contactor  
 R > Resistencia calefactora / Heating resistance

S > Interruptor con piloto neon / Neon light switch  
 T > Termostato de seguridad / Safety thermostat  
 N > Interruptor de nivel / Level switch  
 C1 > Contactor / Contactor  
 R1,2 > Resistencia calefactora / Heating resistance

## 1.10. WARRANTY

1. The product covered by this Certificate is protected for one (1) YEAR after the date of the purchase delivery note of the product.
2. During this warranty period, in case of manufacturing error of the product, the user has the right to its servicing, its replacement (except if one of this options was impossible or disproportionate), to have a reduction in the price and to the resolution of the contract (when the purchaser could not require the servicing or the replacement). The servicing will be free of charge for the user. All the servicing and replacements will be carried out in a reasonable period of time and without major inconvenience for the consumer.
3. The breakdowns, flaws or caused damages in the product will be excluded as consequence of: abuses, improper use, exposure to adverse weather conditions, spilled batteries, blows, falls, bad or inadequate installation, power failure, waste by normal use, servicing carried out by different people indicated by us, as well as modification of this warranty and/or another external causes to the product.
4. Glass parts of our distillers are excluded of this warranty.
5. This warranty will be automatically cancelled in case of mishandling of the product or if it has been used in violation of the standards contained in this manual, regarding use, assemble, operation and service.
6. The resolution of the contract will not proceed when the lack of conformity is of little importance.
7. The warranty does not affect the rights available to the consumer in accordance with the provisions of the law.
8. In case the user needs technical help regarding the guaranteed product, the user will contact the distributor.
9. For any questions that may arise due to the interpretation or execution of this warranty, the provisions of the applicable national legislation will be followed.
10. Please, attach a copy of the delivery note or the purchase invoice with the Warranty.