**DESTILADOR  
DISTILLER**

**D-10 LARGE 4903015**

## 1.1. Tabla de contenidos

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## 1.3. Security

- *Magnetic switch that avoids water still worked in case insufficient water level.*
- *Safety thermostat protects still in the event of water supply failure.*
- *Automatic device to cut water supply in the event of electrical failure.*
- *Device to fill a container without worrying when it finished. It cuts water and electrical supply to avoid useless expenses of water and electricity.*

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#### 1.4. Packing list

- *Water still.*
- *Water inlet tube.*
- *One filter.*
- *One brace 8-16.*
- *Silicone tubes set.*
- *Two nuts Ø 32 with liner*
- *Two nuts Ø 18*
- *Two bent olives*

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## 1.5. Recommendations

### Electric:

- *Electrical wiring must be, at least, of 6 mm<sup>2</sup> section.*
- *It must have ground.*
- *Electrical system must let intensity of 32 A pass.*

### Outlet:


- *The tap water should have screw mouth of 3/4"*

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## 1.6. Assembly

1. For your security, the water still is supplied separately: the metallic base and boiler and condenser, water inlet level and security level.
2. When you assemble the condenser over the boiler, make it carefully fitting on the o-rings. The distilled water outlet should be as you see in the bottom figure.
3. Assemble the water inlet level threading the nut onto the boiler and the silicone tube onto the condenser.
4. Assemble the security level threading the nut onto the boiler on its back side and connect the wire to the base.
5. Close the draining tap if it was opened - it is located on the back side of the water still -.
6. Connect the apparatus to the water mains with the supplied tube.
7. Connect the apparatus to the mains with the supplied electrical wire.

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## 1.7. Controls

1. Revolving switch, 3 positions, located on the right side of the apparatus. Functions: distilling, stop and cleaning.
2. Green pilots: distilling or cleaning.
3. Red pilots: FULL TANK OF DISTILLED WATER and NOT ENOUGH WATER

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## 1.8. Operating

1. *Connect the apparatus to the suitable power supply.*
2. *Put the selector in Cleaning position.*
3. *Turn on the tap water. The filling of the apparatus starts.*
4. *When the water reaches the heaters, put the in Distilling and if the water flow is enough the apparatus starts distilling when the water reaches the security switch.*
5. *The apparatus can work continuous during all the day.*

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## 1.9. Maintenance

- You should clean the still regularly. So, in this way, the results will be better and you can get a still more durable.

- When you finished using the water still, it is recommended to empty the used water and fill again with clean water. You will save many cleanings.

- The periodicity of the cleanings depends on salts and hardness of the supplied water.

- It is strongly recommended to clean the filters placed in the water tube.

Steps to follow:

1. Put the selector in *Cleaning position*.
2. Turn on the *draining tap* on the back of the apparatus. The water keeps going in and out at the same time.
3. Or, If you want to empty completely the apparatus, turn on the *draining tap* with the selector in *OFF* position and the apparatus will be empty.
4. Add a detergent (as *LITOPOL*) and clean the embedded impurities.

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*1. The water stills does not work:*

- Check the plug.*

*2. The red pilot "Not enough water" is light:*

- Check that there is water flow in the tubes.*
- Check that the dirt is not blocked the water flow.*


*3. The apparatus is filled with water and it does not starts working:*

- Check that the water flow is enough.*
- Check that the security switch is working properly.*

*4. The apparatus has water leakages:*

- Check all the external nuts and tighten them, if it is necessary.*
- If the leakages are internal, call manufacturer Technical Service.*

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

1. In order to connect the level sensor, you should connect one tube to the distilled water outlet (9) and another tube to the vacuostat (10). You have to put the level sensor into the bottle of distilled water (11).

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### **ATTENTION**

Since distilled water is in contact with the atmosphere before the conductivity measurement can be carried out, its value is actually much less (up to 10 times less) than what it is actually measured. The measurement of any ultrapure water in contact with atmosphere has a conductivity of 1-2  $\mu\text{S/cm}$  and a pH of 5, due to a small amount of  $\text{CO}_2$  (0.05 ppm) absorbed from the atmosphere.

At a general level, features on distilled water are:

|                   | Conductivity<br>( $\mu\text{S/cm}$ ) | Resistivity<br>( $\text{M}\Omega \cdot \text{cm}$ ) | Silicates<br>( $\text{mg/l}$ ) | Heavy<br>metals<br>( $\text{mg/l}$ ) | Reduced<br>Permanganate<br>(min) | Sodium<br>( $\text{mg/l}$ ) | Calcium<br>hardness<br>( $\text{mg/l}$ ) | Ammonium<br>( $\text{mg/l}$ ) | Bacteria<br>(u.f.c./ml) | pH<br>(at 25°C) |
|-------------------|--------------------------------------|---|--------------------------------|--------------------------------------|----------------------------------|-----------------------------|--|-------------------------------|-------------------------|-----------------|
| Once<br>distilled | 10 – 2                               | 0,1 – 0,5   | 1 – 0,5                        | 1 – 0,5                              | 30                               | 5 – 2                       | 3 – 1                                    | 0,01                          | <10                     | 5 – 6,5         |
| Bidistilled       | 2 – 1                                | 0,5 – 1   | 0,7 – 0,1                      | 0,8 – 0,1                            | 30                               | 1 – 0,5                     | 0,3 – 0,1                                | 0,01                          | <10                     | 5 – 6,5         |