**COMBATEMP S -30°C**  
**COMBATEMP L -30°C**  
**COMBATEMP S -40°C**  
**COMBATEMP L -40°C**

**2101040**  
**2101041**  
**2101181**  
**2101182**

**ARCONES CONSERVADORES**  
**CHEST FREEZERS**

# Indice

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## ENGLISH VERSION

## 1 GENERAL INFORMATION

1) Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributor immediately.

2) Do not install or use the equipment without reading this handbook first.

3) This handbook must always be attached to the equipment and it must be available for all users.

4) If you have any doubts or enquiries, please contact your supplier or J.P. Selecta's technical service.

**5) IMPORTANT! J.P. SELECTA WILL NOT ACCEPT ANY EQUIPMENT TO BE REPAIRED IF IT IS NOT DULY CLEANED.**

6) If any modification, elimination or lacking in maintenance of any device of the equipment by the user transgress the directive 89/655/CEE, the manufacturer is not responsible for the damage that can occur.

7) Do not use the equipment with liquids which can give off vapours capable of making explosive mixtures.

## 2 PACKING LIST

The standard equipment consists of the following components:

| DESCRIPTION                | CODE                                     | QUANTITY         |
|----------------------------|--|------------------|
| COMBATEMP-S<br>COMBATEMP-L | 2101040 / 2101181<br>2101041 / 2101182   | 1                |
| Basket                     | 2101040<br>2101041<br>2101181<br>2101182 | 1<br>2<br>3<br>2 |
| Key                        | ---                                      | 1                |
| Instruction manual         | 80079                                    | 1                |

### 3 TECHNICAL FEATURES

Voltage supply 115-230V 50/60 Hz according to the characteristics plate indications.

| CODE    | Capacity<br>litres | Freezing cap.<br>in 24h/L | Temp.<br>up to °C | Chamber dimensions<br>depth/length/width | Overall dimensions<br>depth/length/width | Power<br>W | Weight<br>Kg |
|---------|--------------------|---------------------------|-------------------|--|--|------------|--------------|
| 2101040 | 203                | 28                        | -30               | 55 / 71 / 50                             | 86 / 83 / 73                             | 160        | 65           |
| 2101041 | 352                | 38                        | -30               | 70 / 100 / 52                            | 91 / 113 / 76                            | 200        | 75           |
| 2101181 | 203                | 30                        | -40               | 55 / 71 / 50                             | 86 / 83 / 73                             | 230        | 70           |
| 2101182 | 352                | 40                        | -40               | 70 / 100 / 52                            | 91 / 113 / 76                            | 400        | 80           |

### 4 EQUIPMENT DESCRIPTION

COMBATEMP chests are designed for adjustable temperatures from -10°C to -30°C and -40°C. Inner chamber of embossed aluminium with guide for supporting baskets or shelves and outer casing epoxy-coated. Top and double door in ABS plastic easy to clean. Flap door with counterbalanced hinges, handle and key lock. Temperature regulation is done by a digital electronic regulator.

### 5 INSTALLATION

**CAUTION!!! IMPORTANT FOR YOUR SAFETY**

Do not use the equipment if it is not earthed.

Be sure that the voltage supply is the same that the one indicated on the characteristics plate.

If you change the plug, bear in mind the following:

|                     |          |
|---------------------|----------|
| Blue cable:         | Neutral. |
| Brown cable:        | Phase.   |
| Yellow/green cable: | Earth.   |

### 6 OPERATION

STARTING UP:

Before the introduction of the first charge, it is advisable to leave the chest working for about two hours at the preset temperature.

1. Connect the equipment into the mains.
2. Press the main switch (1).
3. Adjustment the temperature.

There are two models of cabinets: -30°C and -40°C. Each control panel works in a different way.

## 6.1 Models 2101040 and 2101041 (-30°C)

### Electrical connection

Only operate the appliance with **alternating current (AC)**.


The permissible voltage and frequency are indicated on the type plate. The position of the type plate is shown in the section entitled **Description of the appliance**.

The socket must be properly earthed and protected by a fuse. The tripping current of the fuse must be between 10 A and 16 A.

The socket must not be situated behind the appliance and must be easily accessible.

Do not connect the appliance using an extension cable or extension socket.

Do not use stand-alone inverters (conversion of direct current to alternating current/three-phase current) or energy-saving plugs. Risk of damage to the electronic control system!


### Saving energy

- Always ensure good ventilation. Do not cover ventilation openings or grille.
- Avoid positioning the appliance in direct sunlight or near cookers, radiators and similar sources of heat.
- The energy consumption depends on the installation conditions, e.g. the ambient temperature.
- Keep the time the appliance is open to a minimum.
- Store food logically.
- Keep all food properly packed and covered so as to avoid its frosting up on the outside.
- Always allow hot food to cool to room temperature before placing in the appliance.
- If there is a thick layer of frost in the appliance: defrost the appliance.

Accumulated dust increases the energy consumption:

- Regularly clean the motor compartment ventilation openings to keep them free from dust.

### Operating and control elements

- ① Temperature control
- ② Operating lamp
- ③ Alarm lamp
- ④ Superfrost button
- ⑤ Superfrost lamp

### Switching the appliance on and off

You are recommended to clean the appliance before switching it on for the first time (see "Cleaning").

#### To switch the appliance on

Connect the appliance to the mains > the appliance is switched on.

- The green operating lamp comes on. → On
- The red alarm lamp comes on. → Alarm

The red alarm lamp goes out as soon as the freezer compartment is cold enough.

#### To switch the appliance off

Pull out the mains plug or set the temperature control to 0. →

#### Important

**When the temperature control is turned to "0", only the refrigeration system is switched off. Before cleaning, always disconnect the appliance from the mains. Remove the plug.**

### Setting the temperature

Using a coin, turn the arrow on the temperature control to a setting between 1 and 4.

When the control is turned to 4, the temperature inside the appliance is at its lowest.

For long-term storage of frozen food a temperature of at least -18°C is recommended. →

If the chest is filled up to the bottom stacking mark (see "Freezing"), the temperature control can be set to the 1 setting. →

### Control elements

The **red alarm lamp** comes on if the temperature in the freezer compartment is not cold enough. → Alarm

#### This can be caused by:

- a long power failure;
- a fault in the appliance.


In each case, make sure food has not defrosted or perished.

The **"Superfrost" lamp** lights up when the "Superfrost" facility is switched on for freezing large quantities of fresh food (see "Freezing").

## 6.2 Models 2101181 and 2101182 (-40°C)

### Operating and control elements

- ① **ON/OFF** button (to switch the appliance on and off)
- ② Selection buttons
- ③ **Set** button (Enter)
- ④ **Audible warning** on/off button
- ⑤ Button for calling up stored alarm events
- ⑥ Temperature display

#### Control elements

- Compressor is running
- LED flashing - refrigeration unit switches on after a delay. The compressor will start automatically after the pressure in the refrigerant circuit has equalised.)

**AUX** Temperature display via product sensor is activated

Alarm function

If appears in the display, the appliance has a fault. Consult your nearest customer service point.

#### HACCP (Hazard Analysis Critical Control Point)

The HACCP display means that the power supply and interior temperature of the appliance are recorded.

If HACCP flashes in the display, there has either been a power failure or the temperature in the appliance exceeded the permissible range.

### Switching the appliance on and off

Connect the appliance to the mains - the display reads **OFF**.

**Switching the appliance on:** Keep the **ON/OFF** button pressed for approx. 5 seconds - the display reads **ON**.

No alarm is displayed or sounded when the appliance is switched on for the first time.

If the appliance is disconnected from the mains for a long time after it has been switched on for the first time and if the temperature inside the appliance rises above the upper alarm limit, this will be detected as a fault by the electronic control system (**HACCP** appears in the display).

When the appliance is switched on again, this display must be reset as shown below.

Press button .

Press + for 5 seconds. The display will read **rES**.

The **HACCP** LED will now light up permanently.

Press for 5 seconds.

The electronic control system will switch back to normal operating mode.

**Switching the appliance off:** Keep the **ON/OFF** button pressed for approx. 5 seconds - the display reads **OFF**.

### Setting the temperature

- Press button for 1 second. The temperature display flashes.
- To increase the temperature (warmer): press button .
- To reduce the temperature (colder): press button .
- Press button again.

The desired temperature setting is saved.

### Audible warning signal

The audible warning signal will sound in certain alarm events.

It can be cancelled by pressing button .

### Lid open alarm


When the lid is opened, the LED lights up and the temperature display begins to flash.

When the lid has been left open for more than 60 seconds, the LED begins to flash, and **door** and the temperature indication flash alternately in the display.

The audible warning signal sounds (unless the audible warning signal function has been deactivated).

If the lid has to stay open for longer in order to insert items to be cooled, cancel the audible warning signal by pressing button .

If the lid is open, the audible warning signal can also be muted by pressing the contact switch on the electronics housing. This will disable the audible warning signal for 60 seconds.

### Setting the delay time for the lid open alarm

The time before the audible warning signal sounds after the lid has been opened can be adjusted.

Press for 5 seconds. Display = **r61**

Display = **rA3**

Display = **rC1**

Display = **rC3**

Display = **d8d**

Display = **!** (minute)

Setting range = 1 - 5 minutes

Use buttons and to select the desired setting.

Display = **d8d**

Press for 5 seconds.

The electronic control system will switch back to normal operating mode.

### Deactivating the audible warning signal function

The audible warning signal function can be completely deactivated if necessary.

#### Note

In this case, the sentence stated in these operating instructions "The audible warning signal sounds" must be skipped when reading the section in question.

Press for 5 seconds. Display = **r61**

Display = **tc**

Display = **H4**

Display = **0**

Use buttons and to select the desired setting.

**0** = audible warning signal function activated

**!** = audible warning signal function deactivated

Display = **H4**

Press for 5 seconds.

The electronic control system will switch back to normal operating mode.

### Audible warning signal settings

The audible warning signal will be muted for the current alarm after the button has been pressed. Complete the following steps if you want the audible warning signal to reactivate automatically.

Press for 5 seconds. Display = **r61**

Display = **tc**

Display = **H4**

Display = **H0**

Display = **ASd**

Display = **ASn**

Display = **0**

Display = **!**

Display = **ASn**

Automatic reactivation of the audible warning signal is now active.

The time before the audible warning signal sounds again must be set.

Display = **ASd**

Display = **!** Time in minutes after which the audible warning signal will sound again after the button has been pressed. Setting range = 1 - 120 minutes.

Use buttons and to select the desired setting.

Display = **ASd**

Press for 5 seconds.

The electronic control system will switch back to normal operating mode.

## Alarm messages

### 1. LED flashes

If  appears in the display, the appliance has a fault. Consult your nearest customer service point.

### 2. LED flashes, the display reads *HI* or *LO*

The interior is too warm (*HI*) or too cold (*LO*).

The audible warning signal sounds (unless the audible warning signal function has been deactivated).

#### Note

The alarm parameters can be adjusted. See **Adjusting the alarm parameters**.

### 3. *HA/HF/HACCP* flashes

There has been a power cut (**HF**) of some length or the interior was too warm or too cold (**HA**) during a certain period of time.

Up to three alarm events can be stored and called up.

## Alarm test

This test checks the function of the internal and any external connected alarm device.

The appliance does not stop its refrigerating function during this test.

### Activating the test

Press  +  for 5 seconds.

- The display will change to a temperature value of 0.2°C below the set upper alarm limit.
- The temperature value will now rise by 0.1°C every 2 seconds.
- When the upper alarm limit is reached, *HI0* will appear in the display. An external alarm unit connected to the floating alarm output will now be activated.
- The temperature value will continue to rise up to 0.2°C above the upper alarm limit.
- The same process will take place automatically for the lower alarm limit. *LI0* will appear in the display.

The LED  will be lit during the test.

The electronic control system will automatically switch back to normal operating mode.

### Cancelling the test prematurely

Press  for 5 seconds.

#### Note

If the values of the upper and lower alarm limit (**AL** and **AH** in the section entitled "Adjusting the alarm parameters") are set to 0, *H--* and *L--* will appear in the display during this test.

## Adjusting the alarm parameters

The alarm limits (difference to the set temperature) and the alarm delay (delay until alarm goes on) can be adjusted.

Press  for 5 seconds. Display = *AL*

 Display = *AL3*

 Display = *AL1*

 Display = *AL3*

 Display = *d8d*

 Display = *AL* Lower alarm limit

 Display = temperature difference in °C

Use buttons  and  to select the desired setting.

### Important note

**Set positive values only.**

 Display = *AL*

 Display = *AH* Upper alarm limit

 Display = temperature difference in °C

Use buttons  and  to select the desired setting.

### Important note

**Set positive values only.**

 Display = *AH*

 Display = *Ad*

 Display = alarm delay in minutes


Use buttons  and  to select the desired setting.

 Display = *Ad*

Press  for 5 seconds.

The electronic control system will switch back to normal operating mode.

### Calling up stored alarm events and reading the temperature progression

Scroll through the list using or .

- HA n** Number of temperature alarms
- HA** Last temperature alarm
- HA !** Last temperature alarm but one
- HA 2** Temperature alarm before **HA !**
- HF n** Number of power cuts
- HF** Last power cut
- HF !** Last power cut but one
- HF 2** Power cut before **HF !**
- r t** Period in hours in which the maximum and minimum interior temperatures were measured
- r H** Maximum (highest) measured temperature
- r L** Lowest measured temperature


Select the required item using the button. Press this button again to return to the list.

Note: You can exit the menu at any time by pressing for 5 seconds.


If no button is pressed within 60 seconds, the electronic control system switches back automatically.

### Resetting the recorded temperature progression r t

Complete the following steps if you wish to reset the value saved for **r t** in the previous section to 0.

Press the button or until **r t** appears in the display.

Press for 5 seconds. Display = **r ES**.

The values for **r H** and **r L** (highest and lowest measured interior temperature) are then reset to the current interior temperature.

Press for 5 seconds.

The electronic control system will switch back to normal operating mode.

### Example of an alarm query

Situation: **HA/HF/HACCP** flashes in the display.

Display = **0** There has not been an alarm status with a too high or too low temperature. You must switch to display **HF n**.

Press this button until **HF n** appears in the display.

Display = **! 1** power failure has occurred.

Display = **HF** Last power failure.

Display = **Y 14** Year 2014

Display = **MO5** Month 05 (May)

Display = **d30** Day 30

Display = **h23** Hour 23

Display = **n 14** Minute 14

Display = **t03** The power failure lasted 3 hours.

Press + for 5 seconds. The display will read **r ES**.

The **HACCP** LED will now light up permanently.

**HA/HF** is cancelled in the display.

The electronic control system is now ready for the next alarm.

Press for 5 seconds.

The electronic control system will switch back to normal operating mode.

**Calibrating the control sensor** (standard sensor for temperature control)

Possible tolerances of the control sensor (the displayed temperature compared to the actual interior temperature) can be offset with this function.

Press for 5 seconds. Display = **rE1**

Display = **rA3**

Display = **rC1**

Display = correction value set at the factory

Use buttons or to increase or decrease the correction value in 0.1°C increments.

Display = actual (corrected) interior temperature

Display = **rC1**

Press for 5 seconds.

The electronic control system will switch back to normal operating mode.

**Product sensor** (available accessory)

The temperature may be measured or recorded at any point in the interior using the product sensor.

- Connect sensor (see section entitled **External alarm**)

**Activating the sensor**

Press for 5 seconds. Display = **rE1**

Display = **rA3**

Display = **0**

Display = **1**

Display = **rA3**

Press for 5 seconds.

The electronic control system will switch back to normal operating mode.

If **---** appears in the display, the product sensor has not been activated.

If **E2** appears in the display, the product sensor has not been connected, or is faulty.

**Calibrating the product sensor**

Possible tolerances of the product sensor (the set temperature compared to the actual interior temperature) can be offset with this function.

Press for 5 seconds. Display = **rE1**

Display = **rA3**

Display = **rC1**

Display = **rC3**

Display = **0.0**

Use buttons or to increase or decrease the correction value in 0.1°C increments.

Display = actual (corrected) product sensor temperature

Press for 5 seconds.

The electronic control system will switch back to normal operating mode.

**Switching the temperature display between control sensor and product sensor**

Press for 5 seconds. Display = **rE1**

Display = **1** (control sensor)

Display = **2** (product sensor)

If the product sensor is activated, AUX appears in the display.

Display = **rE1**

Press for 5 seconds.

The electronic control system will switch back to normal operating mode.

### Changing the network address

When connecting several appliances via the RS485 interface, each appliance must have its own network address.

Press for 5 seconds. Display = **261**

Display = **6c**

Display = **44**

Display = **40**

Display = **1**

Use buttons or to change the network address (**1-207**).

Display = **40**

Press for 5 seconds.

The electronic control system will switch back to normal operating mode.

### Resetting the parameters to factory settings

The **alarm limits** and **sensor calibration values** can be reset to the factory settings using this function.

Pull out the mains plug.

Keep pressed and connect the mains plug.

Display = **bn 1**

Display = **5td**

The electronic control system will automatically switch back to normal operating mode.

### Setting the real time clock

The real time clock is preset (CET). Other time zones or summer/winter time must be adjusted manually.

Press for 5 seconds. Display = **261**

Display = **6c**

Display = **4 14** Year 2014

Display = **10** Set the year by pressing the buttons.

= save new setting

Display = **7 07** Month (1-12)

Display = **7** Set the month by pressing the buttons.

= save new setting

Display = **1 01** Day (1-31)

Display = **1** Set the day by pressing the buttons.

= save new setting

Display = **1 03** Days of the week (1 = Monday, 7 = Sunday)

Display = **3** Set the day of the week by pressing the buttons.

= save new setting

Display = **12** Hour (0-23)

Display = **12** Set the hour by pressing the buttons.

= save new setting

Display = **12 48** Minute (0-59)

Display = **48** Set the minutes by pressing the buttons.

= save new setting

Press for 5 seconds.

The electronic control system will switch back to normal operating mode.

#### Note

When **6c** appears in the display, the real time clock must be reset.

## 7 GUARANTEE

This product is guaranteed for one year. The guarantee does not cover damage caused by incorrect use of causes beyond the control of J.P. SELECTA,S.A.

Any manipulation of the equipment by unauthorized personnel by J.P. SELECTA,S.A. cancels the guarantee automatically.

## 8 MAINTENANCE

Before removing the cover disconnect the equipment from the mains.

**The manipulation of the internal electronic circuits by unauthorized personnel can cause irreparable damage. Take it to one of the J.P.SELECTA, s.a. authorized technical services.**

CLEANING:

For the cleaning of the different parts use the following products:

Cleaning of stainless steel: Alcohol

Cleaning of plastic: Alcohol with cotton duster.