Armario refrigerado
Storage refrigerator

Stocklow LS
Stocklow GS

2101283
2101284

Indice

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Index

1	Preface	9
2	Packing list.....	9
3	Use of the equipment	9
4	Technical features	9
5	Operating principle.....	9
6	Control unit	10
7	Handling	10
8	Installation procedure.....	10
9	Connecting to the main power supply	10
10	Maintenance instructions	11
11	Troubleshooting	11
12	Technical service	11
13	Controller configuration	12
14	Technical parameters.....	14
15	Official approval and standards	14
16	Esquemas eléctricos / Electrical diagrams.....	15

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EN 61000-3-3:2002

1 Preface

This instruction manual provides all the necessary information regarding:

- Use of the refrigerator
- Technical specifications
- Installation and handling
- Operator procedures and instructions
- Maintenance operation

The manual is to be considered an integral part of the refrigerator and should be stored in a safe place for father consult to permit a good working life of the refrigerator.

The manufacturer cannot be held liable in the following cases:

- Improper installation (not in accordance with the guidelines indicated herein)
- Misuse of the refrigerator
- Power supply defects
- Improper or inadequate maintenance
- Manipulation or modification without authorization
- Use of non-original spare parts
- Partial or total failure to comply with the instructions

All electrical equipments can be hazardous to health. Current standards and legal requirements must be complied during the installation and use of any equipment.

2 Packing List

- Refrigerated cabinet
- Set of keys for the lock
- 3 trays with supports

3 Use of the equipment

These refrigerators are for preserving samples or products that require preservation below room temperature, with an in-built refrigerated unit.

The operating temperature for refrigeration is:

- between +1°C and +8°C at room temperature of +25°C and 60% HR.

4 Technical features

The refrigerator is a ventilated system, the evaporator is in a separate insulated box on the top.

All the materials used in the manufacture of this unit are guaranteed to be suitable for use with samples.

The refrigerating circuit are in compliance with the current standards.

5 Operating principle

The gas in the refrigerating circuit is in the first time compressed, liquefied and then evaporated in the ventilated evaporator, situated on the top of the equipment.

This cycle involves the absorption of heat from the air in the cooling compartment to cool it. The heat produced is then dissipated to the outside environment by a condenser unit located on the top of the refrigerator.

6 Control unit

The refrigerator is command from a "digital control unit" and a "main switch pilot light" in the top panel of the refrigerator.

The "main switch pilot light" is for turning on the power supply.

The green pilot light comes on to indicate that the unit is connected to the mains and is going to start working.

The green pilot light comes off to indicate that the unit is disconnected and it is not working. The "digital control unit" is for the regulation of all parameters to provide the correct working of the refrigerator. Please check all parameters in chapter 13 in this manual. Keep the manual always near the equipment.

This manual is part of the equipment and it is very important in case of a service call.

7 Handling


The refrigerator arrives in PET film and packed in a cardboard box on a wood pallet.

The refrigerator must be transported and handled with care to avoid posing a hazard to persons or properties.

Never place a refrigerator with an in-built refrigerated unit on its side or turn it upside down as this may damage or impair operation of the refrigerated unit.

We cannot hold liable for any damage or defects arising directly or indirectly from an improper handling of the equipment or non-compliance with the safeguards illustrated above.

8 Installation procedure

· Place the refrigerator in the coolest and best ventilated part of the room. Do not install the refrigerator near a heat source or direct sunlight.

· Remove the straps securing the cardboard packing

Remove the cardboard. Covering

Remove the PET protection film

Remove the metal transport belt from the compressor (star screwdriver) and the door support (without screws). Marked with a yellow label.

· Clean the refrigerator with mild detergent and then dry it with a soft cloth.

· Once located and before starting up, it is recommended that the equipment will be without starting around 24 hours so that the compressor oil and gas keep in their position.


9 Connecting to the main power supply

This operation must be carried out by professionally and qualified persons.

The refrigerators are complete supplied with a power supply cable for the connection to the main power supply.

A thermomagnetic circuit breaker (not supplied) must be installed between the mains power point and the power supply cable of the refrigerator.


9.1 Before proceeding make sure that:

· The mains voltage corresponds to the refrigerator indicated on the characteristics label; to ensure proper operation, it is essential for the power supply voltage to come within a range of +/- 10% of the unit's rated voltage.

· The electronic system to which the refrigerator is connected is made in compliance with current standard requirements.

· The electric connections and the installation of the thermomagnetic circuit breaker have been done by qualified person.

9.2 Connecting steps

- Install a thermomagnetic circuit breaker suitable to the power input of the unit being installed.
- Connect the refrigerator unit to the thermomagnetic circuit breaker outlet.
- Check that the refrigerator is in order as demonstrated by the pilot light incorporated in the main switch coming on.

10 Maintenance instructions

The smooth operation and life of the equipment are mainly determined by correct and regular maintenance.

10.1 Cleaning

Regular cleaning of the refrigerator unit is strongly recommended each month. Please follow the instructions below.

Disconnect the refrigerator power supply cable from the mains prior to carrying out any type of cleaning operation.

10.1.1 Cleaning the refrigerator surface

Clean the refrigerator with mild detergent and then dry it with a soft cloth. Do not use abrasive detergents!


10.1.2 Cleaning inside of the refrigerator

Clean the inside area minimum each month with a detergent suitable for use with samples.

10.1.3 Cleaning the condenser

For an efficient operation of the refrigerator, it is advisable to clean the condenser regularly approx. every 4 months with a dry brush or vacuum cleaner.

11 Troubleshooting

Refrigerator stops working (light off):

- Power supply failure
- Solution:
 - Check that the plug is inserted properly in the socket
 - Check that the switch on/off
 - Check that the mains voltage powers the plug


Refrigerator temperature goes up:

- Unit near to a heat source
- Condenser dirty or close
- Solution:
 - Move the counter or the heat source further away
 - Clean the condenser

12 Technical service


For technical service, please contact the dealer technical department or JP SELECTA's technical service and provide the serial number, and the date of buy.

13 Controller configuration

- 1. Power switch
- 2. Microcomputer controller
- 3. Compressor
- 4. Evaporator
- 5. Evaporator fan motor
- 6. Condenser fan motor
- 7. Condenser
- 8. Drain case
- 9. Shelf
- 10. Lamp

1. New upright air-cooling refrigerator should be first opened and ventilated before its use.
After that, users should use warm water to clean inside.

2. After connecting the power supply, press the "POWER" switch on the controller keyboard (Green Indicator Light ON), the fridge will come to work. The microcomputer controller, installed in the controller keyboard, could automatically adjust the temperature ranges. This intelligent digital controller works as: if the temperature increases and reaches set point plus differential the compressor is started and then turned off when the temperature reaches the set point value again.

3. **SET** To display target set point, select a parameter in programming mode or confirm an operation.
- ❄ (DEF) To start a manual defrost.
 - ▲ (UP) In programming mode, browse parameter codes or increase the displayed value.
 - ▼_{aux} (DOWN) In programming mode, browse the parameter codes or decrease the displayed value.
 - ▼+▲ To lock or unlock the keyboard.
 - SET+▼ To enter in programming mode.
 - SET+▲ To return to room temperature display.

LED	MODE	SIGNIFICATO
❄	On	Compressor enabled
❄	Flashing	Anti short cycle delay enabled (AC parameter)
❄	On	Defrost in progress
❄	Flashing	Dripping in progress
🌀	On	Fans output enabled
🌀	Flashing	Fans delay after defrost
°C	On	Measurement unit
°C	Flashing	Programming mode
°F	On	Measurement unit
°F	Flashing	Programming mode

4. How to see the set point
Press and immediately release the SET key, the set point will be showed;
Press and immediately release the SET key or wait about 5s to return to normal visualisation.
5. How to change the set point
Press the SET key for more than 2 seconds to change the set point value;
The set point value will be displayed and the "°C" or "°F" LED starts blinking;
To change the set value, press UP or DOWN arrows.
To memorise the new set point value, press SET key again or wait 10s.

6. How to start a manual defrost

Press the DEF key for more than 2 seconds and a manual defrost will start. It is recommended to do it once a year

7. How to change a parameter value

To change the parameter's value, operate as follows:

Enter the Programming mode by pressing SET + DOWN keys for 3 sec. ("°C" or "°F" LED starts blinking).

Select the required parameter. Press the SET key to display its value

Use UP or DOWN keys to change its value.

Press SET to store the new value and move to the following parameter.

To exit: Press SET + UP keys or wait for 5 sec without pressing a key.

NOTE: the set value is stored even when the procedure is exited by waiting the time-out to expire.

8. To lock the keyboard

Keep the UP + DOWN keys pressed for more than 3sec.

The "OF" message will be displayed and the keyboard will be locked.

If a key is pressed more than 3s, the "OF" message will be displayed.

9. To unlock the keyboard

Keep the UP + DOWN keys pressed together for more than 3s. till the "ON" message will be displayed.

10. Alarm signalling

Error	Cause	Results
P1	Room probe failure	Compressor output according to "Cy" and "Cn"
P2	Evaporator probe failure	Defrost end is timed
HA	Maximum temperature alarm	Results unchanged
LA	Minimum temperature alarm	Results unchanged
EA	External alarm	Results unchanged
CA	Serious external alarm	All Results OFF
dA	Door open	Compressor and fan restarts

**OUR CONTROL EQUIPMENTS HAVE BEEN
PRECISELY ADJUSTED
BEFORE LEAVING FACTORY.
SO TO AVOID DAMAGING COMPRESSOR UNIT
OR OTHER MALFUNCTIONS, USER MUST NOT MODIFY
THE CONTROLLER PARAMETERS PRIVATELY.**

14 Technical parameters

Description:	Single-door storage refrigerator
Model:	2101283 / 2101284
Prevention class of getting an electric shock:	I
Power source:	indicated on the characteristics label
Rated current:	2.3 A
Temperature range:	from +1 up to +8 °C
Refrigerant:	R134a
Dimensions:	730x845x2130 mm
Net weight:	120 Kg

15 Official approval and standards

Our products fulfil the present E.U. standards, including the CE mark of the European official approval.

89/336/EEC including amendments-electromagnetic compatibility (EMC)

73/23/EEC including amendments-low voltage (LVD)

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