Armario refrigerado
Storage refrigerator

Stocklow LSC

2101289

Indice

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Index

1	Preface	9
2	Packing list.....	9
3	Use of the equipment	9
4	Technical features	9
5	Operating principle.....	9
6	Control unit	10
7	Handling	10
8	Installation procedure.....	10
9	Connecting to the main power supply	10
10	Maintenance instructions	11
11	Troubleshooting	11
12	Technical service	11
13	Controller configuration	12
14	Technical parameters.....	14
15	Official approval and standards	14
16	Esquemas eléctricos / Electrical diagrams.....	15

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1 Preface

This instruction manual provides all the necessary information regarding:

- Use of the refrigerator
- Technical specifications
- Installation and handling
- Operator procedures and instructions
- Maintenance operation

The manual is to be considered an integral part of the refrigerator and should be stored in a safe place for father consult to permit a good working life of the refrigerator.

The manufacturer cannot be held liable in the following cases:

- Improper installation (not in accordance with the guidelines indicated herein)
- Misuse of the refrigerator
- Power supply defects
- Improper or inadequate maintenance
- Manipulation or modification without authorization
- Use of non-original spare parts
- Partial or total failure to comply with the instructions

All electrical equipments can be hazardous to health. Current standards and legal requirements must be complied during the installation and use of any equipment.

2 Packing List

- Refrigerated cabinet
- Set of keys for the lock
- 3 trays with supports

3 Use of the equipment

These refrigerators are for preserving samples or products that require preservation below room temperature, with an in-built refrigerated unit.

The operating temperature for refrigeration is between +2°C and +10°C at room temperature of +32°C and 60% HR. The operating temperature for freezing is between -6°C and -18°C at room temperature of +32°C and 60% HR.

4 Technical features

The refrigerator is a ventilated system, the evaporator is in a separate insulated box on the top.

All the materials used in the manufacture of this unit are guaranteed to be suitable for use with samples.

The refrigerating circuit are in compliance with the current standards.

5 Operating principle

The gas in the refrigerating circuit is in the first time compressed, liquefied and then evaporated in the ventilated evaporator, situated on the top of the equipment.

This cycle involves the absorption of heat from the air in the cooling compartment to cool it. The heat produced is then dissipated to the outside environment by a condenser unit located on the top of the refrigerator.

6 Control unit

The refrigerator is command from a "digital control unit" and a "main switch pilot light" in the top panel of the refrigerator.

The "main switch pilot light" is for turning on the power supply.

The green pilot light comes on to indicate that the unit is connected to the mains and is going to start working.

The green pilot light comes off to indicate that the unit is disconnected and it is not working. The "digital control unit" is for the regulation of all parameters to provide the correct working of the refrigerator. Please check all parameters in chapter 13 in this manual. Keep the manual always near the equipment.

This manual is part of the equipment and it is very important in case of a service call.

7 Handling


The refrigerator arrives in PET film and packed in a cardboard box on a wood pallet.

The refrigerator must be transported and handled with care to avoid posing a hazard to persons or properties.

Never place a refrigerator with an in-built refrigerated unit on its side or turn it upside down as this may damage or impair operation of the refrigerated unit.

We cannot hold liable for any damage or defects arising directly or indirectly from an improper handling of the equipment or non-compliance with the safeguards illustrated above.

8 Installation procedure

· Place the refrigerator in the coolest and best ventilated part of the room. Do not install the refrigerator near a heat source or direct sunlight.

· Remove the straps securing the cardboard packing

Remove the cardboard. Covering

Remove the PET protection film

Remove the metal transport belt from the compressor (star screwdriver) and the door support (without screws). Marked with a yellow label.

· Clean the refrigerator with mild detergent and then dry it with a soft cloth.

· Once located and before starting up, it is recommended that the equipment will be without starting around 24 hours so that the compressor oil and gas keep in their position.


9 Connecting to the main power supply

This operation must be carried out by professionally and qualified persons.

The refrigerators are complete supplied with a power supply cable for the connection to the main power supply.

A thermomagnetic circuit breaker (not supplied) must be installed between the mains power point and the power supply cable of the refrigerator.


9.1 Before proceeding make sure that:

· The mains voltage corresponds to the refrigerator indicated on the characteristics label; to ensure proper operation, it is essential for the power supply voltage to come within a range of +/- 10% of the unit's rated voltage.

· The electronic system to which the refrigerator is connected is made in compliance with current standard requirements.

· The electric connections and the installation of the thermomagnetic circuit breaker have been done by qualified person.

9.2 Connecting steps

- Install a thermomagnetic circuit breaker suitable to the power input of the unit being installed.
- Connect the refrigerator unit to the thermomagnetic circuit breaker outlet.
- Check that the refrigerator is in order as demonstrated by the pilot light incorporated in the main switch coming on.

10 Maintenance instructions

The smooth operation and life of the equipment are mainly determined by correct and regular maintenance.

10.1 Cleaning

Regular cleaning of the refrigerator unit is strongly recommended each month. Please follow the instructions below.

Disconnect the refrigerator power supply cable from the mains prior to carrying out any type of cleaning operation.

10.1.1 Cleaning the refrigerator surface

Clean the refrigerator with mild detergent and then dry it with a soft cloth. Do not use abrasive detergents!


10.1.2 Cleaning inside of the refrigerator

Clean the inside area minimum each month with a detergent suitable for use with samples.

10.1.3 Cleaning the condenser

For an efficient operation of the refrigerator, it is advisable to clean the condenser regularly approx. every 4 months with a dry brush or vacuum cleaner.

11 Troubleshooting

Refrigerator stops working (light off):

- Power supply failure
- Solution:
 - Check that the plug is inserted properly in the socket
 - Check that the switch on/off
 - Check that the mains voltage powers the plug

Refrigerator temperature goes up:

- Unit near to a heat source
- Condenser dirty or close
- Solution:
 - Move the counter or the heat source further away
 - Clean the condenser

12 Technical service

For technical service, please contact the dealer technical department or JP SELECTA's technical service and provide the serial number, and the date of buy.

13 Controller configuration

Display and functions:

During normal operation, the controller displays the value of the probe set using parameter/4(=1 ambient default probe, =2 second probe, =3 third probe). In addition, the display has LEDs that indicate the activation of the control functions (see Table 1), while the 3 keys can be used to activate/deactivate some of the functions (see table 2).

LEDs and associated functions

Table 1

icon	function	normal operation			start up
		ON	OFF	blink	
	Compressor	on	off	request	ON
	fan	on	off	request	ON
	defrost	on	off	request	ON
AUX	aux	output on	output off	-	ON
	alarm	all	no alarm	-	ON
	clock	RTC fitted and enabled, at least 1 time band set	RTC not fitted or disabled, not even 1 time band set	-	ON if RTC fitted

Table of functions activated by the buttons-models S, X, Y & C

Table 2

button		normal operation		start up	
		pressing the button alone	pressed together		
	up ON/OFF	more than 3s toggle ON/OFF	pressed together start/stop continuous cycle	-	
	down defrost	more than 3s: start/stop defrost		Pressed together start parameter reset procedure	for 1s display firmware vers.code
	set mute	-1s: display/set the set point -more than 3s: access parameter Setting menu(enter password 22) -mute audible alarm (buzzer)			for 1s RESET current EZY set

Setting the set point (desired temperature)

- Press SET for 1 sec the set value will start flashing after a few moments.
- Increase or decrease the value using UP or DOWN.
- Press SET to confirm the new value.

Switching the device ON/OFF

Press UP for more than 3 sec. The control and defrost algorithms are now disabled and the instrument displays the message "OFF" alternating with the temperature read by the set probe.

Manual defrost (models S, X, Y and C only)

Press for DOWN more than 3 sec (the defrost starts only the temperature conditions are valid).

Continuous cycle (models S, X, Y and C only)

Press UP and DOWN together for more than 3 sec.

Access and setting type F (frequent) and type C (configuration) parameters

1. Press SET for 3 sec (the display will show "PS")
2. To access the type F and C parameter menu, enter the password "22" using UP/DOWN;
To access the F parameter menu only, press SET (without entering the password);
3. Scroll inside the parameter menu using UP/DOWN;
4. To display/set the values of the parameter displayed, press SET, then UP/DOWN and finally SET to confirm the changed (returning to the parameter menu).
To save all the new values and exit the parameter menu, press SET for 3 sec.
To exit the menu without saving the changed values (exit by timeout) do not press any button for at least 60 sec

Table of alarms

Alam code	buzzer and alarm relay	LED	Description	Parameters involved
E0	active	ON	probe 1 error=control	-
E1	inactive	ON	probe 2 error=defrost	[d0=0/1]
E2	inactive	ON	probe 3 error=condenser	[A4=10]
IA	active	ON	external alarm	[A4=10][+A7]
dOr	active	ON	open door alarm	[A4=7/8][+A7]
LO	active	ON	low temperature alarm	[AL][Ad]
HI	active	ON	high temperature alarm	[Ah][Ad]
EE	inactive	ON	unit parameter error	-
EF	inactive	ON	operating parameter error	-
Ed	inactive	ON	defrost ended by timeout	[dP][dt][d4][A8]
dF	inactive	ON	defrost running	[d6=0]
cht	inactive	ON	condenser dirty pre-alarm	[A4=10]
CHT	active	ON	condenser dirty alarm	[A4=10]
EtC	inactive	ON	clock alarm	

Our products have been precissely modified before leaving factory. To avoid damaging compressor unit or other malfunctions, users must not modify the microcomputer parameters privately.

14 Technical parameters

Description:	Single-door storage refrigerator
Model:	2101289
Prevention class of getting an electric shock:	I
Power source:	indicated on the characteristics label
Rated current:	2.3 A
Temperature range:	from +1 up to +8 °C
Refrigerant:	R134a
Dimensions:	730x845x2130 mm
Net weight:	120 Kg

15 Official approval and standards

Our products fulfil the present E.U. standards, including the CE mark of the European official approval.

Latest version of following directives and standards

2006/42/EC	2014/30/EU
2012/19/EU	2011/65/EU
EN 60335-1	EN 60335-2-89
EN 62233	EN 55014-1
EN 61000-3-2	EN 61000-3-3
EN 55014-2	

Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.

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