PHARMALOW - L

2101300

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1. SAFETY INSTRUCTIONS

- 1) To prevent injury or damage to the unit, the appliance should be unpacked and set up by two people.
- 2) In the event that the appliance is damaged on delivery, contact the supplier immediately before connecting to the mains.
- 3) To guarantee safe operation, ensure that the appliance is set up and connected as described in these operating instructions.
- 4) Disconnect the appliance from the mains if any fail occurs. Pull out the plug, switch off or remove the fuse.
- 5) When disconnecting the appliance, pull on the plug, not on the cable.
- 6) **Any repairs and work on the appliance should only be carried out by the customer service department, as unauthorized work could prove highly dangerous for the user.** The same applies to changing the mains power cable.
- 7) Do not stand on the plinth, drawers or doors or use them to support anything else.
- 8) **If you have a lockable appliance, do not keep the key near the appliance or within reach of children.**
- 9) Do not store explosives or sprays using combustible propellants such as butane, propane, pentane etc. in the appliance. Electrical components might cause leaking gas to ignite. You may identify such sprays by the printed contents or a flame symbol.

2. SETTING UP

- 1) **Avoid positioning the appliance in direct sunlight or near cookers, radiators and similar sources of heat.**
- 2) The floor on which the appliance stands should be horizontal and level. Compensate for uneven floors with the adjustable feet. Leave a minimum space of 10cm between the back of the equipment and the wall.
- 3) Always ensure that there is a good ventilation and that the outward flowing air is able to escape.

3. FEATURES

- Gross capacity (see plate): 250L
- External dimensions: 168 x 60 x 65cm (H x W x D)
- Maximum load per rack: 45kg

4. CONDITIONS OF USE

- Altitude up to 2000 m (795 mbar)
- Maximum relative humidity 65 %
- Mains voltage fluctuations up to $\pm 10\%$ of rated voltage.
- Equipment for indoor use.

5. POWER CONNECTION

Power supply (a.c.) and voltage at the operating point must comply with the details on the type plate, which is attached at the rear part of the refrigerator. The socket must be fused with a 10 A fuse or higher.

To ensure trouble-free operation, do not connect the appliance to the power supply until 6 to 8 hours after delivery. Connect the appliance through a properly earthed plug and socket.

The socket should be easily accessible and not directly behind the appliance.

The wires in the mains lead are coloured in accordance with the following code: green/yellow = earth, blue = neutral, brown = live.

Warning! This appliance must be earthed.

6. SWITCHING THE APPLIANCE ON AND OFF

You are recommended to clean the appliance before switching it on for the first time (see "Cleaning"). When it is connected up (plugged into an approved wall socket) the appliance is ready for operation.

To switch the appliance off: pull out the mains plug.

7. OPERATING

Our products leave the manufacturing process after through testing procedures to ensure they are in proper working order.



To operate them, simply connect them to a GROUNDED electrical outlet to the proper voltage of the equipment. The unit should operate normally.

We strongly recommend against manipulation of the components of this unit by non-authorized personnel.

If you observe any malfunction, contact technical service immediately.

Operating controls

- (1) Lock
- (2) Temperature display
- (3) Light door contact switch
- (4) Temperature controller
- (5) Fan door contact switch
- (6) On/off light switch

Switching the appliance on and off

Switching on

Plug in the mains plug - the appliance is switched on.

Switching off

Remove mains plug.

Switching the interior lighting on and off


The interior light always switches on when the door of the appliance is opened.

Setting the temperature

Set the temperature in the appliance using the control knob.

The lowest temperature in the interior is reached in the Max position.

Adjust the controller between Min and Max as required.

8. COOLING MATERIALS

This appliance can be used for chilling all types of materials. The grid shelves can be moved to accommodate different height packages or objects.

To remove the shelves the door must be opened wider than 90°. Substances which give off or absorb odours and flavours as well as liquids should always be stored covered or in closed containers. High-percentage alcohols should be tightly sealed and stored upright.

The ventilation slots on appliances with recirculated air fans on the inside must be not covered.

9. SHUTTING YOUR APPLIANCE DOWN

If your appliance is to be shut down for some time, switch it off and disconnect the plug. Clean the appliance and leave the door open in order to prevent unpleasant smells.

10. DEFROSTING

The refrigerator compartment defrosts automatically. The water that forms on the rearwall drains into a reservoir at the back of the appliance and evaporates automatically through the compressor heat. Ensure that the water can flow without hindrance through the drain hole (see "Cleaning").

A layer of frost or ice can form on the cooling elements when the appliance has been operating at the coldest setting for a long time.

If a **layer of ice** has formed, the cooling elements should be defrosted and the inner surfaces cleaned.

Do not use any mechanical devices or other artificial aids for defrosting other than those recommended by the manufacturer.

11. CLEANING

Periodic cleaning of the device

To clean the device, follow the instructions below:

- Disconnect the equipment from the mains and remove all products stored in it.
- Open all doors and allow time for the interior to reach room temperature. Remove all interior accessories and clean them with soap and lukewarm water. Dry all accessories completely with a soft cloth.
- Once the chamber has reached room temperature, clean all interior and exterior surfaces with soapy water.

Rinse thoroughly and dry with a soft cloth. Failure to dry the device properly may result in water stains.

Stainless steel cleaners are also available in the market, which can repair and protect the protective layer of steel surfaces.

- Put the accessories back in place and connect the unit to the mains.
- Pitting corrosion or cracks in steel are signs of material deterioration. In this case, apply stainless steel cleaners capable of repairing the steel passivation.

WARNING: Never use steel scrubbers, wire brushed or spatulas to clean the device.

NOTE: Cleaning products used must be alkaline-based or chlorine-free. Any cleaner containing chlorides will damage the protective layer of the stainless steel.

Cleaning of the condenser


The condenser, located behind the back grille of the device, should be checked regularly. Cleaning frequency will depend on the working environment. Air must flow freely through the condenser, so its surface must be free of dirt and grease. Dirty condensers cause compressor fault and product loss. If the condenser battery is dirty or blocked, follow the steps below:

- Disconnect the device from the mains.
- Remove the back grille from the unit.
- On some models, it will be necessary to remove the screws that affix the condensing unit to the skirting board and remove it in order to clean the condenser.
- If the condenser has a protective case, it must be unscrewed and removed.
- Once the surface of the condenser is accessible, it should be cleaned using a vacuum cleaner or a soft brush. Do not use metal brushes.
- If the dirt is excessive, compressed air can be used for cleaning.
- After cleaning, reinsert the protective case, return the condensing unit to its original position, and re-attach all screws.
- Finally, replace the rear grille and connect the device to the mains.

WARNING!: Do not use water to clean the condenser, as it may damage nearby electrical components.

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Over time and due to normal wear, door hinges may move off place slightly. If you notice that the door is not properly aligned, tightened the screws that attach the hinge brackets to the cabinet.

Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.