PHARMALOW

Cod.2101320 (Pharmalow M Printer)

Cod.2101331 (Pharmalow B)

Cod.2101332 (Pharmalow C)

Cod.2101334 (Pharmalow D)

INDICE

ESPAÑOL

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INDEX

ENGLISH

1. GENERAL INFORMATION ABOUT THE REFRIGERATOR	28
1.1. Safety Instructions	29
1.2. Precautions in use	31
2. EQUIPMENT INSTALLATION	31
2.1. Environmental conditions	31
2.2. Location	32
2.3. Levelling	32
2.4. Preparation before use	32
2.5. First Power-On	32
2.6. Operación after power failure	33
3. EQUIPMENT SPECIFICATIONS	34
3.1. Packing List	34
3.2. Equipment Layout	35
3.3. Nameplate	35
4. EQUIPMENT COMMISSIONING	36
5. OPERATION	36
5.1. Control Panel	36
5.2. Temperature visualization	37
5.3. Parameters:	39
6. MAINTENANCE & CLEANING	41
6.1. Cleaning	41
6.2. Maintenance for doors & hinges	42
6.3. Printer	42
6.4. Battery	42
6.5. Technical Service	43
6.6. Spare Parts	44
6.7. Refrigerating System Diagram	45
6.8. Electric Diagram	45
7. TROUBLESHOOTING	46
8. WARRANTY	48

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ENGLISH

RECYCLING INFO ABOUT THIS PRODUCT:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.

1. GENERAL INFORMATION ABOUT THE REFRIGERATOR

This manual has been made in a simple way so that by reading it you can know the operation and maintenance of our refrigerator. It is recommended that you read it carefully and keep it for any reference.

The manufacturer is not responsible for any damage to people or objects that may be caused by failure to comply with the prescriptions contained in this manual. To know all the advantages of this device, please read carefully before proceeding to install it. Anyone who uses this device is recommended to read this user manual. If in doubt, consult your distributor.

This product has been manufactured under strict quality controls and meets all the requirements established by JP Selecta. Before leaving the factory, each unit has been tested to guarantee its quality. This equipment has been manufactured with recyclable materials, through an environmentally friendly production process.

This furniture complies with directive 2014/30/EC, 2014/35/EU and 2011/65/EU. Furthermore, the harmonized European standards defined in the equipment's Declaration of Conformity have been applied.

WARNING! This device must be used only for the purpose described in this manual.

Electrical/electronic equipment at the end of its useful life must be managed by an authorized waste manager.

IMPORTANT: see the **TECHNICAL LABEL** located on the upper side of the equipment that indicates the **TYPE OF REFRIGERANT**.

FOR MODELS WITH R600a / HYDROCARBONS REFRIGERANT:

CAUTION – RISK OF FIRE OR EXPLOSION. FLAMMABLE REFRIGERANT. SHOULD ONLY BE REPAIRED BY A QUALIFIED TECHNICIAN. DO NOT DRILL THE REFRIGERANT PIPE.

1.1. Safety Instructions

When using this product for the first time, pay attention to the meaning of the following warning signs and read the safety precautions carefully so that you can use this equipment safely and correctly.

	Caution! Failure to follow this warning may result in injury or death.		Caution! Failure to follow this notice may result in personal injury or property damage.
	Protective conductor terminal		Warning! Explosion risk
	Warning! Flammable material or fire hazard		Warning! Crushing Hazard
	Warning! Risk of electric shock		Avoid sun exposure.
	Review instruction manual		Handle with care

Caution: Failure to follow precautions may result in serious personal injury or death.

Do not touch the equipment with wet hands to avoid electric shock; The equipment should not be grounded through gas pipes, water supply pipes, telephone lines or lightning rods, which can easily cause electric shock.

This equipment can only be installed by professional technicians or after-sales maintenance personnel, as it may cause electric shock or fire.

Be sure to install the equipment on solid, flat ground and take due care to prevent it from tipping over. If the ground is not solid enough or the installation location is not suitable, the equipment may fall and cause damage to the equipment or personal injury.

Handle the power cord with care to avoid short circuits or open circuits. Turn off the power before pulling out the plug. Hold the plug carefully and remove it. Do not pull on the cords from the plug. Otherwise, it may cause electric shock or fire due to short circuit. Do not coil the power cord, do not press it under furniture or heavy objects, and do not approach heat sources such as compressors.

Firmly insert the power plug into the outlet to ensure firmness and reliability, so as to avoid electrical leakage; After installation, the plug must be within reach so that it can be unplugged in time in an emergency.

Separate special sockets must be used and reliably grounded. The cross-sectional area of the copper conductor in the wall connected to the socket must have sufficient section. Do not extend the power cord without authorization to avoid overheating or fire.

Do not use a power supply that is not specified in the design of the equipment, to avoid overheating, short circuits and other failures. For example, connecting products with a nominal voltage of 110V to a 220V power source may cause failures such as overheating and short circuits in the equipment. For detailed information on AC input voltage, see

the nameplate specifications (nominal voltage $\pm 10\%$). If the operating voltage is too low or too high, a

suitable automatic voltage stabilizer must be installed.

Place the equipment stably and avoid shocks and sudden movements.

Do not place the equipment in a dangerous area and do not operate the equipment near flammable items to avoid explosion or fire accident.

Do not place the equipment in areas exposed to sun or rain, to avoid dangers such as short circuits or overheating.

Do not tilt or place the equipment on its side and do not hit the body of the equipment; Cooling systems are installed on the equipment, which is easy to damage by tilting or impact.

Place the equipment in a dry and dust-free environment to avoid risks such as overheating and short circuits.

In case of sound, smell, smoke, etc. unexpected when the power is turned on, please unplug it in time and contact the manufacturer or supplier.

Place the equipment in a dry and ventilated environment, and ensure that the ventilation grilles of the equipment and instrument surfaces are not blocked or protected by walls or other objects; Do not use it in a poorly ventilated environment to avoid damage caused by the heat released by the equipment.

It is prohibited to disassemble and modify this equipment without authorization, to avoid possible safety risks. In case of non-compliance, the equipment loses its warranty and the user is responsible for any damage that may arise.

It is prohibited to place flammable and explosive dangerous goods, strong corrosive acids, alkalis and other items unsuitable for the equipment in the equipment.

When storing toxic, harmful materials, etc. use the equipment in safe areas. Improper use may cause harm to human health or the environment.

Metal objects, such as nails or iron wires, should not be inserted into any opening, space or outlet on the equipment; Otherwise, electric shock or injury may occur due to accidental contact between the above objects and moving parts.

To ensure the normal operation, ventilation and heat dissipation of the equipment, the back, left and right sides of the cabinet should be at least 30cm from the wall, and the air inlet and outlet should not be blocked by obstacles. In case of non-compliance, the equipment loses its warranty.

This equipment must be connected to a ground wire.

Note: Failure to follow precautions may result in personal injury or equipment failure and related property losses.

When the equipment is working, do not touch the inner surface of the cabinet without wearing protective equipment.

Hold the handle and close the door to avoid pinching your fingers; When the equipment is not used for a long period, unplug it.

When restarting the computer after a power outage or shutdown, check the computer settings first; Otherwise, the stored items may be damaged due to the setting change.

Store the keys correctly to avoid accidents due to involuntary openings or unauthorized personnel.

When handling the equipment, be careful not to tip it over to avoid damage to the equipment or personal injury.

When handling it, it will be lifted from below, with an inclined plane no greater than 45°, and it will be handled with care. Use the equipment in safe areas. Improper use may cause harm to human health or the environment.

1.2. Precautions in use

- Before placing the items into the equipment, first confirm that the temperature in the freezer storage chamber has reached the set value, and then place the items in batches. Each time you place items, it should not exceed 1/3 of the internal volume of the freezer storage chamber to avoid excessive temperature rise.
- The temperature display value of the equipment is the temperature at the temperature sensor in the freezer storage chamber. There is a certain gap between the displayed temperature and the actual temperature in the center of the equipment when the equipment just starts working, but as the equipment enters a stable state, the displayed temperature will gradually approach the actual temperature.
- Use a diluted neutral cleaner to clean the equipment, and do not use brushes, acid, gasoline, soap powder, polishing agent or hot water to clean the equipment, otherwise the painted surface and rubber plastic parts may be damaged. Be careful not to clean plastic rubber parts with volatile solvents such as gasoline.
- When the equipment is not used for a long period, the power supply should be turned off.

2. EQUIPMENT INSTALLATION

2.1. Environmental conditions

- Ambient temperature: 16~32°C, the ideal temperature range is 18~25°C. The use of air conditioning equipment will be required to ensure the indicated conditions. If this point is met, the performance and operation of the equipment may be affected.
- Relative humidity: ≤80%.
- The base of the equipment will not transmit strong vibrations nor will there be an atmosphere of corrosive gases around it.
- Avoid the existence of atmospheres with a large amount of dust.
- Avoid rocking or shaking the equipment.
- Elevation of the equipment's workplace: less than 2000m.
- Indoor use, pollution degree 2 and overvoltage category II.
- For detailed information on AC input voltage, check the equipment nameplate (nominal voltage ±10%).

- There should be no direct sunlight, other sources of cold and/or heat or strong electromagnetic interference, which could affect the normal operation of the control system and could damage the equipment.

2.2. Location

To operate the equipment normally and obtain the best level of performance, the location of the equipment must meet the following requirements:

- It cannot be installed in a narrow and closed space, consider the size of the equipment in relation to the access and required maintenance area.
- The installation floor must be solid, flat, non-combustible and capable of supporting the weight of the equipment during operation;
- It will have good ventilation and direct sunlight will be avoided;
- Each device needs to use a power outlet independently. Make sure the plug and jack are firmly connected;
- Check the working voltage before use. In areas with unstable voltage, consider using a voltage stabilizer suitable for the equipment load to ensure that the input voltage requirements in the installation environment are met;

Attention: Since the ambient temperature has a great influence on the equipment, if the environmental requirements cannot be met, the equipment may not operate normally. Improve the environment before using the equipment.

2.3. Levelling

It is very important that the appliance is perfectly level for correct operation, the doors are aligned and the unit is not subjected to undue stress.

These models are supplied from the factory with height-adjustable legs (maximum height 180 mm). In this case, you must ensure that the ground where the unit is located is level, so that it can work stably.

As an option, wheels can be supplied for all models. In case you want to install wheels, be sure to lock the front brakes so that the unit is completely stable and level. In the "Wheel Installation" section you will find detailed information about wheel adjustment.

2.4. Preparation before use

- Previously remove the outer packaging of all products (including the protective foam from the packaging box).
- Check all the accessories included in the packing list (see section 3)
- Clean the product before using it (see section 9)

2.5. First Power-On

When you use your computer for the first time, follow these steps:

- After the equipment is placed, leveled and cleaned, it should be left idle for more than 24 hours and then turned on to ensure normal operation of the equipment.
- Without load, connect the power cord to a special socket with appropriate specifications
- After turning it on, turn on the power switch of the computer. Move the battery on/off switch to the ON position

located on the back of the unit.

- Check whether the operating temperature of the equipment reaches the required value, observe the normal starting and stopping of the equipment for more than 24 hours, and place a small amount of items in the refrigerator after confirming normal operation.
- Store items in batches, with items not exceeding 1/3 of the cabinet volume at a time. Please make sure the equipment is working properly after being turned off for more than 12 hours before placing the next batch of items.)
- Try not to open the door during cooling, otherwise the temperature will increase.

2.6. Operación after power failure

The equipment has a memory function for the set value. When power is restored after a power outage, the computer will continue to operate as it did before the power outage occurred.

- Once the equipment is turned off, 5 minutes are necessary before it can be turned on again, to avoid damaging the compressor.
- We guarantee normal operation of this equipment under certain conditions, but we are not responsible for any loss or damage to stored items after a power outage.
- In case of failure or shutdown, the temperature in the Refrigerator will increase. If it cannot be repaired in a short time, please take out the stored items and transfer them to a place that meets the storage temperature requirements to avoid damage to the items.
- Before placing items in the equipment, it should be confirmed in advance whether the temperature range of the equipment meets the temperature requirements of the items, to avoid damage to the stored items due to the difference between the settable temperature of the equipment and the required temperature of the articles. Pay attention not to block the air outlet and inlet when placing items on the equipment.
- Due to the inertia of cooling, there is a certain difference between the actual display temperature and the set temperature of the equipment, which is a normal phenomenon.
- Do not change the set temperature frequently in a short time; Otherwise, the expected setting effect may not be achieved due to large temperature inertia; Please ensure that there is a certain air circulation space around the cabinet when placing items, especially do not block the temperature sensor in the cabinet, otherwise it will affect the stability and precise temperature control
- When placing items, if the moisture content of the items is too high or too low, it may affect the degree of humidity in the cabinet, so it is best to keep the items sealed; Humidity will also be affected if the door is opened too frequently and not closed properly.

3. EQUIPMENT SPECIFICATIONS

PHARMALOW refrigeration cabinets are designed for adjustable temperatures between +2°C and +8°C. They are equipped with a hermetic compressor group mounted on an anti-vibration seat with a forced-type finned evaporator and a ventilated-type finned condenser, an air circulation turbine with a stop when the door is opened, temperature sensors, interior LED light, glass door. double, insulating and fog-proof with lock and printer.

White lacquered outdoor furniture. Reversible double glass door with easy-to-replace lock and gasket, with automatic return mechanism and manual interlock.

Interior air circulation for greater temperature homogeneity, Safety thermostat. Microprocessor temperature control that allows you to configure and view the following parameters through the display:

Security password

Set temperature

Set alarms

Calibration

Delay in starting compressor


Code	Capac. liters	Dimensions inside cm			Dimensions external cm			Load per Shelf (kg)	Refrigerant Gas	Nr.Shelves	Rated Power (W)	Weight (Kg)
		Height	Depth	Wide	Height	Depth	Wide					
2101331	130	58	51	55	81	63	65	23	R600a	2+1	150	51
2101332	315	112	53	58	176	67	65	20	R600a	4+1	245	87
2101320	395	135	53	58	199	67	65	20	R600a	6+1	252	95
2101334	525	126	68	61	196	81	72	30	R290	6	380	148

3.1. Packing List

Standard equipment consists of the following components:

Code	Shelves Nr.	Trays	Keys	User Manual	EU Declaration of Conformity
2101331	2	1	1	1	1
2101332	4	1	1	1	1
2101320	6	1	1	1	1
2101334	6	0	1	1	1

3.2 Equipment Layout

3.3 Nameplate

The nameplate is a label permanently affixed to the side of the equipment, which contains important electrical information, as well as data related to the refrigeration system of each unit. In addition, it includes the model and serial number.


4. EQUIPMENT COMMISSIONING

- All JP Selecta products are tested at the factory, evaluating their quality and performance, and verifying they do not present defects.
- When you receive your appliance, it should be carefully inspected for any possible damage that may have occurred during transportation.
- If any damage is detected to the unit, you must retain all packaging material and report such damage on the carrier's delivery note. A claim must be made immediately to the transport company.
- If damage is noticed during or immediately after installation, contact your dealer immediately.

5. OPERATION

5.1 Control Panel

The control system of the refrigeration equipment is defined below.

Screen symbols description:

Door Switch	WIFI	Defrosting	Key Lock
Power Failure	Print	Door Heating	Refrigeration
Mute	Battery Level	Serial Port	Fan

a) Door switch indicator:

When the cabinet door is opened, the indicator lights up. When the door is closed, the indicator is off.

b) WIFI:

The WIFI lamp turns on if the network is activated and is off if it is not.

c) Defrosting:

When the refrigerator enters defrost mode, it will always be on. When the refrigerator comes out of the defrosting state, the indicator is off.

d) Key Lock:

When the keys are in lock mode, none of the keys respond and the key lock lights up. In this case, after pressing the combination “up key”+“down key” for 3 seconds, you will be prompted for a password (default 005). After entering correctly, press the “Set/Mute Multiplex” key and the key lock is deactivated. At this time, the lock indicator stops lighting. In the unlock state, if no key is pressed for 60 seconds, the key lock starts again and the key lock indicator lights up again. If you want to activate the locking mode, press the combination “up key” + “down key” for 3 seconds.

e) Power Failure:

The refrigerator normally operates at 220V/110V. When the input does not receive power, the buzzer will activate, “PF” (power failure) appears on the display intermittently every 3 seconds and the power failure indicator remains illuminated. Once the power is restored, the power failure indicator turns off.

f) Print Indicator:

When the printer is not working, the indicator is off, when the printer is active, the indicator lights up.

g) Door Heating indicator:

When the door heating is turned off, the indicator is off and the opposite when it is active.

h) Refrigeration indicator:

If the compressor is running, the indicator will be on. The other way around if it is stopped.

i) Alarm mute indicator:

When the alarm notification is silenced, the indicator will be on.

j) Low Level Battery Indicator:

When the battery level is less than 8V, the buzzer is activated, the battery indicator will light up and the low battery message “BL” appears on the screen, intermittently every 3 seconds. When the battery voltage recovers to 12V, the buzzer will turn off, as well as the indicator light and the message on the screen will disappear.

k) Serial Port:

When the RS-485 serial port is not connected to any equipment, the indicator is off. When the refrigerator connects to external equipment through the port and there is communication, the serial port indicator will light.

l) Fan indicator:

When the evaporator fan is running, the indicator will be active

5.2 Temperature visualization

Located in the control panel area, it allows you to view the temperature value:


To visualize environmental temperature:

In the key lock state, press  and the display shows the room temperature, returning to its normal display mode after 5 seconds if no other key is pressed, or by pressing the  and  keys simultaneously. In the key unlock state, press  and the display shows the room temperature, returning to its normal display mode after 5 seconds.

To check the humidity, with the keyboard unlocked, press and hold the  and key  , the display will show the humidity, then, if no key is pressed afterwards, the original display will appear after 5 seconds.

a) **Set / mute key** is used to also enable / disable the sound of the alarms. When the device is not in alarm mode, press  and the temperature will be shown for 5 seconds, then returning to the normal display mode. In keyboard unlocked mode, press  for 3 seconds and you will enter the user menu. When the timer is triggered (including cabinet high temperature alarm, door open alarm, sensor failure alarm, etc.) and in the key unlock state, press  the first time, the buzzer will stop ringing and the room temperature will be displayed for 5 seconds, after which normal display will resume (press the mute button to turn off the buzzer). to alarm this abnormal state, for example troubleshooting, and the buzzer will be activated next time for any abnormality). By pressing again, the alarm is activated, the room temperature is displayed for 5 seconds and the chamber temperature and alarm status are displayed.

In the keyboard unlocked state, the key  is used as the "SET" confirmation key. In the keyboard unlocked state and parameter programming mode, press this key to display the parameter values and their names. If pressed for more than 3 seconds, the changes are saved and you return to the starting point.

b) **Up arrow key**  In parameter programming mode, moves to the next parameter or increments the parameter value. For example, if the temperature value is to be programmed, increase the set value. To assign the programming value, press and hold the key and the parameter value will increase rapidly. On the initial screen, if you press the key for 3 seconds, you can import data using a USB pendrive. (see section f)

c) **Down key**  In parameter programming mode, moves to the previous parameter or decreases the parameter value. For example, if the temperature is being changed, its value will be reduced. If the key is held down, the value decreases quickly.

d) **Print key**,  The system can hold up to 7 days of data for printing. Press the key to print the temperature for the desired period of time.

e) **Lighting key**  When the device is turned on, the lights turn off by default. Press this key to turn the lights on and off. When the key is active, the lighting is activated when the door is opened and goes off when the door is closed. If the key is disabled, the refrigerator interior light will not turn on when you open or close the door. This equipment has two functions. Automatic and manual lighting mode. In automatic lighting mode, the lighting is activated when the door is opened and will turn off after 5 seconds. The user can also activate manual mode. In this case, the lighting is always activated (the indicator above the key is illuminated). In this case, the lighting is always active regardless of whether the door is open or closed. Only when the manual lighting is deactivated using the key, the automatic lighting is activated (and the key light goes out).

f) **USB data export.**

Automatic mode: When the pendrive is connected to the USB socket, the buzzer will emit a sound and the message "ON" will appear on the display. The data is stored on the pendrive in PDF format. Once the data has been entered, the buzzer emits another sound and "END" appears on the display. After 6 seconds, the display returns to its normal state.

Note: If there is no data to upload, the display will not show the "ON" or "END" messages.

Manual mode: In keyboard unlock mode, and with the pendrive connected, press the up key for 3 seconds. The message "d01" will appear on the display. Using the up and down keys, adjust "d00..d12" and press the key to generate the file "d00" or generate the PDF file for the months 01..12

5.3 Parameters:

After turning on the equipment, it enters operation mode;

The parameter programming process is described below:

To unlock: On the home screen and using the “up arrow” and “down arrow” keys, press for 3 seconds. Next, the code “000” will appear on the display, enter the password “005” to unlock (With the code “099”, the user password “005” can be reset if it has been accidentally changed). After unlocking, press the key  for 3 seconds, you will enter the user menu. If the display shows “PS1” it means that it has been done correctly.

Use the “up arrow” and “down arrow” keys to move through the list of parameters, the display follows this order:

PS1→b1→b2→Set→H→L→n→y→r→S→F→Pt→Th1→P1→P2→MAX→MIN→CLR, press the key to confirm the parameter.

Use the “up arrow” and “down arrow” keys to move through the parameters

Press the key to view the parameter

Use the “up arrow” and “down arrow” keys to modify the parameters

Press the key to confirm and store the parameter

If the key is held down for more than 3 seconds, the parameter is stored. Likewise, if no key is pressed, parameter programming is exited and the initial menu is displayed.

Below is the list of programmable parameters.

Nr	Menu Item	Parameter Range	Suggested Value	Comments
1	Set	0,0...10,0	5,0	Temp. Setting
2	H	0,0...10,0	5,0	High temp alarm+H; With H=0, the alarm is deactivated, if the temperature is above the alarm temperature, H1 will appear on the display
3	L	0,0...10,0	4,0	High temp alarm-L; With L=0, the alarm is deactivated, if the temperature is below the alarm temperature, L1 will appear on the display
4	Pt	0..240 min	20	Print interval
5	Th1	20,0..50,0°C	40	Environmental limit temperature alarm
6	P1	1. Automatic heating mode mode 1 2. Automatic heating mode		

Nr	Menu Item	Parameter Range	Suggested Value	Comments
		mode 2 3. Automatic heating mode mode 31 4. Always active 5. Always off	1	-
7	P2	1. Temp. average 2. Higher temperature 3. Lower temperature	1	-
8	PS1	-	005	User Password
9	MAX	-	-	Max. Temp. since last reset
10	MIN	-	-	Min. Temp. since last reset
11	CLR	-	0	Reset of max & min temp
12	B1	V1.0-9.9		Hardware version
13	B2	V1.0-9.9		Software version
14	n	10..50		Register year set
15	y	1..12		Register month set
16	r	1..31		Register day set
17	S	00..23		Register hour set
18	F	00..59		Register minute set
19	SCY	0..240	10	Temperatura data storage (minutes)

Alarm codes:

Code	Alarm Description
H1	High temp alarm
L1	Low temp alarm
H2	Ambient high temp.
H3	Condenser overheating
Do	Door opened alarm
PF	Power Fail
BL	Low Battery
ER	Recorder not connected
LOF	Recorder not startup
EE	Communication Error
E1	Control Sensor Error
E2	Upper Sensor Error

E3	Lower Sensor Error
E4	Ambient Sensor Error
E5	Evaporator Sensor Error
E6	Condenser Sensor Error
E7	Humidity Sensor Error

6. MAINTENANCE & CLEANING

6.1 Cleaning

To clean the device, follow the following instructions:

Disconnect the equipment from the mains and remove all products inside.

Open all doors and drawers, and let the interior reach room temperature. Remove all interior accessories and clean with mild soap and warm water. Dry all accessories completely with a soft cloth.

Once the chamber has reached room temperature, clean all interior and exterior surfaces with soapy water. Rinse thoroughly and dry with a soft cloth. Not drying correctly can cause water stains to appear. Likewise, there are stainless steel cleaners that can repair and protect the protective layer of steel surfaces.

Place the accessories in their original position and connect the unit to the network.

Pitting or cracks in steel are signs of deterioration of the material. In this case, apply stainless steel cleaner capable of repairing the passivity of the steel.

NOTE: Never use steel wool, wire brushes or spatulas to clean the appliance or scouring pads.

NOTE: The cleaning products you use must be alkaline-based or chlorine-free. Any cleaner containing chlorides will damage the protective layer of the stainless steel.

Door gaskets maintenance

Door gaskets require regular cleaning to protect its elasticity, ensuring proper closure and preventing mold growth. Cleaning the weather stripping can be done with soapy water. Avoid using harsh cleaners and sharp utensils.

The door gaskets can be easily removed by pulling it out of the door profile, and replaced by pressing it against it.

NOTE: if the weather stripping becomes contaminated with grease and is not cleaned, it will lose its elasticity and break in a short time.

Condenser Cleaning

The condenser, located behind the rear or front grille of the unit, depending on the model, must be checked periodically; it is recommended to do so at least every 6 months. The frequency of cleaning will depend on

the existing work environment. It must be ensured that air circulates freely through the condenser, so its surface must be free of dirt and grease.

The dirty condensers cause malfunctions and product loss. If the condenser coil is dirty or blocked, follow the steps detailed below.

- Disconnect the appliance from the mains.
- Remove the rear or front protection from the unit, depending on the model.
- In some models it will be necessary to remove the screws that adjust the condensing unit to the baseboard, and extract it in order to clean the condenser.
- If the condenser has a protective casing, you must unscrew it and remove it.
- Once the surface of the condenser is free, it must be cleaned using a vacuum cleaner or soft brush. Never use a metal brush.
- If the dirt is excessive, you can use compressed air to clean.
- Once clean, replace the protective casing, return the condensing unit to its original position and replace all screws.
- Finally, replace the protection and connect the appliance to the mains.

WARNING! : Never use water to clean the condenser as it could damage nearby electrical components.

6.2 Maintenance for doors & hinges

Along the time, and the doors are used, the hinges may move slightly. If you notice that the door is starting to misalign, you should tighten the screws that attach the hinge brackets to the cabinet.

6.3 Printer (Optional for some models)

The printer comes installed at the factory with a roll of paper. If the paper roll runs out after a long period of use, you can purchase the same paper roll (size: thermal paper, paper width: 57.5±0.5mm, reel outer diameter: no more than 40 mm) JP Selecta can supply you with this consumable.

Printer panel description:

Open button. Press to open the lid;

SEL button, indicator light, for factory setting. Don't push it;

Lf button. The green indicator is the power indicator and is normally on when the power is turned on;

Install the paper roll:

Press the open button 1, open the cover, install the paper roll, close the cover and allow the head of the paper roll to extend slightly out of the cover.

6.4 Battery

Battery Location: Top of cabinet and bottom inside electrical box.

Battery maintenance

To prolong battery life and prevent the battery from remaining unused for a long time, the refrigerator should work

for more than 24 hours every month for easy charging;

The equipment needs to be connected to the power supply regularly (usually once a month) and turn on the power switch and battery switch for a period of time to complete charging. Each charging time should not be less than 24 hours;

The power switch and battery switch must be turned off after unplugging the equipment, otherwise it will cause loss of battery capacity, and in severe cases, permanent damage to the battery may occur;

The main power switch must be turned off during prolonged power outage or during transportation; Otherwise, prolonged discharge will cause battery power loss and even permanent damage to the battery.

The battery is consumable and has a lifespan of approximately 2 to 3 years. If the battery is used incorrectly, such as power loss. or reaching the end of battery life, it will cause a low battery alarm (cooling is not affected, but there is an alarm and it may influence the use function of the printer, please contact JP Selecta to request replacement Battery replacement and recovery

Turn off the power switch and unplug the plug from the socket.

Remove the screws from the electrical box with a screwdriver. (Note: There are high voltage electrical components in the electrical box, therefore, before opening, unplug the plug.

The electrical box can only be handled by qualified personnel or maintenance personnel;

Remove the battery connection cable; (Before disconnecting the battery fixing cable, pay attention to the cathode/anode of the battery and the sequence of the connecting cables, to avoid burning the control system due to installing the new battery with reversed cathode/anode Generally, the red wire is connected to the cathode (+) and the black wire is connected to the anode (-);

Use a screwdriver to remove the fixing screws of the battery fixing plate and take out the battery;

To obtain a replaced battery that can be recycled, contact your local battery recycling agency.

6.5 Technical Service

WARNING! Make sure that the appliance is disconnected from the mains before carrying out any maintenance or repair operation. These works must be carried out by qualified personnel.

After carrying out the relevant checks, **DO NOT CARRY OUT ANY REPAIRS YOURSELF**. Contact your Technical Assistance Service, providing the model and serial number of the device (located on the characteristics label),

JP Selecta Technical Assistance Service, SAU


Tel.+34 937700877

Email:serviciotecnico@jpselecta.es


6.6 Spare Parts

Material code	Material description
890392981	Temperature sensor for upper part
890291477	Temperature sensor for lower part
890291542	Temperature sensor for control
890294864	Temperature sensor for condenser
890292582	Temperature sensor for evaporator
890322912	Temperature and humidity sensor
890147520	Condenser
890253363	Fan motor (evaporator)
890196880	Universal wheel
890253565	Control board
890188445	Power supply
890280617	Transformer
890253364	Fan motor (condenser)
890195205	gasket
890326919	Display board
890201348	Door lock
800008811	Power switch
890197125	Compressor
890348235	Filter
890215302	USB temperature recorder

6.7 Refrigerating System Diagram

6.8 Electric Diagram

7 TROUBLESHOOTING

Many operating problems arise from causes that can be easily eliminated without the need to contact the Technical Assistance Service. The following list includes various types of problems and how to solve them.

PROBLEM	POSSIBLE REASON	POSSIBLE SOLUTION
Device not working	<ol style="list-style-type: none"> 1. The plug is not connected. 2. Electrical current is not reaching the plug because the fuse has blown or the automatic power limiter has tripped. 	<ol style="list-style-type: none"> 1. Connect the plug to the outlet current and check that electric current exists. 2. Change the fuse or reconnect the automatic power limiter.
Device is not cooling enough	<ol style="list-style-type: none"> 1. Check cutting temperatura in the controller. 2. Poorly closed door or frequent openings. 3. Obstruction of the ventilation grilles of the appliance. 4. Dirty condenser. 5. The device is exposed to direct sunlight or a heat source. 	<ol style="list-style-type: none"> 1. Reduce the cutting temperature. 2. Make sure the door does not remain open for too long. 3. Keep these areas free as indicated in the "installation" chapter of this manual. 4. Clean with compressed air or a hard hair brush (not steel). 5. Change the location of the refrigerator or protect it from these heat sources.
Noisy operation	<ol style="list-style-type: none"> 1. The appliance has not been leveled correctly. 2. Some of the inner tubes rub. 3. Loose fastening screws on some parts. 4. Fan on condenser or evaporator causing vibrations. 5. Compressor oil charge too low. 6. Loose parts in the condensing unit. 	<ol style="list-style-type: none"> 1. Level as indicated in the chapter "installation" of this manual. 2. Separate the tubes that rub. 3. Tighten the loose screws. 4. Level the equipment and tighten any loose screws. 5. If the equipment has been tipped over at any time, leave it in an upright position without connecting for 24 hours so that the oil returns to the compressor. Check for possible oil leaks.
The device is generating too much ice into the evaporator	<ol style="list-style-type: none"> 1. Poorly closed doors. 2. Excessive opening of doors. 3. Defrosting has not been carried out. 	<ol style="list-style-type: none"> 1,2. Make sure the door is not stay open for a long time. 3. Schedule a defrost.

PROBLEM	POSSIBLE REASON	POSSIBLE SOLUTION
Compressor doesn't start	<ol style="list-style-type: none"> 1. Switch open. 2. Blown fuse. 3. Defective wiring. 4. Clixon open. 5. Controller contacts open. 6. Defective relay. 7. Low gas load in the system. 	<ol style="list-style-type: none"> 1. Close the switch. 2. Replace the blown fuse. 3. Check electrical wiring. 4. Check if there is a voltage abnormally low intake current. 5. Defective controller, or device located in an area that is too cold. 6. Replace the relay. 7. Check for leaks.
Compressor starts, but stops due to overload	<ol style="list-style-type: none"> 1. Low voltage. 2. Defective unit wiring. 3. Defective starting capacitor. 4. Defective compressor. 	<ol style="list-style-type: none"> 1. Check if there is a voltage abnormally low power outlet. 2. Check wiring and electrical installation of the equipment. 3. Replace the starting capacitor. 4. Replace the compressor.
High condensing pressure	<ol style="list-style-type: none"> 1. Unit overloaded with product hot. 2. Air or non-condensable gases in the system. 3. Dirty condenser. 4. Defective condenser fan. 5. Appliance located in an area that is too hot. 6. Obstruction in expansion valve or filter. 7. Discharge valve partially closed. 8. Obstruction in discharge line. 	<ol style="list-style-type: none"> 1. Check temperature of the product, and let it cool outside the equipment if it is very hot. 2. Create a vacuum to extract air or gases. 3. Clean the condenser. 4. Replace the fan. 5. Remove equipment from very hot areas. 6. Adjust expansion valve or filter. 7. Adjust discharge valve. 8. Check and adjust discharge line.
Reduced condensing pressure	<ol style="list-style-type: none"> 1. Insufficient refrigerant charge. 2. Leaks in the system. 3. Appliance located in an area that is too cold. 	<ol style="list-style-type: none"> 1. Check if there are gas leaks refrigerant. 2. Repair system leaks. 3. Remove the equipment from very cold areas.

PROBLEM	POSSIBLE REASON	POSSIBLE SOLUTION
The compressor performs short cycles	<ol style="list-style-type: none"> 1. Differential control adjusted in intervals too small. 2. Low refrigerant charge. 3. Excessive refrigerant charge. 4. Leaking discharge valve. 	<ol style="list-style-type: none"> 1. Adjust differential with controller. 2. Check coolant pressure. 3. Adjust the coolant level. 4. Change the discharge valve. 5. Adjust refrigerant charge to avoid excess overpressures. 6. Clean the condenser.

8. WARRANTY

Dear customer, we inform you of the rules on warranty management that our company grants to its products:

- First of all, we inform you that the products manufactured and sold by JP SELECTA, SAU are capital goods intended for industrial and not domestic use. It is therefore that the guarantee applied is not regulated by the consumer and user law but by the commercial guarantee laws.
- The JP SELECTA, SAU warranty covers any manufacturing defect or any hidden defect in the device for a period of one year.
- The warranty does not cover glass breakage after delivery by JP SELECTA, SAU; nor of parts damaged by misuse or normal wear and tear.
- Any intervention on the device that affects the electrical connections, refrigeration part or electronic microcontroller not authorized by our SAT will result in the loss of the remaining warranty period on the machine.
- Any return authorized through the SAT, either for repair or replacement, is inspected at our facilities. If anomalies other than those claimed are detected that are unrelated to our manufacturing or due to misuse or wear, JP SELECTA, SAU will not be responsible for the costs of repair or replacement, which will be assumed by the customer.
- The conditions of the guarantee will not be modified unless an agreement to modify the conditions of the supply contract is previously established in writing with the client.
- **ATTENTION! NO MACHINE WILL BE ACCEPTED FOR REPAIR THAT IS NOT PROPERLY CLEANED AND DISINFECTED**